

Cellfectin® Reagent

Cat. No. 10362-010

Size: 1 ml

Store at +4°C (do not freeze)

Description

Cellfectin® Reagent is suitable for the transfection of DNA into insect and mammalian cells, and is a 1:1.5 (M/M) liposome formulation of the cationic lipid N, N^I, N^{II}, N^{III}-Tetramethyl- N, N^I, N^{II}, N^{III}-tetrapalmityl-spermine (TM-TPS), and dioleoyl phosphatidylethanolamine (DOPE) in membrane-filtered water.

Cellfectin® has been found to be superior for the transfection of Sf9 cells. Refer to the Cell Lines database at www.invitrogen.com for a list of other cell types successfully transfected.

Important Guidelines for Transfection

1. Form complexes using the amounts of DNA and Cellfectin® recommended (see pages 2-3). Optimize as necessary. **Note:** We recommend diluting DNA and Cellfectin® into Sf-900 II SFM (for insect cells; Cat. No. 10902-096) or Opti-MEM® I Reduced Serum Medium (for mammalian cells; Cat. No. 31985-062) before complexing.
2. Cellfectin® is a lipid suspension that may settle with time. Mix thoroughly by inverting the tube 5-10 times before use.
3. Transfect cells at the confluence or cell density recommended (see pages 2-3). Optimize as necessary. Maintain the same seeding conditions between experiments.
4. **Do not** add antibiotics to media during transfection as this causes cell death.
5. For optimal results, perform transfection in medium without serum. Cells may be transfected in the presence of serum, if desired; however, complexes **must** be formed in serum-free medium.
6. Test serum-free media for compatibility with Cellfectin® since some serum-free formulations may inhibit cationic lipid-mediated transfection.
7. Procedures are provided to transfect cells in a 6-well format (60-mm format for stable mammalian cell transfection). For other formats, vary the amounts of DNA, Cellfectin®, cells, and medium used in proportion to the relative surface area of the tissue culture vessel.

Part No.: 10362.pps

Rev. Date: 16 July 2004

This product is distributed for laboratory research only. CAUTION: Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established.

Transfecting Insect Cells with Baculovirus DNA

Use this procedure to transfect Sf9 insect cells in a **6-well format**. All amounts and volumes are given on a per well basis.

1. Plate 9×10^5 Sf9 cells in 2 ml of Sf-900 II SFM containing antibiotics (*i.e.* penicillin/streptomycin/neomycin; Cat. No. 15640-055) at 0.5X final concentration. Allow cells to attach for at least 1 hour.
2. **For each transfection sample**, prepare complexes as follows:
 - a. Dilute 1-2 μg of baculovirus DNA in 100 μl of Sf-900 II SFM without antibiotics.
 - b. Mix Cellfectin® before use, then dilute 1.5-9 μl in 100 μl of Sf-900 II SFM without antibiotics.
 - c. Combine the diluted DNA with diluted Cellfectin® (total volume = 200 μl). Mix gently and incubate for 15-45 minutes at room temperature (solution may appear cloudy).
3. Remove the growth medium from the cells and wash once with Sf-900 II SFM without antibiotics. Remove the wash medium.
4. Add 0.8 ml of Sf-900 II SFM to the complexes (Step 2c), mix gently and add to the cells. Incubate cells at 27°C for 5 hours.
5. Remove the transfection mixture and replace with 2 ml of Sf-900 II SFM containing antibiotics. Incubate cells at 27°C for 48 hours.
6. Assay for gene activity at 48 hours post-transfection; harvest virus at 72 hours.

Transfecting Adherent Mammalian Cells

Use the following procedure to transiently or stably transfect mammalian cells. All amounts and volumes are given on a per well basis.

1. One day before transfection, plate cells in growth medium without antibiotics such that they will be at the recommended confluence at the time of transfection.

Condition	Cell no.	Growth med. vol.	Format	Confluence at txfn
Transient	$1-2 \times 10^5$	2 ml	6-well	60-80%
Stable	$1-2 \times 10^5$	4 ml	60-mm	30-50%

2. **For each transfection sample**, prepare complexes as follows:
 - a. Dilute 1-2 µg of DNA in 100 µl of Opti-MEM® I Reduced Serum Medium (or other medium) without serum.
 - b. Mix Cellfectin® before use, then dilute 2-15 µl of Cellfectin® in 100 µl of Opti-MEM® I Medium (or other medium) without serum.
 - c. Combine the diluted DNA with diluted Cellfectin® (total volume = 200 µl). Mix gently and incubate for 10-15 minutes (15-45 minutes for stable transfection) at room temperature (solution may appear cloudy).
3. Remove the growth medium from the cells and wash once with 2 ml of growth medium without serum. Remove the wash medium.
4. Add 0.8 ml of Opti-MEM® I Medium without serum (1.8 ml for stable transfection) to the complexes (Step 2c), mix gently and add to cells.
5. Incubate cells at 37°C in a CO₂ incubator for 5-24 hours.
6. Replace the medium with 2 ml of growth medium containing serum (4 ml for stable transfection).
7. **Transient:** Test for transgene expression 24-72 hours post-transfection.
Stable cell lines: Passage cells at a 1:5 (or higher dilution) into selective medium 48-72 hours post-transfection.

Transfecting Suspension Mammalian Cells (6-well format)

1. For each transfection sample, form complexes as in Step 2 of the protocol above, with the following exceptions: (1) Dilute DNA and Cellfectin® into 0.5 ml of Opti-MEM® I Medium each and (2) Incubate complexes for 15-45 minutes at room temperature.
2. For each transfection sample, centrifuge 2×10^6 cells and aspirate medium.
3. Resuspend the cell pellet with the DNA-lipid complex solution (Step 1) and transfer to a well of a 6-well plate. Incubate at 37°C in a CO₂ incubator for 4-5 hours.
4. Add 0.5 ml of growth medium containing 30% fetal bovine serum. **Note:** Add PMA and/or PHA, if desired, to enhance promoter activity and gene expression.
5. The following day, add 1-2 ml of complete growth medium. Assay for transgene expression 24-48 hours post-transfection.

Optimizing Transfections

To obtain the highest transfection efficiency and low non-specific effects, optimize transfection conditions by varying cell density, DNA and Cellfectin® concentrations, and transfection incubation time.

Quality Control

Cellfectin® is tested for the absence of microbial contamination using blood agar plates, Sabaraud dextrose agar plates, and fluid thioglycolate medium, and functionally by transfection of CHO-K1 cells with a reporter plasmid.

Limited Use Label License No. 36: Cellfectin® Reagent

This product is the subject of one or more of U.S. Patent Nos. 5,674,908, 5,834,439, 6,110,916, and foreign equivalents owned by Invitrogen. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Invitrogen Corporation, 1600 Faraday Avenue, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500.