

AcTEV™ Protease


Description:

AcTEV™ Protease is an improved version of Tobacco Etch Virus (TEV) protease that is highly site-specific, highly active, and significantly more stable than native TEV protease, resulting in enhanced long-term activity. AcTEV™ Protease specifically recognizes a seven amino acid sequence (Glu-Asn-Leu-Tyr-Phe-Gln-Gly, cleaving between Gln and Gly), making it useful for removing affinity tags from fusion proteins (1,2). Incubation with AcTEV™ Protease releases the protein of interest from the fusion tag. This is an effective way to remove solubility, secretion, detection, and purification tags from recombinant proteins.

AcTEV™ Protease features:

- Highly specific cleavage activity
- Enhanced enzyme stability for prolonged protease activity (Figure 1)
- Activity over a broad temperature (+ 4°C to 30°C) and pH (6.0 to 8.5) range
- Six-histidine sequence to facilitate its removal from the digested protein sample
- Greater than 85% single-band purity with no non-specific protease contamination

Figure 1 - AcTEV™ and TEV protease stability at + 4°C



AcTEV™ and TEV Protease were incubated at +4°C for 0, 6 h, 24 h, and 168 h (1 week). Samples were analyzed on a 4-20% Novex® Tris-Glycine Gel (Panel A). The amount of degraded product was quantitated and plotted as a function of time (Panel B).

Enzyme specifications:

Purified from *E. coli* expressing the AcTEV™ Protease gene. Unit Definition: One unit of AcTEV™ Protease cleaves 85% of a 3 µg control substrate in 1 h at 30°C. Unit Reaction Conditions: 50 mM Tris-HCl (pH 8.0), 0.5 mM EDTA, 1 mM DTT, 3 µg control substrate, and 1 unit enzyme in 30 µl for 1 h at 30°C. AcTEV™ Protease is functionally tested for the absence of any non-specific protease activity.

Contents and Storage:

AcTEV™ Protease is supplied with a vial of 20X TEV buffer [1 M Tris-HCl (pH 8.0), 10 mM EDTA], and a vial of 100 mM DTT. Store at -20°C. Guaranteed stable for 6 months when properly stored.

Product	Concentration	Quantity	Cat. no.
AcTEV™ Protease	10 units/µl	1,000 units (100 mL)	12575-015
	10 units/µl	10,000 units (1 mL)	12575-023

Reference(s):

1. Polayes, D.A. *et al.* (1994) *Focus*® **16**: 1
2. Parks, T.D. *et al.* (1994) *Anal. Biochem.* **216**: 413.

Recommended vectors containing TEV sites

Highest-level protein production in *E. coli*: Champion™ pET151/D-TOPO (Invitrogen Cat. No. K151-01)
Protein production in baculovirus: pFastBac™ HT (Invitrogen Cat. No. 10584-027)

Insert a TEV cleavage sequence into custom primers

Free, online primer design software with one-step addition of a TEV cleavage sequence:
Visit the OligoPerfect™ Designer under the Custom Primers menu at www.invitrogen.com.



These products may be covered by one or more Limited Use Label Licenses (See the Invitrogen catalog or our web site, www.invitrogen.com). By the use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses.

For research use only. Not intended for any animal or human therapeutic or diagnostic use. Printed in the U.S.A. ©2003 Invitrogen Corporation. All rights reserved. Reproduction forbidden without permission.

Corporate headquarters:

1600 Faraday Avenue • Carlsbad, CA 92008 USA • Tel: 760 603 7200 • Fax: 760 602 6500 • Toll Free Tel: 800 955 6288 • E-mail: tech_service@invitrogen.com • www.invitrogen.com

European headquarters:

Invitrogen Ltd • Inchinnan Business Park • 3 Fountain Drive • Paisley PA4 9RF, UK • Tel: +44 (0) 141 814 6100 • Fax: +44 (0) 141 814 6260 • E-mail: eurotech@invitrogen.com