

## Operating Instructions

### Affinity Chromatography Media

## Cellufine<sup>®</sup> ET clean for endotoxin removal

### Description

The Cellufine ETclean is poly( $\epsilon$ -lysine) immobilized Cellufine<sup>®</sup> (cellulose spherical beads). The beads bind and remove endotoxin from your sample solution. The poly( $\epsilon$ -lysine) is a microbial poly(amino acid) that consist of 30-35 lysine residues produced by *Streptomyces albulus*. The poly( $\epsilon$ -lysine) as ligand and the cellulose beads act as matrix and are products of Chisso Corporation.

### Physical-Chemical Characteristics

Product Name	Supplied	Particle size	Pore size*
Cellufine ET clean S	a slurry in 20 % ethanol	ca. 40 - 130 $\mu\text{m}$	$M_{\text{lim}}$ 2000
Cellufine ET clean L			$>M_{\text{lim}}$ $2 \times 10^6$

\*The pore size (molecular weight exclusion;  $M_{\text{lim}}$ ) of the beads was estimated from calibration curves obtained by size exclusion chromatography. Pullulan and maltose were used for the  $M_{\text{lim}}$  determination.

### Column Packing

1. Calculate volume required for the desired bed dimension.
2. Prepare a 40 – 60 % (v/v) slurry with the appropriate elution buffer (high salt). Allow the gel to equilibrate at ambient temperature for one hour.
3. Gently stir or place under vacuum to degas.
4. With column outlet closed, carefully pour the slurry into column. Depending on the volume, a filler tube may be necessary.
5. With the inlet open to release air, insert and affix the top adjuster assembly at the slurry interface.
6. Open the column outlet and begin pumping elution buffer at rate 10 % – 20 % greater than the operational flow rate.
7. After the bed stabilizes, close the column outlet. Then with the inlet open, reposition the end cell on top of the bed. Equilibrate with 10 column volumes of adsorption buffer before sample loading.

## Operating Guidelines

### General Operation

1. CIP for endotoxin free. Wash the Cellufine ET clean column with 5 column volumes of alkali solution. It leaves as it is required time until endotoxin becomes free. Followed by wash with endotoxin-free water.

➤ Note

Alkali solution and require time for endotoxin free.

Alkali solution	Require time for endotoxin free
0.2mol/l NaOH	16hr or over night
0.2mol/l NaOH in 20% EtOH	3 to 5 hr
0.2mol/l NaOH in 95%EtOH	1hr

2. Equilibrate the column with 5 column volumes of a suitable endotoxin-free buffer.
3. Apply sample through the column at a flow rate of 10 to 50cm/h at 4-25 °C.
4. Collect the effluent and determine the endotoxin content of the effluent, as a sample solution after endotoxin-removing treatment.
5. The column can be reused after washing with a cleanup method of (1) to (2). We had already checked that the beads can be regenerated 5 times.

### Recommended Buffers

Generally, 10mM to 50mM sodium phosphate buffer or Tris-HCl buffer, neutral pH can use it well. As long as a sample is stable, simply endotoxin-water is sufficient to use.

Since acid protein may be adsorbed on ET-clean, salt concentration is raised when protein adsorbs. By buffer pH lower than protein  $pI$ , the protein is hard to adsorb on ET-clean.

### Sample Preparation and Load

Prepare samples at a concentration of 1 – 20 mg/ml, in adsorption buffer. Remove insoluble material by centrifugation or microfiltration. If necessary, exchange sample buffer using dialysis, diafiltration or desalting chromatography.

The protein may not be recovered if protein is added to superfluous ET-clean by very low

concentration and very small amount.

### **Selection guide**

The endotoxin-adsorbing capacity of the Cellufine ET clean beads was strongly depend on “exclusion limit” (abbreviates to  $M_{lim}$ .) the adsorbing capacity increased from 500 to 1000  $\times 10^{-6}$ g (LPS from *E. coli* O111:B4) per ml of wet beads, while the  $M_{lim}$  increased from  $2.0 \times 10^3$  to  $>2 \times 10^6$  at pH 7.0 and NaCl concentration of 0.17M. Although Cellufine ET clean-L, having the large  $M_{lim}$  of  $>2 \times 10^6$ , show the greatest endotoxin-removing activity, ionic binding of components other than endotoxin may occur by entry of the components into the pore of the beads.

The beads must be selected as follows:

- (1) To reduce endotoxin from a sample solution containing acidic protein with  $pI$  4.0-6.5, you can use Cellufine ET clean-S beads with a small pore size at pH 5-7 and NaCl concentration of 0.1-0.4 M.
- (2) To reduce endotoxin from a sample solution containing neutral or basic protein with  $pI$  7.0-10.5, you can use Cellufine ET clean-L beads with a large pore size at pH 7-9 and NaCl concentration of 0.1-0.4 M

### **Storage**

Short term (2 weeks or less), bulk and column can be stored in 1 M NaCl in neutral buffer at 2 – 4°C. Longer term storage can be conducted under identical conditions; however, a preservative (e.g. 0.1 % formalin, 0.05 % chloreton or 0.02 % sodium azide) should be added to the buffer. Store at 4 – 8 °C. Do not freeze.

**Shelf Lifetime:** 5 years

### **Storage**

For storage of opened containers, it is recommended that they be kept in a cold room (2 - 8 °C). Do not freeze.

### **Batchwise Method**

The Cellufine ET clean beads must be endotoxin-free. See “General Operation”.


(Example)

- (1) Place 5 ml of Cellufine ET clean beads into a glass-buchner funnels with fritted disc (pore size: 30 micron-m). Add 25 ml of alkali solution in it, and suspend the mixture with a spatula. Stand the

suspension for a required time until endotoxin becomes free and then remove the solution by vacuum.

Alkali solution	Require time for endotoxin free
0.2mol/l NaOH	16hr or over night
0.2mol/l NaOH in 20% EtOH	3 to 5 hr
0.2mol/l NaOH in 95%EtOH	1hr

- (2) Wash the beads, then, with other cleanup solutions (2 M NaCl e.q., endotoxin-free water, and then equilibrate buffer, respectively) by a similar method of (1).
- (3) Suspend 0.2- to 0.4-g portion of wet adsorbent (after removing equilibrate buffer by vacuum) into a flask with 2 ml of sample solution. Shake the suspension for 2 h at 4-25 °C and filter it through a membrane filter (0.8 micron-m) to remove the beads.
- (4) determine the endotoxin content of the filtrate obtained, as a sample solution after endotoxin-removing treatment.
- (5) The beads can be regenerated before each use, with the washing method of (1) to (2).



supplement

**ET clean L can easy to remove LPS from buffers.**

ET clean L easily reduces LPS to low concentrations compared with polymyxin-immobilized agarose gel or DEAE Cellufine A-500(anion-exchanger).

**Materials & Method**

Column : 9mm I.D. x 100mm

Pump : Peristaltic with silicon tube

Buffer : 1M sodium phosphate, pH 7.0 (spiked 2.6 EU/ml LPS)

LPS Removal (pre-washing)

Column & media: wash with 5CV 0.2 M NaOH; let stand for 16 hours, then wash with endotoxin-free water.

Silicon tubing: wash with 0.5 M NaOH, let stand for 16 hours, then wash with endotoxin-free water.

※0.2 M NaOH-20% EtOH is more effective for endotoxin-free.

Chromatography Flow rate : 30ml/h [Residence time 12.8min; linear velocity 47cm/h]

Fraction 6.36ml (total 128 tubes)

Assay

LAL (rate assay, EndospecyES-50M Set; SEIKAGAKU CORPORATION)




Fig. Comparison of the LPS removal capability from 1M phosphate buffer of ET-clean L, and polymyxin-immobilized agarose gel and DEAE Cellufine A-500.

**Product Ordering Information**

Cellufine ET clean L		Cellufine ET clean S	
Pack Size	Catalogue No.	Pack Size	Catalogue No.
Mini-column 1ml x 5	20051	Mini-column 1ml x 5	20151
10ml	681 984 324	10ml	682 985 324
50ml	681 984 326	100ml	682 985 326
500 ml	681 984 328	500 ml	682 985 328
5 Liters	681 984 330	5 Liters	682 985 330
10 Liters	681 984 335	10 Liters	682 985 335

Cellufine ET clean was developed by the Joint Project of Kumamoto University & Chisso Corporation.

## CHISSO CORPORATION

Fine Chemicals Sales Department

2-1, Otemachi 2-Chome, Chiyoda-ku, Tokyo 100-8105, Japan

Phone +81-3-3243-6150, Fax+81-3-3243-6219

E-mail: [cellufine@chisso.co.jp](mailto:cellufine@chisso.co.jp)

<http://www.chisso.co.jp/fine/en/cellufine>