

# Pierce<sup>®</sup> LAL Chromogenic Endotoxin Quantitation Kit

88282

2445.1

Number	Description
88282	<p><b>Pierce LAL Chromogenic Endotoxin Quantitation Kit</b>, sufficient reagents to perform 50 microplate assays</p> <p><b>Kit Contents:</b></p> <p><b><i>Escherichia coli</i> (<i>E. coli</i>) Endotoxin Standard (011:B4)</b>, lyophilized, 1 vial, 15-40 endotoxin units (EU)/mL upon reconstitution</p> <p><b>Limulus Amebocyte Lysate (LAL)</b>, lyophilized, 2 vials, 1.4mL/vial upon reconstitution</p> <p><b>Chromogenic Substrate</b>, lyophilized, 1 vial, 6.5mL/vial upon reconstitution</p> <p><b>Endotoxin-free Water</b>, 1 vial, 30mL</p>

**Storage:** Upon receipt store at 4°C. Product shipped with an ice pack.

## Introduction

The Thermo Scientific Pierce LAL Chromogenic Endotoxin Quantitation Kit is an efficient, quantitative endpoint assay for the detection of gram-negative bacterial endotoxins. Bacterial endotoxin catalyzes the activation of a proenzyme in the modified Limulus Amebocyte Lysate (LAL). The activated proenzyme then catalyzes the splitting of *p*-Nitroaniline (pNA) from the colorless substrate, Ac-Ile-Glu-Ala-Arg-pNA; the activation rate is proportional to the sample endotoxin concentration. After stopping the reaction, the released pNA is photometrically measured at 405-410nm. The correlation between absorbance and endotoxin concentration is linear in the 0.1-1.0EU/mL range. The developed color intensity is proportional to the amount of endotoxin present in the sample and can be calculated using a standard curve.

## Important Product Information

- Accurate pipetting is critical for maintaining consistent results. A repetitive pipettor can aid in normalizing volumes between samples. Ensure pipetting order and rate of reagent addition remain consistent from well to well and row to row.
- All materials (e.g., pipette tips, glass tubes, microcentrifuge tubes and disposable 96-well microplates) must be endotoxin-free.
- Adjust the sample pH to 6-8 using endotoxin-free 0.1M NaOH or 0.1M HCl. Avoid pH-electrode contamination of the sample by testing the pH of a small sample taken from the bulk sample.
- Store samples to be tested to stop all bacteriological activity. Store at 2-8°C for <24 hours and -20°C for >24 hours.
- Maintaining the correct temperature is critical for reproducibility. Use a proper heating block at 37°C±1°C. Do not use a cabinet-style incubator to perform the assay.
- Endotoxin adheres to glass and plastic surfaces; before pipetting, vortex solutions to ensure the correct endotoxin concentrations are measured.
- Glass tubes are preferred for making standard stock solutions; however, polystyrene or polypropylene microcentrifuge tubes (1.5 mL) may also be used. When using microcentrifuge tubes, dedicate the bag of tubes for the assay and follow aseptic techniques.
- If the test sample endotoxin concentration is >1.0EU/mL, dilute the sample five-fold in endotoxin-free water. Re-test.

- Assay inhibition occurs when substances in the test sample interfere with the LAL reaction. In the chromogenic assay, this inhibition results in a lower final absorbance, indicating lower levels of endotoxin than what may be present in the test sample. Determine the lack of product inhibition for each sample undiluted or at an appropriate dilution. See the Additional Information Section for more details.
- Samples turning yellow after addition of the Stop Reagent (25% acetic acid) or possessing significant initial color may require special attention. See the Additional Information Section for more details.
- The kit reagents are “matched” to comply with Food and Drug Administration (FDA) requirements for endotoxin testing. Each LAL lot is tested for functionality using the United States Reference Standard EC-6. The LAL lot is then matched to a lot of our Control Standard Endotoxin (CSE) by testing in parallel with the Reference Standard Endotoxin (RSE). The RSE/CSE correlation assay determines the potency of each CSE lot when used with each matching LAL lot.

## Additional Materials Required

- Disposable endotoxin-free glass tubes or 1.5mL microcentrifuge tubes
- Disposable endotoxin-free pipette tips
- Disposable endotoxin-free 96-well microplates
- Heating block at 37°C±1°C
- Pipettor
- Repetitive pipettor (optional)
- Microplate reader
- 25% acetic acid (Stop Reagent)

## Material Preparation

**Note:** Equilibrate all reagents to room temperature before use.

### A. Endotoxin Standard Stock Solutions

1. Each *E.coli* Endotoxin Standard vial contains ~15-40EU of lyophilized endotoxin; the actual vial concentration is printed on the label. Reconstitute by adding 1mL of room temperature endotoxin-free water to make Endotoxin Standard Stock. For example, a vial with a concentration of 26EU, when reconstituted with 1.0mL of endotoxin-free water, will yield a concentration of 26EU/mL. Shake the solution vigorously for at least 15 minutes on a vortex mixer before use.

**Note:** Store lyophilized *E.coli* Endotoxin Standard at 2-8°C. Reconstituted stock solution is stable for 4 weeks at 2-8°C. Before use, warm the solution to room temperature and vigorously mix for 15 minutes.

2. Prepare Standard Stock Solutions from the Endotoxin Standard Stock using the dilutions and procedures in Table 1.

**Table 1. Dilutions and procedures for preparing Standard Stock Solutions.**

Vial	Volume of Endotoxin Standard Stock (mL)	Volume of vial A (mL)	Endotoxin-free Water (mL)	Final endotoxin concentration (EU/mL)
A	0.05	–	(X-1)/20*	1.0
B	–	0.25	0.25	0.50
C	–	0.25	0.75	0.25
D	–	0.1	0.9	0.1

\*X = endotoxin concentration of the *E. coli* Endotoxin Standard supplied with the kit; refer to the Certificate of Analysis to get the lot-specific concentration.

- A. Prepare a solution containing 1.0EU/mL of endotoxin standard by diluting 0.05mL of the Endotoxin Standard Stock with [(X-1)/20]mL of endotoxin-free water, where X equals the endotoxin concentration of the vial (e.g., if X= 26EU/mL, then dilute 0.05mL of this stock with (26-1)/20, or 1.25mL, of endotoxin-free water). Vigorously vortex the solution for  $\geq 1$  minute before proceeding.
- B. Transfer 0.25mL of the 1.0EU/mL Standard Vial A into 0.25mL of endotoxin-free water in a suitable container to prepare 0.5EU/mL Standard Stock Solution and **vortex the solution vigorously for 1 minute before use.**
- C. Transfer 0.25mL of 1.0EU/mL Standard Vial A into 0.75mL of endotoxin-free water in a suitable container to prepare 0.25EU/mL Standard Stock Solution and **vortex the solution vigorously for 1 minute before use.**
- D. Transfer 0.1mL of 1.0EU/mL Standard Vial A into 0.90mL of endotoxin-free water in a suitable container to prepare 0.1EU/mL Standard Stock Solution and **vortex the solution vigorously for 1 minute before use.**

### B. Limulus Amebocyte Lysate (LAL)

1. The LAL reagent contains lyophilized lysate prepared from the circulating amebocytes of the horseshoe crab *Limulus polyphemus*. Reconstitute immediately before use with 1.4mL of endotoxin-free water and swirl gently to dissolve the powder. Avoid foaming; do not vortex the solution.
2. Store lyophilized LAL protected from light at 2-8°C. Store reconstituted LAL reagent at -20°C or colder for 1 week if frozen immediately after reconstitution. Thaw and use the reconstituted LAL reagent one time; once thawed, gently swirl the reagent to mix before adding to samples.

### C. Chromogenic Substrate

1. Each vial contains ~7mg of lyophilized substrate. Reconstitute the Chromogenic Substrate by adding 6.5mL of endotoxin-free water to yield a final concentration of ~2mM.
2. Store lyophilized Chromogenic Substrate protected from light at 2-8°C. Reconstituted substrate is stable for 4 weeks when stored at 2-8°C.


## Microplate Assay Procedure

**Note:** Equilibrate all reagents to room temperature before use.

1. Equilibrate the microplate in a heating block for 10 minutes at 37°C.
2. With the microplate maintained at 37°C, carefully dispense 50µL of each standard or unknown sample replicate into the appropriate microplate well, cover the plate with the lid and incubate for 5 minutes at 37°C.

**Note:** Each series of determinations must include duplicate runs of a blank and the four endotoxin standards; the blank contains 50µL of endotoxin-free water. If reaction inhibition is possible, see the Additional Information Section.

3. Add 50µL of LAL to each well using a pipettor, cover the plate with the lid and gently shake the plate on a plate shaker for 10 seconds. Incubate the plate at 37°C for 10 minutes.  
**Note:** Ensure pipetting order and rate of reagent addition remain consistent from well-to-well and row-to-row.
4. After exactly 10 minutes, add 100µL of substrate solution to each well. Pipette the substrate solution in the same manner as in step 3. Maintain a consistent pipetting speed.
5. Cover the plate with lid and gently shake on a plate mixer for 10 seconds. Incubate the plate at 37°C for 6 minutes.
6. Add 50µL of Stop Reagent (25% acetic acid). Maintain the same pipetting order as in steps 3 and 4. Gently shake the plate on a plate mixer for 10 seconds.
7. Measure the absorbance at 405-410nm on a plate reader.
8. Subtract the average absorbance of the blank replicates from the average absorbance of all individual standard and unknown sample replicates.
9. Prepare a standard curve by plotting the average blank-corrected absorbance for each standard versus its concentration in EU/mL. The coefficient of determination,  $r^2$ , must be  $\geq 0.98$ .
10. Use the formulated standard curve to determine the endotoxin concentration of each unknown sample (Figure 1).



**Figure 1. Example standard curve for the quantitation of endotoxin in a chromogenic assay.**

## Troubleshooting

Problem	Possible Cause	Solution
Non-linear standard curve	Standard Stock Solutions were not mixed well	Vortex the Endotoxin Standard Stock for 15 minutes before each use
		Vortex all Standard Stock Solutions for 1 minute before each use
Higher absorbance in blank than standards	Pipetting order and rate of reagent addition was irregular	Ensure pipetting order and rate of reagent addition remain consistent from well to well and row to row
		Use a repetitive pipettor
Higher absorbance in blank than standards	Materials (e.g., tips, vials, microplates) were contaminated	Use endotoxin-free materials

## Additional Information

### A. Product Inhibition/Enhancement

Product inhibition occurs when substances in the test sample interfere with the LAL reaction. In the chromogenic assay, this inhibition results in a lower, final absorbance, indicating lower levels of endotoxin than what may be present in the test sample. Determine the lack of product inhibition for each sample undiluted or at an appropriate dilution.

To verify the lack of product inhibition, spike an aliquot or dilution of a test sample with a known amount of endotoxin (e.g., 0.5EU/mL). Assay the spiked sample and the unspiked samples to determine the respective endotoxin concentrations. The difference between the two calculated endotoxin values should equal the known concentration of the spike  $\pm 25\%$ .

### B. Example for Determination of a Non-inhibitory Dilution

**Table 2. Example study for determining inhibition of a sample.**

Sample dilution	Spiked <sup>§</sup> concentration (EU/mL)	Unspiked concentration (EU/mL)	Difference
1:10	0.28	0.18	0.10 = Inhibitory
1:20	0.36	0.11	0.25 = Inhibitory
1:40	0.50	< 0.1	0.50 = Non-inhibitory

<sup>§</sup> Spiked concentrations should all show a value of 0.50EU/mL. The values of 0.28 and 0.36 are indicative of inhibition at the respective dilutions.

Samples showing inhibition on the LAL reaction may require further dilution to overcome the inhibitory effects. Once the non-inhibitory dilution is determined, the exact dilution can be found by testing two-fold dilutions near that dilution. The degree of inhibition or enhancement is dependent on the product concentration.

Beta glucans are polymers of D-glucose found in fungi and plant cell walls with >1000-fold less LAL reactivity than lipopolysaccharides. Inhibition and false-positive colors can occur in samples contaminated with beta glucans. Use appropriate beta glucan blockers if any contamination is possible.

### C. Colored Samples

In the chromogenic assay, samples turning yellow after addition of the Stop Reagent or possessing significant initial color (e.g., tissue culture media) may require special attention. To determine if a sample's intrinsic color will alter the absorbance readings, construct a mock reaction tube by adding 50 $\mu$ L of sample, 150 $\mu$ L of endotoxin-free water and 50 $\mu$ L of Stop Reagent with no incubation. Read the absorbance at 405-410nm. If the absorbance is significantly greater than the absorbance of endotoxin-free water, then the intrinsic color will alter the correct sample absorbance readings. In such cases, include appropriate controls in the assay.

## Related Thermo Scientific Products

88270	Pierce High Capacity Endotoxin Removal Resin, 10mL
88271	Pierce High Capacity Endotoxin Removal Resin, 100mL
88272	Pierce High Capacity Endotoxin Removal Resin, 250mL
88273	Pierce High Capacity Endotoxin Removal Spin Column, 0.25mL, 5 columns
88274	Pierce High Capacity Endotoxin Removal Spin Column, 0.50mL, 5 columns
88276	Pierce High Capacity Endotoxin Removal Spin Column, 1mL, 5 columns
20339	Detoxi-Gel™ Endotoxin Removing Gel
20344	Detoxi-Gel Endotoxin Removing Columns
89896	Pierce Centrifuge Columns, 2mL, 25/pkg

---

<b>89897</b>	<b>Pierce Centrifuge Columns, 5mL, 25/pkg</b>
<b>89898</b>	<b>Pierce Centrifuge Columns, 10mL, 25/pkg</b>
<b>23225</b>	<b>Pierce BCA Protein Assay Kit</b>
<b>22660</b>	<b>Pierce 660nm Protein Assay Kit</b>

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

**No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).**

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at [www.thermoscientific.com/pierce](http://www.thermoscientific.com/pierce). For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2011 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.