



Novagen®

Merck Biosciences

Calbiochem | Novabiochem | Novagen

# InsectDirect™

Protein Expression & Purification System

## When should I use an insect cell expression system?

Researchers dealing with high-throughput (HT) expression screening face the fundamental challenge of trying to rapidly identify the appropriate expression system for many targets in parallel. The focus has traditionally not involved insect expression systems due to their complexity and time-consuming protocols. However, insect cell systems provide an efficient means for expressing proteins that require specific post-translational modifications to retain solubility and activity.

## What does the InsectDirect System offer?

The InsectDirect System provides a rapid and powerful, virus-free method for heterologous protein expression in insect cells. The InsectDirect System consists of the following:

- The pIEx™ series of plasmid expression vectors for rapid, high-level protein expression in insect cells, which eliminates the need to create recombinant baculovirus for protein expression
- The pBiEx™ series of multisystem expression vectors designed to allow rapid characterization/expression of target genes in both *E. coli* and insect cells
- Insect GeneJuice® Transfection Reagent for high-efficiency transfection of insect cells with minimal toxicity
- Insect RoboPop™ Ni-NTA His•Bind® Purification Kit for automation-compatible His•Tag® fusion protein purification



## How does use of the InsectDirect™ System benefit me?

**48 hours**  
VERSUS  
**21 days**

- Allows protein expression in 48 hours, not 2 to 3 weeks
- Provides an option for HT screening of multiple targets
- Produces up to 80 µg target protein per ml culture
- Includes a protocol that can be scaled for larger culture volumes and higher protein yields
- Offers vector choices for secretion signals, N-terminal and C-terminal fusion tags, and “tagless” options

## Insect Cell Expression Vectors

The pIEx™ and pBiEx™ vectors allow rapid, high-level protein expression in insect cells, omitting the time-consuming production of recombinant baculovirus. Both the pIEx and pBiEx vectors feature an hr5 enhancer and ie1 (immediate early) promoter combination derived from AcNPV baculovirus. The hr5/ie1 combination utilizes endogenous insect cell transcription machinery for direct expression in insect cells, thereby avoiding baculovirus infection and its associated cytopathic effects. The pBiEx vectors also allow for expression in *E. coli* by the tightly controlled T7lac promoter.

Vectors	Promoter(s)	Signal sequence	Fusion Tags		Protease cleavage sites <sup>†</sup>	Ek/LIC vector option	Size	Cat. No.
			N-terminal	C-terminal				
pIEx-1	hr5/ie1	No	His•Tag®/S•Tag™	HSV•Tag®	Tb/Ek	Yes <sup>††</sup>	20 µg	71241-3
pIEx-2	hr5/ie1	No	GST•Tag™/His•Tag/S•Tag	HSV•Tag	Tb/Ek	Yes <sup>††</sup>	20 µg	71238-3
pIEx-3	hr5/ie1	Yes	GST•Tag/His•Tag/S•Tag	HSV•Tag	Tb/Ek	Yes <sup>††</sup>	20 µg	71243-3
pIEx-4	hr5/ie1	No	None	S•Tag/His•Tag	None	No	20 µg	71235-3
pIEx-5	hr5/ie1	Yes	None	S•Tag/His•Tag	None	No	20 µg	71242-3
pIEx-6	hr5/ie1	No	His•Tag	S•Tag	Ek	No	20 µg	71333-3
pIEx-7 Ek/LIC*	hr5/ie1	No	His•Tag	S•Tag	Ek	Yes	20 rxn	71339-3
pBiEx-1	hr5/ie1, T7lac	No	His•Tag/S•Tag	HSV•Tag	Tb/Ek	No	20 µg	71234-3
pBiEx-2	hr5/ie1, T7lac	No	GST•Tag/His•Tag/S•Tag	HSV•Tag	Tb/Ek	No	20 µg	71233-3
pBiEx-3	hr5/ie1, T7lac	No	None	S•Tag/His•Tag	None	No	20 µg	71232-3

\* pIEx-7 is available only as an Ek/LIC vector. The pIEx-7 Ek/LIC Vector Kit includes Linearized Ek/LIC Vector; T4 DNA Polymerase, LIC-qualified; polymerase buffer; DTT, EDTA, and dATP solutions; Nuclease-free Water; competent cells; SOC Medium; and controls.


<sup>†</sup> Tb: thrombin; Ek: enterokinase

<sup>††</sup> pIEx-1, pIEx-2, and pIEx-3 vectors are also available in an Ek/LIC Vector Kit. For more information visit our website address listed below.

## Insect GeneJuice® Transfection Reagent

Insect GeneJuice Transfection Reagent is a liposomal formulation optimized for maximal transfection efficiency that offers extremely low toxicity to insect cells. It can be used for both transient and stable transfections in serum-containing or serum-free media. Insect GeneJuice is ideal for HT or large-scale protein expression when using the pIEx or pBiEx vectors for suspension culture transfection of Sf9 and other insect cells.

Product	Size	Cat. No.
Insect GeneJuice®	0.3 ml	71259-3
Transfection	1 ml	71259-4
Reagent	10 × 1 ml	71259-5



Sf9 cells transfected with pIEx-1/β-gal using Insect GeneJuice Transfection Reagent. Transfected cells were stained for β-galactosidase activity using the X-Gal Solution-based BetaBlue™ Staining Kit.


The Insect RoboPop Ni-NTA His•Bind Purification Kit is designed for filtration-based 96-well purification of His•Tag® fusion proteins directly from transfected insect cell cultures. The kit configuration allows for robotic processing of transfected 10-ml suspension cultures and purification of up to 400 µg His•Tag fusion protein per culture based on binding capacity of the resin.

The kit includes:

- Insect PopCulture® Reagent, for protein extraction directly from total cultures
- Benzonase® Nuclease, for viscosity reduction
- Ni-NTA His•Bind Resin and buffers
- 2-ml 96-well Filter Plate and Collection Plate with Sealer

## Insect RoboPop™ Ni-NTA His•Bind® Purification Kit

Product	Cat. No.
Insect RoboPop™ Ni-NTA His•Bind® Purification Kit	71257-3
Components:	
• 50 ml	Insect PopCulture® Reagent
• 10 KU	Benzonase® Nuclease, Purity > 90%
• 10 ml	Ni-NTA His•Bind Resin
• 125 ml	4X Ni-NTA Bind Buffer
• 2 × 125 ml	4X Ni-NTA Wash Buffer
• 50 ml	4X Ni-NTA Elute Buffer
• 1	2-ml 96-well Filter Plate
• 1	Collection Plate with Aluminum Plate Sealer
Note: 1 KU = 1000 units	




### Overview of reagent requirements for various culture sizes

Culture size	Sf9 cells	Plasmid DNA dilution in BacVector® Medium	Insect GeneJuice Transfection Reagent dilution in BacVector Medium	Insect PopCulture Reagent	Benzonase Nuclease	His•Bind Ni-NTA Resin	Resin pool applied to column	His•Bind Ni-NTA 1X Wash Buffer	His•Bind Ni-NTA 1X Elute Buffer
10 ml	9 ml at $1.0 \times 10^6$ /ml	20 µg + 0.5 ml	60 µl + 0.4 ml	0.5 ml	5 µl	0.25 ml	96-well filter plate	2 × 1.5 ml	0.5 ml
100 ml	90 ml at $1.1 \times 10^6$ /ml	200 µg + 5 ml	1 ml + 4 ml	5 ml	50 µl	2.5 ml	30 mm × 20 cm	2 × 15 ml	7.5 ml
200 ml	180 ml at $1.1 \times 10^6$ /ml	400 µg + 10 ml	2 ml + 8 ml	10 ml	100 µl	10 ml	30 mm × 20 cm	2 × 30 ml	15 ml
1 liter (4 × 250-ml cultures)	225 ml at $2.5 \times 10^6$ /ml per flask	500 µg + 12 ml per flask	2.5 ml + 10 ml per flask	12.5 ml per flask	125 µl per flask	25 ml to pooled lysate	30 mm × 23 cm	2 × 100 ml	75 ml

## Small-scale (10 ml) and large-scale (1 liter) Sf9 suspension culture transfections

### A. Heat shock proteins and phospholipase, 10-ml scale



### B. Protein kinase, 1-liter scale



**Panel A.** Sf9 cells in 10-ml suspension cultures ( $1 \times 10^6$  cells/ml) were transfected with 20 µg pEx™ recombinant plasmid using Insect GeneJuice® Transfection Reagent. Total culture extracts were prepared 48 h later by adding Insect PopCulture® Reagent, followed by the addition of Benzonase® Nuclease. Samples were removed at that point to assess total cell protein. Ni-NTA His•Bind® Resin was then added directly to the extracts. The samples were processed robotically using a MultiPROBE® II HT EX Liquid Handling Station (PerkinElmer). Samples of the crude and purified fractions from each transfection were analyzed by SDS-page (10–20% gradient). Purified protein yields were determined by BCA assay.

**Panel B.** Four 250-ml suspension cultures were seeded with Sf9 cells ( $1 \times 10^6$  cells/ml) in serum-free BacVector® Insect Cell Medium. A pEx/protein kinase isoform construct was added to the diluted Insect GeneJuice Transfection Reagent, incubated for 15 min at room temperature, and then added to the flasks. After 48 h incubation at 28°C, Insect PopCulture Reagent and Benzonase Nuclease were added directly to each flask and incubated for 15 min at room temperature to lyse the cells. The His•Tag® fusion protein was manually purified on Ni-NTA His•Bind Resin according to the standard protocol. Protein yield was quantified by a modified Bradford assay.

Prices and availability are subject to change. Copyright© 2004 Merck Biosciences. An Affiliate of Merck KGaA, Darmstadt, Germany. All Rights Reserved. Each product is sold with a limited warranty which is provided with each purchase. Each product is intended to be used for research purposes only. It is not to be used for drug or diagnostic purposes nor is it intended for human use. Merck Biosciences products may not be resold, modified for resale, or used to manufacture commercial products without written approval of Merck Biosciences.

BacVector®, Calbiochem®, GeneJuice®, His•Bind®, His•Tag®, HSV•Tag®, Novabiochem®, Novagen®, and PopCulture® are registered trademarks of EMD Biosciences, Inc. BetaBlue™, GST•Tag™, InsectDirect™, pBiEx™, Perfect Protein™, pEx™, RoboPop™, and S•Tag™ are trademarks of EMD Biosciences, Inc. Benzonase® is a registered trademark of Merck KGaA, Darmstadt, Germany and MultiPROBE® is a registered trademark of PerkinElmer, Inc.

For more information or to order Novagen products, contact us directly or your local sales office.

#### Orders:

Phone +1 608 238 6110

Fax +1 800 776 0999

#### Technical Support:

Phone +1 608 238 6110

E-mail [novagen@merckbiosciences.com](mailto:novagen@merckbiosciences.com)

#### Visit our website:

[www.novagen.com](http://www.novagen.com)

Calbiochem, Novabiochem, and Novagen are brands of Merck Biosciences, an affiliate of Merck KGaA, Darmstadt, Germany.

Printed in the USA.