

## Table of Contents

I. Introduction.....	3
II. Components .....	4
III. Storage .....	4
IV. Overview of <i>Brevibacillus</i> expression system.....	4
IV-1. Selection of expression vectors (1) pNCMO2 (2) pNY326	
IV-2. Cloning into expression vectors	
IV-3. Transformation of <i>Brevibacillus</i>	
IV-4. Detection of products and scaling-up	
V. Procedure.....	10
V-1. <i>Brevibacillus</i> strains	
V-1-1. Genotypes	
V-1-2. Control	
V-1-3. Storage	
V-2. <i>E. coli</i> host strains	
V-3. Construction of expression plasmid using pNCMO2	
V-3-1. Cloning into pNCMO2	
V-3-2. Analysis of pNCMO2 recombinants	
V-3-3. Sequencing	
V-4. Construction of expression plasmid using pNY326	
V-4-1. Cloning into pNY326	
V-4-2. Transformation of <i>Brevibacillus</i>	
V-4-3. Analysis of pNY326 recombinants	
V-5. Transformation of <i>Brevibacillus</i>	
V-5-1. Equipments and reagents required	
V-5-2. Transformation by electroporation	
V-6. Expression of proteins by <i>Brevibacillus</i> recombinants	
V-6-1. Experimental outline	
V-6-2. Culture media	
V-6-3. Culture for secretory production	
V-7. SDS-PAGE analysis	
V-7-1. Sample preparation	
V-7-2. Control	
V-7-3. Analysis of expressed proteins	

V-8. Optimization of protein expression	
V-8-1. In the case of low expression	
V-8-2. In the case of no expression	
V-9. Purification of expressed proteins	
V-10. Medium components	
VI. Appendix.....	17
VI-1. Transformation of <i>Brevibacillus</i> by Tris-PEG method	
VI-1-1. Preparation of solutions	
VI-1-2. Transformation	
VII. References.....	18
VIII. Related products.....	18
IX. Notice.....	18

### **Purchaser's Agreement**

1. *B. choshinensis* SP3, pNY326, pNCMO2 and pNCMO2-BLA were developed and manufactured by Higeta Shoyu Co., Ltd. and sold by TAKARA BIO INC.
2. The use of this product is allowed for academic research only. You have to conclude a license agreement if you use this product or the products produced by using this product for commercial purpose. Contact us if you would like to ask any question about license etc.  
If any questions are received, either Higeta Shoyu Co., Ltd. or an agent for it will answer.  
TaKaRa Bio Inc.  
Tel +81-77-543-7247 Fax +81-77-543-9254
3. All rights of intellectual property related to this product are owned by Higeta Shoyu Co., Ltd. This product is covered by the following Japanese patents and corresponding international patent :  
Japanese patents, 2103774, 3696322, 3734593, 3753945, PCT/JP2004/016912 (international application based on PCT treaty), 3433807.
4. This product, the components or the derivatives of this product, or the products produced by using these items can not be transferred (distributed freely or sold) to third party. (When you purchase this product, you submit the license agreement to us.)

**I. Introduction**

*Brevibacillus choshinensis*, a gram positive bacterium, has excellent ability to produce many kinds of proteins extracellularly. And the *Brevibacillus* Expression System based on this bacterium is well-suited for secretory production of heterologous proteins with high efficiency. Making use of this characteristic of the host bacterium, many successful records have been accumulated so far in this regard. This system has following merits :

- The host bacterium secretes proteins very efficiently.
- The host bacterium produces negligible amount of extracellular protease so that products remain unscathed in the culture medium.
- Proteins are produced as active forms.
- Easy to culture and sterilize the host bacterium.
- The host bacterium is amenable to genetic engineering.
- The host bacterium is guaranteed to be safe organism.

A part of our successful results are listed in Table 1. They cover enzymes, antigens, cytokines. All of them were produced at very high level and confirmed to have native biological activity. In addition to them, proteins from eubacteria, archaeobacteria, eukaryotes, viruses, were successfully produced. Special merit of this system is that it exhibits very high efficiency in the production of eukaryotic secretory proteins. Usually they have some cysteine residues intramolecularly and disulfide bonds must be formed at exact positions precisely, which makes it difficult to achieve efficient production using other bacterial expression systems based on intracellular expression. Meanwhile, this system is very good at secretory production of proteins with native structure irrespective of whether they have disulfide bonds or not.

Table 1. Some examples of successful production of heterologous proteins using *B. choshinensis* host-vector system.

Proteins	Origins	Production (g/L)	References
<b>Enzymes</b>			
$\alpha$ -amylase	<i>B. licheniformis</i>	3.7	
Sphingomyelinase	<i>B. cereus</i>	3.0	
Xylanase	<i>B. halodurans</i>	0.2	
CGTase	<i>B. macerans</i>	1.5	2)
Chitosanase	<i>B. circulans</i>	1.4	
Hyper thermo-stable protease	<i>A. pernix</i>	0.1	
Hyper thermo-stable nuclease	<i>P. horikoshii</i>	0.7	
PDI	human	1.0	3)
<b>Antigens</b>			
Surface antigen	<i>E. rhusiopathiae</i>	0.9	
Surface antigen	<i>T. pallidum</i>	0.8	
<b>Cytokines</b>			
EGF	human	1.5	4)
IL-2	human	0.6	5)
NGF	mouse	0.2	
IFN- $\gamma$	chicken	0.5	6)
TNF- $\alpha$	bovine	0.4	
GM-CSF	bovine	0.2	
GH	flounder	0.2	

The host bacterium exhibits high transformation efficiency by electroporation so that it is easy to construct expression clones. There is a shuttle vector between *B. choshinensis* and *E. coli* as well as the one maintained only in *B. choshinensis*. Thus, either the 2 step construction method using the shuttle vector or 1 step one using the vector for only *B. choshinensis* can be selected.

Two kinds of culture media are recommended for the culture of transformants (components are shown in V-10.). Cultures are performed using test tubes or flasks fixed on an appropriate shaker. After culture, supernatant is separated from cells by centrifugation from which target protein can be purified. The unnecessary to conduct cell disruption makes the purification procedure simpler and easier.

## II. Components

*Brevibacillus* Expression System (Cat.#HB100)

- Expression vector
  - pNY326 DNA 10  $\mu$ g (Cat.#HB111)
  - pNCMO2 DNA 10  $\mu$ g (Cat.#HB112)
- Control vector
  - pNCMO2BLA DNA 1  $\mu$ g (Cat.#HB113)
- Competent cells (for electroporation)
  - Brevibacillus choshinensis* Electro-Cells 100  $\mu$ l  $\times$  10 tubes (Cat.#HB115)

## III. Storage

Expression vectors pNY326, pNCMO2 DNA:  $-20^{\circ}\text{C}$

Control plasmid pNCMO2BLA DNA:  $-20^{\circ}\text{C}$

Competent cells *Brevibacillus choshinensis* Electro-Cells:  $-80^{\circ}\text{C}$

## IV. Overview of *Brevibacillus* expression system

The flow of experiments to produce the target protein using this kit is described below.

### IV-1. Selection of the expression vectors

(1) pNCMO2

pNCMO2 is a shuttle vector between *B. choshinensis* and *E. coli*. At first, expression plasmids can be constructed by employing *E. coli* as the host, then transferred into *B. choshinensis* cells.

The pNCMO2 contains P2 promoter, which is one of 5 promoters driving transcription of cell wall protein (HWP) gene and does not work in *E. coli*. Therefore, it is suitable for cloning genes of which expressions are stressful for *E. coli*. In *B. choshinensis* cells, P2 promoter works as very strong promoter, enabling efficient protein production.




Fig 1. Map of pNCMO2

Because of the strong promoter activity of P2 promoter, it shows detrimental effect on the growth of the transformants. In that case, the use of pNY326 is recommended.

<Features of pNCMO2>

P2 promoter	A part of 5' upstream region of cell wall protein (HWP) gene, which does not work in <i>E. coli</i> , and works as a strong promoter in <i>Brevibacillus</i> .
Sec signal peptide	Secretion signal sequence of HWP is modified to increase the secretion efficiency.
Multiple cloning site	11 restriction sites allowing insertion of target gene into the vector
Terminator	46 bp nucleotides sequence functioning as a terminator is introduced downstream of multiple cloning sites
rep	Gene involving in plasmid replication (pUB110 derived)
Ori	Replication origin allowing replication and maintenance of the plasmid in <i>Brevibacillus</i> (pUB110 derived)
Nm <sup>r</sup>	Neomycin resistance gene, working as selection marker in <i>Brevibacillus</i> .
ColE1 Ori	Replication origin allowing replication and maintenance of the plasmid in <i>E. coli</i> (pUC derived)
Amp <sup>r</sup>	Ampicillin resistance gene, working as selection marker in <i>E. coli</i>



(2) pNY326

pNY326 can replicate in *B. choshinensis* cells and its size (3.4 kb) is much smaller than that of pNCMO2 (5.2 kb). In addition, the promoter activity of the vector is substantially weaker than that of pNCMO2. Taken together, it can be maintained stably even if it expresses rather toxic proteins to the host cells. In some case, higher efficient production of protein can be achieved using this vector because of the more stable growth of the transformants. Even after repeated inoculation to scale up the culture, it still exhibits stable productivity. Thus, this vector is especially suited to culture on large scale. The vector can be maintained only in *B. choshinensis* so that only *B. choshinensis* can be used as the host. The construction of expression plasmid must be constructed in one step using *B. choshinensis* as the host. Because the transformation efficiency must be sufficiently high, the use of *Brevibacillus choshinensis* Electro-Cells ( $10^5$  transformants/ $\mu$ g DNA) is recommended.




Fig 2. Map of pNY326

<Features of pNY326>

P5 promoter	A part of 5' upstream region of cell wall protein (HWP) gene, which works as a relatively weak promoter in <i>Brevibacillus</i> .
Sec signal peptide	Secretion signal sequence of HWP is modified to increase the secretion efficiency.
Multiple cloning site	11 restriction sites allowing insertion of target gene into the vector
Terminator	26 bp nucleotide sequence functioning as a terminator is introduced downstream of multiple cloning sites
rep	Gene involving in plasmid replication (pUB110 derived)
Ori	Replication origin allowing replication and maintenance of the plasmid in <i>Brevibacillus</i> (pUB110 derived)
NmI	Neomycin resistance gene, working as selection marker in <i>Brevibacillus</i>





**IV-2. Cloning into expression vectors**

The expression vectors have secretion signal sequence which derives from *B. choshinensis* cell wall protein gene. DNA fragment encoding the target protein should be cloned into the MCS which is located at the 3' end of the signal sequence. Using two restriction sites, one is on the signal sequence and the other on MCS, the gene of interest can be cloned into the vector with correct orientation.




Fig. 3. Construction of *B. choshinensis* expression plasmid

**IV-3. Transformation of *Brevibacillus***

The transformation of *B. choshinensis* cells can be done by electroporation. Selection marker is neomycin resistance. When the shuttle vector, pNCMO2, is used in combination with *E. coli* host, ampicillin resistance can be used as the selection marker.

If the electroporation device is not available, another transformation method "Tris-PEG method" is applicable. However, the transformation efficiency of "Tris-PEG method" is very low compared to electroporation so that in that case it is recommended to use the completed plasmid which is constructed using *E. coli* as a host.

**IV-4. Detection of the product and scaling-up**

Pick the transformant harboring the expression plasmid and culture it in the liquid medium. After 48-64 hrs' culture, prepare culture supernatant by centrifugation of the culture broth. An aliquot of the supernatant should be subjected to SDS-PAGE analysis to confirm the production of the target protein. If the production is not sufficient, please try scaling-up. Large scale culture of *B. choshinensis* is not difficult.

## V. Procedure

### V-1. *Brevibacillus* strains

Standard genetic engineering technique is applicable.

#### V-1-1. Genotypes

An essential gene for spore formation cascade is disrupted in *B. choshinensis* SP3 so that sterilization of the transformants is easily done. Furthermore, the trace activity due to intracellular protease gene (*imp*) and extracellular protease gene (*emp*) is nullified by gene disruption to secure the intactness of protein products.

#### V-1-2. Control DNA

pNCMO2-BLA : Positive control plasmid for secretory production which have the gene encoding *Bacillus licheniformis*  $\alpha$ -amylase. The size of the protein is about 55 kDa and secreted to the amount of 0.1 g/l or more.

#### V-1-3. Storage

Short term storage (about 1 week)

Pick a single colony and spread it on MTNm plate (see V-10. Medium components). Place it overnight in the incubator kept at 30°C. Seal the plate and store it at room temperature (about 20°C).

Note) Do not store plates in a refrigerator !

Long term storage (longer than 1 month)

Pick a single colony and inoculate it into 2SYNm medium (see V-10. Medium components). Culture it overnight on a shaker. Transfer the broth to a freezing vial and add equal amount of LB medium containing 40 % glycerol. Store it at - 80°C.

### V-2. *E. coli* host strains

The nucleotide sequence encoding lac operator is inserted in pNCMO2 to weaken the promoter activity in *E. coli*. For that reason, F factor integrated strains like JM109 have to be used as host strains. The genotype of JM109 is shown below for reference.

JM109 : *recA1, endA1, hsdR17 thi, hsdR17* (rk<sup>-</sup> mk<sup>+</sup>), *e14<sup>-</sup>* (*mcrA<sup>-</sup>*), *supE44, relA1, Δ (lac-proAB) /F' [ traD36, proAB<sup>+</sup>, lacIq, lacZΔM15 ]*

### V-3. Construction of expression plasmid using pNCMO2

Procedure and precaution for the construction of expression plasmids using pNCMO2 as the vector.

*E. coli* strains having *lacIq* and *recA<sup>-</sup>* like JM109 are recommended as the host.

Insert DNA should be cloned in the vector in-frame at the downstream of secretion signal. Stop codon must be introduced at the 3' end of the gene.

If target protein is prokaryotic secretory protein, better result may be obtained using the protein's own secretory signal. In that case, pNY326 is selected.

#### V-3-1. Cloning into pNCMO2

<Amplification of the gene by PCR>

PCR primers should be designed to clone target gene into the site downstream of the secretion signal. Select appropriate 2 restriction sites available on the expression vector and introduce them into the primers' sequences. Be sure the orientation is correct. Then perform the amplification of the gene by PCR. Please set up optimal condition for PCR according to the characteristics of the gene and PCR enzymes. It is recommended to use high-fidelity PCR enzymes (PrimeSTAR® HS DNA Polymerase (Cat.#R010A) etc.).

## &lt;Construction of expression plasmids via ligation&gt;

Both the PCR product and the vector (0.5-1.0  $\mu$ g) should be digested with 2 species of restriction enzymes. The DNAs are subjected to electrophoresis on agarose gel. Extract and purify DNAs from the pieces of gel.

For ligation reaction, use 100ng of each DNA. Ligation can be done using, for instance, DNA ligation Kit <MightyMix> (Cat.#6023). Use 1/5 of the ligation mixture for the transformation of *E. coli* and the rest for the transformation of *B. choshinensis* after purifying it by ethanol precipitation. As for the procedure thereafter, refer to the corresponding part of V-4. pNY326.

<Transformation of *E. coli*>

Use *E. coli* host of which transformation efficiency is high. It is recommended to use *E. coli* JM109 Competent Cells (Cat.#9052) or *E. coli* JM109 Electro-Cells (Cat.#9022).

**V-3-2. Analysis of pNCMO2 recombinants**

After spreading 100-200  $\mu$ l of the transformed mixture onto LB plates containing 50-100  $\mu$ g/ml ampicillin, incubate the plates at 37°C for 15-18 hr. Select 10-20 ampicillin-resistant colonies and inoculate them into 2 ml of LB liquid medium containing 50-100  $\mu$ g/ml ampicillin.

Grow for 15-18 hr at 37°C with shaking and isolate plasmid DNA.

Perform restriction enzyme analysis using appropriate amount of DNA. Usually the enzymes which were used for the construction of expression plasmid are employed. Confirm the existence of the insert DNA with agarose electrophoresis. If the existence of the insert could be detected, perform sequencing to confirm that the gene has no error introduced by PCR and is in frame with signal sequence.

**V-3-3. Sequencing**

For sequencing, following sequences are used as the forward and the reverse primers (the primer sequences are common in both pNCMO2 and pNY326).

Forward Sequencing Primer : 5'-CGCTTGCAGGATTCGG-3'

Reverse Sequencing Primer : 5-CAATGTAATTGTTCCCTACCTGC-3'

**V-4. Construction of expression plasmids using pNY326**

Procedure and precaution for the construction of expression plasmids using pNY326 as the vector.

*B. choshinensis* Electro-Cells are recommended to use for the transformation when pNY326 is used as the cloning vector, because competent cells with high transformation efficiency are needed.

Insert DNA should be cloned in frame with the signal sequence on the vector. Stop codon should be introduced at the 3' end of the gene. If target gene is a bacterial gene encoding secretory protein, it is recommended to try the original secretion signal as well as *B. choshinensis* cell wall protein signal. In that case, the restriction site of 5' end for the insert DNA is *Bsp*H I site which covers the initiation codon of cell wall gene. If transformants can not be obtained using pNCMO2 or higher stability is required, pNY326 should be selected.

**V-4-1. Cloning into pNY326**

<Amplification of the gene by PCR>

PCR primers should be designed to clone target gene into the site downstream of the secretion signal. Select appropriate 2 restriction sites available on the vector and introduce them into the 5' ends of the primers. Be sure the orientation is correct. Then perform the amplification of the gene by PCR. Please set up optimal condition for PCR according to the characteristics of the gene and PCR enzymes. It is recommended to use high-fidelity PCR enzymes (PrimeSTAR® HS DNA Polymerase (Cat.#R010A) etc.).

<Construction of expression plasmid via ligation>

Both the PCR product and the vector (0.5-1.0  $\mu$ g) should be digested with 2 species of restriction enzymes. The DNAs are subjected to electrophoresis on agarose gel. Extract and purify DNAs from the pieces of gel. For ligation reaction, use 100 ng of each DNA. Ligation can be done using eg. DNA ligation Kit <MightyMix> (Cat.#6023). Use the entire ligation mixture for the transformation of *B. choshinensis* after ethanol precipitation.

**V-4-2. Transformation of *Brevibacillus***

Competent cell of high transformation efficiency should be used as the host (It is recommended to use *Brevibacillus choshinensis* Electro-Cells).

As for the method of transformation, refer to the section V-5 Transformation of *Brevibacillus*.

**V-4-3. Analysis of pNY326 recombinants**

- (1) Pick 10-20 colonies on the MTNm plates at random and inoculate them into 2 ml TMN liquid medium.
- (2) Culture for 15-18 hr at 37°C with shaking. Isolate plasmid DNA. Usually 1-2  $\mu$ g of DNA can be recovered.
- (3) Carry out restriction enzyme analysis using appropriate amount of DNA. Usually the enzymes which were used in the preparation of the vector and the insert DNA for the construction of expression plasmid are used. Confirm the existence of the insert by electrophoresis on agarose.
- (4) If the existence of the insert was detected, perform sequencing to confirm the gene has no error introduced by PCR and is in frame with signal sequence (refer to the section V-3-3 for the primers of sequence).

**V-5. Transformation of *Brevibacillus*****V-5-1. Equipments and reagents required**

Competent cells for electroporation (*Brevibacillus choshinensis* Electro-Cells) (on ice)  
Expression plasmid to express your protein  
Positive control (pNCMO2-BLA)  
Negative control (pNCMO2 or pNY326)  
MTNm plates  
MT liquid medium  
Cuvette for electroporation  
Culture tube  
Sterilized microtube

**V-5-2. Transformation by electroporation**

- (1) Place the tube containing competent cells on ice
- (2) Add 5  $\mu$ l DNA solution (200-500 ng)\* to the competent cells and mix well by flipping the tube (Never vortex!).
- (3) Transfer the mixture to the electroporation cuvette. Store the cuvette on ice for 10 min.
- (4) Pulse the cells under the following condition :  
Electroporator : Gene Pulser II (Bio-Rad)  
Charging voltage : 7.5 kV/cm  
Capacitance : 25  $\mu$ F  
Resistance : 1000  $\Omega$   
Gap of the cuvette : 2 mm  
Electroporator : Electroporator 2510 (Eppendorf)  
Charging voltage : 14 kV/cm  
Gap of the cuvette : 1 mm
- (5) Immediately add 1 ml MT medium and then transfer the mixture into a culture tube
- (6) Incubate the tube at 30°C with shaking at 120 rpm.
- (7) Spread 100  $\mu$ l of the mixture onto a MTNm plate. Centrifuge the remaining mixture (5000 rpm, 5 min, room temp.) to concentrate the cells. Spread the whole amount of cells onto another MTNm plate.
- (8) Incubate plates for 15-18 hrs at 37°C. If the sizes of the colonies are small, extend the incubation further.
- (9) Obtained colonies can be subjected to plasmid check or protein expression analysis.

\*When constructed plasmids as the control plasmids are used, the amount of DNA should be reduced to 10-100 ng.

**V-6. Expression of proteins by *Brevibacillus* recombinants**

After the construction of expression clone is completed, small-scale expression check can be carried out. The general protocol for confirmation of expressed protein is shown below.

**V-6-1. Experimental outline**

Once a positive clone has been identified, expression check can be carried out using positive and negative controls. Following recombinants are used as the controls.

*B. choshinensis* SP3/pNCMO2-BLA (Positive control for secretory expression)

*B. choshinensis* SP3/pNCMO2 or pNY326 (negative control)

The production of the protein is variable among colonies depending on the nature of proteins. Thus, 6-10 colonies should be picked at random (including both large and small colonies) to carry out expression check. If the plates are left over for several days, the potential of protein production may decrease. In that case, do transformation once more to obtain fresh transformants.

**V-6-2. Culture media**

TM and 2SY are recommended as the basic media for expression check. The protein productivity may vary depending on the medium. Because the balance between growth and protein productivity is an important factor for final production rate, it is recommended to try both the media.

**V-6-3. Culture for secretory production**

Check the expression condition by conducting the expression trial using positive control, *B. choshinensis* SP3/pNCMO2-BLA, which expresses *B. licheniformis*  $\alpha$ -amylase (about 55 kDa). By culturing the transformant and the negative control simultaneously and performing the same processing, the expression of target protein can be confirmed. Expression protocol is as follows.

- (1) Pick single colonies and inoculate them into 3 ml of 2SYNm and TMNm in culture tubes (16 mm in diameter). Grow them for 48-64 hrs at 30-33°C with shaking at 120 rpm. During the culture, check the protein production every 24 hrs.
- (2) After the completion of the culture, centrifuge the broth (5000×g, 5 min) to separate supernatant from cells. Cells should be resuspended in the same volume of PBS.
- (3) Supernatant and cellular fraction should be subjected to SDS-PAGE analysis (CBB staining or western blotting) or activity measurement.

**V-7. SDS-PAGE analysis**

Protocol for SDS-PAGE analysis of the protein is as follows.

**V-7-1. Sample preparation**

Add 10  $\mu$ l of 5XSDS-PAGE loading buffer to 40  $\mu$ l of culture supernatant and the cell suspension.  
Mix and boil for 10 min to make samples for electrophoresis.

**V-7-2. Control**

Use following samples as the control.

- a. molecular size marker
- b. standard of your protein
- c. sample from background control (prepared from the culture of the cells harboring only vector)
- d. sample from positive control (prepared from the culture of *B. choshinensis* SP3/pNCMO2-BLA which secretes *B. licheniformis*  $\alpha$ -amylase)

**V-7-3. Analysis of expressed proteins**

It can be confirmed to produce target protein by comparing the position and the density of the specific band with those of the standard protein, which was electrophoresed in parallel with the protein by SDS-PAGE. If the protein production is low, the product is insoluble, or the product is masked by other protein, the protein detection by SDS-PAGE may be difficult. In those cases, western blot analysis using specific antibody, functional evaluation (e.g. specific activity), or protein purification by special technique may be useful for detection of protein production.

**V-8. Optimization of protein expression**

In this system, protein production level is 100 mg/l or more, though the production exceeded 1 g/l in some of the cases. If the productivity of target protein is substantially low or no product can be detected, refer to the guideline described below.

**V-8-1. In the case of low expression**

- a) Try using both the vectors, pNCMO2 and pNY326. The former has stronger promoter activity but does not always exhibit higher productivity. Sometimes pNY326 is the better choice because of the low promoter activity which reduces the growth stress caused by the extremely efficient production of proteins.
- b) Select another type of culture medium. The productivity changes substantially depending on the composition of culture media.
- c) Compare the plasmid copy number with the control. If remarkable reduction of the copy number is detected, use pNY326 as a vector, change the culture medium, or increase the concentration of the antibiotic (400 mg/l).
- d) There is the case that the protein is not suitable for secretory production. In that case, try intracellular production.  
If the production is little, if any, try to produce the protein intracellularly.

**V-8-2. In the case of no expression**

It is recommended to do the same trial as that described in the V-8-1, "In the case of low expression". If no improvement is obtained, following possibilities are conceivable.

- a) Check the secondary structure of mRNA. A strong palindrome structure in the sequence may occur bad effect on the efficiency of translation. In that case, change the sequence to remove the stacking region.
- b) If the amino acid sequence around the signal cleavage site is improper, it may exert some influence on the secretion of the protein. In that case, try to introduce some tag sequences for purification or detection to the N-terminal of the protein, if the addition of some extra amino acids produces no substantial effect on its activity.

**V-9. Purification of expressed proteins**

Purification methods are variable according to the target proteins. If the protein is secreted, the proteins can be easily obtained from supernatant by removing the bacteria, allowing easier purification. After this step, continue purification with ordinary methods (e.g. ion exchange, hydrophobic, affinity chromatography etc.).

**V-10. Medium components**

## 2SY liquid medium

Components :

Glucose	20.0 g/L
Bacto Soytone (Becton Dickinson)	40.0 g/L
Bacto Yeast Extract (Becton Dickinson)	5.0 g/L
CaCl <sub>2</sub> · 2H <sub>2</sub> O	0.15 g/L

\* Autoclave the mixture of glucose and CaCl<sub>2</sub> separately from the rest of the components. Mix them after autoclaving.

## 2SY Nm liquid medium

Add neomycin solution (stock sol. 50 mg/ml) to 2SY liquid medium to the concentration of 50 μg/ml.

## TM liquid medium

Components :

Glucose	10.0 g/L
Polypeptone	10.0 g/L
Meat extract	5.0 g/L
Yeast extract	2.0 g/L
FeSO <sub>4</sub> · 7H <sub>2</sub> O	10 mg/L
MnSO <sub>4</sub> · 4H <sub>2</sub> O	10 mg/L
ZnSO <sub>4</sub> · 7H <sub>2</sub> O	1 mg/L

Adjust the pH to 7.0

\* Autoclave glucose separately from the rest of the components. Mix them after autoclaving.

## TM Nm liquid medium

Add neomycin solution (stock sol. 50 mg/ml) to TM liquid medium to the concentration of 10 μg/ml.

## MT liquid medium

Components :

Glucose	10.0 g/L
Polypeptone	10.0 g/L
Meat extract	5.0 g/L
Yeast extract	2.0 g/L
FeSO <sub>4</sub> · 7H <sub>2</sub> O	10 mg/L
MnSO <sub>4</sub> · 4H <sub>2</sub> O	10 mg/L
ZnSO <sub>4</sub> · 7H <sub>2</sub> O	1 mg/L
MgCl <sub>2</sub>	4.1 g/L

Adjust the pH to 7.0

\* Autoclave glucose separately from the rest of the ingredients. Mix them after autoclaving.

## MTNm plates

Add 7.5 g of agar to 500 ml MT liquid medium and autoclave. Allow the medium to cool to about 50°C before adding neomycin solution (stock sol. 50 mg/ml) to the concentration of 10 μg/ml. Mix the medium by swirling and dispense it to plates.



**VI. Appendix****VI-1. Transformation of *Brevibacillus* by Tris-PEG method**

Highest efficiency of transformation can be obtained by electroporation but if the necessary device is not available, transformation can be done by Tris-PEG method. However, the transformation efficiency is as low as  $10^1$ - $10^2$  transformants/ $\mu$ g plasmid DNA so that the use of the method should be limited to the case that premade pNCMO2-based expression plasmid is used.

**VI-1-1. Preparation of solutions**

Phosphate buffer : dissolve 1.905 g of  $\text{KH}_2\text{PO}_4$  and 0.852 g of  $\text{Na}_2\text{HPO}_4$  into distilled water to make 100 ml solution.

PEG solution : dissolve 40 g of PEG6000 into 50 ml of phosphate buffer. Add  $\text{H}_2\text{O}$  up to 100 ml.

TP medium : Mix phosphate buffer with the same volume of concentrated ( $\times 2$ ) TM medium.

MT medium : TM medium containing 20 mM  $\text{MgCl}_2$ .

50 mM Tris-HCl buffer (pH7.5)

50 mM Tris-HCl buffer (pH8.5)

TE buffer : 10 mM Tris-HCl (pH7.5), 1 mM EDTA

**VI-1-2. Transformation**

- (1) Inoculate *Brevibacillus* cultured overnight cultured into 5 ml TM medium and continue the culture on a shaker.
- (2) Collect the cells in log-phase ( $\text{OD}_{660} = 2.0$ - $3.0$ ) by centrifugation ( $5,000 \times g$ , 5 min., at room temp.) using 50 ml tube.
- (3) Rinse the cells with 5 ml, 50 mM Tris-HCl buffer (pH7.5) and collect them by centrifugation.
- (4) Suspend the cells in 5 ml, 50 mM Tris-HCl buffer (pH8.5) and shake them on a shaker for 60 min. at  $37^\circ\text{C}$ .
- (5) Collect the cells by centrifugation ( $5,000 \times g$ , 5 min, at room temp.) and suspend them in 0.5 ml TP medium.
- (6) Prepare plasmid solution containing about 1  $\mu$ g DNA in 50  $\mu$ l TE buffer plus 50  $\mu$ l TP medium. Add the solution to the cell suspension and mix them. Add 1.5 ml PEG solution and mix them quickly.
- (7) Shake it for 10 min at  $37^\circ\text{C}$  and then collect the cells by centrifugation ( $3000 \times g$ , 10 min, at room temp).
- (8) Add 1 ml MT medium to suspend the cells and shake it for 1 hr at  $37^\circ\text{C}$ .
- (9) Spread 0.1 ml/plate of the suspension onto MTNm plate and incubate them overnight at  $37^\circ\text{C}$ . Pick the colonies resistant to Nm and analyze them afterward. (as for the procedure after this step, refer to "V-4-3. Analysis of pNY326 recombinants").

**VII. References**

- 1) H. Takagi, K. Kadowaki and S. Udaka., Screening and Characterization of Protein-Hyperproducing Bacteria without Detectable Exoprotease Activity. *Agric. Biol. Chem.*, **53** (3), 691-699, 1989
- 2) T. Takano, A. Miyauchi, H. Takagi, K. Kadowaki, K. Yamane, S. Kobayashi., Expression of the Cyclodextrin Glucanotransferase Gene of *Bacillus macerans* in *Bacillus brevis*. *Biosci. Biotech. Biochem.*, **56** (5), 808-809, 1992
- 3) H. Tojo, T. Asano, K. Kato, S. Udaka, R. Horinouchi, A. Kakinuma., Production of Human Protein Disulfide Isomerase by *Bacillus brevis*. *J Biotechnol.* **33** (1), 55-62, 1994
- 4) H. Yamagata, K. Nakahama, Y. Suzuki, A. Kakinuma, N. Tsukakoshi, S. Udaka., Use of *Bacillus brevis* for efficient synthesis and secretion of human epidermal growth factor. *Proc. Natl. Acad. Sci. USA.*, **86**, 3589-3593, 1989
- 5) Y. Takimura, M. Kato, T. Ohta, H. Yamagata, S. Udaka., Secretion of Human Interleukin-2 in Biologically Active Form by *Bacillus brevis* Directly into Culture Medium. *Biosci. Biotechnol. Biochem.*, **61** (11), 1858-1861, 1997
- 6) K. Yashiro, J. W. Lowenthal, T. E. O'Neil, S. Ebisu, H. Takagi, High-Level Production of Recombinant Chicken Interferon-g by *Brevibacillus choshinensis*. Expression and Purification, **23**, 113-120, 2001

**VIII. Related products**

*E. coli* JM109 Competent Cells (Cat.#9052)  
*E. coli* JM109 Electro-Cells (Cat.#9022)  
DNA Ligation Kit <Mighty Mix> (Cat.#6023)  
PrimeSTAR® HS DNA Polymerase (Cat.#R010A)  
pNC-HisT DNA (Cat.#HB121)  
pNC-HisF DNA (Cat.#HB122)  
pNC-HisE DNA (Cat.#HB123)  
pNI DNA (Cat.#HB131)  
pNI-His DNA (Cat.#HB132)

**IX. Notice : Living Modified Organism**

These products (Cat.#HB100 and Cat.#HB115) include a genetically "Living Modified Organism (LMO)" defined in "The Cartagena Protocol on Biosafety". The supplied *Brevibacillus choshinensis* Electro Cells in these kits contain a part of 2  $\mu$ m plasmid DNA derived from *Saccharomyces cerevisiae*.

Please confirm the guidelines or the laws and regulations that you should obey in your country and pay attention for safe handling, storage, transport and disposal.

---

NOTICE TO PURCHASER: LIMITED LICENSE

**[L39] Brevibacillus Expression System**

1. *B. choshinensis* SP3, pNY326, pNCMO2, pNCMO2-BLA, pNC-HisT, pNC-HisF, pNC-HisE, pNI and pNI-His were developed and manufactured by Higeta Shoyu Co., Ltd. and sold by TAKARA BIO INC.
2. The use of this product is allowed for academic research only. You have to conclude a license agreement if you use this product or the products produced by using this product for commercial purpose.

Contact us if you would like to ask any question about license etc. If any questions are received, either Higeta Shoyu Co., Ltd. or an agent for it will answer.

TaKaRa Bio Inc.

Tel +81-77-543-7247 Fax +81-77-543-9254

2. All rights of intellectual property related to this product are owned by Higeta Shoyu Co., Ltd. This product is covered by the following Japanese patents and corresponding international patent: Japanese patents, 2103774, 3696322, 3734593, 3753945, PCT/JP2004/016912 (international application based on PCT treaty), 3433807.
3. This product, the components or the derivatives of this product, or the products produced by using these items cannot be transferred (distributed freely or sold) to a third party. (When you purchase this product, you submit the license agreement to us.)

---

**NOTE :** This product is intended to be used for research purpose only. They are not to be used for drug or diagnostic purposes, nor are they intended for human use. They shall not be used in products as food, cosmetics, or utensils, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.

If you require licenses for other use, please call at +81 77 543 7247 or contact us via our website at [www.takara-bio.com](http://www.takara-bio.com).