

# Size Exclusion Chromatography

TSKgel SW-type

SW

SW<sub>XL</sub>

SuperSW

TSKgel PW-type

PW

PW<sub>XL</sub>

TSKgel Alpha-type

Alpha

SuperAW


TSKgel H-type

H<sub>XL</sub>

H<sub>HR</sub>

SuperH

SuperHZ



# Introduction to TSK-GEL Size Exclusion Chromatography Columns



Tosoh Bioscience provides TSK-GEL columns for both modes of Size Exclusion Chromatography (SEC), which are Gel Filtration Chromatography (GFC) and Gel Permeation Chromatography (GPC). GFC refers to the SEC separation of water soluble polymers in aqueous mobile phases, while GPC refers to the SEC separation of organic soluble polymers using an organic solvent as mobile phase. From a sample perspective, the TSK-GEL SW-type column line, which is based on spherical porous silica particles, is most suitable for analyzing proteins and peptides by gel filtration. Polymer-based TSK-GEL PW columns are best used for GFC analysis of other water soluble polymers, such as oligosaccharides, acrylic acid, etc. Polymer-based TSK-GEL Alpha and SuperAW columns can be used in aqueous solvents, and also with polar organic solvents, thus bridging the gap between GFC and GPC. Organic soluble polymers are best separated by GPC using polystyrene/divinylbenzene-based TSK-GEL H<sub>XL</sub>, H<sub>HR</sub>, SuperH or SuperHZ columns. *Table I* compares the characteristics of the various TSK-GEL column lines for SEC.

Tosoh Corporation has a proud history of innovation in size exclusion chromatography. The porous silica-based SW columns were originally introduced in 1978 as TSKgel G2000SW, G3000SW and G4000SW, using 10 and 13 μm spherical particles. The second generation of SEC columns for biopolymer analysis was introduced in 1987. By reducing the particle size from 10 to 5 μm for G2000SW<sub>XL</sub> and G3000SW<sub>XL</sub>, and from 13 to 8 μm for G4000SW<sub>XL</sub>, analysis times were cut in half without sacrificing resolution. The third generation of SW-type columns was introduced in 1998 under the name SuperSW, a column line that features 4 μm spherical particles packed in 4.6 mm ID columns rather than the wider (7.8 mm ID) SW<sub>XL</sub> columns. TSK-GEL SuperSW columns are available in two pore sizes that mimic the pore structure of the G3000SW<sub>XL</sub> and G2000SW<sub>XL</sub> columns. Due to their higher column efficiency and smaller column volume, these columns provide better resolution and higher sensitivity in sample limited cases. This trend towards smaller, more efficient particles packed into narrower bore columns has been expanded to include columns for the analysis of non-biological aqueous and non aqueous polymers. SuperAW, SuperH and SuperHZ columns represent the state of the art column technology for high throughput SEC analysis.

Bulk polymeric Toyopearl GFC resins for process scale separations are available in convenient LABPAK samplers (up to 150 mL) and in bulk quantities of up to 500 mL (see the Bulk

Resins section of this catalog). For larger volumes of Toyopearl GFC media, please request a copy of the process media catalog.

## Column Selection

The complete TSK-GEL SW, PW, Alpha and SuperAW column lines for GFC consist of many packings, available in various pore and particle sizes. The main criterion in choosing between the TSK-GEL SW, PW, Alpha and SuperAW SEC columns is the molecular weight of the sample and its solubility. The fact that the TSK-GEL SW columns are based on silica and the TSK-GEL PW, Alpha and SuperAW columns are derived from a hydrophilic polymer network has less impact on the separation than the particle and pore size differences. See the Molecular Weight Range Tables in each column type section to determine the best column choice based on pore size.

Due to their higher resolving power, the TSK-GEL SW columns are suitable for the separation of monodisperse biopolymers such as proteins and nucleic acids. TSK-GEL PW columns are commonly used for the separation of synthetic water-soluble polymers because they exhibit a much larger separation range, better linearity of calibration curves, and less adsorption than the TSK-GEL SW columns. While a TSK-GEL SW column is typically the first column to try for biopolymers, TSK-GEL PW columns have demonstrated good results for smaller peptides (<1,000 Da), protein aggregates, DNA fragments, and viruses.

The TSK-GEL Alpha Series columns offer a new alternative for performing SEC. Their compatibility with a wide range of solvents makes them useful for both gel filtration chromatography (GFC) and gel permeation chromatography (GPC). TSK-GEL SuperAW columns are based upon the same chemistry as Alpha columns but have smaller particle sizes and shorter, narrower column dimensions for high-throughput applications. Unlike the Alpha columns, mixed bed formats are included in the SuperAW product offerings for samples with wide ranges of molecular weights/hydrodynamic radius. See *Figure 1* for more information on compatibility of all TSK-GEL SEC columns and solvents.

Non-aqueous GPC separations are performed with TSK-GEL H-type pre-packed columns. Our TSKgel MultiporeH<sub>XL</sub>-M column contains porous polystyrene divinylbenzene packings with various sized pores on a single bead. TSK-GEL SuperH and SuperHZ columns are offered for high throughput applications where reductions in run time and solvent consumption are critical. *Table XII (see p. 61)* provides a list of possible shipping solvents along with compatible solvents for each H-type packing.

See *Table II* for suggestions on how to choose between TSK-GEL SW, TSK-GEL PW and TSK-GEL Alpha column types.

Features	Benefits
Rigid hydrophilic packings	Minimal swelling and excellent physical strength Low adsorption resulting in high mass recovery
Four series of SEC columns with different ranges of solvent compatibility ( <i>Figure 1</i> )	Suitable for aqueous GFC and non-aqueous GPC
Easy scale up	Analytical and preparative pre-packed SEC columns

## Characteristics of TSK-GEL Size Exclusion Column Lines


Column Line	TSK-GEL SW / SW <sub>XL</sub> / SuperSW	TSK-GEL PW / PW <sub>XL</sub>	TSK-GEL Alpha / SuperAW	TSK-GEL H <sub>XL</sub> / H <sub>HR</sub> / SuperH / SuperHZ
Resin Type	Silica	Methacrylate	Methacrylate	PS-DVB
No. of Available Pore Sizes	3	7	5	8
pH stability	2.5 - 7.5	2.0 - 12.0	2.0 - 12.0	1.0 - 14.0
Solvent Compatibility	100% polar	50% polar	100% polar and nonpolar	100% nonpolar, limited polar
Max. Temperature	30°C	80°C*	80°C	60°C (1000H-3000H) 80°C (4000H-GMH) 140°C (H <sub>HR</sub> , H- HT)
Max. Flow Rate (mL/min)	0.4 (SuperSW) 1.2 (SW and SW <sub>XL</sub> )	1.0 (PW <sub>XL</sub> ) 1.2 (PW)	0.6 (SuperAW) 1.0 (Alpha)	1.0 (H <sub>HR</sub> ) 1.2 (H <sub>XL</sub> ) 0.8 (SuperH) 0.7 (SuperHZ 6.0 mmID) 0.4 (SuperHZ 4.6 mmID)
Pressure** (kg/cm <sup>2</sup> )	8-120	10-40	20-40	10-60
Application Focus	proteins	water-soluble polymers	intermediate polar polymers	organic-soluble polymers

\* Except for the TSKgel G-DNA-PW, which can be operated up to 50°C and the 55 mm ID TSK-GEL PW-type columns, which can be operated up to 60°C. When operating below 10°C, reduce the flow rate to ensure that the maximum pressure is not exceeded.

\*\* Depends on column dimensions and particle size

Note: The operating conditions and specifications for each column are listed on the Operating Conditions and Specifications sheet (OCS) shipped with the column and in the Ordering Information section at the end of each section.

Figure 1



## Column selection guide for high performance Gel Filtration Chromatography


Sample		Column selection		Selection criteria	
		First choice	Alternative		
Carbo-hydrates	polysaccharides		TSKgel GMPW <sub>XL</sub>	G5000PW <sub>XL</sub> and G3000PW <sub>XL</sub>	large pore size, linear calibration curve, small particles, high resolving power
	oligosaccharides		TSKgel G-Oligo-PW or TSKgel G2000PW	G2500PW <sub>XL</sub>	small particles, high resolving power
Nucleic acids	DNA fragments	large	TSKgel G-DNA-PW or TSKgel G5000PW <sub>XL</sub>	—	large pore size, small particles, high resolving power
		medium and small	TSKgel G4000SW <sub>XL</sub> / SW/ TSKgel BioAssist G4SW <sub>XL</sub> , TSKgel SuperSW3000, or G3000SW <sub>XL</sub> / SW/ TSKgel BioAssist G3SW <sub>XL</sub>	—	suitable pore sizes
	RNA		TSKgel G4000SW <sub>XL</sub> / SW/ TSKgel BioAssist G4SW <sub>XL</sub> , TSKgel SuperSW3000, or G3000SW <sub>XL</sub> / SW/ TSKgel BioAssist G3SW <sub>XL</sub>		suitable pore sizes
	oligonucleotides		TSKgel G2500PW <sub>XL</sub>	—	small pore size, ionic interaction
Proteins	normal size small-medium proteins		TSKgel SuperSW3000, G3000SW <sub>XL</sub> / SW/ TSKgel BioAssist G3SW <sub>XL</sub> , TSKgel G4000SW <sub>XL</sub> / SW, TSKgel BioAssist G4SW <sub>XL</sub> , TSKgel SuperSW2000, or G2000SW <sub>XL</sub> / SW/ TSKgel BioAssist G2SW <sub>XL</sub>	G3000PW <sub>XL</sub> or G4000PW <sub>XL</sub>	small particles small to medium range pore sizes
	large proteins	low density lipoprotein	TSKgel G6000PW <sub>XL</sub> or TSKgel G5000PW <sub>XL</sub>	—	large pore sizes
		gelatin	TSKgel GMPW <sub>XL</sub>	G5000PW <sub>XL</sub> and G3000PW <sub>XL</sub>	large pore size, linear calibration curve
Peptides	large		TSKgel SuperSW3000, G3000SW <sub>XL</sub> / SW/ TSKgel BioAssist G3SW <sub>XL</sub> or G2000SW <sub>XL</sub> / SW/ TSKgel BioAssist G2SW <sub>XL</sub>	SuperSW2000 or G3000PW <sub>XL</sub>	small to medium range pore size, versatile
	small		TSKgel G2500PW <sub>XL</sub>	SuperSW2000 or G2000SW <sub>XL</sub> / SW	linear calibration curve, high resolving power
Viruses			TSKgel G6000PW <sub>XL</sub> or TSKgel G5000PW <sub>XL</sub>	—	large pore size, high resolving power
Synthetic polymers			TSKgel GMPW <sub>XL</sub> or TSKgel Alpha-M	G5000PW and G3000PW <sub>XL</sub> or Alpha-5000 and Alpha-3000	large pore size, low adsorption, linear calibration curve
Synthetic oligomers	nonionic and cationic		TSKgel G-Oligo-PW, TSKgel G2500PW <sub>XL</sub> or TSKgel Alpha-2500	G2500PW or SuperAW2500	small pore size, high resolving power
	anionic		TSKgel G2500PW <sub>XL</sub> or TSKgel Alpha-2500	G2500PW or SuperAW2500	small pore size, ionic interaction

**Polynucleotides**

TSKgel G2000SW, G3000SW, G4000SW, and G5000PW columns are effective in separating double-stranded DNA fragments and ribosomal and transfer RNA. The choice of column is dependent on sample molecular weight. Small nucleic acids are adequately analyzed by using TSK-GEL SW columns. Larger nucleic acids should be analyzed with TSK-GEL PW columns of larger pore size, such as the TSKgel G-DNA-PW and TSKgel G5000PW columns. Calibration curves for double-stranded DNA fragments on TSK-GEL SW type columns and a TSKgel G5000PW column are shown in *Figure 2*; *Table III* lists the recommended TSK-GEL SW and TSK-GEL PW columns for separating double-stranded DNA and RNA fragments.

Separation of four *E. coli* RNAs, shown in *Figure 3*, confirms the better performance of TSK-GEL SW columns for samples with a wide molecular weight range. The sample consists of 4S tRNA (25,000 Da), 5S rRNA (39,000 Da), 16S rRNA (560,000 Da), and 23S rRNA (1,100,000 Da). All four polynucleotides are within the molecular weight range recommended for TSK-GEL SW type columns. The two chromatograms demonstrate a superior separation with the TSKgel G4000SW column.

**Figure 2**



**Recommended TSK-GEL SW and TSK-GEL PW columns for separating double-stranded DNA and RNA fragments**


Base pairs of DNA	Recommended column
< 55	TSKgel G2000SW <sub>XL</sub> /SW or G3000SW <sub>XL</sub> /SW or SuperSW2000 or 3000
55 – 110	TSKgel G3000SW <sub>XL</sub> /SW or SuperSW3000
110 – 375	TSKgel G4000SW <sub>XL</sub> /SW
375 – 1500	TSKgel G5000PW <sub>XL</sub> /PW
1000 – 7000	TSKgel G-DNA-PW

RNA (Da) Molecular Weight	Recommended column
< 60,000	TSKgel G2000SW <sub>XL</sub> /SW or G3000SW <sub>XL</sub> /SW or SuperSW2000 or 3000
60,000 – 120,000	TSKgel G3000SW <sub>XL</sub> /SW or SuperSW3000
120,000 – 1,200,000	TSKgel G4000SW <sub>XL</sub> /SW
1,200,000 – 10,000,000	TSKgel G5000PW <sub>XL</sub>

Note: To determine the approximate molecular weight (Da) of a DNA fragment, multiply the number of base pairs by 650.

**Figure 3**



## Peptides

As Figure 4 shows, the calibration curves for TSKgel G2500PW<sub>XL</sub> and TSKgel G2000SW<sub>XL</sub> are very similar for samples below 3,000 Da. The data was generated using 17 samples ranging in size from myoglobin (17,800 Da) to glycine (75 Da). While the curves are similar in shape through this range of sample sizes, each sample molecule behaved differently on the two columns, indicating additional sample-resin interaction. For example, although an organic solvent was used to reduce hydrophobic effects, the elution of hydrophobic peptide leu-enkephalin was delayed on the TSKgel G2500PW<sub>XL</sub> column.

Small peptides may be difficult to chromatograph by aqueous GFC, due to complex non-size effects such as ionic and hydrophobic interactions. The addition of organic solvents and

buffered salt solutions overcome these effects. Figure 5 compares the separation of two mixtures of peptides on both TSKgel G2500PW<sub>XL</sub> and TSKgel G2000SW<sub>XL</sub> columns to demonstrate which might be superior for a particular type of peptide. The first group of peptides had molecular weights ranging from 6,500 Da to 555 Da. In the second group, the range was from 75 Da to 17,800 Da. The chromatograms confirm that TSKgel G2000SW<sub>XL</sub> columns give higher resolution for most peptide mixtures, but do not perform as well as TSKgel G2500PW<sub>XL</sub> columns at peptide molecular weights lower than 1,000 Da. For very small peptides, the TSK-GEL PW<sub>XL</sub> column type is preferable.

Figure 4





Figure 5



## Proteins

In general, TSK-GEL SW columns have a fractionation range that matches the molecular weight range of most proteins and peptides. However, analytical methods that utilize combinations of TSK-GEL SW and TSK-GEL PW columns in series for isolating large lipoproteins (high and low density), chylomicron, etc. are illustrated in *Figure 6*.


As shown in *Figure 7*, successful elution of the hydrophobic amphoteric polymer collagen (a connective tissue protein) and gelatin on the TSKgel GMPW column requires the addition of 20% CH<sub>3</sub>CN to 0.1M NaNO<sub>3</sub>. Peak areas are reduced and elution is not reproducible when the organic solvent is omitted from the elution buffer for hydrophobic samples.

## IgG

The most suitable column and mobile phase will depend on the particular components that need to be measured. TSK-GEL SW<sub>XL</sub> columns can provide fast, simplified quality control analyses of proteins, peptides and other large-size biopharmaceuticals. For example, a therapeutic solution of intravenous IgG may contain albumin as a stabilizer, and both proteins must be quantified following manufacture. Although literature reports describe the separation of these two proteins by many other chromatographic methods, long analysis times and complex gradient elutions are required. A method developed on TSKgel G3000SW<sub>XL</sub> provides quantitative separation of the two proteins in 15 minutes with a simple, isocratic elution system. As shown in *Figure 8*, human albumin can be separated from a 20-fold excess of IgG, and quantified using an optimized elution buffer. This simple separation method can be applied to the isolation of other IgGs, such as monoclonal antibodies in ascites fluid.

Figure 7


## Use of mobile phase additives to improve the elution profile of amphoteric polymers



Column: TSKgel GMPW, two 17  $\mu$ m, 7.5 mm ID x 60 cm L columns in series  
 Sample: 0.5 mL of 0.05-0.1% collagen, Type VI, or gelatin, Type I  
 Elution: A. distilled water, B. 0.1 M NaNO<sub>3</sub> in water, C. 20% CH<sub>3</sub>CN in 0.1 M NaNO<sub>3</sub>  
 Flow Rate: 0.5 mL/min  
 Detection: RI

Figure 6


## The relation between molecular weight of lipoproteins and elution volume for the combination GFC columns



Column: 7.5 mm ID x 60 cm L  
 Sample:  $\nabla$  chylomicron  
 $\blacktriangledown$  VLDL  
 $\bullet$  LDL  
 $\blacksquare$  HDL<sub>2</sub>  
 $\circ$  HDL<sub>3</sub>  
 $\square$  albumin  
 $\blacktriangle$  ovalbumin  
 Elution: 0.1 M Tris-HCl buffer, pH 7.4  
 Flow Rate: 1.0 mL/min  
 Detection: on-line post column reaction

Figure 8

## QC test for albumin in intravenous human IgG



Column: TSKgel G3000SW<sub>XL</sub>, 5  $\mu$ m, 7.8 mm ID x 30 cm L  
 Sample: 5  $\mu$ L of Venilon<sup>®</sup>, containing 237.5 mg of IgG and 12.5 mg of albumin  
 Elution: 0.1 M Na<sub>2</sub>SO<sub>4</sub> in 0.05 M sodium phosphate buffer, pH 5.0  
 Flow Rate: 1.0 mL/min  
 Detection: UV @ 280 nm

## Cleaning Procedures

Please see Appendix A for instructions on how to maintain and clean TSK-GEL Size Exclusion columns.

### Improving resolution with SEC columns

- Verify that the column is not overloaded (refer to sample load information in the SW column section).
- Decrease the dead volume in the HPLC system by using the shortest tubing lengths and the smallest tubing ID possible without exceeding the maximum pressure for the column.
- Decrease the flow rate, but not lower than 0.3 mL/min for TSK-GEL SW and SW<sub>XL</sub> columns because increased diffusion will occur.

*Note: TSK-GEL SuperSW columns show optimum resolution at flow rates below 0.3mL/min because of their small diameter.*

- If using a 30 cm column, add an additional 30 cm column or switch to a TSK-GEL<sub>XL</sub>-type column that has a smaller particle size. Resolution will improve more than two-fold. Alternatively, if using the SW<sub>XL</sub>-type columns, switch to the TSK-GEL SuperSW that offer increased sensitivity and resolution due to 50% more theoretical plates.

### Literature

For additional information describing applications of Size Exclusion Chromatography columns, or how to select the optimal Size Exclusion Chromatography column, please contact our Technical Service specialists or refer to Tosoh HPLC database on our website: [www.tosohbioscience.com](http://www.tosohbioscience.com).

# Silica-based TSK-GEL SW, TSK-GEL SW<sub>XL</sub> and TSK-GEL SuperSW columns for aqueous Gel Filtration Chromatography of proteins and peptides

## Highlights

- **-NEW-** PEEK column hardware available for SW<sub>XL</sub> packings.
- TSK-GEL SW-type packings are comprised of rigid spherical silica gel chemically bonded with hydrophilic compounds.
- Low adsorption and well defined pore size distribution, which are required for high performance SEC.
- Particles having three different pore sizes are available packed as the standard SW or SW<sub>XL</sub> (higher resolution) columns.
- TSK-GEL SuperSW columns feature a 4 µm particle size packing in a narrower 4.6 mm ID column for the highest resolution and sensitivity in the SW series.
- The TSK-GEL SuperSW as well as the larger particle size TSK-GEL QC-PAK (fast analysis) columns are available in two pore sizes. (The properties and molecular weight ranges of these packings are summarized in *Table IV*).
- Semi-preparative (21.5 mm ID) and preparative (55 mm, 108 mm ID) stainless steel columns are available for precise scale-up to commercial production.

## Column selection

The three different pore sizes for the TSK-GEL SW packings results in different exclusion limits for several sample types, as shown by the calibration curves in *Figure 9*. From this data, recommended separation ranges for different sample types can be made for each column (*Table IV*). The ranges for polyethylene glycol (PEG) should be consulted when choosing a column for straight chain molecules, dextran for branched molecules, and globular proteins for globular molecules.

Standard TSK-GEL SW columns are packed with 10 µm particles (13 µm particles in G4000SW columns). The higher performance TSK-GEL SW<sub>XL</sub> versions are packed with 5 µm particles (8 µm particles in G4000SW<sub>XL</sub> columns). *Figure 10* shows the improved resolution and reduced analysis time with a protein standard mixture on the TSK-GEL SW<sub>XL</sub> compared with SW analytical columns. The calibration curves for protein standards for the three types of SW columns are compared on the TSK-GEL SW<sub>XL</sub> packings in *Figure 11*.

For preliminary research or reducing quality control testing time, the 15 cm long QC-PAK columns provide analysis times half as long as those on standard 30 cm columns, while retaining baseline resolution of many protein mixtures. They are available in glass or stainless steel packed with the high performance 5 µm TSK-GEL SW<sub>XL</sub> packings. TSKgel QC-PAK GFC 200 columns contain the G2000SW<sub>XL</sub> packing, while TSKgel QC-PAK GFC 300 columns contain the G3000SW<sub>XL</sub> packing.

*Table IV*

### Properties and separation ranges for TSK-GEL SW-type packings

TSK-GEL packing	Particle Size (µm)	Pore Size (Å)	Molecular weight of sample (Da)		
			Globular proteins	Dextrans	Polyethylene glycols and oxides
SuperSW2000	4	125	5,000–1.5 x 10 <sup>5</sup>	1,000–3 x 10 <sup>4</sup>	500–15,000
G2000SW <sub>XL</sub>	5	125	5,000–1.5 x 10 <sup>5</sup>	1,000–3 x 10 <sup>4</sup>	500–15,000
BioAssist G2SW <sub>XL</sub>	5	125	5,000–1.5 x 10 <sup>5</sup>	1,000–3 x 10 <sup>4</sup>	500–15,000
QC-PAK GFC 200	5	125	5,000–1.5 x 10 <sup>5</sup>	1,000–3 x 10 <sup>4</sup>	500–15,000
G2000SW	10, 13	125	5,000–1 x 10 <sup>5</sup>	1,000–3 x 10 <sup>4</sup>	500–15,000
SuperSW3000	4	250	1 x 10 <sup>4</sup> –5 x 10 <sup>5</sup>	2,000–7 x 10 <sup>4</sup>	1,000–3.5 x 10 <sup>4</sup>
G3000SW <sub>XL</sub>	5	250	1 x 10 <sup>4</sup> –5 x 10 <sup>5</sup>	2,000–7 x 10 <sup>4</sup>	1,000–3.5 x 10 <sup>4</sup>
BioAssist G3SW <sub>XL</sub>	5	250	1 x 10 <sup>4</sup> –5 x 10 <sup>5</sup>	2,000–7 x 10 <sup>4</sup>	1,000–3.5 x 10 <sup>4</sup>
QC-PAK GFC 300	5	250	1 x 10 <sup>4</sup> –5 x 10 <sup>5</sup>	2,000–7 x 10 <sup>4</sup>	1,000–3.5 x 10 <sup>4</sup>
G3000SW	10	250	1 x 10 <sup>4</sup> –5 x 10 <sup>5</sup>	2,000–7 x 10 <sup>4</sup>	1,000–3.5 x 10 <sup>4</sup>
G4000SW <sub>XL</sub>	8	450	2 x 10 <sup>4</sup> –7 x 10 <sup>6</sup>	4,000–5 x 10 <sup>5</sup>	2,000–2.5 x 10 <sup>5</sup>
BioAssist G4SW <sub>XL</sub>	8	450	2 x 10 <sup>4</sup> –7 x 10 <sup>6</sup>	4,000–5 x 10 <sup>5</sup>	2,000–2.5 x 10 <sup>5</sup>
G4000SW	13, 17	450	2 x 10 <sup>4</sup> –7 x 10 <sup>6</sup>	4,000–5 x 10 <sup>5</sup>	2,000–2.5 x 10 <sup>5</sup>

Data generated using the following conditions:

**Columns:** Two 4 µm, 4.6 mm ID x 30 cm L TSK-GEL SuperSW columns in series; two 5 µm, 7.8mm ID x 30 cm L TSK-GEL SW<sub>XL</sub> columns in series; two 10 µm, 7.5 mm ID x 60 cm L TSK-GEL SW columns in series

**Elution:** Globular proteins: 0.3 M NaCl in 0.1 M (0.05 M for SW<sub>XL</sub> columns) phosphate buffer, pH 7  
Dextrans and polyethylene glycols and oxides (PEOs): distilled water

### Polyethylene oxide, dextran and protein calibration curves for TSK-GEL SW columns





Figure 10

### Higher resolution with 5 $\mu\text{m}$ TSK-GEL SW<sub>XL</sub> compared with 10 $\mu\text{m}$ TSK-GEL SW columns




Column: Left: TSK-GEL SW, two 10  $\mu\text{m}$  7.5 mm ID x 30 cm L columns in series  
 Right: TSK-GEL SW<sub>XL</sub>, one 5  $\mu\text{m}$ , 7.8 mm ID x 30 cm L column  
 Sample: 1. glutamate dehydrogenase (55,000 Da)  
 2. lactate dehydrogenase (36,500 Da)  
 3. enolase (67,000 Da)  
 4. adenylate kinase (6,000 Da)  
 5. cytochrome C (12,400 Da)  
 Elution: 0.3 M NaCl in 0.05 M phosphate buffer, pH 6.9  
 Flow Rate: 1.0 mL/min  
 Detection: UV @ 220 nm

To further improve efficiency and resolution, TSK-GEL SuperSW columns are filled with 4  $\mu\text{m}$  particles. These columns feature a 4.6 mm ID to provide better sensitivity in sample-limited cases. A minimum of 30,000 plates/column is guaranteed. Two pore sizes are available, 125 Å and 250 Å, which provide molecular weight ranges as shown in *Table IV*. *Figure 12* shows similarities of the calibration curves of the TSK-GEL SuperSW and SW<sub>XL</sub>. Chromatograms of a standard protein mixture in *Figure 13* illustrate the increased sensitivity and resolution of the TSKgel SuperSW3000 compared with the TSKgel G3000SW<sub>XL</sub>. *Figure 14* depicts the same benefits on the SuperSW2000.


The eluent plays an important role in determining the appropriate column for a separation. Since the structure of globular proteins is compact, the molecular weight exclusion limit for proteins is higher than for linear polymers such as polyethylene oxides, dextrans, or double-stranded DNA. However, when denaturing eluents are used, the exclusion limit for proteins becomes more like the exclusion limit for linear polymers. *Table V* lists the suggested TSK-GEL SW and TSK-GEL PW columns for the separation of proteins and peptides with different eluents. *Figure 15* shows a calibration curve for proteins denatured with 0.1% SDS on the SuperSW3000.

Figure 11

Protein calibration curves for TSK-GEL SW<sub>XL</sub> columns


Column: TSK-GEL SW<sub>XL</sub> columns, 5 or 8  $\mu$ m, 7.8 mm ID x 30 cm L  
 Sample: 1. thyroglobulin (660,000 Da); 2. IgG (156,000 Da);  
 3. bovine serum albumin (67,000 Da);  
 4. ovalbumin (43,000 Da); 5. peroxidase (40,200 Da); 6.  $\beta$ -lactoglobulin (35,000 Da);  
 7. myoglobin (16,900 Da);  
 8. ribonuclease A (13,700 Da); 9. cytochrome C (12,400 Da); 10. glycine tetramer (246 Da)  
 Elution: 0.3 M NaCl in 0.1 M sodium phosphate buffer, pH 7  
 Detection: UV @ 220 nm

Figure 12

Calibration curves for TSK-GEL SuperSW and SW<sub>XL</sub>


Sample: proteins: 1. thyroglobulin (670,000 Da);  
 2.  $\gamma$ -globulin (155,000 Da); 3. BSA (66,300 Da);  
 4.  $\beta$ -lactoglobulin (18,400 Da); 5. lysozyme (14,300 Da);  
 6. cytochrome C (12,400 Da); 7. triglycine (189 Da)  
 Elution: 0.15 M phosphate buffer (pH 6.8)  
 Flow Rate: 0.35 mL/min for SuperSW; 1.0 mL/min for SW<sub>XL</sub>  
 Temperature: 25°C  
 Detection: UV @ 280 nm (220 nm for triglycine)

Figure 13

Comparison of TSKgel SuperSW3000 and TSKgel G3000SW<sub>XL</sub> for the separation of proteins

Column: A. TSKgel G3000SW<sub>XL</sub>, 7.8 mm ID x 30 cm L; B. TSKgel Super SW3000, 4.6 mm ID x 30 cm L  
 Sample: 5  $\mu$ L of a mixture of 1. thyroglobulin, 0.5 mg/mL (660,000 Da); 2.  $\gamma$ -globulin, 1.0 mg/mL (150,000 Da);  
 3. ovalbumin, 1.0 mg/mL (43,000 Da); 4. ribonuclease A, 1.5 mg/mL (12,600 Da);  
 5. *p*-aminobenzoic acid, 0.01 mg/mL (137 Da)  
 Elution: 0.1 M Na<sub>2</sub>SO<sub>4</sub> in 0.1 M phosphate buffer with 0.05% NaN<sub>3</sub>, pH 6.7  
 Flow Rate: 1.0 mL/min for G3000SW<sub>XL</sub>; 0.35 mL/min for SuperSW3000  
 Temperature: 25°C  
 Detection: UV @ 220 nm


### Comparison of chromatograms of proteins on TSKgel SuperSW2000 and TSKgel G2000SW<sub>XL</sub>



Column: A. TSKgel G2000SW<sub>XL</sub>, 7.8 mm ID x 30 cm L; B. TSKgel SuperSW2000, 4.6 mm ID x 30 cm L  
 Sample: 1. thyroglobulin (0.2 mg/mL); 2. albumin (1.0 mg/mL); 3. ribonuclease A (1.0 mg/mL); 4. *p*-aminobenzoic acid (0.01 mg/mL)  
 Injection Volume: 5  $\mu$ L  
 Elution: 0.1 M phosphate buffer + 0.1 M Na<sub>2</sub>SO<sub>4</sub> + 0.05% NaN<sub>3</sub> (pH 6.7)  
 Flow Rate: 0.35 mL/min for SuperSW2000; 1.0 mL/min for G2000SW<sub>XL</sub>  
 Temperature: 25°C  
 Detection: UV @ 280 nm

Figure 15

### Effect of buffer concentration on calibration curves of proteins on SuperSW3000 in 0.1% SDS eluent



Sample: proteins\*: 1. phosphorylase (94,000 MW);  
 2. BSA (67,000 MW); 3. ovalbumin (43,000 MW);  
 4. carbonic anhydrase (30,000 MW);  
 5. soybean trypsin inhibitor (20,100 MW);  
 6.  $\alpha$ -lactalbumin (14,400 MW)  
 Eluent: 50 mM-200 mM phosphate buffer (pH 6.8) containing 0.1% SDS  
 Flow Rate: 0.35 mL/min  
 Temperature: 25°C  
 Detection: UV @ 280 nm

\* Proteins are denatured in phosphate buffer containing SDS and DTT at 40°C during 15 min.

Table V

### Recommended TSK-GEL SW and TSK-GEL PW columns for separating peptides and proteins

Sample (Da)	Recommended column
<b>Peptides</b>	
with 45% CH <sub>3</sub> CN in 0.1% TFA <sup>a</sup>	
500 – 2.5 x 10 <sup>4</sup>	TSKgel G2000SW <sub>XL</sub> <sup>b</sup> /BioAssist G2SW <sub>XL</sub> <sup>b</sup> /SuperSW2000 <sup>b</sup>
100 – 5 x 10 <sup>3</sup>	TSKgel G2500PW <sub>XL</sub> <sup>c</sup>
1,000 – 5 x 10 <sup>4</sup>	TSKgel G3000PW <sub>XL</sub> <sup>d</sup>
<b>Peptides and proteins</b>	
with 0.1% SDS in 0.1 M phosphate, pH 7.0	
1,000 – 2.5 x 10 <sup>4</sup>	TSKgel G2000SW <sub>XL</sub> /SW/BioAssist G2SW <sub>XL</sub> or SuperSW2000
10,000 – 7 x 10 <sup>4</sup>	TSKgel G3000SW <sub>XL</sub> /SW/BioAssist G3SW <sub>XL</sub> or SuperSW3000
15,000 – 3 x 10 <sup>5</sup>	TSKgel G4000SW <sub>XL</sub> or SW
with 6M guanidine-HCl in 0.1M phosphate, pH 6.0	
1,000 – 2.5 x 10 <sup>3</sup>	TSKgel G2000SW <sub>XL</sub> /SW/BioAssist G2SW <sub>XL</sub> or SuperSW2000
2,000 – 7 x 10 <sup>4</sup>	TSKgel G3000SW <sub>XL</sub> /SW <sup>e</sup> /BioAssist G3SW <sub>XL</sub> or SuperSW3000
3,000 – 4 x 10 <sup>5</sup>	TSKgel G4000SW <sub>XL</sub> /SW
<b>Proteins</b>	
with 0.3 M NaCl in 0.05 ~ 0.1 M phosphate, pH 7.0	
5 x 10 <sup>3</sup> – 1 x 10 <sup>5</sup>	TSKgel G2000SW <sub>XL</sub> /SW/BioAssist G2SW <sub>XL</sub> or SuperSW2000
1 x 10 <sup>4</sup> – 5 x 10 <sup>5</sup>	TSKgel G3000SW <sub>XL</sub> /SW/BioAssist G3SW <sub>XL</sub> or SuperSW3000
2 x 10 <sup>4</sup> – 7 x 10 <sup>6</sup>	TSKgel G4000SW <sub>XL</sub> or SW

#### Notes:


- Proteins may denature in this eluent, increasing the apparent molecular weight.
- Recommended as first column for peptide separation since it offers high resolution across a wide fractionation range.
- Recommended for small peptides. Lower overall resolution than the SW column, but durable to alkaline solutions.
- Recommended for larger peptides. TSKgel G3000SW<sub>XL</sub> is also effective with mild buffer for peptide fragments with molecular weight >10,000 Da.
- TSKgel G4000SW or SW<sub>XL</sub> is recommended for polypeptides above 60,000 Da.

## Applications


### Proteins

The effect of different concentrations of surfactant on the separation of membrane proteins is seen in *Figure 16*. As the concentration of octaethyleneglycol dodecylether increases to 0.05%, the main peak becomes sharper and recovery increases. The TSKgel SuperSW3000 provides an excellent high resolution separation of IgG<sub>1</sub> from mouse ascites fluid as can be seen in *Figure 17*.

*Figure 16*




*Figure 17*



### Enzymes

Mobile phase conditions in GFC are optimized to ensure none or minimal interaction of the sample with the packing material. This gentle technique allows for high recovery of enzymatic activity. For example, crude samples of peroxidase and glutathione S-transferase were separated in only 15 minutes on a TSKgel G3000SW<sub>XL</sub> column, and activity recovery was 98% and 89% respectively. The elution profiles of the separations in *Figure 18* show that all of the activity eluted in a narrow band of about 1.5 mL.

*Figure 18*



## Optimizing GFC with TSK-GEL SW, SW<sub>XL</sub> and SuperSW Columns

### Sample Load

In all modes of chromatography, high sample loads distort peak shapes and cause an overall decrease in efficiency due to column overload. Sample loads may be increased by using organic solvents to enhance the solubility of the sample or by using higher column temperatures to lower the viscosity of the mobile phase. When recovering protein activity, however, these techniques are usually not recommended.


The sample load limitations of TSK-GEL SW and TSK-GEL SW<sub>XL</sub> analytical and semi-preparative columns are shown in *Figure 19*. At protein injections higher than 0.8 mg, HETP increases rapidly and column efficiency decreases with TSK-GEL SW and TSK-GEL SW<sub>XL</sub> analytical columns. Sample concentrations in the

range of 1-20 mg/mL are recommended, but proteins can be loaded at higher concentrations and higher total loads than synthetic macromolecules. For example, a TSKgel G3000SW semi-preparative column can provide high efficiency for 100 mg of bovine serum albumin (BSA) shown in *Figure 19*. For a synthetic, linear polyethylene glycol (7,500 Da), sample loads higher than 20 mg cause a rapid increase in HETP.


SuperSW columns provide excellent resolution at a sample load of less than 0.3 mg of BSA as seen in *Figure 20*. A sample volume less than 20  $\mu$ L is recommended for the SuperSW2000 column, and less than 80  $\mu$ L for the SuperSW3000 column, illustrated in *Figure 21*. However, a 5  $\mu$ L injection ensures optimal results.

If a larger sample load is necessary, a 7.5 mm ID x 60 cm L SW column has twice the capacity of a 30 cm SW<sub>XL</sub> column with very little loss of resolution.


*Figure 19*



*Figure 20*



### Comparison of sample volume on TSKgel SuperSW and TSKgel SW<sub>XL</sub>



Column: A. TSKgel SuperSW2000, 4.6 mm ID x 30 cm L and TSKgel G2000SW<sub>XL</sub>, 7.8 mm ID x 30 cm L  
 B. TSKgel SuperSW3000, 4.6 mm ID x 30 cm L and TSKgel G3000SW<sub>XL</sub>, 7.8 mm ID x 30 cm L  
 Sample: 1 mg/mL of bovine serum albumin  
 Elution: 0.1 M Na<sub>2</sub>SO<sub>4</sub> in 50 mM phosphate buffer (pH 6.7)  
 Flow Rate: 0.20 mL/min for Super SW; 0.57 mL/min for SW<sub>XL</sub>  
 Detection: UV @ 280 nm

### Selecting mobile phase buffers

Proper selection of elution conditions is necessary to maximize molecular sieving mechanisms and to minimize secondary effects, such as ionic and hydrophobic interactions between the sample and the column packing material. Under conditions of high ionic strength (>1.0 M), hydrophobic interactions may occur. Under low ionic strength (<0.1 M), ionic interactions are more likely to occur. Secondary interactions take place more often with small solutes since the residual silanol sites of TSK-GEL SW packings are mainly located in small pores, which are accessible to small solutes. In general, the use of relatively high ionic strength buffers can be used to overcome secondary interactions. A neutral salt, such as sodium sulfate, is often added to increase buffer ionic strength.

While the extent of secondary interactions is not as great for proteins as for smaller molecules, complex interactions can occur. For each protein, there will be an optimum buffer type and concentration that results in the highest resolution and recovery. As shown in *Table VI*, resolution among protein samples varies according to buffer type and strength. Although each salt used has the same molar concentration in the mobile phase, each provides a different ionic strength, from the lowest in strength, sodium chloride, to the highest, ammonium phosphate.

The ionic species of the mobile phase will also affect the separation, shown in *Table VI*, by the difference in resolution values for the different salts. In addition, calibration curves for proteins in potassium phosphate buffers are more shallow than those generated in sodium phosphate buffers, while the slope of

the curve in Sorenson buffer (containing both Na<sup>+</sup> and K<sup>+</sup>) is midway between the slopes generated with either cation alone. *Table VII* illustrates the impact of different buffer conditions on mass recovery for six sample proteins. In this case, the mass recovery of proteins is higher with sodium or potassium phosphate buffers (pH 6.9) than with Tris-HCl buffers (pH 7.8).

If hydrophobic interaction occurs between the sample and the column matrix, up to 100% water soluble organic, such as acetonitrile, acetone, methanol or ethanol, can be added to the mobile phase. Alternatively, 1% nonionic detergent (0.1% SDS), 8M urea, or 6 M guanidine can be used. Viscous mobile phases may require a flow rate reduction if the maximum pressure is exceeded (listed on the Operating Condition and Specification [OCS] sheet or the Ordering Information section at the end of this catalog section).

### Choice of mobile phase pH

The optimal buffer pH for a separation depends on the protein's isoelectric point (pI) and the stability of the packing. In general, Tosoh Bioscience recommends the use of a buffer with a pH that is different from the pI of the protein so the protein will carry an overall net positive or negative charge. This generally increases the protein's stability. The pH stability of TSK-GEL SW packings ranges from 2.5 to 7.5. If the separation must be performed under basic conditions, the pH stable polymer-based TSK-GEL PW or TSK-GEL Alpha columns should be used.

Table VI

**Resolution factors ( $R_s$ ) of proteins with various ionic species added to the elution buffer**

	<i>Da</i>	<i>0.2 M (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub></i>	<i>0.2 M Na<sub>2</sub>SO<sub>4</sub></i>	<i>0.2 M (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub></i>	<i>0.2 M MgCl<sub>2</sub></i>	<i>0.2 M NaCl</i>
Glutamate dehydrogenase	280,000	}	1.4	1.1	1.5	1.3
Alcohol dehydrogenase	150,000		0.4	0.5	0.2	0.5
Glutathione reductase	113,000		0.9	1.1	0.9	0.9
Enolase	67,000		1.1	0.8	1.3	0.8
Ovalbumin	43,000		1.3	1.3	1.4	1.9
Trypsinogen	24,000		1.0	0.5	0.2	0.8
Cytochrome C	12,400					

Column: TSKgel G3000SW, 7.5 mm ID x 30 cm L

Elution: 0.05 M GTA buffer, pH 7.0 (G: 3,3-dimethylglutaric acid, T: trishydroxyaminomethane, A: 2-amino-2-methyl-1,3-propanediol) plus indicated salt

Detection: UV @ 220 nm

Table VII

**Effect of buffer composition on mass recovery of proteins on TSKgel G3000SW**

<i>Protein</i>	<i>Recovery (%)</i>		
	<i>A. Sodium phosphate buffer</i>	<i>B. Potassium phosphate buffer</i>	<i>C. Tris-HCl buffer</i>
Cytochrome C	98	101	92
Lysozyme	92	96	75
$\alpha$ -Chymotrypsinogen	95	98	90
IgG	95	98	88
Thyroglobulin	94	94	85
Ovalbumin	96	92	66

Column: TSKgel G3000SW, 7.5 mm ID x 60 cm L

Elution: A. 0.2 M NaH<sub>2</sub>PO<sub>4</sub> and 0.2 M Na<sub>2</sub>HPO<sub>4</sub>, pH 6.9;

B. 0.2 M KH<sub>2</sub>PO<sub>4</sub> and 0.2 M K<sub>2</sub>HPO<sub>4</sub>, pH 6.9;

C. 0.2 M NaCl and 0.05 M Tris-HCl, pH 7.8

Flow Rate: 1.0 mL/min

Detection: UV @ 220 nm

## Ordering Information

### Analytical and preparative TSK-GEL Size Exclusion silica-based column products: typical properties

See Appendix A for total column volumes and void volumes for each column size

Part #	Description	ID (mm)	Length (cm)	Particle Size (µm)	Min. Number Theoretical Plates	Flow Rate (mL/min)		Maximum Pressure Drop (kg/cm <sup>2</sup> )
						Range	Max.	
<b>Glass columns</b>								
16214	QC-PAK GFC 200, Glass, 125 Å	8.0	15	5	10,000	0.5 – 1.0	1.2	40
16216	QC-PAK GFC 300, Glass, 250 Å	8.0	15	5	10,000	0.5 – 1.0	1.2	40
08799	G2000SW, Glass, 125 Å	8.0	30	10	10,000	0.4 – 0.8	0.8	20
08800	G3000SW, Glass, 250 Å	8.0	30	10	10,000	0.4 – 0.8	0.8	20
08801	G4000SW, Glass, 450 Å	8.0	30	13	8,000	0.4 – 0.8	0.8	20
14464	G3000SW, Glass, 250 Å	20.0	30	13	6,000	3.0 – 6.0	8.0	8
<b>Stainless steel columns</b>								
18674	SuperSW2000, 125 Å	4.6	30	4	30,000	0.1 – 0.35	0.4	120
18675	SuperSW3000, 125 Å	4.6	30	4	30,000	0.1 – 0.35	0.4	120
08540	G2000SW <sub>XL</sub> , 125 Å	7.8	30	5	20,000	0.5 – 1.0	1.2	70
08541	G3000SW <sub>XL</sub> , 250 Å	7.8	30	5	20,000	0.5 – 1.0	1.2	70
08542	G4000SW <sub>XL</sub> , 450 Å	7.8	30	8	16,000	0.5 – 1.0	1.2	35
16215	QC-PAK GFC 200, 125 Å	7.8	15	5	10,000	0.5 – 1.0	1.2	40
16049	QC-PAK GFC 300, 250 Å	7.8	15	5	10,000	0.5 – 1.0	1.2	40
05788	G2000SW, 125 Å	7.5	30	10	10,000	0.5 – 1.0	1.2	20
05789	G3000SW, 250 Å	7.5	30	10	10,000	0.5 – 1.0	1.2	25
05790	G4000SW, 450 Å	7.5	30	13	8,000	0.5 – 1.0	1.2	15
05102	G2000SW, 125 Å	7.5	60	10	20,000	0.5 – 1.0	1.2	40
05103	G3000SW, 250 Å	7.5	60	10	20,000	0.5 – 1.0	1.2	50
05104	G4000SW, 450 Å	7.5	60	13	16,000	0.5 – 1.0	1.2	30
06727	G2000SW, 125 Å	21.5	30	13	10,000	3.0 – 6.0	8.0	10
06728	G3000SW, 250 Å	21.5	30	13	10,000	3.0 – 6.0	8.0	15
06729	G4000SW, 450 Å	21.5	30	17	8,000	3.0 – 6.0	8.0	10
05146	G2000SW, 125 Å	21.5	60	13	20,000	3.0 – 6.0	8.0	20
05147	G3000SW, 250 Å	21.5	60	13	20,000	3.0 – 6.0	8.0	30
05148	G4000SW, 450 Å	21.5	60	17	16,000	3.0 – 6.0	8.0	20
07428	G2000SW, 125 Å	55.0	30	20	4,000	15.0 – 25.0	50.0	10
07481	G3000SW, 250 Å	55.0	30	20	4,000	15.0 – 25.0	50.0	10
07429	G2000SW, 125 Å	55.0	60	20	9,000	15.0 – 25.0	50.0	15
07482	G3000SW, 250 Å	55.0	60	20	9,000	15.0 – 25.0	50.0	15
<b>PEEK Columns</b>								
20027	BioAssist G2SW <sub>XL</sub> -NEW-	7.8	30	5	20,000	0.5 – 1.0	1.2	70
20026	BioAssist G3SW <sub>XL</sub> -NEW-	7.8	30	5	20,000	0.5 – 1.0	1.2	70
20027	BioAssist G4SW <sub>XL</sub> -NEW-	7.8	30	8	16,000	0.5 – 1.0	1.2	35

Pore Sizes: SuperSW2000, QC-PAK GFC 200, G2000SW, G2000SW<sub>XL</sub> and BioAssist G2SW<sub>XL</sub> = 125 Å; SuperSW3000, QC-PAK GFC 300, G3000SW, G3000SW<sub>XL</sub> and BioAssist G3SW<sub>XL</sub> = 250 Å; G4000SW, G4000SW<sub>XL</sub> and BioAssist G4SW<sub>XL</sub> = 450 Å.

## Ordering Information

**Analytical and preparative TSK-GEL Size Exclusion silica-based column products: typical properties** (continued)

See Appendix A for total column volumes and void volumes for each column size

Part #	Description	ID (mm)	Length (cm)	Particle Size (µm)	
<b>Guard columns</b>					
08805	SW Guard column, Glass	8.0	4.0	10	For all 8 mm ID SW and QC-PAK glass
14465	SW Guard column, Glass	20.0	4.0	20	For P/N 14464 SW glass
18762	SuperSW Guard column	4.6	3.5	4	For 4.6 mm ID SuperSW (contains SuperSW3000 packing)
08543	SW <sub>XL</sub> Guard column	6.0	4.0	7	For all SW <sub>XL</sub> and P/Ns 16215 and 16049 (contains 3000SW <sub>XL</sub> packing)
18008	SW <sub>XL</sub> Guard column, PEEK	6.0	4.0	7	For all BioAssist SW <sub>XL</sub> , PEEK columns
05371	SW Guard column	7.5	7.5	10	For all 7.5 mm ID SW (contains 3000SW packing)
05758	SW Guard column	21.5	7.5	13	For all 21.5 mm ID SW
07427	SW Guard column	45.0	5.0	20	For 55 mm ID SW
<b>Bulk packing</b>					
08544	SW <sub>XL</sub> Top-Off, 1g wet gel			5	For SW <sub>XL</sub> and QC-PAK
06819	SW Top-Off, 1g wet gel			10	For all 7.5 mm ID SW

# Polymer-based TSK-GEL PW and TSK-GEL PW<sub>XL</sub> columns for Gel Filtration Chromatography of water soluble polymers

## Highlights

- **-NEW-** PEEK column hardware available for G6000PW packings for ultra-low sample adsorption during virus analysis.
- Hydrophilic, rigid, spherical, porous methacrylate beads.
- Excellent chemical and mechanical stability.
- pH range of 2 to 12, with up to 50% organic solvent
- Temperatures up to 80°C (50°C for TSKgel G-DNA-PW).
- Wide molecular weight separation range up to  $8 \times 10^6$  Da for linear polymers
- Available in analytical (7.5 mm & 7.8 mm ID), semi-preparative (21.5 mm ID) and preparative (55 mm and 108 mm ID) columns.

Polymeric TSK-GEL PW and TSK-GEL PW<sub>XL</sub> columns are designed for Gel Filtration Chromatography (GFC) of water soluble organic polymers, proteins, peptides, polysaccharides, oligosaccharides, DNA, and RNA.

## Column Selection

The properties and molecular weight separation ranges for all TSK-GEL PW type columns are summarized in *Table VIII*. Resins with six pore sizes are available. Particle sizes of resins packed in TSK-GEL PW columns are 10 μm, 17 μm, 20 μm, or 22 μm. The resins in TSK-GEL PW<sub>XL</sub> columns are 6 μm, 10 μm, or 13 μm. Thus, PX<sub>XL</sub> resins are higher performance versions of the TSK-GEL PW resins due to the smaller particle sizes. Specialty resin-based columns

include the mixed-bed TSKgel GMPW and TSKgel GMPW<sub>XL</sub> for samples with a broad molecular weight range. They also include TSKgel G-Oligo-PW and TSKgel G-DNA-PW columns for oligosaccharides and for DNA or RNA respectively.

The eluent plays an important role in determining the appropriate column for a separation. Since the structure of globular proteins is compact, the molecular weight exclusion limit for proteins is higher than for linear polymers such as polyethylene oxides, dextrans or double-stranded DNA. However, when denaturing eluents are used, the exclusion limit for proteins is comparable to the exclusion limit for linear polymers. Therefore, consult *Table V* in the SW section for help in choosing a column. Use the molecular weight ranges for PEO in *Table VIII* when choosing a column for linear molecules; dextran for branched molecules, and globular proteins for globular molecules. Calibration curves for polyethylene glycols chromatographed on TSK-GEL PW and TSK-GEL PW<sub>XL</sub> columns are shown in *Figure 22*. Protein calibration curves on TSKgel PW<sub>XL</sub> columns are presented in *Figure 23*. Although many methods for polymer analysis have been satisfactorily developed on TSK-GEL PW columns, higher resolution can be achieved with a TSK-GEL PW<sub>XL</sub> column. The smaller particle sizes of the TSK-GEL PW<sub>XL</sub> columns provide almost 2.5 times the resolution of their TSK-GEL PW counterparts. In addition, with shorter TSK-GEL PW<sub>XL</sub> columns, higher resolution separations are possible in less than half the time, as shown in *Figure 24* with polyethylene glycol and polyethylene oxide standards.

For analytical purposes, the TSK-GEL PW<sub>XL</sub> columns are preferred. For preparative work, or for other cases in which large amounts of sample must be used, the 60 cm TSK-GEL PW columns are recommended because of their increased loading capacity.

*Table VIII*

**Properties and molecular weight separation ranges for TSK-GEL PW packings**

TSKgel column	Particle Size* (μm)	Average pore Size (Å)	Molecular weight of sample (Da)		
			Polyethylene glycols & oxides	Dextrans**	Globular proteins**
G1000PW	10	<100	up to 1,000	—	<2,000
G2000PW	10, 17, 20	125	up to 2,000	—	<5,000
G2500PW <sub>XL</sub>	6	<200	up to 3,000	—	<8,000
G2500PW	10, 17, 20				
G3000PW <sub>XL</sub>	6	200	up to $5 \times 10^4$	up to $6 \times 10^4$	$500-8 \times 10^5$
G3000PW	10, 17, 20				
G4000PW <sub>XL</sub>	10	500	$2,000-3 \times 10^4$	$1,000-7 \times 10^5$	$1 \times 10^4-1.5 \times 10^6$
G4000PW	17, 22				
G5000PW <sub>XL</sub>	10	1000	$4,000-1 \times 10^6$	$5 \times 10^4-2.5 \times 10^6$	$<1 \times 10^7$
G5000PW	17, 20, 22				
G6000PW <sub>XL</sub> /BioAssist G6PW <sub>XL</sub>	13	>1000	$4 \times 10^4-8 \times 10^6$	$5 \times 10^5-5 \times 10^7$	$<2 \times 10^8$
G6000PW	17, 25				
GMPW <sub>XL</sub>	13	<100-1000	$500-8 \times 10^6$	$<5 \times 10^7$	$<2 \times 10^8$
GMPW	17				
G-Oligo-PW	6	125	up to 3,000	—	<3,000
G-DNA-PW	10	>1000	$4 \times 10^4-8 \times 10^6$	$<5 \times 10^7$	$<2 \times 10^8$

Column: TSK-GEL PW columns, 7.5 mm ID x 60 cm L; TSKgel PW<sub>XL</sub>, G-Oligo-PW & G-DNA-PW, 7.8 mm ID x 30 cm L

Elution: Polyethylene glycols and oxides: distilled water; dextrans and proteins: 0.2 M phosphate buffer, pH 6.8

Flow Rate: 1.0 mL/min

Note: \*Larger particle sizes of each group are for 21.5 mm ID x 60 cm L semi-preparative and 55 mm or 108 mm ID x 60 cm L preparative columns.

\*\*Maximum separation range determined from estimated exclusion limits.

### Polyethylene glycol and oxide calibration curves on TSK-GEL PW and TSK-GEL PW<sub>XL</sub> columns





Figure 23

### Protein calibration curves on TSKgel PW<sub>XL</sub> columns




Column: 1. G3000PW<sub>XL</sub>  
2. G4000PW<sub>XL</sub>  
3. G5000PW<sub>XL</sub>  
4. G6000PW<sub>XL</sub>  
5. GMPW<sub>XL</sub>

Sample: a. thyroglobulin (660,000 Da)  
b.  $\gamma$ -globulin (150,000 Da)  
c. albumin (67,000 Da)  
d. ovalbumin (43,000 Da)  
e.  $\beta$ -lactoglobulin (36,000 Da)  
f. myoglobin (16,900 Da)  
g. cytochrome C (12,400 Da)

Elution: 0.2 M phosphate buffer (pH 6.8)  
Flow Rate: 1.0 mL/min  
Detection: UV @ 280 nm

### Faster analysis and higher resolution with TSK-GEL PW<sub>XL</sub> columns



Column: A. TSKgel G2500PW, two 10  $\mu$ m, 7.5 mm ID x 60 cm L columns in series  
B. TSKgel G2500PW<sub>XL</sub>, two 6  $\mu$ m, 7.8 mm ID x 30 cm L columns in series  
C. TSKgel G4000PW, 17  $\mu$ m, 7.5 mm ID x 60 cm L  
D. TSKgel G4000PW<sub>XL</sub>, 10  $\mu$ m, 7.8 mm ID x 30 cm L

Sample: A. & B.: polyethylene glycol 200  
C. & D.: polyethylene oxide standards: SE-150, SE-15 and SE-2 in 100  $\mu$ L

Elution: A. & B.: distilled water; C. & D.: 0.1 M NaCl  
Flow Rate: 1.0 mL/min  
Temperature: A. & B.: 25°C; C. & D.: 50°C  
Detection: RI

### TSKgel Oligo PW

The specialty column TSKgel Oligo PW is designed for high resolution separations of nonionic and cationic oligomers. Because of the presence of residual cationic groups, this column is not recommended for separating anionic materials. The polyethylene glycol and polyethylene oxide calibration curves for TSKgel G-Oligo-PW (not shown) are identical to TSKgel G2500PW<sub>XL</sub> shown in Figure 22.

### TSKgel G-DNA-PW

The TSKgel G-DNA-PW column is dedicated to the separation of large polynucleotides, such as DNA and RNA fragments of 500 to 5,000 base pairs. Figure 25 shows a calibration curve for double-stranded DNA on the G-DNA-PW column. The exclusion limits for double-stranded DNA fragments are lower than those for rRNAs, indicating that double-stranded DNA fragments have a larger effective molecular weight in solution than rRNAs of the same molecular weight. The packing has very large pores (>1000 Å) and a small particle size (10 µm).


For the separation of large DNA fragments greater than 1,000 base pairs, a four-column system is typically required. Baseline resolution of DNA fragments up to 7,000 base pairs can be achieved, provided there is a two-fold difference in the chain length of the fragments. Figure 26 shows the elution of double-stranded DNA fragments, obtained from pBR322 DNA cleaved by both *Eco*RI and *Bst*NI, on four TSKgel G-DNA-PW columns in series. The eluted peaks were collected and subjected to polyacrylamide gel electrophoresis, which showed almost complete separation of the 1060, 1857, and 4362 base pair fragments. Although lower flow rates typically yield better separations of most fragments, the resolution of the 1857 and 4362 base pair fragments was slightly greater at the higher flow rate, as shown in Figure 26B.

### TSKgel GMPW and TSKgel GMPW<sub>XL</sub>

When the molecular weight range of the sample is broad or unknown, Tosoh Bioscience offers two mixed-bed columns, TSKgel GMPW and TSKgel GMPW<sub>XL</sub>, for analysis. The TSKgel GMPW column and its high resolution counterpart, TSKgel GMPW<sub>XL</sub>, are packed with the G2500, G3000 and G6000 PW or corresponding PW<sub>XL</sub> resins. They offer a broad molecular weight separation range. As shown in Figure 22, the calibration curve for polyethylene glycols and oxides on these mixed-bed columns

Figure 25 covers the range of 100-1,000,000 Da.


### Calibration curve for double-stranded DNA



Column: TSKgel G-DNA-PW, four 10 µm, 7.8 mm ID x 30 cm L columns in series  
Sample: *Eco*RI and *Bst*NI-cleaved pBR322 DNA, void volume determined with λ-DNA  
Elution: 0.3 M NaCl in 0.1 M Tris-HCl, pH 7.5, plus 1 mM EDTA  
Flow Rate: 0.15 mL/min  
Detection: UV @ 260 nm

Figure 26

### Separation of large DNA fragments on a TSKgel G-DNA-PW column



Column: TSKgel G-DNA-PW, four 10 µm, 7.8 mm ID x 30 cm L columns in series  
Sample: 60 µL of *Eco*RI and *Bst*NI - cleaved pBR322 DNA, base pairs: a. 4362, b. 1857, c. 1060 & 928, d. 383, e. 121, f. 13  
Elution: 0.3 M NaCl in 0.1 M Tris-HCl, pH 7.5, plus 1mM EDTA  
Flow Rate: A. 0.15 mL/min, B. 0.5 mL/min  
Detection: UV @ 260 nm

## Applications

### Nucleic acids

Desalting of nucleosides can be accomplished using the TSKgel G2500PW<sub>XL</sub> as depicted in *Figure 27*. In this case secondary interaction effects induce adenosine and uridine to elute after the void volume in the unbuffered water mobile phase.

### Polysaccharides


TSK-GEL PW columns are recommended for polysaccharide analysis, due to their ability to separate a wide molecular weight distribution. Nonionic polysaccharides are the least complicated molecules to analyze by size exclusion chromatography because they seldom exhibit secondary interactions with solid support. Two TSK-GEL PW columns (TSKgel G5000PW and TSKgel G3000PW) in series are effective for the characterization of clinical dextran.

Cationic samples can be adsorbed on the resin by electrostatic interaction. If the polymer is strongly cationic, a fairly high salt concentration is required to prevent ionic interactions. *Figure 28* demonstrates the effect of increasing sodium nitrate concentration on peak shapes for a cationic polymer, DEAE-dextran. A mobile phase of 0.5 M acetic acid with 0.3 M Na<sub>2</sub>SO<sub>4</sub> can also be used.

An effective separation of the anionic hydrophilic glucosaminoglycan, hyaluronic acid, is shown in *Figure 29* on a TSKgel G6000PW and TSKgel G4000PW column in series with 0.2 M sodium chloride.

Figure 28


### Use of sodium nitrate to elute cationic polymers



Column: TSKgel GMPW, two 17  $\mu$ m, 7.5 mm ID x 60 cm L columns in series  
Sample: 0.5 mL of 0.05-0.1% of the DEAE-dextran cationic polymer  
Elution: 0.1 M, 0.2 M, 0.4 M or 0.8 M NaNO<sub>3</sub> in water  
Flow Rate: 0.5 mL/min  
Detection: RI

Figure 27


### Desalting of nucleosides



Column: TSKgel G2500PW<sub>XL</sub>, 7.8 mm ID x 30 cm L  
Sample: A. 0.5 NaCl; B. uridine; C. adenosine  
Eluent: distilled water  
Flow Rate: 1.0 mL/min  
Detection: UV @ 260 nm

Figure 29

### Separation of hyaluronic acid




Column: TSKgel G6000PW + G4000PW, two 7.5 mm ID x 60 cm L columns in series  
Sample: hyaluronic acid  
Elution: 0.2 M NaCl  
Flow Rate: 0.9 mL/min  
Temperature: 40°C

## Oligomers

The TSKgel G-Oligo-PW column is designed for high resolution separations of nonionic and cationic oligomers. *Figure 30* demonstrates excellent resolution of chito-oligosaccharides obtained by using the smaller, 6  $\mu\text{m}$  particle size packing in TSKgel G-Oligo-PW columns as compared with the resolution obtained with a TSKgel G2000PW column. The pore sizes in both TSKgel G-Oligo-PW and TSKgel G2000PW columns are about 125  $\text{\AA}$  and both resins bear approximately 0.2  $\mu\text{eq/mL}$  of cationic groups. Because of the presence of cationic groups, neither column is recommended for separating anionic materials. However, for nonionic oligomers, TSKgel G-Oligo-PW columns provide higher resolution than TSKgel G2500PW<sub>XL</sub> columns. The calibration curve for hydrolyzed  $\beta$ -cyclodextrin on two TSKgel G-Oligo-PW columns in series is shown in *Figure 31*.

Figure 30






## Polymers

Sodium polyacrylate, an anionic polymer, is effectively separated on two TSKgel GMPW columns in *Figure 32*. The addition of 0.01 M  $\text{NaNO}_3$  results in normal elution and peak shape overcoming the ionic repulsion between the anionic sample and the resin.

Polyethyleneimine, a cationic polymer, is separated in *Figure 33* on two TSKgel GMPW<sub>XL</sub> columns with a 0.5 M acetic acid and 0.5 M sodium acetate as mobile phase.

Figure 31





## Optimizing GFC with TSK-GEL PW and PW<sub>XL</sub>

### Selecting mobile phase buffers

In an ideal SEC separation, the mechanism is pure sieving, with no chemical interaction between the column matrix and the sample molecules. In practice, however, a small number of weakly charged groups on the surface of all TSK-GEL PW type packings can cause changes in elution order from that of an ideal system. Fortunately, the eluent composition can vary greatly with TSK-GEL PW columns, to be compatible with a wide range of neutral, polar, anionic and cationic samples. *Table IX* lists appropriate eluents for GFC of all polymer types on TSK-GEL PW type columns.

For some nonionic, nonpolar polymers, such as polyethylene glycols, normal chromatograms can be obtained by using distilled water. More polar ionic polymers may exhibit abnormal peak shapes or minor peaks near the void volume when eluted with distilled water, due to ionic interactions between the sample and residual charged groups on the resin surface. To eliminate ionic interactions, a neutral salt such as sodium nitrate or sodium sulfate is added to the aqueous eluent. Generally, a salt concentration of 0.1 M to 0.5 M is sufficient to overcome undesirable ionic interactions.

### Hydrophobic samples

TSK-GEL PW-type resins are more hydrophobic than polysaccharide gels such as cross-linked dextran. The extent of hydrophobic interaction increases as the salt concentration of the eluent increases, but it can be reduced by the addition of an organic solvent modifier such as acetonitrile. Water-soluble organic solvents are frequently used as modifiers to suppress hydrophobic interactions between the sample and the resin surface.

Modifiers are also used for optimizing the elution of both charged and neutral hydrophobic polymers. Typical examples for a variety of sample types are given in *Table IX*. All TSK-GEL PW-type column packings are compatible with 20% aqueous solutions of methanol, ethanol, propanol, acetonitrile, dimethyl formamide, dimethyl sulfoxide, formic acid, and acetic acid. In addition these columns can be run in 50% aqueous acetone.

**Recommended eluents for GFC of water-soluble polymers on TSK-GEL PW type columns**

<i>Type of polymer</i>	<i>Typical sample</i>	<i>Suitable eluent</i>
Nonionic hydrophilic	polyethylene glycol soluble starch, methyl cellulose, pullulan dextran, hydroxyethyl cellulose, polyvinyl alcohol, polyacrylamide	distilled water 0.01 N NaOH 20% DMSO Buffer or salt solution (e.g., 0.1–0.5 M NaNO <sub>3</sub> )
Nonionic hydrophobic	polyvinylpyrrolidone	Buffer or salt solution with organic solvent (e.g., 20% CH <sub>3</sub> CN in 0.1 M NaNO <sub>3</sub> )
Anionic hydrophilic	sodium chondroitin sulfate, sodium alginate, carboxymethyl cellulose, sodium polyacrylate, sodium hyaluronate	Buffer or salt solution (e.g., 0.1 M NaNO <sub>3</sub> )
Anionic hydrophobic	sulfonated lignin sodium salt, sodium polystyrenesulfonate	Buffer or salt solution with organic solvent (e.g., 20% CH <sub>3</sub> CN in 0.1 M NaNO <sub>3</sub> )
Cationic hydrophilic	glycol chitosan, DEAE-dextran, poly(ethyleneimine), poly(trimethylaminoethyl methacrylate) iodide salt	0.5 M acetic acid with 0.3 M Na <sub>2</sub> SO <sub>4</sub> , or 0.8 M NaNO <sub>3</sub>
Cationic hydrophobic	poly(4-vinylbenzyltrimethylammonium chloride), poly(N-methyl-2-vinylpyridinium) iodide salt	0.5 M acetic acid with 0.3 M Na <sub>2</sub> SO <sub>4</sub>
Amphoteric hydrophilic	peptides, proteins, poly-and oligosaccharides, DNA, RNA	Buffer or salt solution (e.g., 0.1 M NaNO <sub>3</sub> )
Amphoteric hydrophobic	blue dextran, collagen, gelatin, hydrophobic proteins hydrophobic peptides	Buffer or salt solution with organic solvent (e.g., 20% CH <sub>3</sub> CN in 0.1 M NaNO <sub>3</sub> or 35–45% CH <sub>3</sub> CN in 0.1% TFA)

## Ordering Information

### Analytical and preparative TSK-GEL Size Exclusion polymer-based column products: typical properties

Part #	Description	ID (mm)	Length (cm)	Particle Size ( $\mu\text{m}$ )	Min. Number Theoretical Plates	Flow Rate (mL/min)		Maximum Pressure Drop (kg/cm <sup>2</sup> )
						Range	Max.	
<b>Stainless steel columns</b>								
08031	G-Oligo-PW, 125 Å	7.8	30	6	14,000	0.5 – 0.8	1.0	40
08032	G-DNA-PW, >1.000 Å	7.8	30	10	10,000	0.2 – 0.5	0.6	20
08020	G2500PW <sub>XL</sub> , <200 Å	7.8	30	6	14,000	0.5 – 0.8	1.0	40
08021	G3000PW <sub>XL</sub> , 200 Å	7.8	30	6	14,000	0.5 – 0.8	1.0	40
08022	G4000PW <sub>XL</sub> , 500 Å	7.8	30	10	10,000	0.3 – 0.8	1.0	20
08023	G5000PW <sub>XL</sub> , 1.000 Å	7.8	30	10	10,000	0.3 – 0.8	1.0	20
08024	G6000PW <sub>XL</sub> , >1.000 Å	7.8	30	13	7,000	0.3 – 0.8	1.0	20
08025	GMPW <sub>XL</sub> , 100-1.000 Å	7.8	30	13	7,000	0.3 – 0.8	1.0	20
05760	G1000PW, <100 Å	7.5	30	10	5,000	0.5 – 1.0	1.2	20
05761	G2000PW, 125 Å	7.5	30	10	5,000	0.5 – 1.0	1.2	20
08028	G2500PW, <200 Å	7.5	30	10	5,000	0.5 – 1.0	1.2	20
05762	G3000PW, 200 Å	7.5	30	10	5,000	0.5 – 1.0	1.2	20
05763	G4000PW, 500 Å	7.5	30	17	3,000	0.5 – 1.0	1.2	10
05764	G5000PW, 1.000 Å	7.5	30	17	3,000	0.5 – 1.0	1.2	10
05765	G6000PW, >1.000 Å	7.5	30	17	3,000	0.5 – 1.0	1.2	10
08026	GMPW, 100-1.000 Å	7.5	30	17	3,000	0.5 – 1.0	1.2	10
05105	G2000PW, 125 Å	7.5	60	10	10,000	0.5 – 1.0	1.2	40
08029	G2500PW, <200 Å	7.5	60	10	10,000	0.5 – 1.0	1.2	40
05106	G3000PW, 200 Å	7.5	60	10	10,000	0.5 – 1.0	1.2	40
05107	G4000PW, 500 Å	7.5	60	17	6,000	0.5 – 1.0	1.2	20
05108	G5000PW, 1.000 Å	7.5	60	17	6,000	0.5 – 1.0	1.2	20
05109	G6000PW, >1.000 Å	7.5	60	17	6,000	0.5 – 1.0	1.2	20
08027	GMPW, 100-1.000 Å	7.5	60	17	6,000	0.5 – 1.0	1.2	20
05150	G2000PW, 125 Å	21.5	60	17	10,000	1.0 – 6.0	8.0	20
08030	G2500PW, <200 Å	21.5	60	17	10,000	1.0 – 6.0	8.0	20
05151	G3000PW, 200 Å	21.5	60	17	10,000	1.0 – 6.0	8.0	20
05152	G4000PW, 500 Å	21.5	60	22	6,000	1.0 – 6.0	8.0	20
05153	G5000PW, 1.000 Å	21.5	60	22	6,000	1.0 – 6.0	8.0	20
05154	G6000PW, >1.000 Å	21.5	60	25	6,000	1.0 – 6.0	8.0	20
07926	G3000PW, 200 Å	55.0	60	20	4,500	15.0 – 25.0	30.0	15
07927	G5000PW, 1.000 Å	55.0	60	20	4,500	15.0 – 25.0	30.0	15
<b>Guard columns</b>								
08034	Oligo Guard column	6.0	4.0	12.0	For G-Oligo-PW			
08033	PW <sub>XL</sub> Guard column	6.0	4.0	12.0	For 7.8 mm ID PW <sub>XL</sub> and G-DNA-PW (contains 3000PW packing)			
06763	PW-L Guard column	7.5	7.5	12.0	For 7.5 mm ID G1000PW and G2000PW (contains 2000PW packing)			
06762	PW-H Guard column	7.5	7.5	12.0	For 7.5 mm ID G2500PW - G6000PW + GMPW (contains 3000PW packing)			
06757	PW-L Guard column	21.5	7.5	17.0	For 21.5 mm ID G2000PW			
06758	PW-H Guard column	21.5	7.5	17.0	For 21.5 mm ID G2500PW through G6000PW			
07924	PW Guard column	45.0	5.0	20.0	For 55 mm ID G3000PW + G5000PW			
<b>Bulk packing</b>								
08035	PW <sub>XL</sub> Top-Off, 1g wet resin			10.0	For all PW <sub>XL</sub> and G-DNA-PW			


# Polymer based TSK-GEL Alpha and SuperAW columns for Gel Filtration and Gel Permeation Chromatography of water soluble and organic soluble polymers

## Highlights

- A unique hydrophilic, polyvinyl resin.
- Exhibits strong mechanical stability and minimal swelling characteristics.
- A wide range of solvent compatibility, from 100% water to 100% non-polar organic solvents (Figure 34 and Figure 35).
- The reduced particle size and shorter column length of TSK-GEL SuperAW columns provides equivalent resolution in ½ the time for high throughput applications.
- Unlike polystyrene-divinylbenzene (PS-DVB) resins that may adsorb polymers due to hydrophobic interaction, both the TSK-GEL Alpha and SuperAW columns allow for the separation of polymers soluble in methanol.
- Provide accurate molecular weight (MW) determination of samples in dimethyl formamide and exhibit normal retention of polystyrene polymers.
- System peaks from salts in the eluent elute away from the oligomer of interest (Figure 36), providing accurate MW determinations.

Figure 35


## Solvent Compatibility of TSK-GEL SuperAW Series



Column: TSK-GEL SuperAW Series (6.0 mm ID x 15 cm L)  
 Eluent: Water  
 Flow rate: 0.6 mL/min  
 Temperature: 25°C  
 Detection: Refractive index detector  
 Sample: Ethylene glycol  
 Inj. volume: 5 µL (2.5 g/L)

Figure 34

## Solvent compatibility for TSKgel Alpha-3000 for organic solvents




Conditions of solvent change  
 Flow Rate: 1.0 mL/min  
 Temperature: 25°C  
 Time for purge: 8h

Conditions for TP measurement  
 Sample: ethylene glycol  
 Flow Rate: 1.0 mL/min  
 Temperature: 25°C  
 Detection: RI

Figure 36

## Effect of salt peak on chromatogram



Column: TSKgel Alpha-M, two 7.8 mm ID x 30 cm L columns in series  
 Sample: Polyimide  
 Elution: 60 mM LiBr in DMF  
 Flow Rate: 1.0 mL/min  
 Temperature: 40°C  
 Detection: RI

## Column Selection

The TSK-GEL Alpha Series consists of six columns with three particle sizes: 7, 10, and 13  $\mu\text{m}$ . These columns span a wide MW separation range from 100 to more than  $1 \times 10^6$  Da when using polyethylene oxide (PEO) as a MW standard. Exclusion limits for the TSK-GEL Alpha columns for PEO in water and polystyrene (PS) in dimethylformamide (DMF) are shown in Table X. Calibration curves for the TSK-GEL Alpha Series columns are shown in Figure 37 for polyethylene oxide, polyethylene glycol and polystyrene standards.

The TSK-GEL SuperAW series contains a similar chemistry as the Alpha series but offers the benefit of smaller particle sizes (4  $\mu\text{m}$  to 9  $\mu\text{m}$ ) and smaller column dimensions. Reductions in analysis time and mobile phase consumption result making these columns ideal for high throughput applications. TSK-GEL SuperAW columns are offered in 5 discrete exclusion ranges and 1 mixed bed, polymers up accommodate PEO polymers standards over several million Dalton. (see Figure 38)

Table X

## Exclusion limits for TSK-GEL Alpha Series and SuperAW Series columns


Column	Particle Size ( $\mu\text{m}$ )	Exclusion limit (Da) for various standards and eluents		
		PEO <sup>a</sup> /H <sub>2</sub> O	PS <sup>b</sup> /10mM LiBr in DMF	PEG <sup>c</sup> /10mM LiBr in MeOH
Alpha-2500	7	$5 \times 10^3$	$1 \times 10^4$	$1 \times 10^4$
Alpha-3000	7	$9 \times 10^4$	$1 \times 10^5$	$6 \times 10^4$
Alpha-4000	10	$4 \times 10^5$	$1 \times 10^6$	$3 \times 10^6$
Alpha-5000	10	$1 \times 10^6$	$7 \times 10^6$	N.D.
Alpha-6000	13	$>1 \times 10^7$	$>1 \times 10^7$	N.D.
Alpha-M	13	$>1 \times 10^7$	$>1 \times 10^7$	N.D.
SuperAW2500	4	$5 \times 10^3$	$8 \times 10^3$	$1 \times 10^4$
SuperAW3000	4	$9 \times 10^4$	$8 \times 10^4$	$1 \times 10^5$
SuperAW4000	6	$1 \times 10^6$	$6 \times 10^5$	$6 \times 10^6$
SuperAW5000	7	$1 \times 10^{6*}$	N.D.	N.D.
SuperAW6000	9	$1 \times 10^{7*}$	N.D.	N.D.
SuperAWM-H	9	$1 \times 10^{7*}$	N.D.	N.D.

N.D.= not determined a Polyethylene oxide b Polystyrene Divinyl Benzene c Polyethylene glycol

\* Exclusion limit for SuperAW5000, SuperAW6000, and SuperAWM-H are estimated, respectively

Figure 37


**Polyethylene oxide (PEO), polyethylene glycol (PEG) and polystyrene (PS) calibration curves for TSK-GEL Alpha columns**



Column: Alpha Series, 7.8 mm ID x 30 cm L  
 Eluent: A. H<sub>2</sub>O; B. 10 mM LiBr in Methanol; C. 10 mM LiBr in DMF  
 Flow Rate: 1.0 mL/min  
 Temperature: A. 25°C; B. 25°C; C. 40°C  
 Detection: RI

Figure 38

**Polyethylene oxide (PEO), polyethylene glycol (PEG) and ethylene glycol (EG) calibration curves for TSK-GEL SuperAW Series**



Column: TSK-GEL SuperAW Series (6.0 mm ID x 15 cm L)  
 Eluent: A. Water; B. MeOH containing 10mM LiBr; C. DMF containing 10 mM LiBr  
 Flow rate: 0.6 mL/min  
 Temperature: 25°C  
 Detection: Refractive index detector  
 Samples: Standard polyethylene oxide, polyethylene glycol, ethylene glycol

## Applications

The versatility of using TSK-GEL Alpha columns with various polar solvents is illustrated in *Figure 39* for the analysis of cellulose derivatives. A TSKgel Alpha-M column was used to separate ethylcellulose with the polar solvent DMF and ethylhydroxyethyl cellulose with methanol.

The separation of polyvinylalcohol with different degrees of saponification is shown in *Figure 40*. This separation was performed with a TSKgel Alpha-5000 and a TSKgel Alpha-3000 column in series using a hexafluoroisopropanol mobile phase.


Separation of the components in cleansing gel are shown in *Figure 41*.

A separation of sodium dodecyl sulfate is shown in *Figure 43* on three Alpha columns in series: the TSKgel Alpha-4000, 3000 and 2500. This separation is performed in 50 mM lithium bromide in DMF.

The separation of dextran in *Figure 42* indicates the decreased run time that results when using TSKgel SuperAW2500 in comparison to a conventional 6  $\mu$ m TSKgel G2500 PW<sub>XL</sub>.

*Figure 40*


### Polyvinylalcohol characterization using TSKgel Alpha-5000 and Alpha-3000 columns in series



Column: TSKgel Alpha-5000 and Alpha-3000, 7.8 mm ID x 30 cm L in series  
Sample: degree of saponification of polyvinyl alcohol: A. 75%; B. 88%; C. 100%  
Eluent: hexafluoroisopropanol (HFIP)  
Flow Rate: 0.5 mL/min  
Temperature: 40°C  
Detection: RI

*Figure 39*


### TSKgel Alpha-M separation of cellulose derivatives



Column: TSKgel Alpha-M, 7.8 mm ID x 30 cm L  
Sample: A. 50  $\mu$ L ethylcellulose, 0.1%; B. 50  $\mu$ L ethylhydroxyethylcellulose, 0.1%  
Elution: A. 10 mM LiBr in DMF; B. 10 mM LiBr in methanol  
Flow Rate: 0.5 mL/min  
Temperature: 40°C  
Detection: RI


*Figure 41*

### TSKgel Alpha-2500 separation of cleansing gel




Column: TSKgel Alpha-2500, 7.8 mm ID x 30 cm L  
Sample: cleansing gel: 1. glyceryl monostearate POE (20); 2. glyceryl tri 2-ethylhexanoate; 3. glycerin; 4. sorbitol  
Eluent: methanol  
Flow Rate: 0.5 mL/min  
Temperature: 40°C  
Detection: RI

### Comparison of TSKgel SuperAW2500 and TSKgel G2500PW<sub>XL</sub>



Column: TSKgel SuperAW2500 (6.0 mm ID x 15 cm L)  
 TSKgel G2500PW<sub>XL</sub> (7.8 mm ID x 30 cm L)  
 Eluent: Water  
 Flow rate: 0.6 mL/min (TSKgel SuperAW2500)  
 1.0 mL/min (TSKgel G2500PW<sub>XL</sub>)  
 Temperature: 25°C  
 Detection: Refractive index detector  
 Sample: Dextran T-40 hydrolysate

### Separation of sodium dodecyl sulfate



Columns: (1) TSKgel Alpha-4000, 7.8 mm ID x 30 cm L  
 (2) TSKgel Alpha-3000, 7.8 mm ID x 30 cm L  
 (3) TSKgel Alpha-2500, 7.8 mm ID x 30 cm L  
 in series  
 Sample: 200  $\mu$ L of 1.7% sodium dodecyl sulfate  
 Elution: 50 mM lithium bromide in DMF  
 Flow Rate: 1.0 mL/min  
 Detection: RI

## Ordering Information

### TSK-GEL Alpha Series and TSK-GEL SuperAW Series

Part #	Description	ID (mm)	Length (cm)	Particle Size ( $\mu$ m)	Min. Number Theoretical Plates	Flow Rate (mL/min)		Maximum Pressure Drop (kg/cm <sup>2</sup> )
						Range	Max.	
<b>Stainless steel columns</b>								
18339	Alpha-2500	7.8	30	7	16,000	0.5 – 0.8	1.0	40
18340	Alpha-3000	7.8	30	7	16,000	0.5 – 0.8	1.0	40
18341	Alpha-4000	7.8	30	10	10,000	0.3 – 0.6	1.0	30
18342	Alpha-5000	7.8	30	10	10,000	0.3 – 0.6	1.0	30
18343	Alpha-6000	7.8	30	13	7,000	0.3 – 0.6	1.0	20
18344	Alpha-M (mixed bed)	7.8	30	13	7,000	0.3 – 0.6	1.0	20
<b>Guard columns</b>								
18345	Alpha Guard column	6	4	13				
<b>Stainless steel columns</b>								
19315	SuperAW2500 <b>-NEW-</b>	6.0	15	4	16,000	0.3 – 0.6	0.6	60
19316	SuperAW3000 <b>-NEW-</b>	6.0	15	4	16,000	0.3 – 0.6	0.6	60
19317	SuperAW4000 <b>-NEW-</b>	6.0	15	6	10,000	0.3 – 0.6	0.6	40
19318	SuperAW5000 <b>-NEW-</b>	6.0	15	7	10,000	0.3 – 0.6	0.6	30
19319	SuperAW6000 <b>-NEW-</b>	6.0	15	9	7,000	0.3 – 0.6	0.6	20
19320	SuperAWM-H <b>-NEW-</b>	6.0	15	9	7,000	0.3 – 0.6	0.6	20
<b>Guard columns</b>								
19321	SuperAW-L Guard Column <b>-NEW-</b>		4.6	3.5	7			
19322	SuperAW-H Guard Column <b>-NEW-</b>		4.6	3.5	23			

# Polymer-based TSK-GEL H<sub>XL</sub>, TSK-GEL H<sub>HR</sub>, TSK-GEL SuperH, TSK-GEL SuperHZ columns for Gel Permeation Chromatography of organic soluble polymers

## Highlights

- **-NEW-** Reduce solvent consumption during high throughput analysis with narrow bore SuperH (6.0 mm ID) and SuperHZ (4.6 and 6.0 mm ID). (Figure 45)
- Porous, highly cross-linked, spherical polystyrene divinylbenzene (PS-DVB) resin.
- Four different TSK-GEL H-type columns are available. Each of these are packed with different particle sizes (Table XI).
- H-type packings are available in eight pore sizes.
- Expanded molecular weight ranges with exclusion limits from 1,000 Da to an estimated  $4 \times 10^6$  Da.
- Minimal shrinking and swelling of the column bed.
- Chemically and thermally stable.
- Novel multi-pore distribution in the TSKgel MultiporeH<sub>XL</sub>-M column provides linear calibration curves over a wider MW range.
- Mixed bed columns with optimized particle and pore sizes to prevent polymer sheering.

## Column Selection

TSK-GEL H-type packings are available in eight pore sizes and span four different chemistries. The exclusion limits, pore sizes and application areas for the four chemistries of TSK-GEL H-type columns are listed in Table XI. For polymer samples with a broad molecular weight range, packings of several pore sizes are provided in the mixed-bed columns: TSK-GEL SuperHZM series, TSK-GEL SuperHM series, TSKgel GMH<sub>XL</sub>, TSKgel GMH<sub>HR</sub> and selected high temperature versions provide linear calibration curves up to several million Daltons. Additionally, TSK-GEL SuperHZM and SuperHM series are offered in several particle sizes appropriately matched to the molecular weight of the polymer sample in order to prevent sheering during the analysis. The "Super" prefix refers to the efficiency of the column. The super series columns contain ultra efficient particles as low as 3 μm, housed in 15 cm length columns. The smaller particle allows for equivalent resolution to conventional H<sub>XL</sub> columns, with 50% less run time due to the shorter column length. The super series columns are an excellent choice for high throughput polymer analysis. Calibration curves for all four H-type columns are provided in Figures 46, 47, 48, 49. A comparison of the 3μm TSKgel SuperHZ columns to conventional larger particle TSKgel H<sub>XL</sub> is shown in Figure 44.

Table XI

### Comparison of TSK-GEL H-type resins

Series Type	SuperHZ	H <sub>XL</sub>	SuperH	H <sub>HR</sub>
Application focus	High-throughput polymer analysis with ultra low polymer adsorption, limited solvent compatibility range.	Conventional polymer analysis with ultra low polymer adsorption, limited solvent compatibility range.	High-throughput polymer analysis with expanded solvent compatibility.	Conventional polymer analysis with expanded solvent compatibility range.
Particle size	3, 5 and 10 μm, depending on pore size	5 and 9 μm, depending on pore size	3 and 5 μm, depending on pore size	5 μm
Theoretical plates <sup>1</sup>	16,000/15 cm column	16,000/30 cm column	16,000/15 cm column	16,000/30 cm column
Maximum temperature	G1000 - G4000 60°C G5000 - mixed 80°C	G1000 - G4000 60°C G5000 - mixed 80°C	140°C	140°C
Standard shipping solvent	THF	THF <sup>2</sup>	THF <sup>2</sup>	THF <sup>2</sup>
THF can be switched to	benzene, chloroform, toluene, xylene, dichloromethane <sup>3</sup> and dichloroethane <sup>3</sup>		see Table XII (page 61)	
Other shipping solvents available?	yes <sup>4</sup>		no	
Number of solvent substitutions	One time only		Several <sup>5</sup>	
Solvent exchange instructions	Linear gradient with a 2%/min rate of change with a flow rate <0.25 mL/min.	Linear gradient with a 2%/min rate of change with a flow rate <0.5 mL/min.	Linear gradient with a 2%/min rate of change according to flow rates listed in Table XIII	

1) Theoretical plates listed are based on smallest particle size listed

2) High-temperature columns (HT) are shipped with OCDB (Orthochlorodivinylbenzene) as standard shipping solvent.


3) Dichloromethane and Dichloroethane are not available for G1000 pore size columns.

4) See Table XII for available shipping solvents

5) After switching to a very polar solvent such as acetone, switching back to a nonpolar solvent is not recommended.

Figure 44


Comparison of resolution on TSKgel H<sub>HR</sub> and H<sub>XL</sub> columns



Column: A: TSKgel GMH<sub>HR</sub>-H, 7.8 mm ID x 30 cm L  
 B: TSKgel GMH<sub>XL</sub>, 7.8 mm ID x 30 cm L  
 Sample: polystyrene standards 5 µL  
 Elution: THF  
 Temperature: 25°C  
 Detection: UV @ 254 nm

Figure 45


Comparison between SuperHZ and GMH<sub>XL</sub>



Columns: (A) TSKgel SuperHZ (4000, 3000, 2500) (4.6 mm ID x 15cm L)  
 (B) TSKgel GMH<sub>XL</sub> (4000, 3000, 2500) (7.8 mm ID x 30 cm L)  
 Eluent: THF  
 Flow rate: (A) 0.35 mL/min  
 (B) 1.0 mL/min  
 Temperature: 40°C  
 Detection: RI  
 Sample: Phenolic resin  
 Inj. volume: (A) 5 µL, (B) 30 µL

Figure 46

Calibration curves for TSKgel H<sub>XL</sub> columns with polystyrene standards



Column size: 7.8 mm ID x 30 cm L  
 Sample: polystyrene standards  
 Eluent: THF  
 Flow Rate: 1.0 mL/min  
 Temperature: 25°C  
 Detection: UV @ 254 nm

Figure 47




Figure 48





Figure 49



**Solvent compatibility for TSK-GEL H-type columns**

The column shipping solvent limits the types of solvents that can be used with each column type.

Column type	Shipping solvent *	Can be replaced with
SuperHZ and H <sub>XL</sub> <sup>1</sup>	Tetrahydrofuran <sup>3,4</sup>	benzene, chloroform, toluene, xylene, dichloromethane, dichloroethane
	Acetone**	carbon tetrachloride <sup>5</sup> , <i>o</i> -dichlorobenzene, dimethylformamide, dodecane, dimethyl sulfoxide, dioxane, ethylacetate, FC-113, hexane, pyridine, hexafluoroisopropanol/chloroform, methyl ethyl ketone, quinoline, cyclohexane,
	Chloroform**	<i>m</i> -cresol in chloroform, up to 10% hexafluoroisopropanol/chloroform
	Dimethylformamide**	dimethyl sulfoxide, dioxane, tetrahydrofuran, toluene
	<i>o</i> -Dichlorobenzene**	1-chloronaphthalene, trichlorobenzene
SuperH and H <sub>HR</sub> <sup>2</sup>	Tetrahydrofuran <sup>3</sup>	acetone, decahydronaphthalene, ethanol, quinoline, benzene, <i>o</i> -dichlorobenzene, ethyl acetate, tetrafluoroethylene, carbon tetrachloride <sup>5</sup> , dichloromethane, n-hexane, tetrahydrofuran, chloroform, 1,4-dioxane, hexafluoroisopropanol, toluene, 1-chloronaphthalene, dimethylacetacetamide, methyl ethyl ketone, trichlorobenzene, <i>m</i> -cresol, dimethylformamide, -methylpyrrolidone, <i>o</i> -chlorophenol/chloroform, methyl sulfoxide, pyridine

**Notes:**

1 In case of SuperHZ and H<sub>XL</sub>, keep flow rate below <0.5 mL/min during solvent change. Solvent can be changed one way/one time only.

2 In case of SuperH and H<sub>HR</sub>, see table XIII for appropriate flow rates for solvent exchange. After switching to a very polar solvent, switching to a nonpolar is not recommended.

3 All TSK-GEL H<sub>XL</sub>, H<sub>HR</sub>, SuperHZ, SuperH and GMH analytical columns are shipped containing tetrahydrofuran (THF), except GMH-HT columns, which contain *o*-dichlorobenzene.

4 THF in G1000H<sub>XL</sub> columns cannot be replaced with dichloromethane or dichloroethane.

5 Prolonged exposure to carbon tetrachloride can corrode the stainless steel parts of a column and an HPLC system.

\* 100% methanol cannot be used with H-type columns; use this solvent with TSK-GEL SW-type columns or Alpha-type columns.

\*\* TSK-GEL H columns may be specially ordered with this shipping solvent.

Table XIII

**Recommended flow rates for solvent exchange in SuperH and H<sub>HR</sub> columns**


Solvent	Flow Rate (mL/min)	
	SuperH	H <sub>HR</sub>
	6.0 mm ID x 15 cm L	7.8 mm ID x 30 cm L
<i>n</i> -Hexane	0.5	0.9
methyl ethyl ketone	0.4	0.7
dichloromethane, ethyl acetate	0.35	0.6
toluene, chloroform	0.3	0.5
dimethylformamide	0.2	0.4
carbon tetrachloride, pyridine	0.15	0.3
dimethyl sulfoxide, dioxane, ethanol, N-methylpyrrolidone, <i>o</i> -dichlorobenzene	0.1	0.2
quinoline, hexafluoroisopropanol, 1-chloronaphthalene	0.05	0.1

The TSKgel MultiporeH<sub>XL</sub>-M column offers a unique packing material and new strategy for the precise analysis of polymers by Gel Permeation Chromatography (GPC). Until now, the GPC separation of a sample containing a wide range of molecular weight polymers was performed by one of two strategies. One strategy combines columns with different pore sizes of packing material in series. The other strategy employs a single column with a blend of different pore sizes packing materials, commonly referred to as a mixed bed. Mixed bed columns do not always provide linear calibration curves, which may result in broad or split peaks. With the introduction of the TSKgel MultiporeH<sub>XL</sub>-M column, a new strategy has been developed using a single column containing a novel polystyrene packing material with a multi-pore size distribution. *Figure 50* gives an illustration of the three strategies. A scanning electron microscope photograph of the TSKgel Multipore column in *Figure 51* shows several

different pore sizes with continuous distribution in every bead. This results in sharper peaks without inflections that may be observed using mixed-bed columns.

The pore size distributions of the TSKgel MultiporeH<sub>XL</sub>-M column and a mixed-bed column are shown in *Figure 52*. The mixed-bed column shows a sharp maximum for pores with a diameter of 0.08  $\mu\text{m}$ , though the overall pore size distribution ranges from 0.006 to 0.6  $\mu\text{m}$  in diameter. In the case of the TSKgel MultiporeH<sub>XL</sub>-M column, the pore size distribution exhibits a wider maximum range from 0.02 to 0.1  $\mu\text{m}$  in diameter. This difference in pore size distribution may explain the reason for the inflection phenomenon. A comparison of calibration curves for polystyrene standards on the TSKgel MultiporeH<sub>XL</sub>-M column, the TSKgel GMH<sub>HR</sub>-H and PLgel Mixed-C, are shown in *Figure 53*. Both the TSKgel GMH<sub>HR</sub>-H and PLgel Mixed-C columns are mixed-bed columns.

*Figure 50*



*Figure 51*  
SEM photograph of the TSKgel multipore packing



*Figure 52*





Figure 53

### Calibration curves for Multipore and mixed-bed columns



Column: TSKgel Multipore H<sub>XL</sub>-M, 7.8 mm ID x 30 cm L;  
 TSKgel GMH<sub>HR</sub>-H, 7.8 mm ID x 30 cm L;  
 PLgel Mixed-C, 7.5 mm ID x 30 cm L

Sample: polystyrene standards

Elution: THF

Flow Rate: 1.0 mL/min

Temperature: 40°C

Detection: UV @ 254 nm

\*PLgel Mixed-C is a product from Polymer Laboratories Ltd.

## Optimizing GPC with TSK-GEL H-type columns


### Sample loading

In general, column efficiency improves with decreasing sample concentration and decreasing sample volume. When particle size decreases, all other conditions constant, these trends are more pronounced. Similarly, reducing column volume by either reducing column length or column diameter also results in reduced sample volumes and concentrations at which efficiency starts to deteriorate.

Efficiency decreases rapidly as the injection volume increases from 100  $\mu$ L to 300  $\mu$ L for polystyrene with MW <110,000 Da as seen with the mixed bed TSKgel GMH<sub>XL</sub> in Figure 54. Sample concentrations from 0.001-0.5 mg can be loaded on the 7.5 or 7.8 mm ID x 30 cm L column for an analytical separation (Figure 55). The smaller, more efficient particles in the Super series columns are more sensitive to injection volumes and sample concentrations. Figure 56 shows that injection volumes of 1  $\mu$ L to 10  $\mu$ L are recommended for TSKgel SuperHZ2500. Recommended sample concentrations vary dependant on sample molecular weight. For polymers below  $1 \times 10^6$  Da concentrations can range up to 1 g/L. Larger polymers should not exceed 0.1 g/L to 0.5 g/L according to Figure 57.

Figure 54

### Dependence of the number of theoretical plates (N) on injection volume in the separation of polystyrene on TSKgel GMH<sub>XL</sub>



Column: TSKgel GMH<sub>XL</sub>, 7.8 mm ID x 30 cm L

Sample: polystyrene, F-1, 10, 80, 450 (0.125 mg/mL)

Elution: THF


Flow Rate: 1.0 mL/min

Detection: UV @ 254 nm

Temperature: 40°C

Figure 55

### Dependence of the number of theoretical plates (N) on sample concentration in the separation of polystyrene on TSKgel GMH<sub>XL</sub>



Column: TSKgel GMH<sub>XL</sub>, 7.8 mm ID x 30 cm L

Sample: polystyrene, F-1, 10, 80, 450 (0.125 mg/mL)

Elution: THF


Flow Rate: 1.0 mL/min

Detection: UV @ 254 nm

Temperature: 40°C

Figure 56


### Relationship between HETP and Sample injection volume in SuperHZ2500



Column: TSKgel SuperHZ2500  
 6.0 mm ID x 15 cm L, 4.6 mm ID x 15 cm L  
 Eluent: THF  
 Flow rate: 0.6 mL/min (6.0 mm ID x 15 cm L),  
 0.35 mL/min (4.6 mm ID x 15 cm L)  
 Temperature: 25°C  
 Sample: DCHP

Figure 57

### Effect of sample concentration on elution time with SuperHZ-M



Column: TSKgel SuperHZM-M (4.6 mm ID x 15 cm L x 2)  
 Eluent: THF  
 Flow rate: 0.35 mL/min  
 Temperature: 40°C  
 Detection: RI  
 Sample: Standard polystyrene  
 Injection vol: 10 µL

## Applications for TSK-GEL H-type columns

Many industrial applications detailing analyses of organic polymers have been developed with TSK-GEL H-type columns. *Table XV* gives examples of the recommended solvent for each application.

*Table XV*


Recommended solvent	Application
THF	polystyrene, epoxy resin, phenoxy resin, polycarbonate, polyisoprene, polyvinyl acetate, polyvinyl chloride, monoglycerides, fatty acids, poly (butadiene-as) {as resin}, polybutadiene, poly (methyl methacrylate), poly (styrene-butadiene), poly (styrene-acrylonitrile) {as resin}
N,N-Dimethylformamide (DMF) +/- 5mM LiBr	polyvinyl chloride, polyvinyl fluoride, urea resins, polyurethane, polystyrene, polyester, polyimido ether, polyimido ester, polyphenol (aqueous solution), polyacrylonitrile
<i>o</i> -Dichlorobenzene (ODCB)	polyethylene, polypropylene
Chloroform	polycarboxylic ether, acrylic resin, epoxy resin, polystyrene
<i>m</i> -Cresol/Chloroform	nylon, polyester, polyamide, poly (ethylene terephthalate)
Toluene	polybutadiene, polysiloxane

## Phthalate esters

*Figure 58* demonstrates the high resolving power of a TSKgel G1000H<sub>XL</sub> column for low molecular weight phthalate esters. Resolution was close to baseline, even though the molecular weights of the esters differed by less than 50 Da.

*Figure 58*

### High resolution of phthalate esters on TSKgel G1000H<sub>XL</sub>



Column: TSKgel G1000H<sub>XL</sub>, 7.8 mm ID x 30 cm L

Sample: 1. polystyrene (10,200 Da), 2. dioctylphthalate (391 Da), 3. dibutylphthalate (278 Da), 4. dipropylphthalate (250 Da), 5. diethylphthalate (222 Da), 6. dimethylphthalate (194 Da), 7. n-propylbenzene (120 Da), 8. ethylbenzene (116 Da), 9. toluene (92 Da), 10. benzene (78 Da)

Elution: THF


Flow Rate: 1.0 mL/min

Detection: UV @ 254 nm

## Phenol resin

The TSKgel GMH<sub>XL</sub>-L column has been designed to provide a complete profile for high molecular weight samples that contain low molecular weight additives. The calibration curve for this mixed-bed column is shallow in the low molecular weight range of oligomers. Sample adsorption is not observed. For example, the complete profile of a phenol resin, with high resolution of the low molecular weight components, is shown in *Figure 59*. Other applications for the TSKgel GMH<sub>XL</sub>-L column include analyses of paint materials, bond and adhesive components, and synthetic polymer additives.

*Figure 59*



## Polymethylmethacrylate

The effect of different pore size distributions in the mixed beds of TSKgel GMH<sub>HR</sub>-H and TSKgel GMH<sub>HR</sub>-M is illustrated in *Figure 60*. The TSKgel GMH<sub>HR</sub>-M produces sharper peaks in the  $8 \times 10^5$  to  $1 \times 10^6$  Da range.


*Figure 61* shows two TSKgel Super HZM columns in series for the separation of polymethyl methacrylate. A low dead volume HPLC system as well as small sample volume and concentration are required to achieve optimum efficiency using the 4.6 mm ID columns.

*Figure 60*

### Comparison of TSKgel GMH<sub>HR</sub>-H and -M columns with polymethylmethacrylate standards

A. TSKgel GMH<sub>HR</sub>-H


B. TSKgel GMH<sub>HR</sub>-M



Column: A. TSKgel GMH<sub>HR</sub>-H, 7.8 mm ID x 30 cm L;  
B. TSKgel GMH<sub>HR</sub>-M, 7.8 mm ID x 30 cm L  
Sample: polymethylmethacrylate: 1. 820,000 Da, 2. 67,000 Da, 3. 10,200 Da, 4. 1,950 Da  
Solvent: 5 mM sodium trifluoroacetate in hexafluoroisopropanol  
Flow Rate: 1.0 mL/min  
Detection: UV @ 220 nm  
Temperature: 40°C

*Figure 61*

### Chromatogram of polymethyl methacrylate




Column: TSKgel SuperHZM-M x 2 (15 cm L each)  
Eluent: THF  
Flow rate: 0.35 mL/min (4.6 mm ID)  
0.6 mL/min (6.0 mm ID)  
Temperature: 40°C  
Detection: RI  
Sample: Polymethyl methacrylate (1 g/L)  
Inj. volume: 5  $\mu$ L (4.6 mm ID)  
9  $\mu$ L (6.0 mm ID)

## Fatty acids

In *Figure 62*, two TSKgel G2000H<sub>XL</sub> columns in series separate a mixture of fatty acids ranging from C<sub>4</sub> to C<sub>30</sub>.


*Figure 62*



## Acrylic polymer

*Figure 63* shows the separation of an acrylic polymer on the TSKgel MultiporeH<sub>XL</sub>-M column compared with two commercially available mixed-bed columns. The arrows illustrate the inflections seen in the chromatographic separations on mixed-bed columns and the improvement achieved when using the TSKgel MultiporeH<sub>XL</sub>-M column. For faster run times, the TSKgel SuperH<sub>ZM</sub> can perform the separation in under 10 minutes (*Figure 65*).


*Figure 63*




## Epoxy resin

*Figure 64* compares the separation of epoxy resin on the TSKgel MultiporeH<sub>XL</sub>-M and TSK-GEL H<sub>XL</sub> columns. The multi-pore size distribution of the TSKgel MultiporeH<sub>XL</sub>-M column gives a broad range of separation and linearity. Four TSK-GEL H<sub>XL</sub> columns with different pore sizes are required to achieve the same broad range of separation. Note the appearance of the inflections in the chromatogram when using the TSK-GEL H<sub>XL</sub> columns.

*Figure 64*



*Figure 65*



## Ordering Information

### Analytical and preparative TSK-GEL GPC column products: typical properties

Part #	Description	ID (mm)	Length (cm)	Particle Size ( $\mu\text{m}$ )	Min. Number Theoretical Plates	Flow Rate (mL/min)		Maximum Pressure Drop (kg/cm <sup>2</sup> )
						Range	Max.	
<b>Stainless steel columns</b>								
17352	G1000H <sub>HR</sub>	7.8	30	5	16,000	0.5 – 1.0	2.0	50
17353	G2000H <sub>HR</sub>	7.8	30	5	16,000	0.5 – 1.0	1.0	50
17354	G2500H <sub>HR</sub>	7.8	30	5	16,000	0.5 – 1.0	1.0	50
17355	G3000H <sub>HR</sub>	7.8	30	5	16,000	0.5 – 1.0	1.0	50
17356	G4000H <sub>HR</sub>	7.8	30	5	16,000	0.5 – 1.0	1.0	50
17357	G5000H <sub>HR</sub>	7.8	30	5	16,000	0.5 – 1.0	1.0	50
17358	G6000H <sub>HR</sub>	7.8	30	5	16,000	0.5 – 1.0	1.0	50
17359	G7000H <sub>HR</sub>	7.8	30	5	16,000	0.5 – 1.0	1.0	50
17362	GMH <sub>HR</sub> -L mixed-bed	7.8	30	5	16,000	0.5 – 1.0	1.0	50
18055	GMH <sub>HR</sub> -N mixed-bed	7.8	30	5	16,000	0.5 – 1.0	1.0	50
17392	GMH <sub>HR</sub> -M mixed-bed	7.8	30	5	16,000	0.5 – 1.0	1.0	50
17360	GMH <sub>HR</sub> -H mixed-bed	7.8	30	5	16,000	0.5 – 1.0	1.0	50
16131	G1000H <sub>XL</sub>	7.8	30	5	16,000	0.5 – 1.0	1.0	50
16134	G2000H <sub>XL</sub>	7.8	30	5	16,000	0.5 – 1.0	1.2	50
16135	G2500H <sub>XL</sub>	7.8	30	5	16,000	0.5 – 1.0	1.2	50
16136	G3000H <sub>XL</sub>	7.8	30	6	16,000	0.5 – 1.0	1.2	35
16137	G4000H <sub>XL</sub>	7.8	30	6	16,000	0.5 – 1.0	1.2	35
16652	GMH <sub>XL</sub> -L mixed-bed	7.8	30	6	16,000	0.5 – 1.0	1.2	35
16138	G5000H <sub>XL</sub>	7.8	30	9	14,000	0.5 – 1.0	1.2	15
16139	G6000H <sub>XL</sub>	7.8	30	9	14,000	0.5 – 1.0	1.2	15
16140	G7000H <sub>XL</sub>	7.8	30	9	14,000	0.5 – 1.0	1.2	15
16141	GMH <sub>XL</sub> mixed-bed	7.8	30	9	14,000	0.5 – 1.0	1.2	15
07112	GMH <sub>XL</sub> -HT	7.8	30	13	5,500	0.5 – 1.0	1.2	15
18403	Multipore H <sub>XL</sub> -M	7.8	30	5	16,000	0.5 – 1.0	1.0	35
17990	TSKgel SuperH1000 <b>-NEW-</b>	6.0	15	3	16,000	0.3 – 0.6	0.8	70
17991	TSKgel SuperH2000 <b>-NEW-</b>	6.0	15	3	16,000	0.3 – 0.6	0.8	60
17992	TSKgel SuperH2500 <b>-NEW-</b>	6.0	15	3	16,000	0.3 – 0.6	0.8	60
17993	TSKgel SuperH3000 <b>-NEW-</b>	6.0	15	3	16,000	0.3 – 0.6	0.8	40
17994	TSKgel SuperH4000 <b>-NEW-</b>	6.0	15	3	16,000	0.3 – 0.6	0.8	40
17995	TSKgel SuperH5000 <b>-NEW-</b>	6.0	15	3	16,000	0.3 – 0.6	0.8	40
17996	TSKgel SuperH6000 <b>-NEW-</b>	6.0	15	5	10,000	0.3 – 0.6	0.8	30
17997	TSKgel SuperH7000 <b>-NEW-</b>	6.0	15	5	10,000	0.3 – 0.6	0.8	30
17998	TSKgel SuperHM-L <b>-NEW-</b>	6.0	15	3	16,000	0.3 – 0.6	0.8	60
17999	TSKgel SuperHM-N <b>-NEW-</b>	6.0	15	3	16,000	0.3 – 0.6	0.8	40
18000	TSKgel SuperHM-M <b>-NEW-</b>	6.0	15	3	16,000	0.3 – 0.6	0.8	40
18001	TSKgel SuperHM-H <b>-NEW-</b>	6.0	15	3 and 5	16,000	0.3 – 0.6	0.8	40

## Ordering Information

### Analytical and preparative TSK-GEL GPC column products: typical properties (continued)

Part #	Description	ID (mm)	Length (cm)	Particle Size ( $\mu\text{m}$ )	Min. Number Theoretical Plates	Flow Rate (mL/min)		Maximum Pressure Drop (kg/cm <sup>2</sup> )
						Range	Max.	
<b>Stainless steel columns</b>								
19309	TSKgel SuperHZ1000 <b>-NEW-</b>	4.6	15	3	16,000	0.15 – 0.35	0.4	56
19302	TSKgel SuperHZ1000 <b>-NEW-</b>	6.0	15	3	16,000	0.25 – 0.6	0.7	56
19310	TSKgel SuperHZ2000 <b>-NEW-</b>	4.6	15	3	16,000	0.15 – 0.35	0.4	50
19303	TSKgel SuperHZ2000 <b>-NEW-</b>	6.0	15	3	16,000	0.25 – 0.6	0.7	50
19311	TSKgel SuperHZ2500 <b>-NEW-</b>	4.6	15	3	16,000	0.15 – 0.35	0.4	40
19304	TSKgel SuperHZ2500 <b>-NEW-</b>	6.0	15	3	16,000	0.25 – 0.6	0.7	40
19312	TSKgel SuperHZ3000 <b>-NEW-</b>	4.6	15	3	16,000	0.15 – 0.35	0.4	30
19305	TSKgel SuperHZ3000 <b>-NEW-</b>	6.0	15	3	16,000	0.25 – 0.6	0.7	30
19313	TSKgel SuperHZ4000 <b>-NEW-</b>	4.6	15	3	16,000	0.15 – 0.35	0.4	35
19306	TSKgel SuperHZ4000 <b>-NEW-</b>	6.0	15	3	16,000	0.25 – 0.6	0.7	35
19960	TSKgel SuperHZM-N <b>-NEW-</b>	4.6	15	3	16,000	0.15 – 0.35	0.4	35
19961	TSKgel SuperHZM-N <b>-NEW-</b>	6.0	15	3	16,000	0.25 – 0.6	0.7	35
19962	TSKgel SuperHZM-M <b>-NEW-</b>	4.6	15	3 and 5	16,000	0.15 – 0.35	0.4	20
19963	TSKgel SuperHZM-M <b>-NEW-</b>	6.0	15	3 and 5	16,000	0.25 – 0.6	0.7	20
19964	TSKgel SuperHZM-H <b>-NEW-</b>	4.6	15	10	9,000	0.15 – 0.35	0.4	10
19965	TSKgel SuperHZM-H <b>-NEW-</b>	6.0	15	10	9,000	0.25 – 0.6	0.7	10
<b>Guard columns</b>								
18404	MultiporeH <sub>XL</sub> -M Guard column	6.0	4.0	5	For P/N 18403			
07113	H <sub>XL</sub> -L Guard column	6.0	4.0	7	For G1000H <sub>XL</sub> - G4000H <sub>XL</sub> , GMH <sub>XL</sub> -L			
13727	H <sub>XL</sub> -H Guard column	6.0	4.0	13	For G5000H <sub>XL</sub> - G7000H <sub>XL</sub> + GMH <sub>XL</sub>			
17368	Guard column H <sub>HR</sub> -L and H <sub>HR</sub>	6.0	4.0	5	For G1000-4000H <sub>HR</sub> and P/N 17362			
17369	Guard column H <sub>HR</sub> -H and H <sub>HR</sub>	6.0	4.0	5	For G5000-7000H <sub>HR</sub> and P/Ns 18055, 17392 and 17360			
19314	SuperHZ Guard Column <b>-NEW-</b>	4.6	2.0	4	For 4.6 mm ID columns of SuperHZ1000-4000 and HZM-N &-M			
19668	SuperHZ Guard Column <b>-NEW-</b>	4.6	2.0	10	For 4.6 mm ID columns of SuperHZM-H			
19666	SuperHZ Guard Column <b>-NEW-</b>	4.6	3.5	4	For 6.0 mm ID columns of SuperHZ1000-4000 and HZM-N &-M			
19667	SuperHZ Guard Column <b>-NEW-</b>	4.6	3.5	10	For 6.0 mm ID columns of SuperHZM-H			