

# HiTrap HIC Columns

**HiTrap Phenyl FF (high sub), 1 ml and 5 ml**

**HiTrap Phenyl FF (low sub), 1 ml and 5 ml**

**HiTrap Phenyl HP, 1 ml and 5 ml**

**HiTrap Butyl FF, 1 ml and 5 ml**

**HiTrap Butyl-S FF, 1 ml and 5 ml**

**HiTrap Octyl FF, 1 ml and 5 ml**

**HiTrap HIC Selection Kit**

HiTrap™ Phenyl FF (high sub), HiTrap Phenyl FF (low sub), HiTrap Phenyl HP, HiTrap Butyl FF, HiTrap Butyl-S FF, and HiTrap Octyl FF are prepacked 1 ml and 5 ml, ready to use hydrophobic interaction chromatography (HIC) columns. The column together with modern chromatography media, provides fast, reproducible and easy separation in a convenient format.

HiTrap™ HIC Selection Kit consists of six Hydrophobic Interaction Chromatography media (HIC) with different hydrophobic characteristics. The kit provides the possibility to screen for the most appropriate HIC medium for specific applications. The recommended test procedures help optimize salt concentration, sample loading, resolution and other chromatographic parameters. The six different media are prepacked in ready to use 1 ml HiTrap columns.

Separations are easily performed with a syringe, a pump or a chromatography system such as ÄKTAdesign™ or FPLC™ System.



<b>Code No.</b>	<b>Product</b>	<b>No. supplied</b>
17-1355-01	HiTrap Phenyl FF (high sub)	5 x 1 ml
17-5193-01	HiTrap Phenyl FF (high sub)	5 x 5 ml
17-1353-01	HiTrap Phenyl FF (low sub)	5 x 1 ml
17-5194-01	HiTrap Phenyl FF (low sub)	5 x 5 ml
17-1351-01	HiTrap Phenyl HP	5 x 1 ml
17-5195-01	HiTrap Phenyl HP	5 x 5 ml
17-1357-01	HiTrap Butyl FF	5 x 1 ml
17-5197-01	HiTrap Butyl FF	5 x 5 ml
17-0978-13	HiTrap Butyl-S FF	5 x 1 ml
17-0978-14	HiTrap Butyl-S FF	5 x 5 ml
17-1359-01	HiTrap Octyl FF	5 x 1 ml
17-5196-01	HiTrap Octyl FF	5 x 5 ml

<b>Connectorkit</b>		
<b>Connectors supplied</b>	<b>Usage</b>	<b>No. supplied</b>
1/16" male/luer female	Connection of syringe to top of HiTrap column	1
Tubing connector flangeless/M6 female	Connection of tubing (e.g. Peristaltic Pump P1) to bottom of HiTrap column*	1
Tubing connector flangeless/M6 male	Connection of tubing (e.g. Peristaltic Pump P1) to top of HiTrap column**	1
Union 1/16" female/M6 male	Connection to original FPLC System through bottom of HiTrap column	1
Union M6 female/1/16" male	Connection to original FPLC System through top of HiTrap column	1
Stop plug female, 1/16"	Sealing bottom of HiTrap column	2, 5 or 7

\* Union 1/16" female/M6 male is also needed.

\*\* Union M6 female/1/16" male is also needed.

<b>Code No.</b>	<b>Product</b>	<b>No. supplied</b>
11-1034-53	HiTrap HIC Selection Kit	
	<b>Including:</b>	
	HiTrap Phenyl FF (high sub)	1 x 1 ml
	HiTrap Phenyl FF (low sub)	1 x 1 ml
	HiTrap Phenyl HP	1 x 1 ml
	HiTrap Butyl FF	1 x 1 ml
	HiTrap Butyl-S FF	1 x 1 ml
	HiTrap Octyl FF	1 x 1 ml

### **Connectorkit**

<b>Connectors supplied</b>	<b>Usage</b>	<b>No. supplied</b>
1/16" male/luer female	Connection of syringe to top of HiTrap column	1
Tubing connector flangeless/M6 female	Connection of tubing (e.g. Peristaltic Pump P1) to bottom of HiTrap column*	1
Tubing connector flangeless/M6 male	Connection of tubing (e.g. Peristaltic Pump P1) to top of HiTrap column**	1
Union 1/16" female/M6 male	Connection to original FPLC System through bottom of HiTrap column	1
Union M6 female/1/16" male	Connection to original FPLC System through top of HiTrap column	1
Stop plug female, 1/16"	Sealing bottom of HiTrap column	2, 5 or 7

\* Union 1/16" female/M6 male is also needed.

\*\* Union M6 female/1/16" male is also needed.

# Table of contents

1. Media properties	5
2. Column	7
3. Hydrophobic Interaction Chromatography	9
4. Operation	19
5. Determination of binding capacity	21
6. Clearing and regeneration	22
7. Scaling up	22
8. Storage	23
9. Further information	23
10. Ordering Information	24

# 1. Media properties

The HIC media are based on the highly cross-linked beaded agarose matrices, Sepharose™ Fast Flow and Sepharose High Performance.

The media have excellent flow properties with high physical and chemical stabilities. All Sepharose matrices show virtually no non-specific adsorption and are resistant to microbial degradation due to the presence of the unusual sugar, 3,6-anhydro-L-galactose. The hydrophobic ligands are coupled to the monosaccharide units via glycidylethers. The resulting ether bonds are both stable and uncharged. Characteristics of the different HiTrap HIC media are listed in Table 1 and their chemical stability is shown in Table 3.

**Phenyl Sepharose High Performance** is based on a 34 µm matrix and is ideal for laboratory and intermediate process scale separations and for final step purifications where high resolution is needed. The degree of substitution gives Phenyl Sepharose High Performance a selectivity similar to that of Phenyl Sepharose 6 Fast Flow (low sub).

**Phenyl Sepharose 6 Fast Flow (low sub) and Phenyl Sepharose 6 Fast Flow (high sub)** are based on a 90 µm matrix. They are ideal for initial and intermediate step purifications requiring a matrix with medium to high hydrophobicity. The availability of two degrees of substitution grades increases the possibility of finding the best selectivity and capacity for a given application.

**Butyl Sepharose 4 Fast Flow** is based on a 90 µm matrix. It is intended for initial and intermediate step purifications requiring a matrix with low to medium hydrophobicity. Butyl Sepharose 4 Fast Flow often works efficiently with rather low salt concentrations. The mechanism of adsorption and desorption onto the butyl ligand is different from that onto the phenyl ligand, giving a different selectivity.

**Butyl-S Sepharose 6 Fast Flow** is based on a 90 µm matrix. The main differences between Butyl-S Sepharose 6 Fast Flow and Butyl Sepharose 4 Fast Flow lie in the length of their spacer arms, the concentration of the immobilized ligands, and the type of connector atom (O-ether or S-ether) linking each ligand to the Sepharose base matrix. Butyl-S Sepharose 6 Fast Flow contains a sulfur atom as a linker between the spacer arm and

**Table 1.** HiTrap HIC media characteristics.

	Phenyl Sepharose High Performance	Phenyl Sepharose 6 FastFlow (low sub)	Phenyl Sepharose 6 FastFlow (high sub)	Butyl Sepharose 4 FastFlow	Butyl-S Sepharose 6 FastFlow	Octyl Sepharose 4 FastFlow
Hydrophobic ligand	Phenyl	Phenyl	Phenyl	n-Butyl	Butyl-S	n-Octyl
Ligand density	25 µmol/ml medium	25 µmol/ml medium	40 µmol/ml medium	40 µmol/ml medium	10 µmol/ml medium	5 µmol/ml medium
Mean bead size	34 µm	90 µm	90 µm	90 µm	90 µm	90 µm
Bead size range	24-44 µm	45-165 µm	45-165 µm	45-165 µm	45-165 µm	45-165 µm
Bead structure	6% crosslinked agarose, spherical	6% crosslinked agarose, spherical	6% crosslinked agarose, spherical	4% crosslinked agarose, spherical	6% crosslinked agarose spherical	4% crosslinked agarose, spherical
pH stability						
Short term	2-14	2-14	2-14	2-14	2-14	2-14
Working	3-13	3-13	3-13	3-13	3-13	3-13
Long term	3-13	3-13	3-13	3-13	3-13	3-13
Storage	0.01 M NaOH or 20% ethanol	0.01 M NaOH or 20% ethanol	0.01 M NaOH or 20% ethanol	0.01 M NaOH or 20% ethanol	0.01 M NaOH or 20% ethanol	0.01 M NaOH or 20% ethanol

Please note the following:

**pH stability, long term:** pH interval where the medium can be operated without significant change in function.

**pH stability, working:** pH interval where the medium binds protein as intended or is needed for elution, without adverse long-term effect.

**pH stability, short term:** pH interval to which the medium can be subjected for cleaning- or sanitization-in-place (accumulated 90-100 hours at room temperature) without significant change in function.

the butyl ligand. It is the least hydrophobic medium in the GE Healthcare portfolio and is intended for purification or removal of strongly hydrophobic biomolecules at low salt concentrations, with high recovery and low risk of denaturation.

**Octyl Sepharose 4 Fast Flow** is based on a 90 µm matrix. It has a different hydrophobic character from the phenyl and butyl ligands and is an important complement to the other hydrophobic matrices.

## 2. Column

The column is made of polypropylene, which is biocompatible and non-interactive with biomolecules. The top and bottom frits are manufactured from porous polyethylene. It is delivered with a stopper on the inlet and a snap-off end on the outlet.

The separation can be easily achieved using a syringe together with the supplied luer adaptor, a peristaltic pump, or in a chromatography system such as ÄKTA™ or FPLC System.

**Note:** To prevent leakage it is essential to ensure that the adaptor is tight.

**Note:** The column cannot be opened or refilled.

Characteristics of HiTrap column are listed in Table 2.

**Table 2.** Characteristics of HiTrap column.

Column volume	1 ml or 5 ml
Column dimension	0.7 × 2.5 cm (1 ml) or 1.6 × 2.5 cm (5 ml)
Recommended flow rate	HiTrap 1 ml: 1 ml/min HiTrap 5 ml: 5 ml/min
Maximum flow rate*	HiTrap 1 ml: 4 ml/min HiTrap 5 ml: 20 ml/min
Maximum back pressure	0.3 MPa, 3 bar

\* Room temperature, aqueous buffers. Stay below maximum pressure.

**Table 3.** Chemical stability.

	Phenyl Sephacrose High Performance	Phenyl Sephacrose 6 FastFlow (low sub)	Phenyl Sephacrose 6 FastFlow (high sub)	Butyl Sephacrose 4 FastFlow	Butyl-S Sephacrose 6 FastFlow	Octyl Sephacrose 4 FastFlow
1 M NaOH	X	X	X	X	X	X
1 M acetic acid	X	n.d.	n.d.	n.d.	n.d.	n.d.
1 mM HCl	n.d.	n.d.	n.d.	X	X	X
3 M $(\text{NH}_4)_2\text{SO}_4$	n.d.	X	X	n.d.	X	n.d.
70 % ethanol	X	X	X	X	X	X
30 % isopropanol	X	X	X	X	X	X
6 M GuHCl	X	X	X	X	X	X
8 M Urea	X	X	X	X	X	X

X = Functionally stable for 7 days at +40 °C


n.d. = not determined



# 3. Hydrophobic Interaction Chromatography

## Overview

Separation of biomolecules on HIC media is based on an interplay between the hydrophobicity of the medium, the nature and composition of the sample, the prevalence and distribution of surface-exposed hydrophobic amino acid residues, and the type and concentration of salt in the binding buffer. Unlike reversed phase chromatography (RPC), which is a separation method closely related to HIC, the adsorption of biological solutes to HIC adsorbents is promoted, or otherwise modulated, by the presence of relatively high concentrations of anti-chaotropic salts such as ammonium sulfate, sodium sulfate, etc (Fig 1). Desorption of bound solutes is achieved simply by stepwise or gradient elution with buffers of low salt content.



**Fig 1.** The Hofmeister series of some anions and cations arranged according to their effects on the solubility of protein in aqueous solutions.

HIC media available from GE Healthcare are produced as a graded series of hydrophobic adsorbents based on alkyl or aryl ligands attached to a hydrophilic base matrix, e.g. Sepharose. In each instance, the type and concentration of ligand has been optimized to cover the range of hydrophobicities of the proteins in a biological extract, varying from weak to moderate to strong hydrophobic proteins. This strategy results in HIC adsorbents for “all occasions” where the emphasis is on high recovery, purity, and reduced risk for denaturation of the target proteins in a biological extract.

## Factors affecting HIC

The main parameters to consider when selecting a HIC medium and optimizing its chromatographic performance are:

- The nature of the base matrix (e.g. agarose, organic co-polymers, etc.)
- Structure of the ligand
- Concentration of the ligand
- Characteristics of the target protein and other sample components
- Type of salt
- Concentration of salt
- Temperature
- pH

Of these parameters, **the type and concentration of ligand as well as the type and concentration of salt added during the adsorption step** are of paramount importance in determining the outcome of a HIC event. In general, the type of immobilized ligand determines its adsorption selectivity toward the proteins in a sample while its concentration determines its adsorption capacity.

In general, HIC media fall into two groups, depending on their interactions with sample components. Straight alkyl chains (butyl, octyl) show a “pure” hydrophobic character, while aryl ligands (phenyl) show a mixed mode behaviour, where both aromatic and hydrophobic interactions as well as lack of charge play simultaneous roles. The choice of ligand must be

determined empirically through screening experiments for each individual separation problem.

**The characteristics of the target protein** (in a HIC context) are usually not known since minimal data are available in this respect. There are some published data regarding the hydrophobicity indices for a number of purified proteins based on amino acid composition, the number and distribution of surface-exposed hydrophobic amino acids, and the order of their elution from RPC columns but few, if any, have proved to be useful when purifying a protein in a real biological sample. For this and other reasons, the adsorption behavior of a protein exposed to a HIC medium has to be determined on a case-by-case basis.

**The solvent** is one of the most important parameters which influences capacity and selectivity in HIC. In general, the adsorption process is more selective than the desorption process. It is therefore important to optimize the start buffer with respect to pH, type of solvent, type of salt and concentration of salt. The addition of various "salting-out" salts to the sample promotes ligand-protein interactions in HIC. As the concentration of salt is increased, the amount of bound protein increases up to the precipitation point for the protein. Each type of salt differs in its ability to promote hydrophobic interactions and it may be worthwhile testing several salts. The most commonly used salts are  $(\text{NH}_4)_2\text{SO}_4$ ,  $\text{Na}_2\text{SO}_4$ ,  $\text{NaCl}$ ,  $\text{KCl}$  and  $\text{CH}_3\text{COONH}_4$ . At a given concentration, ammonium sulphate often gives the best resolution of a mixture of standard proteins compared to other salts. If sodium chloride is used, a concentration of up to 3-4 M is usually needed. Due to instability, ammonium sulphate is not suitable when working at pH values above 8.0. Sodium sulphate is also a very good salting-out agent but protein solubility problems may exclude its use at high concentrations.

**The effect of pH** is not well established. In general, an increase in pH above 8.5 weakens hydrophobic interactions whereas a decrease in pH below 5.0 results in an apparent increase in the retention of proteins on HIC adsorbents. In the range of pH 5-8.5, the effect seems to be minimal or insignificant.

In regard to **temperature**, it is generally accepted that the binding of proteins to HIC adsorbents is entropy driven, which implies that the solute-adsorbent interaction increases with increase in temperature. In some instances, the reverse effect has been observed. In practical work, one should

be aware that a downstream purification process that is developed at room temperature might not be reproduced in the cold room, or vice versa. In other instances, temperature control is mandatory in order to obtain reproducible results from run to run.

Sometimes it is necessary to weaken the protein-ligand interactions by including different additives. Commonly used are water-miscible alcohols (propanol, ethylene glycol), detergents (SDS) and solutions of chaotropic salts (lithium perchlorate, urea, guanidine hydrochloride).

Further information about hydrophobic interaction chromatography can be found in the Handbook "HIC & RPC, Principles and Methods", see Ordering Information. Check also [www.amershambiosciences.com/hitrap](http://www.amershambiosciences.com/hitrap) and [www.chromatography.amershambiosciences.com](http://www.chromatography.amershambiosciences.com) for further information. For example can reference lists be found on this site (Code No. 18-1156-78, 18-1156-79, 18-1156-80, 18-1156-81).

## Screening experiments

Protein binding to HIC adsorbents is promoted by moderate to high concentrations of "salting-out" salts, which also have a stabilizing influence on protein structure. Elution is achieved by a linear or step-wise decrease in concentration of the salt. The HIC medium should bind the protein of interest at a reasonably low concentration of salt. As mentioned before, binding conditions are dependent on the salt chosen. The salt concentration should be below that which causes precipitation of proteins in the sample. We recommend 1 M ammonium sulphate as a good starting buffer for screening experiments. If the substance does not bind, a more hydrophobic medium should be chosen. If the substance binds so strongly that non-polar additives are required for elution, a less hydrophobic medium should be tried. The bound protein should be eluted from the column with high recovery.

Figures 2 to 7 show some typical elution profiles from screening experiments and recommendations for further experimental work. The shaded area shows the elution position of the protein of interest.





Fig 2.

**Result:** The protein is eluted early in the gradient. Resolution is not satisfactory.

**Discussion:** Not much can be gained by changing the salt concentration. A decrease in the salt concentration of the start buffer will decrease binding capacity and might even lead to coelution with the unbound fraction. An increase in the salt concentration might lead to the co-adsorption of unwanted impurities leading to a decrease in selectivity. Change of pH in the start buffer might result in stronger binding and higher selectivity. Trying another salt may also improve performance.

**Next step:** Repeat the experiment at a lower or higher pH or with a salt of higher salting-out strength. If no improvement in selectivity is obtained, try a medium with a different ligand or with a higher degree of ligand substitution.




**Fig 3.**

**Result:** The protein is eluted near the end of the gradient. Resolution is not satisfactory.

**Discussion:** A decrease in the initial salt concentration will weaken the binding, resulting in earlier elution of the protein. It may, however, not have a positive effect on selectivity, since the contaminants are eluted very close to, both before and after, the protein of interest. Changing the pH of the start buffer or changing to another salt may have more impact on resolution and should be tried.

**Next step:** Repeat the experiment at a lower or higher pH value or with a salt with lower salting-out strength. If no improvement in selectivity is obtained, try a medium with a different ligand or with a lower degree of ligand substitution.




**Fig 4.**

**Result:** The protein is eluted in the middle of the gradient. Resolution is not satisfactory.

**Discussion:** Changing the concentration of salt in the start buffer will have a limited effect since the contaminants are eluted very close to, both before and after, the protein of interest. Changing the pH of the start buffer or changing to a different salt may have impact on resolution and should be tried. Trying a different gradient slope may also be an effective way to increase the resolution.

**Next step:** Repeat the experiment at a lower or higher pH or with a salt with higher salting-out effect. If no improvement in selectivity is obtained, try a medium with a different ligand or with a higher degree of ligand substitution.




**Fig 5.**

**Result:** The protein is eluted early in the gradient. Resolution is satisfactory.

**Discussion:** This can be a good choice of medium. The fact that the protein is eluted early in the gradient indicates that the binding capacity might be low. An increase in salt concentration to compensate for this may lead to a decrease in selectivity since some of the unbound proteins might be adsorbed together with the protein of interest. Not much can be gained by changing the pH of the start buffer since the resolution was considered satisfactory.

**Next step:** If low capacity is a problem and a moderate increase in salt concentration leads to loss in resolution, try a medium with a different ligand or a medium with a higher degree of ligand substitution.






**Fig 6.**

**Result:** The product is eluted near the end of the gradient. Resolution is satisfactory.

**Discussion:** This can be a good choice of medium. A disadvantage might be that some of the most hydrophobic substances bind so strongly that they are difficult to remove from the media. Decreasing the concentration of salt in the start buffer will give earlier elution of the protein and reduce the risk that strongly bound proteins are difficult to elute. It also leads to a reduction in cycle time.

**Next step:** If problems with strong binding of contaminants are discovered, try a medium with a different ligand or with a lower degree of ligand substitution.



**Fig 7.**

**Result:** The protein is eluted in the middle of the gradient. Resolution is satisfactory.

**Discussion:** The choice of ligand is very good and there is little risk for strong binding of the most hydrophobic contaminants.

## 4. Operation

### Buffer preparation

Water and chemicals should be of high purity. When using high salt concentration buffers, especially ammonium sulphate, it is possible to treat the buffer with activated charcoal for approx. 2 hours and filter through 0.45  $\mu\text{m}$  filter to avoid UV baseline drift.

### Sample preparation

Since adsorption is carried out at high salt concentration, the composition of the sample should be adjusted to the pH and ionic strength of the start buffer (high salt buffer) for consistent and reproducible results. When possible, dissolve the sample in start buffer. Buffer exchange can be carried out using HiPrep™ 26/10 Desalting, HiTrap Desalting or PD-10 columns for ionic strengths up to ~1.5 M, (at higher ionic strengths there is a risk that the medium will shrink).

Another way to increase the ionic strength is by addition of solid salt (note: precipitation may occur due to high local salt concentrations), by addition of the salt as a high concentration stock solution or by dilution of the sample with start buffer followed by pH adjustment.

The sample should be fully solubilized. We recommend centrifugation or filtration immediately before loading on the column to remove particulate material (0.45  $\mu\text{m}$  filter). Never apply turbid solutions to the column. Turbidity indicates sample insolubility which may be due to incorrect ionic strength.

High sample viscosity causes high back pressure, instability of the sample zone and gives irregular flow pattern with decreased resolution. High back pressure can also damage the column packing. Recommended maximum sample viscosity corresponds to a protein concentration of ~50 mg/ml in aqueous solution. If lipids or other very hydrophobic substances are present in the sample, they may interact very strongly with the HIC column, diminishing capacity and being very difficult to remove. In such cases, using a slightly less hydrophobic column as a pre-column can be very efficient. The pre-column should be chosen to bind the most hydrophobic material and allow the substance of interest to pass through under starting conditions.

## Column equilibration

1. Fill the syringe or pump tubing with elution buffer (low salt buffer). Remove the stopper. To avoid introducing air into the column, connect the column "drop to drop" to either the syringe (via the adaptor) or to the pump tubing.
2. Remove the snap-off end at the column outlet.
3. Wash the column with 5 ml elution buffer at 1 ml/min.
4. Wash with 5-10 ml start buffer (high salt buffer). If air is trapped in the column, wash with buffer until the air disappears.

**Note:** If a P1-pump is used a max flow rate of 1-3 ml/min can be run on a HiTrap 1 ml column packed with Sepharose High Performance media.

## Elution with linear descending gradients

A linear decrease of the salt concentration is the most frequently used type of elution in hydrophobic interaction chromatography. Recommended buffers are 50 mM sodium phosphate, 1.0 M ammonium sulphate, pH 7.0 as start buffer and 50 mM sodium phosphate, pH 7.0 as elution buffer. Continuous gradients can be prepared in different ways depending on available equipment:

- A peristaltic pump and a gradient mixer  
e.g. pump P-1, gradient mixer GM-1
  - A one pump system, e.g. ÄKTAprime™ plus
  - A two pump system, e.g. ÄKTAdesign or FPLC
1. Equilibrate the column (see column equilibration).
  2. Adjust the sample to the chosen starting pH and ionic strength (see sample preparation).
  3. Apply the sample.
  4. Wash with 5-10 ml of the start buffer until the UV trace of the effluent returns to near baseline.
  5. Start the gradient elution. A gradient volume of 10-20 ml is usually enough.

6. Regenerate the column by washing with 5 ml of distilled water followed by 5 ml of start buffer. The column is now ready for a new sample. Avoid storage of the column or the chromatography system in high salt buffer to prevent crystal build-up in the equipment.

## **Elution with stepwise descending gradients**

Stepwise elution is the sequential use of the same buffer at different ionic strengths. It is technically simple and fast and suitable for syringe operation. It is often used for sample concentration and sample clean up. Stepwise elution gives small peak volumes and the resolution depends on the difference in elution power between each step. When stepwise elution is applied, one has to keep in mind the danger of artifactual peak when a subsequent step is executed too early after a tailing peak. For this reason it is recommended to start with a continuous gradient to characterise the sample and its chromatographic behaviour.

## **5. Determination of binding capacity**

The amount of sample which can be applied to a column depends on the capacity of the column and the degree of resolution required. The capacity is dependent on the sample composition, chosen starting conditions such as pH, ionic strength, buffer salts and the flow rate at which the adsorption is done. The dynamic capacity can be determined by frontal analysis using real sample:

1. Equilibrate the column (see column equilibration).
2. Adjust the sample to the chosen starting pH and ionic strength (see sample preparation).
3. Determine the concentration of the actual protein in the sample by UV, PAGE, ELISA or other appropriate techniques.
4. Apply the sample solution to the column with a pump or a syringe at the flow rate to be used in the purification method. Collect fractions and continue to apply sample until the column is saturated.
5. Wash the column with 5-10 ml start buffer until the baseline is stable.

6. Elute bound proteins with 2-5 ml of elution buffer (low salt buffer) and collect the eluate.
7. Analyse fractions and eluates from steps 4 and 6 for the protein in question and determine the breakthrough profile. The practical capacity is the amount that can be applied without any breakthrough and the total capacity available is determined by analysing eluate from step 6.

## 6. Cleaning and regeneration

HIC adsorbents can normally be regenerated by washing with distilled water. To prevent gradual build-up of contaminants on the column, regular cleaning is recommended. Precipitated proteins can be removed by washing with 5 ml of 0.5-1.0 M NaOH followed by 5 ml of water at a flow rate of 1 ml/min. Strongly bound substances can be removed by washing with 5-10 ml of up to 70% ethanol or 30% isopropanol.

## 7. Scaling up

Columns and media for scale-up are available. For quick scale up of purifications, two or three HiTrap HIC columns of the same type can be connected in series (backpressure will increase). All HiTrap HIC columns included in HiTrap HIC Selection Kit are available as individual packages of 5 x 1 ml and 5 x 5 ml. Further scale up can be achieved using the prepacked columns HiPrep 16/10 Phenyl (high sub), HiPrep 16/10 Phenyl (low sub), HiPrep 16/10 Butyl, HiPrep 16/10 Octyl or HiLoad™ Phenyl Sepharose High Performance or bulk media packs. See ordering information.

## 8. Storage

Store the HiTrap HIC columns equilibrated with 5–10 ml 20% ethanol or 0.01 M NaOH. The recommended storage temperature is +4 to +30 °C.

## 9. Further information

Visit [www.amershambiosciences.com/hitrap](http://www.amershambiosciences.com/hitrap) and [www.chromatography.amershambiosciences.com](http://www.chromatography.amershambiosciences.com) for further information. Useful handbooks are also available, see ordering information.

## 10. Ordering Information

<b>Product</b>	<b>No. Supplied</b>	<b>Code No.</b>
<b>Pre-packed columns</b>		
HiTrap HIC Selection Kit, 6 different HIC media	6 × 1 ml	11-0034-53
HiTrap Phenyl FF (high sub)	5 × 1 ml	17-1355-01
HiTrap Phenyl FF (high sub)	5 × 5 ml	17-5193-01
HiTrap Phenyl FF (low sub)	5 × 1 ml	17-1353-01
HiTrap Phenyl FF (low sub)	5 × 5 ml	17-5194-01
HiTrap Phenyl HP	5 × 1 ml	17-1351-01
HiTrap Phenyl HP	5 × 5 ml	17-5195-01
HiTrap Butyl FF	5 × 1 ml	17-1357-01
HiTrap Butyl FF	5 × 5 ml	17-5197-01
HiTrap Butyl-S FF	5 × 1 ml	17-0978-13
HiTrap Butyl-S FF	5 × 5 ml	17-0978-14
HiTrap Octyl FF	5 × 1 ml	17-1359-01
HiTrap Octyl FF	5 × 5 ml	17-5196-01
HiPrep 16/10 Phenyl FF (high sub)	1 × 20 ml	17-5095-01
HiPrep 16/10 Phenyl FF (low sub)	1 × 20 ml	17-5094-01
HiPrep 16/10 Butyl FF	1 × 20 ml	17-5096-01
HiPrep 16/10 Octyl FF	1 × 20 ml	17-5097-01
HiLoad 16/10 Phenyl Sepharose High Performance	1 × 20 ml	17-1085-01
HiLoad 26/10 Phenyl Sepharose High Performance	1 × 53 ml	17-1086-01



<b>Related products</b>	<b>Quantity</b>	<b>Code No.</b>
<b>Bulk media</b>		
Phenyl Sepharose High Performance	75 ml <sup>1</sup>	17-1082-01
Phenyl Sepharose 6 Fast Flow (low sub)	25 ml	17-0965-10
	200 ml <sup>1</sup>	17-0965-05
Phenyl Sepharose 6 Fast Flow (high sub)	25 ml	17-0973-10
	200 ml <sup>1</sup>	17-0973-05
Butyl Sepharose 4 Fast Flow	25 ml	17-0980-10
	200 ml <sup>1</sup>	17-0980-01
Butyl-S Sepharose 6 Fast Flow	25 ml	17-0978-10
	200 ml <sup>1</sup>	17-0978-02
Octyl Sepharose 4 Fast Flow	25 ml	17-0946-10
	200 ml <sup>1</sup>	17-0946-02

<sup>1</sup> Larger quantities are available. Please contact your local representative for further information.

<b>Accessories</b>	<b>No. Supplied</b>	<b>Code No.</b>
1/16" male/luer female*	2	18-1112-51
Tubing connector flangeless/M6 female*	2	18-1003-68
Tubing connector flangeless/M6 male*	2	18-1017-98
Union 1/16" female/M6 male*	6	18-1112-57
Union M6 female /1/16" male*	5	18-3858-01
Union luerlock female/M6 female	2	18-1027-12
HiTrap/HiPrep, 1/16" male connector for ÄKTAdesign	8	28-4010-81
Stop plug female, 1/16" <sup>†</sup>	5	11-0004-64
Fingertight stop plug, 1/16" <sup>‡</sup>	5	11-0003-55

\* One connector included in each HiTrap package.

† Two, five, or seven stop plugs female included in HiTrap packages depending on the product.

‡ One fingertight stop plug is connected to the top of each HiTrap column at delivery.

## Literature

---

Hydrophobic Interaction Chromatography & Reversed Phase Chromatography, Principles and Methods, Handbook	11-0012-69
Convenient Protein Purification, HiTrap Column Guide	18-1129-81

---



[www.amershambiosciences.com/hitrap](http://www.amershambiosciences.com/hitrap)  
[www.gehealthcare.com](http://www.gehealthcare.com)

GE Healthcare  
Amersham Biosciences AB  
Björkgatan 30  
751 84 Uppsala  
Sweden

GE Healthcare  
Amersham Biosciences Europe GmbH  
Munzinger Strasse 9  
D-79111 Freiburg  
Germany

GE Healthcare  
Amersham Biosciences UK Ltd  
Amersham Place  
Little Chalfont  
Buckinghamshire, HP7 9NA  
UK

GE Healthcare  
Amersham Biosciences Corp  
800 Centennial Avenue  
P.O. Box 1327  
Piscataway, NJ 08855-1327  
USA

GE Healthcare  
Amersham Biosciences KK  
Sanken Bldg.  
3-25-1, Hyakunincho  
Shinjuku-ku, Tokyo 169-0073  
Japan

ÄKTA, ÄKTAdesign, ÄKTAprime, FPLC, HiLoad, HiPrep, HiTrap and Sepharose are trademarks of GE Healthcare Ltd., a General Electric Company. GE tagline and GE monogram are trademarks of General Electric Company.

Separating Miraculin with Butyl-S Sepharose 6 Fast Flow is subject to US patent No 5,886,155. Licenses are available from BioResources International, Inc., of Somerset, N.J., U.S.A.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. GE Healthcare reserves the right, subject to any regulatory and contractual approval, if required, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your local GE Healthcare representative for the most current information.

© 2005 General Electric Company – All rights reserved.

Amersham Biosciences AB, a General Electric Company going to market as GE Healthcare.



imagination at work

11-0034-52 AB