

# Characterizing the Separation Properties of Two New Mixed-Mode Sorbents

A. Topol<sup>a</sup>, J. Pezzini<sup>b</sup>, X. Santarelli<sup>b</sup>, M.P. Barouillet<sup>c</sup>, G. Clotent-Sanchez<sup>c</sup>, M. Cerruti<sup>d</sup>, G. Magistrelli<sup>e</sup>, J. Jenco<sup>a</sup>, V. Brochier<sup>a</sup> and R. Gantier<sup>a</sup>

<sup>a</sup>Pall Life Sciences. <sup>b</sup>ESTBB (EA 4135). <sup>c</sup>RMSB (UMR5536), Bordeaux, France. <sup>d</sup>UPS3044, Saint Christol, France. <sup>e</sup>NovImmune, Plan les Ouates, Switzerland.

## INTRODUCTION

Mixed mode chromatography offers new selectivities compared to **conventional techniques** such as ion exchange and hydrophobic interaction. New mixed-mode chromatography sorbents combine in a single tool the ability to exploit both the ionic and hydrophobic characteristics of a protein, overcoming the salt issue.

HEA and PPA HyperCel™ sorbents carry mixed-mode synthetic ligands, immobilized on a **robust and scalable** cross-linked cellulose matrix (HyperCel) that confers high porosity, chemical stability and low non-specific binding. The ligands include aliphatic (HEA – hexylamine) and aromatic (PPA – phenylpropylamine) amines. The aromatic group on PPA HyperCel sorbent confers an enhanced hydrophobicity compared to HEA HyperCel sorbent.

A panel of proteins of different isoelectric points (pI) and hydrophobicity (GRAVY<sup>1</sup> index) was used at different pH and ionic strength to characterize the interaction mechanisms with HEA and PPA HyperCel sorbents. The mixed-mode sorbents were also used to purify industrially relevant proteins. Their unique selectivity, illustrated by the separation of protein or peptide isoforms, was demonstrated.

## CHARACTERIZATION OF INTERACTIONS WITH MODEL PROTEINS

Four different proteins with different pI and hydrophobicity were used to characterize the interactions between the mixed-mode sorbents and proteins:

**Acidic proteins:** bovine serum albumin (BSA, pI = 4.7, GRAVY<sup>1</sup> = -0.433) and ovalbumin (pI = 4.6, GRAVY = -0.008). **Basic proteins:** α-chymotrypsinogen (α-chymo, pI = 8.5, GRAVY = +0.051) and lysozyme (pI = 9.1, GRAVY = -0.314).

(1) GRAVY = GRand AVerage of hYdrophobicity - Gasteiger, E., et al, Protein Identification and Analysis Tools on the ExPASy Server; John M. Walker Ed.: The Proteomics Protocols Handbook, Humana Press (2005) 571-607.

### Method used

- Salt and pH influence: α-chymo in different pH and salt conditions (Figure 1).
- Protein separation: capture of a four-protein mixture in PBS, pH 7.4. Elution by pH step and gradient (Figure 2).

### Results

- **Influence of salt and pH on binding of α-chymo on HEA and PPA HyperCel sorbents** (Figure 1):
  - ▶ Limited or no salt (presence or not) influence on binding
  - ▶ Limited binding at acidic pHs
- **Separation of protein mixture** (Figure 2):
  - ▶ Only acidic (BSA and ovalbumin on both sorbents) or hydrophobic basic (α-chymo on PPA HyperCel sorbent) protein bound in PBS, pH 7.4.
  - ▶ Basic proteins (α-chymo on PPA HyperCel sorbent) are eluted first, followed by acidic ones (BSA and ovalbumin).




Figure 1. Influence of salt and pH on binding of α-chymo on HEA and PPA HyperCel sorbents.




Figure 2. Separation of a protein mixture on HEA and PPA HyperCel sorbents. Load: PBS, pH 7.4, step elution at pH 5.4, then pH gradient from 5.4 to 2.6.

- Binding relies more on hydrophobic than on electrostatic interactions.
- Different interaction compared to standard IEX or HIC sorbents.
- Electrostatic repulsion between positive charges on protein and sorbents for elution (Figure 3).




Figure 3. Adsorption and desorption mechanism on HEA and PPA HyperCel mixed-mode sorbents.

## SEPARATION OF RECOMBINANT CHEMOKINE ISOFORMS

PPA HyperCel sorbent was used for the purification of a chemokine (cytokine) from a PEAK cell culture supernatant with 10% fetal calf serum (courtesy of NovImmune).




Figure 4. A = Separation of chemokine isoforms (MW 8.0–8.2 kDa and 8.5 kDa) by a pH gradient on PPA HyperCel sorbent. B = SELDI analysis of the fractions eluted at 115 mL (blue) and 125 mL (red).

### Method used

- Load in PBS, pH 7.4 on PPA HyperCel sorbent. Elution at pH 5.0, 4.0 and 2.6.
- Analysis by ELISA (test developed at NovImmune) and SELDI-MS.

### Results

- Fractions eluted at 95, 115 and 125 mL (Figure 4A) were analysed on ELISA (not shown) and SELDI-MS (Figure 4B).
- ELISA: Fractions 115 mL and 125 mL positive (higher quantity for 125 mL).
- SELDI-MS analysis: Chemokine in fractions 115 mL and 125 mL but presence of molecules of different molecular weight: 8.0–8.2 kDa in fraction 115; 8.5 kDa in fraction 125.

- PPA HyperCel sorbent can separate chemokine isoforms of minor difference.
- Mixed-mode sorbents can differentiate closely related molecules, which may not be done by ion exchange (IEX) or hydrophobic interaction chromatography (HIC).

## PURIFICATION OF RECOMBINANT F(ab')<sub>2</sub> FRAGMENT

HEA HyperCel sorbent was used as a capture chromatography step to purify a recombinant F(ab')<sub>2</sub> fragment obtained through baculovirus expression in insect cells (SF9).

### Method used

- Load on HEA HyperCel sorbent at pH 6, pH 4 and pH 2.
- Analysis by SDS-PAGE, ELISA and BCA assays.

### Results

- No protein in the flowthrough or elution at pH 6.
- F(ab')<sub>2</sub> fragment eluted at pH 4 while HCPs eluted at pH 2.
- F(ab')<sub>2</sub> fragment recovered at 82% and with a 39-fold purification factor (ELISA analysis, data not shown).




Figure 5. Purification of F(ab')<sub>2</sub> fragment using HEA HyperCel sorbent. The fractions identified on the chromatogram during load and elution were analysed using ELISA (not shown) and SDS-PAGE to identify and quantify the F(ab')<sub>2</sub>.

HEA HyperCel sorbent is very efficient to selectively isolate industrially relevant proteins such as recombinant F(ab')<sub>2</sub> fragments expressed in insect cells, while maintaining a good yield.

## CONCLUSION

- HEA and PPA HyperCel mixed-mode sorbents provide different selectivities compared to existing chromatography techniques and offer new options to purify industrially relevant proteins.
- Proteins and peptides with only minor differences can be discriminated (e.g., isoforms).