

Optimizing elution conditions on Capto MMC using Design of Experiments

Abstract

The multimodal ligand of Capto™ MMC generally requires elution conditions that differ from those used with traditional ion exchangers. This application note describes how Design of Experiments (DoE) can be used to optimize elution conditions and thus recovery of a target protein. The results are summarized as follows:

- The effect of quantitative design factors, such as concentration of eluting salt, pH and ionic strength of the buffer, correlated positively with the recovery of the target protein. The qualitative factor, type of salt, also affected recovery.
- The use of DoE proved to be an efficient way of investigating the effect of several parameters on target protein recovery. Furthermore, DoE revealed a curvature effect that could have been overlooked if multivariate analysis had not been applied.

Introduction

Capto MMC has a ligand with multimodal functionality (Figure 1). This means that different types of interactions, such as ionic interactions, hydrogen bonding and hydrophobic interactions, are possible. The design of the ligand enables binding of proteins at high conductivity, however, it does limit elution using only increasing salt concentration. Thus a different elution protocol has to be used compared to traditional ion exchangers.




Fig 1. The multimodal ligand of Capto MMC.

Generally, Capto MMC requires an increase in both pH and salt concentration for effective elution of bound protein. Additionally, the ionic strength of the buffer and the type of salt will also affect elution. Elution conditions for Capto MMC were investigated by measuring the recovery of bovine serum albumin (BSA) in an elution volume corresponding to 3 column volumes as a function of four different parameters: salt concentration, type of salt, buffer ionic strength and pH.



Material and methods

Frontal analysis

Chromatography system:	ÅKTAexplorer™ 10
Column:	Tricorn™ 5/100
Protein:	BSA
Loading buffer:	50 mM sodium acetate, 250 mM NaCl, pH 4.75 (30 mS/cm)
Loading concentration:	4 mg/ml
Amount loaded:	40 mg BSA
Elution buffer:	see DoE below
Flow velocity:	600 cm/h (1 min residence time)

Design of Experiments

The DoE setup was as follows:

Factors:		Responses:
Salt concentration	0.25, 1.00, 1.75 M	Recovery in elution volumes corresponding to 3 column volumes (CVs)
Salt type	NaCl or NH ₄ Cl	
Buffer ionic strength (BIS)*	0.026 – 0.300 M	
pH [†]	5.75, 6.25, 6.75	

* Calculated from the Henderson-Hasselbalch equation without correction for pKa shifts at different ionic strengths.

[†] The buffers used were sodium malonate (pH 5.75 and pH 6.25) and sodium phosphate (pH 6.75).




Fig 2. Coefficient plot for the model of recovery.

The design was a full factorial design in multiple levels in order to resolve main, interaction and an observed curvature effect [1]. The software used for the design and evaluation was MODDE 7.0 (Umetrics AB).

Results and discussion

All of the tested factors affected recovery of BSA as shown by the coefficient plot (Figure 2).

All of the main factors correlated positively with the recovery of BSA, indicating that an increase in any or all of the factors - pH, salt concentration or buffer ionic strength (BIS) - improves recovery. Although no interaction effects were identified a curvature effect was observed, suggesting that recovery will decrease if the ionic strength of the buffer is above 0.23 M.

In the DoE, the effect of salt type was also investigated. NaCl and NH₄Cl were studied, but the Hofmeister series can be consulted for choice of eluting salt [2]. The response surface plots for the recovery of BSA illustrate the results (Figure 3). Full recovery was obtained under different conditions. The highest pH (6.75) allowed for the greatest variation in salt concentration and buffer ionic strength. Under these conditions NH₄Cl gave better recovery than NaCl.




Fig 3. Response surface plots for the recovery of BSA. Recovery is plotted as a function of salt concentration and buffer ionic strength (BIS) for two different salts. Labels correspond to recovery in percent.

Figure 3 illustrates the optimum values for the selected parameters in the DoE required to achieve full recovery. Buffer ionic strength should be about 0.23 M at pH 6.75, which was the highest pH tested in this study. Salt concentration is dependent on the type of salt used; in this DoE, NH_4Cl gave a better recovery than NaCl at lower salt concentrations. For NH_4Cl , full protein recovery was achieved when the salt concentration ranged from about 0.5 M and upwards and the buffer ionic strength was 0.23 M at pH 6.75. It is possible to obtain full recovery at a lower pH with NH_4Cl , but it should be noted that the salt concentration range narrows with decreasing pH.

The response surface plots also show that full protein recovery is possible when elution is performed with NaCl, but only at the highest pH tested (pH 6.75). However, the salt concentration range is narrower for NaCl than for NH_4Cl . For NaCl, full protein recovery is achievable at a salt concentration of 1.3 M and above when buffer ionic strength is 0.23 M at pH 6.75.

In this DoE, the effects of selected parameters on the recovery of a pure protein were investigated. In a real feed application, other factors such as purity and removal of certain contaminants will need to be considered. As proteins/contaminants have different elution characteristics under different conditions, DoE may be used to optimize the purity of the target protein or removal of certain contaminants.

Conclusions

Both a change in pH and an increase in salt concentration were needed to maximize recovery of the target protein. The lyotropic [2] character of the eluting salt influenced recovery – elution with NH_4Cl gave a higher recovery than elution with NaCl. A high buffer ionic strength also facilitated recovery through a more rapid change of pH during elution and possibly through the lyotropic effect of the buffer salt.

In this study we have shown that DoE is a useful tool for scouting the best elution conditions on Capto MMC. DoE is beneficial in that it can reveal interactions and quadratic effects* that may otherwise go unnoticed. DoE can be preferentially applied for optimizing elution conditions in real feed applications where more responses, e. g., purity and concentration of major impurities are relevant for the application.

The DoE experiments were performed on Tricorn 5/100 columns. Scouting is most conveniently performed on prepacked HiTrap™ columns as they save both sample and time.

* The quadratic effect/term is used to quantify deviations from linearity in the data.

References

1. Box, G. E. P., Hunter, W. G., and Hunter, S. *Statistics for Experimenters*, John Wiley & Sons, Inc., New York (1978).
2. Melander, W., and Horváth, C. Salt effects on hydrophobic interactions in precipitation and chromatography of proteins: an interpretation of the lyotropic series. *Arch. Biochem. Biophys.* **183** 200-215 (1977).

www.gehealthcare.com
www.bioprocess.amershambiosciences.com

GE Healthcare
Amersham Biosciences AB
Björkgatan 30
751 84 Uppsala
Sweden

ÅKTExplorer, Capto, HiTrap, Tricorn are trademarks of GE Healthcare Ltd, a General Electric Company. GE tagline and GE Monogram are trademarks of General Electric Company.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. GE Healthcare reserves the right, subject to any regulatory and contractual approval, if required, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation.

© 2005 General Electric Company – All rights reserved.

Amersham Biosciences AB, a General Electric Company going to market as GE Healthcare.

GE Healthcare Amersham Biosciences AB
Björkgatan 30, 751 84 Uppsala, Sweden

GE Healthcare Amersham Biosciences Europe GmbH
Munzinger Strasse 9, D-79111 Freiburg, Germany

GE Healthcare Amersham Biosciences UK Ltd
Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Amersham Biosciences Corp
800 Centennial Avenue, P.O. Box 1327, Piscataway,
NJ 08855-1327, USA

GE Healthcare Amersham Biosciences KK
Sanken Bldg, 3-25-1, Hyakunincho, Shinjuku-ku,
Tokyo 169-0073, Japan

Asia Pacific Tel: +852 2811 8693 Fax: +852 2811 5251 • Australasia Tel: + 61 2 9899 0999 Fax: +61 2 9899 7511 • Austria Tel: 01/57606-1619 Fax: 01/57606-1627 • Belgium Tel: 0800 73 888 Fax: 03 272 1637 • Canada Tel: 800 463 5800 Fax: 800 567 1008 • Central, East, & South East Europe Tel: +43 1 982 3826 Fax: +43 1 985 8327 • Denmark Tel: 45 16 2400 Fax: 45 16 2424 • Finland & Baltics Tel: +358 (0)9 512 39 40 Fax: +358 (0)9 512 39 439 • France Tel: 01 69 35 67 00 Fax: 01 69 41 96 77 • Germany Tel: 0761/4903-490 Fax: 0761/4903-405 • Italy Tel: 02 27322 1 Fax: 02 27302 212 • Japan Tel: +81 3 5331 9336 Fax: +81 3 5331 9370 • Latin America Tel: +55 11 3933 7300 Fax: +55 11 3933 7304 • Middle East & Africa Tel: +30 210 9600 687 Fax: +30 210 9600 693 • Netherlands Tel: 0165 580 410 Fax: 0165 580 401 • Norway Tel: 815 65 555 Fax: 815 65 666 • Portugal Tel: 21 417 7035 Fax: 21 417 3184 • Russia & other C.I.S. & N.I.S Tel: +7 (095) 232 0250 Fax: +7 (095) 230 6377 • South East Asia Tel: 60 3 8024 2080 Fax: 60 3 8024 2090 • Spain Tel: 93 594 49 50 Fax: 93 594 49 55 • Sweden Tel: 018 612 1900 Fax: 018 612 1910 • Switzerland Tel: 0848 8028 12 Fax: 0848 8028 13 • UK Tel: 0800 616928 Fax: 0800 616927 • USA Tel: 800 526 3593 Fax: 877 295 8102

