

Capto Q

Capto™ Q is a strong anion exchange medium for packed bed chromatography that allows increased speed and throughput in capture and intermediate purification (Fig 1). It combines high capacity with high flow velocity and low backpressure, at the same time reducing process cycle times and increasing productivity. Capto Q is a BioProcess™ medium that meets the demands of large-scale biopharmaceutical manufacturers for fast, efficient and cost-effective protein purification.

- Raised productivity with high dynamic binding capacity at high flow
- Increased yield with rapid mass transfer
- Reduced process time with high volume throughput
- Cost-effective processing with smaller unit operations

Media characteristics

High throughput in downstream purification requires separation media that combine mechanical strength of the matrix with a pore structure that allows fast mass transfer of target molecules. Capto Q is based on a highly rigid, chemically modified agarose. The production methodology offers outstanding pressure/flow properties with full control over the pore structure. Capto Q represents the first product in a line of new media for general-use high-flow, low-pressure solutions for bioprocessing. Table 1 summarizes the basic characteristics of Capto Q.

High flow rates and low backpressure in large-scale columns

High flow velocities increase the productivity of large-scale bioprocessing operations and allow processing of large volumes in one working shift. High flow velocities also reduce exposure of the target protein to proteases. Due to its mechanical strength and the low backpressure generated with Capto Q, columns with this medium can be



Figure 1. Capto Q increases throughput and productivity in large-scale bioprocessing operations

operated at high flow velocities with a wide range of bed heights also in large scale. Typical flow velocities for Capto Q in a 1 m diameter column with 20 cm bed height are over 700 cm/h, with a backpressure below 3 bar.

It is important to prove pressure/flow performance in larger-scale columns where no wall support is present. Figure 2 compares pressure/flow curves in a BPG™ 300 column for the Capto base matrix with that of Sepharose™ 6 Fast Flow. Although the bead and pore sizes are similar between the two matrices, the pressure/flow properties of Capto Q are significantly better. This is a result of the exceptional mechanical stability of the Capto base matrix.

Strong anion exchanger with fast mass transfer and high dynamic binding capacity

The functional group on Capto Q is a quaternary amine (Q) group for anion-exchange. This group is linked to the matrix via a dextran surface extender.



Table 1. Characteristics of Capto Q

| | |
|---------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Matrix | highly cross-linked agarose with dextran surface extender |
| Ion exchange type | strong anion, Q |
| Charged group | -N ⁺ (CH ₃) ₃ |
| Total ionic capacity | 0.16-0.22 mmol Cl ⁻ /ml medium |
| Particle size* | 90 µm (d _{50v}) |
| Flow velocity | at least 700 cm/h in columns up to 1 m in diameter with 20 cm bed height at 20 °C using process buffers with the same viscosity as water at < 3 bar (0.3 MPa) |
| Dynamic binding capacity [†] | > 100 mg BSA/ml medium |
| pH stability [‡] | |
| short term | 2 – 14 |
| working | 2 – 12 |
| long term | 2 – 12 |
| Working temperature | +4 to +30 °C |
| Chemical stability | all commonly used aqueous buffers 1 M acetic acid, 1 M NaOH [§] , 8 M urea, 6 M guanidine hydrochloride, and 70% ethanol |
| Avoid | oxidizing agents, anionic detergents |

* d_{50v} is the median particle size of the cumulative volume distribution.

[†] Dynamic binding capacity at 10% breakthrough as measured at a residence time of 1 minute, 600 cm/h in a Tricorn™ 5/100 column with 10 cm bed height in a 50 mM Tris-HCl buffer, pH 8.0.

[‡] Short term pH: pH interval that the medium can be subjected to, for cleaning- or sanitization-in-place (accumulated 90–300 h at room temperature) without significant change in function.

Working pH: pH interval where the medium binds protein as intended or is needed for elution, without adverse long-term effect.

Long term pH: pH interval where the medium can be operated without significant change in function.

[§] No significant change in ionic binding capacity and carbon content after 1 week storage in 1 M NaOH at 40°C.




Figure 2. Pressure-flow properties of Capto Q compared to Sepharose 6 Fast Flow. Running conditions: BPG 300 (30 cm I.D.), open bed at settled bed height equal to 23 cm, with water at 20 °C.

The agarose matrix has an open pore structure that allows fast mass transfer and gives high dynamic binding capacity. The surface extender further increases capacity and mass transfer, resulting in high dynamic binding capacity also at high flow velocities. High dynamic binding capacity is important for obtaining high yield of the target protein. It also contributes to shortening the overall processing time as the total number of cycles may be reduced. Dynamic binding capacities at different residence times are shown in figure 3.




Figure 3. The agarose matrix has an open pore structure that allows fast mass transfer. Dynamic binding capacity for bovine serum albumin (BSA) is shown as a function of residence time. For Q Sepharose Fast Flow, residence times below 2 minutes can not be achieved in large-scale columns.

Agarose assures protein friendliness and high chemical stability

Agarose is a natural polymer with hydrophilic properties that minimizes structural changes of the target molecule and unspecific adsorption to the matrix. It is compatible to the broad range of buffer systems commonly used in bioprocess chromatography for protein purification. Therefore optimized conditions for binding and elution can be applied, assuring mild treatment of the target molecules. The chemical stability of agarose prolongs media lifetime even if harsh cleaning-in-place procedures are applied. Characteristics such as capacity, elution behavior and pressure/flow properties are not affected by common cleaning solutions in process chromatography. These properties make agarose media well accepted by regulatory bodies and have resulted in proven track records in biotherapeutic manufacturing.

Rigid medium for cost-effective purification

The rigidity of the Capto Q base matrix translates into improved process economics. Its characteristics allow a wider working range of flow velocities, bed heights and sample viscosities, all of which affect processing costs in a positive way. High flow velocities increase volume throughput and reduce process time, longer bed heights means smaller equipment

and reduced foot-print, and high-flow processing of viscous samples means less dilution and shorter cycle times.

The available degree of freedom in process design for a medium can be illustrated as its “window of operation”. Figure 4 shows schematically the ranges for key operating variables for Capto Q and Sepharose 6 Fast Flow. Given a maximum allowed pressure, it predicts the allowable combinations of column bed heights and operating velocities. The pressure limits, shown as blue and red curves, are based on a 1 meter diameter column with 20 cm bed height and maximum operating velocities of 700 and 250 cm/h, respectively. At this point, the pressure is 3 bar for Capto Q and 1 bar for Sepharose 6 Fast Flow. For Sepharose 6 Fast Flow, 1 bar represents the highest recommended pressure. For Capto Q, 3 bar corresponds to the maximum pressure for many low-pressure systems; the medium as such can normally be run to the maximum pressure rating of low and medium pressure columns.




Figure 4. The highly rigid Capto base matrix allows a much larger window of operation (area below the curves) at large-scale than Sepharose 6 Fast Flow. This is particularly true at bed heights of 25-30 cm and above. Data correspond to a 1 m diameter column, 20 cm bed height, 20°C and viscosity of water. Red and blue curves correspond to pressure limits of 1 and 3 bar, respectively. Green contours give the residence time in the column in minutes.

The size of the area below the pressure limit curves represents the window of operation, or the available operating range for the respective medium. As shown in Figure 4, this is significantly larger for Capto Q than for Sepharose 6 Fast Flow based media, especially when bed heights increase to 25-30 cm or more. At these bed heights, Capto Q can still be run at flow velocities of 300-400 cm/h or more. Thus, the high mechanical stability of Capto Q allows practical and cost-effective use of smaller diameter columns.

A large window of operation also allows flexibility even if the viscosity of the feed is high. Doubling viscosity halves the operational velocity. For a feedstock with a viscosity of 2 cP at a bed height of 30 cm, the flow velocity of Capto Q is 235 cm/h compared to 80 cm/h for Sepharose 6 Fast Flow.

Figure 4 also shows contours of the residence time in the column. Residence time is important because it affects dynamic binding capacity; the longer target molecules reside in the bed, the more of them adsorb to the medium. It is possible to increase the residence time by either decreasing the flow velocity, or increasing the column bed height. For Capto Q, increasing bed heights assures longer residence times also when high flow velocities are applied.

HiTrap columns for fast screening

Time and sample can be saved in the early stages of development by using small, prepacked HiTrap™ columns to screen for the most suitable media and to develop the basic separation method (Fig 5). In combination with an automated system such as ÄKTAexplorer™ it is easy and fast to scout for optimum selectivity, pH, ionic strength etc. Capto Q is available as HiTrap prepacked 1 ml and 5 ml columns. Basic characteristics of HiTrap Capto Q columns are summarized in Table 2.

Table 2. Characteristics of HiTrap columns

| | |
|------------------------|--------------------------------------------------|
| Column volumes | 1 ml and 5 ml |
| Column dimensions | 0.7 x 2.5 cm (1 ml); 1.6 x 2.5 cm (5 ml) |
| Maximum flow rates | HiTrap 1 ml: 4 ml/min, HiTrap 5 ml: 20 ml/min |
| Recommended flow rates | HiTrap 1 ml: 1ml/min, HiTrap 5 ml: 5 ml/min |
| Maximum backpressure | 0.3 MPa, (3 bar) |

Application

Recent developments in upstream processing result in larger feed volumes and increased protein expression levels. The combination of high volume throughput and high capacity makes Capto Q the optimal choice for processing large amounts of protein in a fast and efficient way. As a strong anion exchanger, its behavior can be easily controlled and its application area predicted by buffer choice and pI of the target protein.

The Capto Q ligand is a quaternary amine identical to the Q ligand used in other media. Despite this, minor differences in selectivity can occur between the media as illustrated in Fig 5.

Improved productivity based on high flow features

Scale-up modeling and productivity calculations based on experimental data at small and pilot scale indicate that it is possible to capture and recover 103 kg of Green Fluorescent Protein (GFP) from an *Escherichia coli* homogenate in 24 h by using Capto Q in a 1.6 m i.d. column at 20 cm bed height. Assuming the same process cycle conditions, the corresponding amount for Q Sepharose Fast Flow is 30 kg. This example supports the argument that Capto Q as being particularly suitable for high throughput and high productivity capture purification. For further details see Application Note 11-0026-20.

Column: HiTrap, 1 ml
 Sample: GFP in *E. coli* homogenate
 Buffer A: 50 mM Tris/HCl, pH 8.2
 Buffer B: 50 mM Tris/HCl, pH 8.2 + 1 M NaCl
 Flow: 1 ml/min (156 cm/h)
 Gradient: 0 – 100% B, 15 column volumes
 System: ÄKTAexplorer 100




Figure 5. Five different anion exchange media available in prepacked 1 ml HiTrap columns were fast and easily screened on ÄKTAexplorer. The selectivity and separation ability for GFP, expressed in *E. coli*, on the different media were compared. GFP absorbs specifically at 490 nm and thus can be easily assayed for throughout purification. Based on the combination of the UV curves 280 nm and 490 nm Capto Q gave optimal selectivity.

Another example of how the high flow agarose technology improves throughput in capture processes is demonstrated by MabSelect™ - a protein A medium that has allowed a five-fold increase in process throughput for capture of monoclonal antibodies. The Capto Q matrix is manufactured using the same technology as MabSelect. However, the Capto matrix is even more rigid than the MabSelect matrix, to allow even higher flow velocities.

Operation

Fast method development

Fast method development can easily be achieved on ÄKTAexplorer with quick initial screening on selectivity and process conditions using prepacked HiTrap columns. Further optimization and method development using Tricorn or XK columns allows straightforward scale-up. The UNICORN™ software on ÄKTAdesign™ systems makes it simple to transfer the optimized method to a production scale process system. For more information about method development and optimization, consult the handbook, "Ion Exchange Chromatography & Chromatofocusing: Principles and Methods", (11-0004-21).

Fully scalable

Capto Q belongs to the BioProcess range of media that are developed and supported for production-scale chromatography. This includes validated manufacturing methods, secure supply and Regulatory Support Files (RSF) to assist process validation and submission to regulatory authorities.

Scale-up from laboratory purification to pilot and production scales is performed by keeping the residence time constant. The residence time is equal to the bed height (cm) divided by the flow velocity (cm/h) applied during sample load. Thus it is possible to scale-up to higher bed heights, provided that the different scales are operated with the same residence time.


This flexibility during scale-up/scale-down is a clear advantage, as lab-scale experiments are often performed with short bed heights to reduce time and sample volume. To utilize the full potential of Capto Q, we recommend bed heights of 20 cm and higher in large scale.

Figure 6 shows scale-up of capture purification of recombinant GFP, expressed in *E. coli* on Capto Q. Process development and optimization was performed in a Tricorn 5/100 column with 10 cm bed height. The process was scaled-up 400 times to pilot scale in a FineLINE™ 70 column with 20 cm bed height. The two-fold increase in bed height was accompanied by a two-fold increase in flow velocity, keeping the residence time constant between the scales. GFP was eluted in ~2.5 CV in FineLINE-scale, where the recovery exceeded 90%.

Cleaning and sanitization

Cleaning-in-place (CIP) is a cleaning procedure that removes contaminants such as lipids, precipitates, or denatured proteins that may remain in the packed column after regeneration. Regular CIP also prevents the build-up of these contaminants in the media bed and helps to maintain the capacity, flow properties and general performance of Capto Q.

A Column: Tricorn 5/100, 10 cm bed height (column volume 2 ml)
 Medium: Capto Q
 Sample: GFP in *E. coli* homogenate, 12 ml
 Buffer A: 50 mM Tris/HCl, pH 8.2
 Buffer B: 50 mM Tris/HCl, pH 8.2 + 1 M NaCl
 Flow rate: 300 cm/h
 Gradient: 10% B (6 CV), 33% B (6 CV), 100% B (4 CV)
 System: ÄKTAexplorer 100



B Column: FineLINE 70, 20 cm bed height (column volume 808 ml)
 Medium: Capto Q
 Sample: GFP in *E. coli* homogenate, 4835 ml
 Buffer A: 50 mM Tris/HCl, pH 8.2
 Buffer B: 50 mM Tris/HCl, pH 8.2 + 1 M NaCl
 Flow rate: 600 cm/h
 Gradient: 10% B (6 CV), 33% B (6 CV), 100% B (4 CV)
 System: ÄKTApilot™




Figure 6. Straightforward scale-up between A. Tricorn 5/100 and B. FineLINE 70 (400 times).

A specific CIP protocol should be designed for each process according to the type of contaminants present. The frequency of CIP depends of the nature and the condition of the starting material, but one CIP cycle is generally recommended every 1–5 separation cycles.

Capto Q withstands all standard CIP procedures, e. g. 1 M NaOH, 2 M NaCl or 70% ethanol.

Equipment

Capto Q can be used together with most equipment available for chromatography from lab scale to production scale. Due to the high rigidity of the medium, packing procedures are slightly different compared to Sepharose 6 Fast Flow based media. For details on packing lab scale columns, see Instruction manual, for process scale columns see Application note 11-0026-21. Appropriate columns from GE Healthcare are shown in Table 3.

Table 3: Appropriate columns

| Column family | Range (bed diameters) |
|-----------------|-----------------------|
| Tricorn | 5 mm, 10 mm |
| XK | 16 mm, 26 mm |
| FineLINE | 35 mm – 350 mm |
| BPG | 100 mm – 300 mm* |
| BioProcess LPLC | 100 mm – 1200 mm |
| Chromaflow™ | 400 mm – 2000 mm |

* The pressure rating of BPG 450 is too low to use it with Capto Q.

Storage

Store unused media and prepacked columns at +4 to +30°C in 20% in ethanol.

Ordering information

Capto Q is also available as prepacked column in the HiTrap format.

| Product | Pack size | Code No |
|----------------|-----------|------------|
| Capto Q | 25 ml | 17-5316-10 |
| Capto Q | 500 ml | 17-5316-01 |
| Capto Q | 5 L | 17-5316-04 |
| Capto Q | 10 L | 17-5316-05 |
| Capto Q | 60 L | 17-5316-60 |
| HiTrap Capto Q | 5x1ml | 11-0013-02 |
| HiTrap Capto Q | 5x5ml | 11-0013-03 |

All bulk media and prepacked column products are supplied in 20% ethanol. For more information, contact your local GE Healthcare representative.

Literature

| | |
|----------------------------------------------------------------------------|------------|
| Ion Exchange Chromatography & Chromatofocusing: Principles and Methods | 11-0004-21 |
| Capture of Green Fluorescent Protein on Capto Q – application note | 11-0026-20 |
| Methods for packing Capto Q in production-scale columns - application note | 11-0026-21 |

Further information

For the latest news, product information or handbooks, visit:

www.amershambiosciences.com

www.chromatography.amershambiosciences.com

www.bioprocess.amershambiosciences.com

or contact your local GE Healthcare representative.

www.gehealthcare.com

GE Healthcare

Björkgatan 30

751 84 Uppsala

Sweden

ÄKTAdesign, ÄKTaexplorer, ÄKTApilot, BioProcess, BPG, Capto, Chromaflow, Dropdesign, FineLINE, HiTrap, MabSelect, Sepharose, Tricorn, and UNICORN are trademarks of Amersham plc, a General Electric company going to market as GE Healthcare. GE tagline and GE Monogram are trademarks of General Electric Company.

Separating viral particles with this product may require a license under United States pat 6,537,793 B2 and foreign equivalents owned by Centelion SAS.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. GE Healthcare reserves the right, subject to any regulatory and contractual approval, if required, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation.

© 2005 General Electric Company – All rights reserved.

Amersham Biosciences AB, a General Electric Company going to market as GE Healthcare.

GE Healthcare Amersham Biosciences AB
Björkgatan 30, 751 84 Uppsala, Sweden

GE Healthcare Amersham Biosciences Europe GmbH
Munzinger Strasse 9, D-79111 Freiburg, Germany

GE Healthcare Amersham Biosciences UK Ltd
Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Amersham Biosciences Corp
800 Centennial Avenue, P.O. Box 1327, Piscataway,
NJ 08855-1327, USA

GE Healthcare Amersham Biosciences KK
Sanken Bldg. 3-25-1, Hyakunincho, Shinjuku-ku,
Tokyo 169-0073, Japan

Asia Pacific Tel: +852 2811 8693 Fax: +852 2811 5251 • Australasia Tel: +61 2 9899 0999 Fax: +61 2 9899 7511 • Austria Tel: 01/57606-1619 Fax: 01/57606-1627 • Belgium Tel: 0800 73 888 Fax: 03 272 1637 • Canada Tel: 800 463 5800 Fax: 800 567 1008 • Central, East, & South East Europe Tel: +43 1 982 3826 Fax: +43 1 985 8327 • Denmark Tel: 45 16 2400 Fax: 45 16 2424 • Finland & Baltics Tel: +358-(0)9-512 39 40 Fax: +358 (0)9 512 39 439 • France Tel: 01 69 35 67 00 Fax: 01 69 41 96 77 • Germany Tel: 0761/4903-490 Fax: 0761/4903-405 • Italy Tel: 02 27322 1 Fax: 02 27302 212 • Japan Tel: +81 3 5331 9336 Fax: +81 3 5331 9370 • Latin America Tel: +55 11 3933 7300 Fax: +55 11 3933 7304 • Middle East & Africa Tel: +30 210 9600 687 Fax: +30 210 9600 693 • Netherlands Tel: 0165 580 410 Fax: 0165 580 401 • Norway Tel: 815 65 555 Fax: 815 65 666 • Portugal Tel: 21 417 7035 Fax: 21 417 3184 • Russia & other C.I.S. & N.I.S Tel: +7 (095) 232 0250, 956 1137 Fax: +7 (095) 230 6377 • South East Asia Tel: 60 3 8024 2080 Fax: 60 3 8024 2090 • Spain Tel: 93 594 49 50 Fax: 93 594 49 55 • Sweden Tel: 018 612 1900 Fax: 018 612 1910 • Switzerland Tel: 0848 8028 12 Fax: 0848 8028 13 • UK Tel: 0800 616928 Fax: 0800 616927 • USA Tel: 800 526 3593 Fax: 877 295 8102



imagination at work