

UltraLink™ Immobilized Jacalin

3747 N. Meridian Road
P.O. Box 117
Rockford, IL 61105

53145

0510

Introduction

Jacalin is an α -D-galactose binding lectin extracted from jack-fruit seeds (*Artocarpus integrifolia*).¹ This lectin is a glycoprotein of approximately 40,000 MW composed of four identical subunits. Jacalin immobilized on supports such as agarose has been useful for the purification of human serum or secretory IgA.²⁻¹² IgA can be separated from human IgG and IgM in human serum or colostrum.² IgD is reported to bind to jacalin.⁴ This support is also useful for removing contaminating IgA from IgG samples. Immobilized jacalin is preferable to protein A because protein A will bind some IgA along with IgG and is not effective for IgA removal.³

UltraLink™ Immobilized Jacalin is jacalin on the azlactone-activated support 3M Emphaze™ Biosupport Medium AB 1. This activated matrix has been coupled to a variety of ligands. 3M Emphaze™ Biosupport Medium AB 1 is a hydrophilic, high capacity, highly crosslinked, copolymeric, porous support. The reactive groups couple to amines yielding a stable bond and minimal ligand leakage. UltraLink™ supports are more sturdy than agarose allowing for their use in FPLC and large-scale applications.

Characteristics of 3M Emphaze™ Biosupport Medium AB 1

Capacity:	See label for specific capacity
pH stability of matrix:	1-13
Particle Size (average)	50-80 microns (typically 60 μ m)
Exclusion limit (proteins):	>2,000,000 daltons
Surface area (average):	>250 m ² /g of beads
Pore volume (average):	>1.2 ml/g of beads (>60% of bead volume)
Pore size:	1000 Å
Maximum pressure:	100 psi (6.9 bar)*
Maximum linear velocity:	3000 cm/hour

***NOTE:** The indicated maximum pressure of 100 p.s.i. refers to the maximum pressure drop across a column the support can withstand. It does not necessarily refer to the indicated system pressure shown on a liquid chromatography apparatus, since the system may not actually be measuring the pressure drop across the column.

Product Description

NUMBER

53145

DESCRIPTION

UltraLink™ Immobilized Jacalin, 5 ml gel

See product label for lot-specific binding capacity.

Loading: Approximately 2-3 mg Jacalin per ml of gel.

Store at 4°C. Do not freeze.

Telephone 800-8-PIERCE or 815-968-0747

Fax 815-968-7316 or 800-842-5007

Sample Protocol for Medium Pressure Chromatography

UltraLink™ supports are extremely valuable for use in medium pressure chromatography applications. When packed into a 3 mm (diameter) x 14 cm (height) glass column, UltraLink™ supports have been run up to approximately 400 p.s.i. (system pressure) with no visual compression of the gel or adverse effects on the chromatography. Typically these columns were run at linear velocities of 85-3000 cm/hour with excellent separation characteristics. Purification time can be reduced from 8-12 hours to several minutes when an FPLC procedure is used instead of a gravity flow method.

Medium Pressure Chromatographic Applications

The 3M Emphaze™ Biosupport Medium can be used in Fast Performance Liquid Chromatography (FPLC) applications. With MPC (medium pressure chromatographic) applications the critical factor for success is limiting the pressure drop across the column that the support experiences. Due to variations with the instrument, the indicated gauge pressure of an MPC apparatus may not actually measure the pressure drop across the column.

A more universal criteria for MPC applications, therefore, would be to measure the linear velocity of the buffers through the column, a pressure-independent measurement. The linear velocity is defined as the velocity of the buffer front passing through the gel bed. Linear velocity is usually expressed in cm/hour, and in the case of the 3M Emphaze™ Biosupport Medium, the maximum linear velocity that the support can tolerate in a 1 cm diameter by 10 cm high bed is 3000 cm/hr.

The linear flow rate through a cylindrical column can be calculated if the inside diameter of the column is known and the column effluent is collected and measured for a given time to determine the flow rate. The calculations for determining linear velocity are shown on the following page:

Radius measured in cm = r

Column cross-sectional area in square centimeters = πr^2

1 ml = 1 cubic centimeter

Measured flow rate (ml/minute) = Column effluent collected in ml/minute

$$\text{Linear velocity per minute} = \frac{(\text{measured flow rate in ml/minute or cm}^3/\text{minute})}{\text{cross-sectional area in cm}^2}$$

Linear velocity per hour = (linear velocity per minute)(60 minutes/hour)

$$\text{Measured flow rate per hour} = \frac{(\text{cm}^3/\text{minute})(60 \text{ minute/hour})}{\pi r^2} = \frac{\text{Linear velocity expressed as cm/hour}}{\text{Column cross-sectional area}}$$

Materials

- A. UltraLink™ Immobilized Jacalin, 5 ml
- B. Human serum
- C. Binding buffer: Phosphate buffered saline (PBS), 0.1 M sodium phosphate, 150 mM NaCl, pH 7.4
- D. Elution buffer: 0.1 M melibiose or 0.1 M α -D-galactose in PBS
- E. Column

Method

1. Gravity pack under flow UltraLink™ Jacalin into a medium pressure chromatography column.
2. Store the column in 0.02% sodium azide if column will not be used immediately.
3. Equilibrate the column with at least 5 column volumes of binding buffer.
4. Add human serum to binding buffer in equal volumes, centrifuge if necessary and add or inject supernatant onto column.
5. Wash the column with an additional 8 column volumes of binding buffer or until baseline is reached using absorbance at 280 nm to measure effluent.

Telephone 800-8-PIERCE or 815-968-0747

Fax 815-968-7316 or 800-842-5007

6. Elute the bound IgA1 using 2 ml aliquots of elution buffer until baseline absorbance is reached. Collect the eluate from the aliquots in separate test tubes and monitor absorbance at 280 nm to determine baseline.
7. Buffer exchange the sample into binding buffer using a gel filtration/desalting column to remove the melibiose or galactose.
8. Regenerate the column by washing with 20 column volumes of PBS.

References

1. Kumar, G.S., Appukuttan, P.S. and Basu, D. (1982). α -D-Galctose-specific lectin from jack fruit (*Artocarpus integra*) seed. *J. Biosci.* **4(3)**, 257-261.
2. Roque-Barreira, M.C. and Campos-Neto, A. (1985). Jacalin, an IgA-binding lectin. *J. Immunol.* **134(30)**, 1740-1743.
3. Van Kamp, G.J. (1979). IgA contamination of IgG prepared on a protein A column. *J. Immunol. Meth.* **27**, 301-305.
4. Aucouturier, P., Mihaesco, E., Mihaesco, C. and Preud'Homme, J.-L. (1987). Characterization of jacalin, the human IgA and IgD binding lectin from jackfruit. *Mol. Immunol.* **24(5)**, 503-511.
5. Biewenga, J., Steneker, I. and Hameleers, D.M.H. (1988). Effect of serum albumin on the recovery of human IgA1 from immobilized jacalin. *J. Immunol. Meth.* **115**, 199-207.
6. Bunn-Moreno, M.M. and Campos-Neto, A. (1981). Lectins(s) extracted from seeds of artocarpus integrifolia (Jackfruit): potent and selective stimulators(s) of distinct human T and B cell functions. *J. Immunol.* **127(2)**, 427-429.
7. Gregory, R.L., Rundegren, J. and Arnold, R.R. (1987). Separation of human IgA1 and IgA2 using jacalin-agarose chromatography. *J. Immunol. Meth.* **99**, 101-106.
8. Jarvis, G.A. and Griffiss, J.M. (1989). Human IgA1 initiates complement-mediated killing of *Neisseria Meningitidis*. *J. Immunol.* **143(5)**, 1703-1709.
9. MacDermott, R.P., Nash, G.S., Bertovich, M.J., Mohrman, R.F., Kodner, I.J., Delacroix, D.L. and Vaerman, J.-P. (1986). Altered patterns of secretion of monomeric IgA and IgA subclass 1 by intestinal mononuclear cells in inflammatory bowel disease. *Gastroent.* **91(2)**, 379-385.
10. Saxon, A., Tsui, F. and Martinez-Maza, O. (1987). Jacalin, an IgA-binding lectin, inhibits differentiation of human B cells by both a direct effect and by activating T-suppressor cells. *Cell. Immunol.* **104**, 134-141.
11. Kondoh, H., Kobayashi, K., Hagiwara, K. and Kajii, T. (1986). Jacalin, a jackfruit lectin, precipitates IgA1 but not IgA2 subclass on gel diffusion reaction. *J. Immunol. Meth.* **88**, 171-173.
12. Wood, G.M., Trejdosiewicz, L.K. and Losowsky, M.S. (1987). ELISA for measurement of secretory IgA distinct from monomeric IgA. *J. Immunol. Meth.* **97**, 269-274.
13. Coleman, P.L., Walker, M.M., Milbrath, D.S. and Stauffer, D.S. (1990). Immobilization of protein A at high density on azlactone-functional polymeric beads and their use in affinity chromatography. *J. Chrom.* **512**, 345-363.
14. Coleman, P.L., Walker, M.M., Heilmann, S.M., Krepski, L.R., Rasmussen, J.K. and Jensen, K.M. (1988). Affinity chromatography on a novel support: azlactone-acrylamide copolymer beads. *FASEB J.* **2: A1770** (#8563).
15. Coleman, P.L., Milbrath, D.S., Walker, M.M., Heilmann, S.M., Rasmussen, J.K. and Krepski, L.R. (1990). Azlactone copolymer beads: applications in bioseparations. *J. Cell. Biochem.* **44**, 19 (S14D).
16. Coleman, P.L., Walker, M.M., Reese, C.L. and Milbrath, D.S. (1991). Effect of polyanionic salts on immobilization of protein A and antibody on azlactone-functional beads. *FASEB J.* **A805** (#2528).
17. Hermanson, G.T., Mallia, A.K. and Smith, P.K. in *Immobilized Affinity Ligand Techniques*. 1992, *Academic Press*, 28-30, 90-95.
18. Milbrath, D.S., Coleman, P.L., Walker, M.M. and Stauffer, D.S. (1990). Azlactone-functional supports useful in affinity chromatography and other bioseparations. *AIChE Extended Abstracts* #104E.
19. Milbrath, D.S., Coleman, P.L., Walker, M.M., Heilmann, S.M., Rasmussen, J.K. and Krepski, L.R. (1989). Azlactone polymer supports for bioseparations. *ACS Abstracts*.
20. Rasmussen, J.K., Heilmann, S.M., Krepski, L.R., Jensen, K.M., Mickelson, J. and Johnson, K. (1991/1992). Crosslinked, hydrophilic, azlactone-functional polymeric beads: a two-step approach. *Reactive Polymers* **16**, 199-212.
21. Rasmussen, J.K., Hembre, J.I., Koski, N.I., Milbrath, D.S., *et al.* (1992). Mechanistic studies in reverse-phase suspension copolymerization of vinyl dimethylazlactone methylenebis (acrylamide). *Makromol. Chem., Macromol. Symp.* **54/55**, 535-550.
22. Rasmussen, J.K., Heilmann, S.M., Krepski, L.R., Smith II, H.K., *et al.* (1990). Hydrophilic, crosslinked, azlactone-functional beads- a new reactive support. *Polymer Reprints* **31(2)**, 442-443.

3M Emphaze Biosupport Medium AB 1 is covered by US Patent No. 4,871,824 (Heilmann *et al*) and European Patent Publication O 392,735 A2 Emphaze is a trademark of 3M Corporation.

©Copyright Pierce Chemical Company, 1995. Printed in U.S.A.

Telephone 800-8-PIERCE or 815-968-0747

Fax 815-968-7316 or 800-842-5007