

July 2005

---

# **Qproteome™ Mammalian Protein Preparation Handbook**

For preparation of total protein from  
mammalian cells



---

[WWW.QIAGEN.COM](http://WWW.QIAGEN.COM)

Trademarks: QIAGEN®, Qproteome™, (QIAGEN Group) Benzonase® (Merck KGaA, Germany), Coomassie® (ICI [Imperial Chemical Industries] Organics Inc.)

Benzonase endonuclease is covered by several patents, i.e. the US patent No. 5.173,418 and EP No. 0229866. Nycomed Pharma A/S, Denmark claims worldwide patent rights to Benzonase endonuclease which are licensed exclusively to MERCK KGaA, Darmstadt (Germany).

© QIAGEN 2005. All rights reserved.

# Contents

<b>Kit Contents</b>	<b>4</b>
<b>Storage</b>	<b>4</b>
<b>Product Use Limitations</b>	<b>4</b>
<b>Product Warranty and Satisfaction Guarantee</b>	<b>4</b>
<b>Technical Assistance</b>	<b>5</b>
<b>Quality Control</b>	<b>5</b>
<b>Safety Information</b>	<b>5</b>
<b>Introduction</b>	<b>6</b>
<b>Protocols</b>	<b>7</b>
■ <b>Isolation of Total Proteins from Mammalian Cells</b>	<b>7</b>
■ <b>Acetone Precipitation of Protein Fractions</b>	<b>8</b>
<b>Troubleshooting Guide</b>	<b>8</b>
<b>Ordering Information</b>	<b>9</b>
<b>QIAGEN Distributors</b>	<b>11</b>

## Kit Contents

<b>Qproteome Mammalian Protein Prep Kit</b>	
<b>Catalog no.</b>	<b>37901</b>
<b>Number of preps</b>	<b>100</b>
Mammalian Cell Lysis Buffer	100 ml
Benzonase®	2000 Units (25 U/ $\mu$ l)
Protease Inhibitor Solution (100x)	4 x 300 $\mu$ l
Handbook	1

## Storage

All kit components should be stored at 2–8°C. For longer storage, Benzonase® should be stored at –20°C.

## Product Use Limitations

Qproteome Kits are developed, designed, and sold for research purposes only. They are not to be used for human diagnostic or drug purposes or to be administered to humans unless expressly cleared for that purpose by the Food and Drug Administration in the USA or the appropriate regulatory authorities in the country of use. All due care and attention should be exercised in the handling of many of the materials described in this text.

## Product Warranty and Satisfaction Guarantee

QIAGEN guarantees the performance of all products in the manner described in our product literature. The purchaser must determine the suitability of the product for its particular use. Should any product fail to perform satisfactorily due to any reason other than misuse, QIAGEN will replace it free of charge or refund the purchase price. We reserve the right to change, alter, or modify any product to enhance its performance and design. If a QIAGEN product does not meet your expectations, simply call your local Technical Service Department or distributor. We will credit your account or exchange the product — as you wish.

A copy of QIAGEN terms and conditions can be obtained on request, and is also provided on the back of our invoices. If you have questions about product specifications or performance, please call QIAGEN Technical Services or your local distributor (see inside back cover).

## Technical Assistance

At QIAGEN we pride ourselves on the quality and availability of our technical support. Our Technical Service Departments are staffed by experienced scientists with extensive practical and theoretical expertise in molecular biology and the use of QIAGEN® products. If you have any questions or experience any difficulties regarding Qproteome Kits or QIAGEN products in general, please do not hesitate to contact us.

QIAGEN customers are a major source of information regarding advanced or specialized uses of our products. This information is helpful to other scientists as well as to the researchers at QIAGEN. We therefore encourage you to contact us if you have any suggestions about product performance or new applications and techniques.

For technical assistance and more information please call one of the QIAGEN Technical Service Departments or local distributors (see inside back cover).

## Quality Control

In accordance with QIAGEN's ISO-certified Total Quality Management System, each lot of Qproteome Kits is tested against predetermined specifications to ensure consistent product quality.

## Safety Information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate material safety data sheets (MSDSs). These are available online in convenient and compact PDF format at [www.qiagen.com/ts/msds.asp](http://www.qiagen.com/ts/msds.asp) where you can find, view, and print the MSDS for each QIAGEN kit and kit component.

### 24-hour emergency information

Emergency medical information in English, French, and German can be obtained 24 hours a day from:

Poison Information Center Mainz, Germany

Tel: +49-6131-19240

## Introduction

In order to be able to efficiently purify and analyze proteins, they must first be released from their host cell in a soluble form. Therefore, the first step in protein purification and analysis is cell lysis. Mammalian cells are bounded by a plasma membrane, which is composed of phospholipids and proteins. In contrast to bacterial or yeast cells, this membrane can be easily disrupted, but any lysis method should be gentle enough to deliver active, native-conformation proteins.

The Qproteome Mammalian Protein Prep Kit provides gentle but efficient detergent-based lysis of mammalian cells to deliver a total protein fraction suitable for any downstream application. Cell pellets are incubated with Mammalian Cell Lysis Buffer that is supplemented with Benzonase<sup>®</sup> — an enzyme that digests nucleic acids — and a protease inhibitor, which protect the released proteins from enzymatic degradation. No specialized equipment is required. Cells are simply incubated with lysis buffer and centrifuged to separate cell debris. The supernatant contains the total protein fraction, which can be further purified or directly analyzed, for example by SDS-PAGE.

**Table 1. Protein Yields and Protein Solubilization Using the Qproteome Mammalian Protein Prep Kit**

Cell line	Protein yield ( $\mu\text{g}$ )*	Solubilization (%)*
Cos-7	900	93
HeLa	1700	97
CHO	800	98
NIH	940	97
Jurkat	780	96

\* Average of 3 preparations, each of  $5 \times 10^6$  cells.

# Protocol: Isolation of Total Proteins from Mammalian Cells

The volumes given in this protocol are suitable for processing of 5–10 x 10<sup>6</sup> mammalian cells. When processing larger or smaller cultures, adjust the volume of buffer used accordingly.

## Equipment and reagents be supplied by the user

- Ice-cold PBS
- Cell scraper
- 15 ml conical tube
- Microcentrifuge tubes
- Microcentrifuge
- Rotary shaker

## Important notes before starting

- It is possible that the Protease Inhibitor does not thaw at room temperature. If this is the case, heat the protease inhibitor solution to 37°C with agitation and cool to room temperature.
- When processing small volumes, it may be necessary to dilute Benzonase<sup>®</sup> using Mammalian Lysis Buffer.

## Procedure

### Cell collection

1. **Aspirate cell-culture medium from culture plate.**
2. **Wash cells twice with 5 ml ice-cold PBS.**
3. **Add 10 ml ice-cold PBS.**
4. **Remove cells from culture plate by gentle scraping with cell-scraper and transfer cells to a pre-chilled 15 ml conical tube.**
5. **Centrifuge cell suspension for 5 min at 450 x g in a centrifuge pre-cooled to 4°C. Discard supernatant. Keep cell pellet on ice.**

## Cell lysis

6. **Add 1 U Benzonase<sup>®</sup> and 10  $\mu$ l of Protease Inhibitor (100x) to 1 ml Mammalian Lysis Buffer.**
7. **Resuspend 5–10 x 10<sup>6</sup> cells in the Mammalian Lysis Buffer prepared in step 6 and incubate on a rotary shaker for 5 minutes at 4°C.**
8. **Centrifuge the suspension for 10 minutes at 14,000 x g in a microcentrifuge tube pre-cooled to 4°C.**
9. **Transfer the supernatant into a new microcentrifuge tube.**

The supernatant contains the total protein fraction. For some applications, such as 2-D gel analysis, the protein fraction may need to be concentrated. This can be achieved by acetone precipitation (see below).

## Protocol: Acetone Precipitation of Protein Fractions

This protocol is suitable for concentrating and desalting protein samples for downstream applications such as 2D-PAGE.

1. **Add four volumes of ice-cold acetone to the protein fraction and incubate for 15 min on ice.**
2. **Centrifuge for 10 min at 12,000 x g in a pre-cooled microcentrifuge at 4°C. Discard the supernatant and air dry the pellet.**

Do not overdry the pellet as this may make it difficult to resuspend.

3. **Depending on the application, resuspend the pellet in the required sample buffer.**

## Troubleshooting Guide

This troubleshooting guide may be helpful in solving any problems that may arise. The scientists in QIAGEN Technical Services are always happy to answer any questions you may have about either the information and protocol in this handbook or molecular biology applications (see back page for contact information).

### Comments and Suggestions

---

#### **Inconsistent results in protein quantification assays**

The extraction buffer contains components that may interfere with protein quantification assays

Perform a precipitation step (e.g., using acetone, see above) to remove interfering substances.

## Ordering Information

Product	Contents	Cat. no.
Qprotome Mammalian Protein Prep Kit	For approximately 100 protein preparations from cultured mammalian cells: Buffer, Reagents, Protease Inhibitor Solution, Benzonase®	37901
<b>Related products</b>		
Qproteome Bacterial Protein Prep Kit	For soluble protein preparations from up to 4.5 liters of bacterial culture: Lysis Buffer, Lysozyme, Benzonase®	37900
Qproteome Nuclear Protein Kit	For 6 nuclear protein preparations: Buffers, Reagents, Protease Inhibitor Solution, Benzonase®	37582
Qproteome Nuclear Subfractionation Kit	For 6 nuclear protein preparations: Buffers, Reagents, Nuclear protein Fractionation Columns (6), Nuclear Protein Fractionation Resin, Protease Inhibitor Solution, Benzonase®	37531
Qproteome Albumin/IgG Depletion Kit	For albumin/IgG depletion of 6 serum or plasma samples: Albumin/IgG Depletion Spin Columns (6)	37521
Qproteome Murine Albumin Depletion Kit	For albumin depletion of 6 murine serum or plasma samples: Murine Albumin Depletion Spin Columns (6)	37591
Qproteome Total Glycoprotein Kit	For 6 total glycoprotein preps: Buffers, Lectin Spin Columns (6), Detergent Solution, Protease Inhibitor Solution, Collection Tubes (6 x 2 ml)	37541
Qproteome Mannose Glycoprotein Kit	For 6 mannose glycoprotein preps: ConA, GNA, and LCH Lectin Spin Columns (2 each); Buffers; Detergent Solution; Protease Inhibitor Solution; Collection Tubes (6 x 2 ml)	37551

<b>Product</b>	<b>Contents</b>	<b>Cat. no.</b>
Qproteome Sialic Glycoprotein Kit	For 6 sialic acid glycoprotein preps: WGA, SNA, and MAL Lectin Spin Columns (2 each); Buffers; Detergent Solution; Protease Inhibitor Solution; Collection Tubes (6 x 2 ml)	37561
Qproteome O-Glycan Glycoprotein Kit	For 6 O-glycan glycoprotein preps: ALL and PNA Lectin Spin Columns (3 each); Buffers; Protease Inhibitor Solution; Collection Tubes (6 x 2 ml)	37571
Qproteome Soluble Protein Separation Kit	For 10 soluble protein fractionations: Fractionation Buffer, Precipitation Reagents, Protease Inhibitor Solution, Benzonase <sup>®</sup>	37512
Qproteome Cell Compartment Kit	For 10 subcellular fractionations: Extraction buffers, Protease Inhibitor Solution, Benzonase <sup>®</sup>	37502

# QIAGEN Companies

Please see the back cover for contact information for your local QIAGEN office.

## QIAGEN Distributors

### Argentina

Tecnolab S.A.  
Tel: (011) 4555 0010  
Fax: (011) 4553 3331  
E-mail: info@tecnolab.com.ar  
Web site: www.tecnolab.com.ar

### Brazil

Uniscience do Brasil  
Tel: 011 3622 2320  
Fax: 011 3622 2323  
E-mail: info@uniscience.com  
Web site: www.uniscience.com

### Chile

Biosonda SA  
Tel: 562 209 6770  
Fax: 562 274 5462  
E-mail: ventas@biosonda.cl  
Web site: www.biosonda.cl

### Croatia

INEL Medicinska Tehnika d.o.o.  
Tel: (01) 2984-898  
Fax: (01) 6520-966  
E-mail: inel-medicinska-tehnika@zg.htnet.hr

### Cyprus

Scientronics Ltd  
Tel: 02-357 22 765416  
Fax: 02-357 22 764614  
E-mail: a.sarpetsas@biotronics.com.cy

### Czech Republic

BIO-CONSULT spol. s.r.o.  
Tel/Fax: (420) 2 417 29 792  
E-mail: info@bioconsult.cz  
Web site: www.bioconsult.cz

### Egypt

Clinilab  
Tel: 52 57 212  
Fax: 52 57 210  
E-mail: Clinilab@link.net

### Estonia

PEAL-Est OÜ  
Tel: (051) 65 830  
Fax: (07) 383 360  
E-mail: langel@ut.ee

### Greece

BioAnalytica S.A.  
Tel: (210)-640 03 18  
Fax: (210)-646 27 48  
E-mail: bioanalyt@hol.gr  
Web site: www.bioanalytica.gr

### Hungary

Kasztel-Med Co. Ltd.  
Tel: (01) 385 3887  
Fax: (01) 381 0695  
E-mail: info@kasztel.hu  
Web site: www.kasztel.hu

### India

Genetix  
Tel: (011)-2542 1714  
or (011)-2515 9346  
Fax: (011)-2546 7637  
E-mail: genetix@nda.vsnl.net.in

### Israel

Westburg (Israel) Ltd.  
Tel: 08-6900655  
or 1-800 20 22 20 (toll free)  
Fax: 08-6900650  
E-mail: info@westburg.co.il  
Web site: www.westburg.co.il

### Korea

LRS Laboratories, Inc.  
Tel: (02) 924-86 97  
Fax: (02) 924-86 96  
E-mail: webmaster@lrslab.co.kr  
Web site: www.lrslab.co.kr

### Malaysia

RESEARCH BIOLABS SDN. BHD.  
Tel: (603)-8070 3101  
Fax: (603)-8070 5101  
E-mail: biolabs@tm.net.my  
Web site: www.researchbiolabs.com

### Mexico

Quimica Valaner S.A. de C.V.  
Tel: (55) 55 25 57 25  
Fax: (55) 55 25 56 25  
E-mail: ventas@valaner.com  
Web site: www.valaner.com

### New Zealand

Biolab Ltd  
Tel: (09) 980 6700  
or 0800 933 966  
Fax: (09) 980 6788  
E-mail: biosciences@nzl.biolabgroup.com  
Web site: www.biolabgroup.com/nzl

### Poland

Syngen Biotech Sp.z.o.o.  
Tel: (071) 798 58 50 - 52  
Fax: (071) 798 58 53  
E-mail: info@syngen.pl  
Web site: www.syngen.pl

### Portugal

IZASA PORTUGAL, LDA  
Tel: (21) 424 7312  
Fax: (21) 417 2674  
E-mail: consultasbiotec@izasa.es

### Saudi Arabia

Abdulla Fouad Co. Ltd.  
Tel: (03) 8324400  
Fax: (03) 8346174  
E-mail: sadiq.omar@abdulla-fouad.com

### Singapore

Research Biolabs Pte Ltd  
Tel: 6777 5366  
Fax: 6778 5177  
E-mail: sales@researchbiolabs.com  
Web site: www.researchbiolabs.com

### Slovak Republic

BIO-CONSULT Slovakia spol. s.r.o.  
Tel/Fax: (02) 5022 1336  
E-mail: bio-cons@cdicon.sk  
Web site: www.bioconsult.cz

### Slovenia

MEDILINE d.o.o.  
Tel: (01) 830-80-40  
Fax: (01) 830-80-70  
or (01) 830-80-63  
E-mail: mediline@siol.net

### South Africa

Southern Cross Biotechnology (Pty) Ltd  
Tel: (021) 671 5166  
Fax: (021) 671 7734  
E-mail: info@scb.co.za

### Spain

IZASA, S.A.  
Tel: (93) 902.20.30.90  
Fax: (93) 902.22.33.66  
E-mail: consultasbiotec@izasa.es

### Taiwan

TAIGEN Bioscience Corporation  
Tel: (02) 2880 2913  
Fax: (02) 2880 2916  
E-mail: taigen@ms10.hinet.net

### Thailand

Theera Trading Co. Ltd.  
Tel: (02) 412-5672  
Fax: (02) 412-3244  
E-mail: theetrad@samart.co.th

### Turkey

Medek Medikal Ürünler ve Sağlık Hizmetleri A. S.  
Tel: (216) 302 15 80  
Fax: (216) 302 15 88  
E-mail: makalp@med-ek.com

### All other countries

QIAGEN GmbH, Germany

**Australia** ■ **QIAGEN Pty Ltd** ■ PO Box 641 ■ DONCASTER VIC 3108

Orders 03-9840-9800 ■ Fax 03-9840-9888 ■ Technical 1-800-243-066

**Austria** ■ **QIAGEN Vertriebs GmbH** ■ Löwengasse 47 / Top 6 ■ A-1030 Wien

Orders: 0800/28-10-10 ■ Fax: 0800/28-10-19 ■ Technical: 0800/28-10-11

**Belgium** ■ **QIAGEN Benelux B.V.** ■ Spoorstraat 50 ■ 5911 KJ Venlo ■ The Netherlands

Orders 0800-79612 ■ Fax 0800-79611 ■ Technical 0800-79556

**Canada** ■ **QIAGEN Inc.** ■ 2800 Argentia Road ■ Unit 7 ■ Mississauga ■ Ontario ■ L5N 8L2

Orders 800-572-9613 ■ Fax 800-713-5951 ■ Technical 800-DNA-PREP (800-362-7737)

**China** ■ **QIAGEN China Representative Office** ■ Room 2213 ■ POS PLAZA ■ Century Avenue 1600 ■

Pudong ■ Shanghai 200122 ■ Telephone 021-51162555 ■ Fax 021-51162500

**France** ■ **QIAGEN S.A.** ■ 3 avenue du Canada ■ LP 809 ■ 91974 COURTABOEUF CEDEX

Orders 01-60-920-920 ■ Fax 01-60-920-925 ■ Technical 01-60-920-930

**Germany** ■ **QIAGEN GmbH** ■ QIAGEN Strasse 1 ■ 40724 Hilden

Orders 02103-29-12000 ■ Fax 02103-29-22000 ■ Technical 02103-29-12400

**Italy** ■ **QIAGEN S.p.A.** ■ Via Grosio, 10/10 ■ 20151 Milano

Orders 02-33430411 ■ Fax 02-33430426 ■ Technical 02-33430414

**Japan** ■ **QIAGEN K.K.** ■ Forefront Tower II ■ 13-1, Kachidoki 3 Chome ■ Chuo-ku, Tokyo 104-0054

Telephone 03-5547-0811 ■ Fax 03-5547-0818 ■ Technical 03-5547-0811

**Luxembourg** ■ **QIAGEN Benelux B.V.** ■ Spoorstraat 50 ■ 5911 KJ Venlo ■ The Netherlands

Orders 8002-2076 ■ Fax 8002-2073 ■ Technical 8002-2067

**The Netherlands** ■ **QIAGEN Benelux B.V.** ■ Spoorstraat 50 ■ 5911 KJ Venlo

Orders 0800-0229592 ■ Fax 0800-0229593 ■ Technical 0800-0229602

**Switzerland** ■ **QIAGEN AG** ■ Garstligweg 8 ■ 8634 Hombrechtikon

Orders 055-254-22-11 ■ Fax 055-254-22-13 ■ Technical 055-254-22-12

**UK and Ireland** ■ **QIAGEN Ltd.** ■ QIAGEN House ■ Fleming Way ■ Crawley ■ West Sussex, RH10 9NQ

Orders 01293-422-911 ■ Fax 01293-422-922 ■ Technical 01293-422-999

**USA** ■ **QIAGEN Inc.** ■ 27220 Turnberry Lane ■ Suite 200 ■ Valencia ■ CA 91355

Orders 800-426-8157 ■ Fax 800-718-2056 ■ Technical 800-DNA-PREP (800-362-7737)

