

January 2005

Qproteome Nuclear Subfractionation Handbook

For fractionation of nuclear and nucleic-acid
binding proteins from eukaryotic cell lysates



WWW.QIAGEN.COM

Trademarks: QIAGEN® (QIAGEN Group); Benzonase® (Merck KGaA, Germany); Parafilm® (American Can Company)

Benzonase endonuclease is covered by several patents, i.e. the US patent No. 5,173,418 and EP No. 0229866. Nycomed Pharma A/S, Denmark claims worldwide patent rights to Benzonase endonuclease which are licensed exclusively to MERCK KGaA, Darmstadt (Germany).

© 2005 QIAGEN, all rights reserved.

Contents

Kit Contents	4
Storage	4
Product Use Limitations	5
Product Warranty and Satisfaction Guarantee	5
Technical Assistance	5
Safety Information	6
Introduction	7
Principle and procedure	9
Preparation of the cytosolic fraction	9
Preparation of the nucleic-acid binding protein fraction	9
Further fractionation of nucleic-acid-binding proteins	10
Preparation of the “insoluble” nuclear protein fraction	10
Protocols	11
■ Fractionation of Nuclear Proteins from Mammalian-Cell Lysates	11
■ Acetone Precipitation of Protein Fractions	16
Troubleshooting Guide	16
Ordering Information	17
QIAGEN Companies	19
QIAGEN Distributors	19

Kit Contents

Qproteome Nuclear Subfractionation Kit	
Catalog no.	37531
Number of preps	6
Lysis Buffer NL	10 ml
Extraction Buffer NX1	1 ml
Extraction Buffer NX2	1 ml
Nuclear Protein Fractionation Resin	4 ml
Nuclear Protein Fractionation Columns	6
Dilution Buffer ND	5 ml
Elution Buffer NE1	10 ml
Elution Buffer NE2	5 ml
Elution Buffer NE3	5 ml
Detergent Solution NP	0.25 ml
DTT Stock Solution	0.25 ml of a 1 M solution
Benzonase [®]	2000 Units (25 U/ μ l)
Protease Inhibitor Solution (100x)	300 μ l

Storage

Buffers NL, NX1, NX2, ND, NE1, NE2, and NE3 and Nuclear Protein Fractionation Columns should be stored at room temperature (15–25°C).

Nuclear Protein Fractionation Resin, Detergent Solution NP, and Protease Inhibitor Solution (100x) should be stored at 4°C.

Benzonase[®] and DTT Stock Solution should be stored at –20°C.

Product Use Limitations

Qproteome Nuclear Subfractionation Kits are developed, designed, and sold for research purposes only. They are not to be used for human diagnostic or drug purposes or to be administered to humans unless expressly cleared for that purpose by the Food and Drug Administration in the USA or the appropriate regulatory authorities in the country of use. All due care and attention should be exercised in the handling of many of the materials described in this text.

Product Warranty and Satisfaction Guarantee

QIAGEN guarantees the performance of all products in the manner described in our product literature. The purchaser must determine the suitability of the product for its particular use. Should any product fail to perform satisfactorily due to any reason other than misuse, QIAGEN will replace it free of charge or refund the purchase price. We reserve the right to change, alter, or modify any product to enhance its performance and design. If a QIAGEN product does not meet your expectations, simply call your local Technical Service Department or distributor. We will credit your account or exchange the product — as you wish.

A copy of QIAGEN terms and conditions can be obtained on request, and is also provided on the back of our invoices. If you have questions about product specifications or performance, please call QIAGEN Technical Services or your local distributor (see inside back cover).

Technical Assistance

At QIAGEN we pride ourselves on the quality and availability of our technical support. Our Technical Service Departments are staffed by experienced scientists with extensive practical and theoretical expertise in molecular biology and the use of QIAGEN® products. If you have any questions or experience any difficulties regarding Qproteome Nuclear Subfractionation Kits or QIAGEN products in general, please do not hesitate to contact us.

QIAGEN customers are a major source of information regarding advanced or specialized uses of our products. This information is helpful to other scientists as well as to the researchers at QIAGEN. We therefore encourage you to contact us if you have any suggestions about product performance or new applications and techniques.

For technical assistance and more information please call one of the QIAGEN Technical Service Departments or local distributors (see inside back cover).

Safety Information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate material safety data sheets (MSDSs). These are available online in convenient and compact PDF format at www.qiagen.com/ts/msds.asp where you can find, view, and print the MSDS for each QIAGEN kit and kit component.

24-hour emergency information

Emergency medical information in English, French, and German can be obtained 24 hours a day from:

Poison Information Center Mainz, Germany

Tel: +49-6131-19240

Introduction

A powerful strategy for proteomics is the combination of classical biochemical fractionation techniques — which enrich particular subcellular structures or classes of protein with distinct functions or characteristics — with mass spectrometric identification of individual proteins.

The functional architecture of the eukaryotic cell nucleus is of great interest to cell biologists. The identification of nuclear proteins — especially nucleic-acid-binding proteins (e.g., transcription factors) — is important for an understanding of genome regulation and function, and provides clues about the molecular function of novel proteins.

The nucleus contains a cell's genetic information and is the site of gene expression. Biological processes involving nucleic acids, such as transcription, replication, recombination, and DNA repair all involve the action of proteins that bind nucleic acids in a sequence-specific manner. These proteins interact with other proteins, which may or may not bind directly to nucleic acids. As a result of these interactions, large functional complexes are formed and anchored to specific nucleic acid sites. Many auxiliary proteins, enzymes, and complexes in the nucleus have important, general functions and are present in small to moderate amounts. However, proteins that bind at selected sites, such as transcription factors bound to specific promoters, represent only 0.01–0.001% of total cellular protein. A fractionation of nuclear proteins, of which over 1000 have been identified, is a prerequisite for analysis of such low-abundance proteins using current methods.

Fractionation of Cell Lysate and Subsequent Subfractionation





Figure 1 **A** Fractions obtained after processing cell lysates using the Qproteome Nuclear Subfractionation Kit. **C**: Cytosolic fraction; **D**: nucleic-acid-binding protein fraction; **I**: “Insoluble” fraction. **B** Eluates 1–3 after subfractionation of fraction N.

Nuclear Protein Fractionation Procedure



Principle and procedure

The Nuclear Subfractionation Kit is designed for specific enrichment of nuclear proteins and provides efficient subfractionation of nucleic-acid-binding proteins from cultured mammalian cells. The Nuclear Subfractionation Kit provides standardized sample preparation for reliable quantitative and qualitative analysis of proteins for targeted proteomics. Nuclear proteins are separated from complex protein mixtures, allowing detection of low-abundance proteins (e.g., most transcription factors). Starting material for one fractionation procedure is 5×10^6 – 1×10^7 cells. The procedure has been used successfully with several different mammalian cell lines including HeLa, Jurkat, NIH3T3, HEK293, and Cos (Table 1).

Table 1. Typical Total Protein Yields in Cytosolic, DNA Binding Protein, and Histone Fractions

Cell type	Number of cells processed	Cytosolic fraction	Nucleic-acid-binding proteins	Histones
HeLa	5.4×10^6	1.9 mg	0.34 mg	0.40 mg
HEK293	4.0×10^7	9.4 mg	0.55 mg	0.89 mg
Jurkat	1.5×10^7	2.6 mg	0.23 mg	0.53 mg
Cos	3.7×10^6	2.9 mg	0.14 mg	0.18 mg
NIH3T3	7.9×10^6	1.9 mg	0.28 mg	0.34 mg

Preparation of the cytosolic fraction

Cells are incubated in hypotonic buffer, causing them to swell. Detergent added to the lysis buffer ruptures the plasma membrane and centrifugation is used to separate the cytosolic fraction (supernatant) from the cell nuclei (pellet).

Preparation of the nucleic-acid binding protein fraction

The nuclear pellet is washed to remove cytosolic contaminants from the isolated cell nuclei. Washed cell nuclei are incubated in a buffer containing a high salt concentration. During this incubation step, the nuclei shrink and nucleic-acid-binding proteins (e.g., transcription factors) separate from nucleic acids and diffuse through the nuclear pore to the exterior of the nucleus. A centrifugation step separates the nucleic-acid-binding proteins (which are found in the supernatant) from the nuclear debris, which includes genomic DNA (the pellet). The expected yield from one fractionation procedure is 200–300 μg specific nucleic-acid-binding proteins, at a concentration of approximately 5 $\mu\text{g}/\mu\text{l}$. This nucleic-acid-binding protein fraction can be used for activity tests (e.g., gel-shift assays,

transcription factor activity assays) and characterization of function without further processing. Alternatively, this fraction can be further fractionated using phosphocellulose columns.

Further fractionation of nucleic-acid-binding proteins

For 2-D PAGE analysis — where a reduction in sample complexity is often desired — fractionation of the nucleic-acid-binding protein sample is recommended. nucleic-acid-binding proteins are fractionated by diluting the nucleic-acid-binding protein fraction and loading it onto a phosphocellulose column, and eluting proteins with increasing salt concentrations (Elution Buffer NE1, 0.1 M NaCl; Elution Buffer NE2, 0.35 M NaCl; Elution Buffer NE3, 1 M NaCl). Subfractionation using stabilized phosphocellulose gravity flow columns reduces the complexity of the nuclear-acid binding protein fraction and allows more sensitive detection and analysis of proteins expressed at low levels (Table 2). For most downstream applications, concentration of the phosphocellulose fractions is required. A protocol for concentration using acetone precipitation is provided on page 16.

Table 2. Typical Protein Yields in Nucleic-Acid-Binding Protein Subfractions

Cell type	Number of cells processed	Subfraction NE1	Subfraction NE2	Subfraction NE3
		Yield (% total)	Yield (% total)	Yield (% total)
HeLa	4.0 x 10 ⁶	70 µg (41%)	50 µg (29%)	50 µg (29%)
HEK293	1.0 x 10 ⁷	130 µg (52%)	60 µg (24%)	60 µg (24%)
Jurkat	1.3 x 10 ⁷	60 µg (38%)	50 µg (32%)	45 µg (29%)
NIH3T3	2.6 x 10 ⁶	50 µg (29%)	60 µg (35%)	60 µg (35%)

Preparation of the “insoluble” nuclear protein fraction

Extraction of “insoluble” nuclear proteins is performed by incubation of the nuclear debris with a buffer containing Benzonase[®], a DNase/RNase which digests genomic DNA and releases nuclear proteins intimately associated with DNA (e.g., histones).

Protocol: Fractionation of Nuclear Proteins from Mammalian-Cell Lysates

Important notes before starting

- All steps are performed at 4°C. Use pre-cooled buffers and equipment. Separated protein fractions should be stored at –80°C.
- For downstream applications, such as 2-D gel analysis, phosphocellulose column fractions must be concentrated. This can be achieved by acetone precipitation (see page 16).
- Starting material for one fractionation procedure using the Nuclear Subfractionation Kit protocol is 5×10^6 – 1×10^7 cells.
- Dilute 10 μ l of the 1 M DTT Stock solution with 90 μ l deionized, sterile water to give a concentration of 0.1 M. Store 0.1 M DTT Solution at –20°C.
- Immediately before starting the protocol, supplement the buffers used in the protocol with Protease Inhibitor Solution, DTT Stock Solutions, and Benzonase[®] as shown in Table 3.

To be supplied by the user

- PBS
- 15 ml conical tube
- Microcentrifuge tubes
- Parafilm[®]
- Optional: Acetone stored at –20°C

Table 3. Preparing Buffers for the Nuclear Subfractionation Kit Protocol.

Buffer	Required volume per prep.	Protease Inhibitor Solution	1 M DTT	0.1 M DTT	Benzonase®
Lysis Buffer NL (for protocol steps 6 and 11)	1000 μ l	10 μ l	–	5 μ l	–
Extraction Buffer NX1 (step 14)	50 μ l	0.5 μ l	–	–	–
Dilution Buffer ND (step 17)	150 μ l	1.5 μ l	–	3 μ l	–
Elution Buffer NE1 (steps 23 and 26)	1300 μ l	13 μ l	1.5 μ l	–	–
Elution Buffer NE2 (step 24)	500 μ l	5 μ l	4 μ l	–	–
Elution Buffer NE3 (step 25)	500 μ l	5 μ l	1.5 μ l	–	–
Extraction Buffer NX2 (step 26)	100 μ l	1 μ l	–	1 μ l	1 μ l

Procedure

Cell collection

- 1. Aspirate cell-culture medium from culture plate.**
- 2. Wash cells twice with 5 ml ice-cold PBS.**
- 3. Add 10 ml ice-cold PBS.**
- 4. Remove cells from culture plate by gentle scraping with cell-scraper and transfer cells to a pre-chilled 15 ml conical tube.**
- 5. Centrifuge cell suspension for 5 min at 450 x g in a centrifuge pre-cooled to 4°C. Discard supernatant. Keep cell pellet on ice.**

Cell lysis and isolation of cell nuclei

6. Gently resuspend the cells collected in step 5 in 500 μ l Lysis Buffer NL (supplemented with Protease Inhibitor Solution and 0.1 M DTT, see Table 3) by pipetting up and down several times. Incubate for 15 min on ice.
7. Transfer the resuspended cells to a clean pre-chilled microcentrifuge tube.
8. Add 25 μ l Detergent Solution NP to the cell suspension and vortex for 10 s at maximum speed.
9. Centrifuge cell suspension for 5 min at 10,000 x g in a microcentrifuge pre-cooled to 4°C.
10. Transfer the supernatant (cytosolic fraction) into a new microcentrifuge tube and store at -80°C.
11. Resuspend the pellet (which contains cell nuclei) in 500 μ l Nuclear Protein Lysis Buffer NL (supplemented with Protease Inhibitor Solution and 0.1 M DTT, see Table 3) by vortexing for 5 s at maximum speed.
12. Centrifuge the suspension of nuclei for 5 min at 10,000 x g in a pre-cooled microcentrifuge at 4°C.
13. Discard the supernatant and save the nuclear pellet.

Extraction of nucleic-acid-binding proteins

14. Resuspend the nuclear pellet from step 13 in 50 μ l Extraction Buffer NX1 (supplemented with Protease Inhibitor Solution, see Table 3) by pipetting up and down. Incubate suspension for 30 min with gentle agitation (e.g., 750 rpm in a Thermomixer) at 4°C.
15. Centrifuge the suspension of nuclei for 10 min at 12,000 x g in a microcentrifuge pre-cooled to 4°C.
16. Transfer the supernatant (which contains nucleic-acid binding proteins) into a new microcentrifuge tube.

Store the supernatant and the pellet at -80°C until you are ready to proceed with further purification steps. Alternatively, this fraction can be used directly for activity assays (e.g., gel-shift assays).

Continue with step 17 to use the supernatant for fractionation with Nuclear Protein Fractionation Resin. This fractionation is recommended to reduce the complexity of samples for 2-D PAGE analysis.

Continue with step 26 to use the pellet for extraction of "insoluble" nuclear proteins (e.g., histones).

Subfractionation of nucleic-acid-binding proteins with Nuclear Protein Fractionation Columns

17. Add 150 μ l Dilution Buffer ND (supplemented with Protease Inhibitor Solution and 0.1 M DTT, see Table 3) to a 50 μ l aliquot of nucleic-acid-binding protein supernatant from step 16.

IMPORTANT: Carry out this and all following steps at 4°C.

18. Prepare a Nuclear Protein Fractionation Column by detaching the top cap and breaking off the bottom closure.

19. Resuspend Nuclear Protein Fractionation Resin by gently vortexing. Pipet 400 μ l Nuclear Protein Fractionation Resin into a Nuclear Protein Fractionation Column and allow the resin to settle.

20. Pipet 800 μ l Elution Buffer NE1 (supplemented with Protease Inhibitor Solution and 1 M DTT, see Table 3) into the Nuclear Protein Fractionation Column.

21. Cover the opening of the column with a small square of Parafilm. Press gently on the center of the Parafilm until buffer begins to flow through the column. Remove the Parafilm and allow the buffer to run through the column.

22. Pipet the 200 μ l diluted nucleic-acid-binding protein supernatant from step 17 onto the equilibrated Nuclear Protein Fractionation Column. Discard the first 200 μ l void volume.

23. Pipet 500 μ l Elution Buffer NE1 (supplemented with Protease Inhibitor Solution and 1 M DTT, see Table 3) onto the column and save the flow-through. Label this fraction as Fraction NE1 and store at -80°C.

24. Pipet 500 μ l Elution Buffer NE2 (supplemented with Protease Inhibitor Solution and 1 M DTT, see Table 3) onto the column and save the flow-through. Label this fraction as Fraction NE2 and store at -80°C.

25. Pipet 500 μ l Elution Buffer NE3 (supplemented with Protease Inhibitor Solution and 1 M DTT, see Table 3) onto the column and save the flow-through. Label this fraction as Fraction NE3 and store at -80°C.

Note: before 2-D PAGE analysis, elution fractions must be concentrated and desalted. This can be achieved by acetone precipitation (see page 16).

Extraction of “insoluble” nuclear proteins (e.g., histones)

- 26. Resuspend the pellet from step 16 in 100 μ l Extraction Buffer NX2 (supplemented with Benzonase[®], Protease Inhibitor Solution, and 0.1 M DTT Stock Solution, see Table 3) by pipetting up and down. Incubate suspension for 1 h with gentle agitation (e.g., 750 rpm in a Thermomixer) at 4°C.**
- 27. Centrifuge pellet suspension for 10 min at 12,000 x g in a pre-cooled microcentrifuge at 4°C. Transfer the supernatant to a new microcentrifuge tube and store at -80°C.**

Protocol: Acetone Precipitation of Protein Fractions

This protocol is suitable for concentrating and desalting protein samples for downstream applications such as SDS-PAGE.

1. **Add four volumes of ice-cold acetone to the protein fraction and incubate for 15 min on ice.**
2. **Centrifuge for 10 min at 12,000 x g in a pre-cooled microcentrifuge at 4°C. Discard the supernatant and air dry the pellet.**

Do not overdry the pellet as this may make it difficult to resuspend.

3. **Depending on the application, resuspend the pellet in the required sample buffer.**

Troubleshooting Guide

	Comments and suggestions
Low protein concentration in protein fractions	
Too few cells in starting material	Use a minimum of 5×10^6 cells per preparation (see Tables 1 and 2).
Poor protein compartmentalization	
a) Incomplete removal of cytosolic fraction	Ensure that all traces of supernatant are removed after the nuclear pellet wash step.
b) Nuclei are disrupted during nucleic-acid-binding protein extraction	Gently resuspend the nuclear pellet in extraction buffer NX1 by pipetting up and down. Do not vortex.
No or low protein activity	
Fractionated proteins are degraded	Ensure all buffers and equipment are cooled to 4°C during the entire procedure. Ensure that Protease Inhibitor Solution has been added to buffers. Snap-freeze eluted proteins using liquid nitrogen.

Ordering Information

Product	Contents	Cat. no.
Qproteome Nuclear Subfractionation Kit	For 6 nuclear protein preparations: Buffers, Reagents, Nuclear protein Fractionation Columns (6), Nuclear Protein Fractionation Resin, Protease Inhibitor Solution, Benzonase®	37531
Related products		
Qproteome Total Glycoprotein Kit	For 6 total glycoprotein preps: Buffers, Lectin Spin Columns (6), Detergent Solution, Protease Inhibitor Solution, Collection Tubes (6 x 2 ml)	37541
Qproteome Mannose Glycoprotein Kit	For 6 mannose glycoprotein preps: ConA, GNA, and LCH Lectin Spin Columns (2 each); Buffers; Detergent Solution; Protease Inhibitor Solution; Collection Tubes (6 x 2 ml)	37551
Qproteome Sialic Glycoprotein Kit	For 6 sialic acid glycoprotein preps: WGA, SNA, and MAL Lectin Spin Columns (2 each); Buffers; Detergent Solution; Protease Inhibitor Solution; Collection Tubes (6 x 2 ml)	37561
Qproteome O-Glycan Glycoprotein Kit	For 6 O-glycan glycoprotein preps: ALL and PNA Lectin Spin Columns (3 each); Buffers; Protease Inhibitor Solution; Collection Tubes (6 x 2 ml)	37571
Qproteome Albumin/IgG Depletion Kit	For albumin/IgG depletion of 6 serum or plasma samples: Albumin/IgG Depletion Spin Columns (6)	37521
Qproteome Soluble Protein Separation Kit	For 10 soluble protein fractionations: Fractionation Buffer, Precipitation Reagents, Protease Inhibitor Solution, Benzonase®	37512
Qproteome Cell Compartment Kit	For 10 subcellular fractionations: Extraction buffers, Protease Inhibitor Solution, Benzonase®	37502

Notes

QIAGEN Companies

Please see the back cover for contact information for your local QIAGEN office.

QIAGEN Distributors

Argentina

Tecnolab S.A.
Tel: (011) 4555 0010
Fax: (011) 4553 3331
E-mail: info@tecnolab.com.ar
Web site: www.tecnolab.com.ar

Brazil

Uniscience do Brasil
Tel: 011 3622 2320
Fax: 011 3622 2323
E-mail: info@uniscience.com
Web site: www.uniscience.com

Chile

Biosonda SA
Tel: 562 209 6770
Fax: 562 274 5462
E-mail: ventas@biosonda.cl
Web site: www.biosonda.cl

China

Gene Company Limited
Tel: (852)2896-6283
Fax: (852)2515-9371
E-mail:
Hong Kong:
info@genehk.com
Beijing:
info_bj@genecompany.com
Shanghai:
info_sh@genecompany.com
Chengdu:
gene@public.cd.sc.cn
Guangzhou:
info_gz@genecompany.com

Croatia

INEL Medicinska Tehnika d.o.o.
Tel: (01) 2984-898
Fax: (01) 6520-966
E-mail:
inel-medicinska-tehnika@zg.htnet.hr

Cyprus

Scientronics Ltd
Tel: 02-357 22 765416
Fax: 02-357 22 764614
E-mail:
a.sarpetsas@biotronics.com.cy

Czech Republic

BIO-CONSULT spol. s.r.o.
Tel/Fax: (420) 2 417 29 792
E-mail: info@bioconsult.cz
Web site: www.bioconsult.cz

Egypt

Clinilab
Tel: 52 57 212
Fax: 52 57 210
E-mail: Clinilab@link.net

Estonia

PEAI-Est OÜ
Tel: (051) 65 830
Fax: (07) 383 360
E-mail: langel@ut.ee

Greece

BioAnalytica S.A.
Tel: (210)-640 03 18
Fax: (210)-646 27 48
E-mail: bioanalyt@hol.gr
Web site: www.bioanalytica.gr

Hungary

Kasztel-Med Co. Ltd.
Tel: (01) 385 3887
Fax: (01) 381 0695
E-mail: info@kasztel.hu
Web site: www.kasztel.hu

India

Genetix
Tel: (011)-2542 1714
or (011)-2515 9346
Fax: (011)-2546 7637
E-mail: genetix@nda.vsnl.net.in

Israel

Westburg (Israel) Ltd.
Tel: 08-6900655
or 1-800 20 22 20 (toll free)
Fax: 08-6900650
E-mail: info@westburg.co.il
Web site: www.westburg.co.il

Korea

LRS Laboratories, Inc.
Tel: (02) 924-86 97
Fax: (02) 924-86 96
E-mail: webmaster@lrslab.co.kr
Web site: www.lrslab.co.kr

Malaysia

RESEARCH BIOLABS SDN. BHD.
Tel: (603)-8070 3101
Fax: (603)-8070 5101
E-mail: biolabs@tm.net.my
Web site: www.researchbiolabs.com

Mexico

Química Valaner S.A. de C.V.
Tel: (55) 55 25 57 25
Fax: (55) 55 25 56 25
E-mail: ventas@valaner.com
Web site: www.valaner.com

Poland

Syngen Biotech Sp.z.o.o.
Tel: (071) 798 58 50 - 52
Fax: (071) 798 58 53
E-mail: info@syngen.pl
Web site: www.syngen.pl

Portugal

IZASA PORTUGAL, LDA
Tel: (21) 424 7312
Fax: (21) 417 2674
E-mail: consultasbiotec@izasa.es

Saudi Arabia

Abdulla Fouad Co. Ltd.
Tel: (03) 8324400
Fax: (03) 8346174
E-mail:
sadiq.omar@abdulla-fouad.com

Singapore

Research Biolabs Pte Ltd
Tel: 6777 5366
Fax: 6778 5177
E-mail: sales@researchbiolabs.com
Web site: www.researchbiolabs.com

Slovak Republic

BIO-CONSULT Slovakia spol. s.r.o.
Tel/Fax: (02) 5022 1336
E-mail: bio-cons@cdicon.sk
Web site: www.bioconsult.cz

Slovenia

MEDILINE d.o.o.
Tel: (01) 830-80-40
Fax: (01) 830-80-70
or (01) 830-80-63
E-mail: mediline@siol.net

South Africa

Southern Cross Biotechnology
(Pty) Ltd
Tel: (021) 671 5166
Fax: (021) 671 7734
E-mail: info@scb.co.za

Spain

IZASA, S.A.
Tel: (93) 902.20.30.90
Fax: (93) 902.22.33.66
E-mail: consultasbiotec@izasa.es

Taiwan

TAIGEN Bioscience Corporation
Tel: (02) 2880 2913
Fax: (02) 2880 2916
E-mail: taigen@ms10.hinet.net

Thailand

Theera Trading Co. Ltd.
Tel: (02) 412-5672
Fax: (02) 412-3244
E-mail: theetrad@samart.co.th

Turkey

Medek Medikal Ürünler
ve Sağlık Hizmetleri A. S.
Tel: (216) 302 15 80
Fax: (216) 302 15 88
E-mail: makialp@med-ek.com

All other countries

QIAGEN GmbH, Germany

Australia ■ **QIAGEN Pty Ltd** ■ PO Box 641 ■ DONCASTER VIC 3108

Orders 03-9840-9800 ■ Fax 03-9840-9888 ■ Technical 1-800-243-066

Belgium ■ **QIAGEN Benelux B.V.** ■ Spoorstraat 50 ■ 5911 KJ Venlo ■ The Netherlands

Orders 0800-79612 ■ Fax 0800-79611 ■ Technical 0800-79556

Canada ■ **QIAGEN Inc.** ■ 2800 Argentia Road ■ Unit 7 ■ Mississauga ■ Ontario ■ L5N 8L2

Orders 800-572-9613 ■ Fax 800-713-5951 ■ Technical 800-DNA-PREP (800-362-7737)

France ■ **QIAGEN S.A.** ■ 3 avenue du Canada ■ LP 809 ■ 91974 COURTABOEUF CEDEX

Orders 01-60-920-920 ■ Fax 01-60-920-925 ■ Technical 01-60-920-930

Germany ■ **QIAGEN GmbH** ■ QIAGEN Strasse 1 ■ 40724 Hilden

Orders 02103-29-12000 ■ Fax 02103-29-22000 ■ Technical 02103-29-12400

Italy ■ **QIAGEN S.p.A.** ■ Via Grosio, 10/10 ■ 20151 Milano

Orders 02-33430411 ■ Fax 02-33430426 ■ Technical 02-33430414

Japan ■ **QIAGEN K.K.** ■ Forefront Tower II ■ 13-1, Kachidoki 3 Chome ■ Chuo-ku, Tokyo 104-0054

Telephone 03-5547-0811 ■ Fax 03-5547-0818 ■ Technical 03-5547-0811

Luxembourg ■ **QIAGEN Benelux B.V.** ■ Spoorstraat 50 ■ 5911 KJ Venlo ■ The Netherlands

Orders 8002-2076 ■ Fax 8002-2073 ■ Technical 8002-2067

The Netherlands ■ **QIAGEN Benelux B.V.** ■ Spoorstraat 50 ■ 5911 KJ Venlo

Orders 0800-0229592 ■ Fax 0800-0229593 ■ Technical 0800-0229602

Switzerland ■ **QIAGEN AG** ■ Garstligweg 8 ■ 8634 Hombrechtikon

Orders 055-254-22-11 ■ Fax 055-254-22-13 ■ Technical 055-254-22-12

UK and Ireland ■ **QIAGEN Ltd.** ■ QIAGEN House ■ Fleming Way ■ Crawley ■ West Sussex, RH10 9NQ

Orders 01293-422-911 ■ Fax 01293-422-922 ■ Technical 01293-422-999

USA ■ **QIAGEN Inc.** ■ 27220 Turnberry Lane ■ Valencia ■ CA 91355

Orders 800-426-8157 ■ Fax 800-718-2056 ■ Technical 800-DNA-PREP (800-362-7737)

