



For life science research only. Not for use in diagnostic procedures. FOR *IN VITRO* USE ONLY.

cOmplete Lysis-M

Reagent set for highly efficient protein extraction from mammalian cells by rapid lysis and simultaneous protection of extracted proteins against a multitude of proteases

Cat. No. 04 719 956 001

Version September 2005

Store cOmplete tablets at +2 to +8°C
Store Lysis-M Reagent at +15 to +25°C

1. What this Product Does

Number of Reactions

The set is designed for the lysis of approximately 100 g of mammalian cells.

Kit Contents

Label	Contents
Lysis-M Reagent	200 ml
cOmplete, Mini Protease Inhibitor Cocktail Tablets	<ul style="list-style-type: none"> • 20 tablets, supplied in <i>EASYpacks</i> (foil blisters) • Each tablet is sufficient for a volume of 10 ml solution. • Each tablet contains 3.7 mg EDTA, resulting in a 1 mM EDTA solution in 10 ml.

Storage and Stability

- If stored at room temperature the Lysis-M Reagent is stable through the expiration date printed on the label.
- If stored dry at +2 to +8°C the cOmplete, Mini Protease Inhibitor Cocktail Tablets are stable through the expiration date printed on the label

Application

cOmplete Lysis-M allows very efficient and gentle extraction of proteins from both the cytoplasm and the nucleus of cultured mammalian cells. Efficient lysis of mammalian cells occurs in only 5 minutes at room temperature, eliminating the need for scraping, sonication or freeze-thaw cycles. The Lysis-M protein extraction reagent for mammalian cells contains a mild detergent in 25 mM bicine buffer (pH 7.6). The protein yields obtained with this kit are significantly higher compared to those obtained using sonification.

Lysis-M Reagent is compatible with many different applications, including reporter assays (*e.g.*, β -galactosidase, luciferase, chloramphenicol acetyltransferase), immunoassays (*e.g.*, Western blots, ELISAs, RIAs), protein assays (*e.g.*, protein kinase A, protein kinase C, and tyrosine kinase), and protein purification. Furthermore, the cell lysate is compatible with protein assays such as Coomassie staining and BCA (2'-Benzoyloxycinnamaldehyde) protein assays, and the reagent can be removed by dialysis.

- Ⓢ cOmplete, Mini tablets contain EDTA. If the protein of interest is to be purified by IMAC (immobilized metal-chelate affinity chromatography), *e.g.*, Poly-His tagged recombinant proteins, EDTA has to be eliminated (*e.g.*, by dialysis) prior to the chromatography. Alternatively, the product cOmplete Lysis-M, EDTA-free can be used (see table "Ordering Information").

2. How To Use this Product

2.1 Before You Begin

General Remarks

For adherent mammalian cells the maximum cell lysis without cell scraping can be obtained by using the volumes of Lysis-M Reagent specified in table 1. Estimate the volume of cells (if unknown) to calculate the required volume of Lysis-M. For example, 2×10^6 of HeLa cells is equivalent to 20 mg of cells (~10 μ l of a packed cell volume) and requires 200 μ l of Lysis-M Reagent. A smaller volume of Lysis-M Reagent may be used if more concentrated cell extracts are preferred. In this case the cells must be scraped for maximum recovery.

Safety precautions

Observe the usual precautions to be taken when handling chemicals.
⚠ Consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation of Working Solutions

The tablets can be added directly to the provided Lysis-M Reagent. One cOmplete, Mini tablet is sufficient for the inhibition of the proteolytic activity in 10 ml Lysis-M Reagent. Dissolve the tablet in 10 ml of the provided Lysis-M Reagent by incubating for 2 min at RT, afterwards vortex shortly.

2.2 Protocol for Lysis of Adherent Mammalian Cells

- 1 • Remove (decant) culture media from the adherent cells grown in monolayer culture.
 - Optional: Wash cells once in washing buffer (*e.g.* PBS*) if the culture medium contains reagents that could interfere with subsequent protein analysis.
- 2 • Add the appropriate amount of Lysis-M Reagent containing cOmplete to each plate or well (see Table 1 below).
 - Incubate for 5 min at RT with gentle shaking.
- 3 • Collect the cell lysate.
 - Ⓢ The lysate can be used directly for analysis in the presence of the cell debris.
 - Transfer the lysate to a microcentrifuge tube.
 - Centrifuge the lysate at $\sim 14,000 \times g$ for 5 – 10 min. The soluble proteins are separated from the insoluble fraction and the cell debris during centrifugation.
- 4 Transfer the supernatant containing soluble protein to a new reaction tube and proceed with further analysis.

Table 1: Suggested volumes of Lysis-M Reagent containing cOmplete to use for different sizes of standard culture plates.

Plate Size/Surface Area	Volume of Lysis-M Reagent + cOmplete
100 mm ¹⁾	500 – 1,000 μ l
60 mm	250 – 500 μ l
6-well plate	200 – 400 μ l per well
24-well plate	100 – 200 μ l per well
96-well plate	50 – 100 μ l per well

¹⁾ Cells grown in 100 mm plates typically contain 10^7 cells (50 mg). The typical yield resulting from the extraction of 10^7 cells is approx. 3 mg of total protein.

2.3 Protocol for Lysis of Mammalian Cells in Suspension

- 1 • Collect cells by centrifugation at $2,500 \times g$ for 10 min.
• Decant the supernatant.
Optional: Wash cells once in washing buffer (e.g. PBS*) if the culture medium contains reagents that could interfere with subsequent protein analysis. Centrifuge the cells at $2,500 \times g$ for 10 min after washing.
- 2 • Add at least 1 ml of Lysis-M Reagent containing cOmplete for each 100 mg (~100 μ l) of wet cell pellet.
 - ⓐ First, add 1/10 of the final recommended volume of Lysis-M Reagent containing cOmplete to the cells if large amounts of cells are used. Resuspend the pellet by pipetting up and down. Then add the rest of the Lysis-M Reagent containing cOmplete to the cell suspension.
 - ⓑ Expect to obtain approximately 6 mg of total protein from 100 mg of wet cell pellet depending on cell type.
- 3 Incubate the lysate for 10 min with gentle shaking.
Pellet cell debris by centrifugation at $\sim 14,000 \times g$ for 15 min.
- 4 Transfer the supernatant containing soluble protein to a new reaction tube and proceed with further analysis.

3. Typical Result

Cos-7 cells at confluency were harvested in Lysis-M Reagent containing cOmplete. The extracted proteins were analysed by SDS-PAGE (5 μ l/lane).
M: marker
W : whole protein fraction
S: supernatant fraction
P: pellet fraction




Fig. 1: SDS-PAGE analysis and Coomassie blue staining of proteins extracted from Cos-7 cells.

When directly dissolved in Lysis-M Reagent, the cOmplete, Mini Protease Inhibitor Cocktail maintained full functionality for efficient inhibition of a multitude of proteases (including serine-, cysteine- and metalloproteases).

Typical values for the inhibition of different proteases and protease mixtures by cOmplete, Mini in Lysis-M Reagent are shown in table 2.

Table 2: Inhibition of different proteases by cOmplete Protease Inhibitor Tablets.

Protease or protease mixture	Enzyme concentration (μ g/ml)	% inhibition after immediate addition to the protease
Pancreatic extract	20	90%
Thermolysin	0.5	96%
Trypsin	0.2	76%
Papain	330	96%

One cOmplete, Mini tablet was added per 10 ml Lysis-M Reagent. Proteolytic activity was determined with the Roche Applied Science Universal Protease Substrate* (casein, resorufin-labeled). When extractions or single-step isolations are necessary in the acid pH range, simply include pepstatin* along with cOmplete, Mini tablets to ensure aspartic (acid) protease inhibition. All experiments were performed at room temperature.

*available from Roche Applied Science

4. Troubleshooting

Observation	Possible Cause	Recommendation
Low protein yield	Protein expression is low	Optimize the transfection procedure.
	Insufficient amounts of Lysis-M Reagent	Add more Lysis-M Reagent.
	Lysis-M Reagent was unable to penetrate the cell membrane	Increase incubation time and shake more vigorously during incubation.

5. Additional Information on this Product

Product Description

cOmplete Lysis-M protein extraction reagent for mammalian cells contains a mild detergent in 25 mM bicine buffer (pH7.6). This simple extraction method allows very efficient and gentle extraction of proteins from both the cytoplasm and the nucleus of cultured mammalian cells. Efficient lysis of mammalian cells occurs in only 5 minutes at room temperature, eliminating the need for scraping, sonication or freeze-thaw cycles. The protein yields obtained with this kit are 20 to 25% higher compared to three cycles of freeze-thaw and approximately 20% higher than 2 minutes of sonication (with 50% pulse).

Proteases are ubiquitous in all living cells. As soon as cells are disrupted, proteases are released and can quickly degrade any protein (1). This can drastically reduce the yield of protein during isolation and purification. The cOmplete, Mini tablets provided with this kit allow the inhibition of a broad spectrum of serine, cysteine and metalloproteases as well as calpains. Due to the optimized composition of the tablets they show excellent protease-inhibiting effects and are therefore very well suited for the protection of proteins isolated from mammalian cells. cOmplete, Mini contains both irreversible and reversible protease inhibitors. A significant advantage is that the protease inhibitor tablets can be directly dissolved in the Lysis-M protein extraction reagent of the kit. The extracted proteins can be further purified or analyzed in downstream applications. cOmplete, Mini tablets eliminate the time-consuming search for the right protease inhibitor. The ready-to-use water-soluble, non-toxic tablets work optimally in combination with the Lysis-M Reagent.

References

- 1 Beynon RJ, Bond JS. (1986) Catabolism of intracellular protein: molecular aspects. *Am J Physiol.*, **251** (2 Pt 1),141-52.

Quality Control

The inhibitory power of cOmplete, Mini has been demonstrated with many proteases and protease mixtures. In these experiments substantially higher concentrations of proteases were used compared to the concentration usually present in extracts. The inhibitory activity of each lot is tested with a concentrated pancreas extract and a concentrated pronase solution. The proteolytic activities are thereby typically inhibited by 95% after one hour (detection with Universal Protease Substrate, casein, resorufin-labeled*).

The efficiency of cell lysis using Lysis-M Reagent is determined for each lot by functional testing.

6. Supplementary Information

6.1 Text Conventions

To make information consistent and memorable, the following text conventions are used in this package insert:

Text Convention	Use
Numbered Instructions labeled ①, ②, etc.	Steps in a procedure that must be performed in the order listed
Asterisk *	Denotes a product available from Roche Applied Science

Symbols

In this package insert the following symbols are used to highlight important information:

Symbol	Description
ⓘ	Information Note: Additional information about the current topic or procedure.
⚠	Important Note: Information critical to the success of the procedure or use of the product.

Abbreviations

In this Instruction Manual the following abbreviations are used:

Abbreviation	Meaning
f.c.	final concentration
min	minute(s)
PAGE	polyacrylamide gel electrophoresis
RT	room temperature

6.2 Ordering Information

Roche Applied Science offers a large selection of reagents and systems for life science research. For a complete overview of related products and manuals, please visit and bookmark our home page www.roche-applied-science.com.

For additional information on protease inhibition, please visit our Special Interest Site at: www.roche-applied-science.com/proteaseinhibitor

	Product	Pack Size	Cat. No.	
Complete Lysis	Lysozyme	10 g	10 837 059 001	
	TriPure Isolation Reagent	50 ml 200 ml	11 667 157 001 11 667 165 001	
	DNase I from bovine pancreas	100 ml sterile	11 284 908 001	
	DNase I recombinant	2 × 10,000 U	04 536 282 001	
	cOmplete Protease Inhibitor Cocktail Tablets in EASYpacks	cOmplete	20 tablets in foil blisters (for 50 ml each)	04 693 116 001
cOmplete, Mini		30 tablets in foil blisters (for 10 ml each)	04 693 124 001	
cOmplete, EDTA-free		20 tablets in foil blisters (for 50 ml each)	04 693 132 001	
cOmplete, Mini, EDTA-free		30 tablets in foil blisters (for 10 ml each)	04 693 159 001	
cOmplete Protease Inhibitor Cocktail Tablets in glass vials		cOmplete	20 tablets in a glass vial (for 50 ml each) 3 × 20 tablets in a glass vial (for 50 ml each)	11 697 498 001 11 836 145 001
		cOmplete, Mini	25 tablets in a glass vial (for 10 ml each)	11 836 153 001
	cOmplete, EDTA-free	20 tablets in a glass vial (for 50 ml each)	11 873 580 001	
	cOmplete, Mini, EDTA-free	25 tablets in a glass vial (for 10 ml each)	11 836 170 001	
Kits and Sets	Pefabloc SC PLUS	Set I: contains 100 mg Pefabloc SC and 5 ml PSC protector solution Set II: contains 1g Pefabloc SC and 2 × 25 ml PSC protector solution	11 873 601 001 11 873 628 001	
	Protease Inhibitor Set	Small quantities of 10 most commonly used protease inhibitors	11 206 893 001	
	Universal Protease Substrate (Casein, resorufin-labeled)	15 mg 40 mg	11 080 733 001 11 734 334 001	
Individual Protease Inhibitors	Aprotinin	10 mg 50 mg 100 mg	10 236 624 001 10 981 532 001 11 583 794 001	
	Bestatin	10 mg 50 mg	10 874 515 001 11 359 070 001	
	Calpain Inhibitor I	25 mg	11 086 090 001	
	Calpain Inhibitor II	25 mg	11 086 103 001	
	Chymostatin	10 mg	11 004 638 001	
	E-64	5 mg 10 mg 25 mg	11 585 673 001 10 874 523 001 11 585 681 001	

Product	Pack Size	Cat. No.
Leupeptin	5 mg	11 017 101 001
	25 mg	11 017 128 001
	50 mg	11 034 626 001
	100 mg	11 529 048 001
α ₂ -Macroglobulin	25 inhibitory units	10 602 442 001
Pefabloc SC	100 mg	11 429 868 001
	500 mg	11 585 916 001
	1 g	11 429 876 001
Pepstatin	2 mg	10 253 286 001
	10 mg	11 359 053 001
	50 mg	11 524 488 001
PMSF	1 g	10 236 608 001
	10 g	10 837 091 001
	25 g	11 359 061 001
TLCK - HCl	100 mg	10 874 485 001
	250 mg	10 874 493 001
Trypsin Inhibitor (chicken, egg white)	1 g	10 109 878 001
Trypsin Inhibitor (soybean)	50 mg	10 109 886 001
	500 mg	10 109 894 001
Buffers in a Box, Pre-mixed PBS Buffer, 10×	4 l	11 666 789 001

Buffers

6.3 Notice to Purchaser

Trademarks

COMPLETE, BUFFERS IN A BOX, and TRIPURE are Trademarks of Roche.

PEFABLOC is a trademark of Pentapharm AG, Basel, Switzerland.

Contact and Support

To ask questions, solve problems, suggest enhancements or report new applications, please visit our **Online Technical Support** Site at:

www.roche-applied-science.com/support

To call, write, fax, or email us, visit the Roche Applied Science home page, www.roche-applied-science.com, and select your home country.

Country-specific contact information will be displayed.

Use the Product Search function to find Pack Inserts and Material Safety Data Sheets.



Diagnostics

Roche Diagnostics GmbH
Roche Applied Science
68298 Mannheim
Germany