




Protein Assays From Bio-Rad

Kits for total protein quantitation

www.bio-rad.com/proteinassaykits

Where to Use a Protein Assay?




Setting Up a Standard Curve

Determine protein concentration by plotting the absorbance vs. concentration of known standards. Use the resulting curve to determine the concentration of unknown proteins based on their absorbance.

Note:
The best standard to use is a purified sample of your target protein. If this is not available, use Bio-Rad's bovine serum albumin or bovine γ -globulin to make your standard curve.

- For best results**
- Run a standard curve with each protein assay
 - Run at least 3 replicates of all standards and samples
 - Process the sample and standards the same way to ensure that differences in color intensity are due only to differences in protein concentration
 - Make sure the sample and standard fall within the same concentration range




Standard curve generation using known standards. A, typical standard curve for Lowry-based assays, including DC protein assay and RC DC protein assay; B, typical standard curve for Bradford-based assays, including Bio-Rad protein assay and Quick Start Bradford protein assay.



Quick Start™ Bradford Protein Assay

Quick and easy to use, includes standards and prediluted reagent

Bradford-Based*




Catalog #	Description
500-0201	Quick Start Bradford Protein Assay Kit 1, includes 1 L 1x dye reagent, BSA standard (5 x 2 mg/ml)
500-0202	Quick Start Bradford Protein Assay Kit 2, includes 1 L 1x dye reagent, BSA standard set (14 x 2 ml of 0.125–2.0 mg/ml)
500-0203	Quick Start Bradford Protein Assay Kit 3, includes 1 L 1x dye reagent, bovine γ -globulin standard (5 x 2 mg/ml)
500-0204	Quick Start Bradford Protein Assay Kit 4, includes 1 L 1x dye reagent, bovine γ -globulin standard set (14 x 2 ml of 0.125–2.0 mg/ml)



Bio-Rad Protein Assay

Includes standards and reagent

Bradford-Based*




Catalog #	Description
500-0001	Bio-Rad Protein Assay Kit I, includes dye concentrate and bovine γ -globulin standard
500-0002	Bio-Rad Protein Assay Kit II, includes dye concentrate and bovine serum albumin standard



DC™ Protein Assay

Detergent compatible

Lowry-Based**



Catalog #	Description
500-0111	DC Protein Assay Kit I, includes DC protein assay reagents package and bovine γ -globulin standard
500-0112	DC Protein Assay Kit II, includes DC protein assay reagents package and bovine serum albumin standard

Related Products


Catalog #	Description	Catalog #	Description
168-1130	iMark™ Microplate Absorbance Reader	170-2511	trUView Cuvettes, pack of 100
170-2502	Standard Cuvette, 1–3.5 ml, quartz	170-2525	SmartSpec™ Plus Spectrophotometer
170-2503	Semimicrovolume Cuvette, 0.5–1.4 ml, quartz	223-9950	Standard Disposable Polystyrene Cuvettes, 3.5 ml, pack of 100
170-2504	Microvolume Cuvette, 200–700 μ l, quartz	223-9955	Semimicrovolume Disposable Polystyrene Cuvettes, 1.5 ml, pack of 100
170-2505	Submicrovolume Cuvette, 80–100 μ l, quartz		
170-2510	trUView™ Cuvettes, pack of 50		



RC DC™ Protein Assay

Reductant and detergent compatible

Lowry-Based**



Catalog #	Description
500-0121	RC DC Protein Assay Kit I, includes RC reagents package, DC protein assay reagents package, bovine γ -globulin standard
500-0122	RC DC Protein Assay Kit II, includes RC reagents package, DC protein assay reagents package, bovine serum albumin standard

Reagent Compatibility

	QSBPA	BRPA	DCPA	RDCPA
Acetate	—	0.6 M	—	—
Acetone	10%	✓	—	—
Acetonitrile	10%	—	—	—
Acidic pH	—	✓	—	—
Adenosine	—	0.001 M	—	—
Amino acids	—	—	✓	—
Ammonium sulfate	1 M	1 M	0.5 M	—
Ampholytes pH 3–10	0.5%	0.5%	—	—
ASB-14	0.03%	—	—	—
Ascorbic acid	0.05 M	—	—	—
ATP	—	0.001 M	—	—
Barbital	—	✓	—	—
BES	—	2.5 M	—	—
Bis-Tris, pH 6.5	0.2 M	—	—	—
Boric acid	—	✓	—	—
Brij-35	—	—	1%	—
C ₁₂ E ₈	—	—	0.2%	—
Cacodylate-Tris	—	0.1 M	—	—
Calcium chloride	0.04 M	—	0.05 M	—
CDTA	—	0.05 M	—	—
CHAPS	10%	—	1%	2%
CHAPSO	10%	—	1%	—
Citrate	—	0.05 M	—	—
Deoxycholate	—	0.1%	—	—
Deoxycholic acid	0.2%	—	—	—
Dithioerythritol (DTE)	0.01 M	—	—	—
Dithiothreitol (DTT)	0.01 M	1 M	0.001 M	0.1 M
DMSO	5%	—	—	—
DNA	—	1 mg/ml	—	—
Eagle's MEM	✓	✓	—	—
Earle's salt solution	✓	✓	—	—
EDTA	0.2 M	0.1 M	0.025 M	0.1 M
EGTA	0.2 M	0.05 M	—	—
Ethanol	10%	✓	—	—
Formic acid	—	1.0 M	—	—
Fructose	—	✓	—	—
Glucose	20%	✓	—	—
Glutathione	—	✓	—	—
Glycerol	5%	99%	—	—
Glycine	0.1 M	0.1 M	—	—
Guanidine HCl	2 M	✓	0.4 M	—
Hank's salt solution	✓	✓	—	—
HCl	0.1 M	—	0.5 M	—
HEPES	0.1 M	0.1 M	—	—
Imidazole	0.2 M	—	—	0.5 M
Laemmli buffer	—	—	—	✓
Magnesium chloride	1 M	1 M	—	—
Malic acid	—	0.2 M	—	—
2-Mercaptoethanol	1 M	1 M	X	5%
MES	0.1 M	0.7 M	—	—
Methanol	10%	✓	—	—
Modified Dulbecco's PBS	✓	—	—	—
MOPS	0.1 M	0.2 M	—	—
NAD	0.002 M	0.001 M	—	—
NaSCN	—	3 M	—	—
Nonidet P-40	0.25%	—	2%	—
Octyl β -glucoside	0.5%	—	1%	—
Octyl β -thioglucoyanoside	1%	—	—	—
PBS	✓	—	—	—
Peptone	—	✓	—	—
Phenol Red	0.5 mg/ml	—	—	—
Phenol	—	5%	—	—
Phosphate	—	1 M	—	—
PIPES	0.2 M	0.5 M	—	—
PMSF	0.002 M	—	—	—
Polyadenylic acid	—	0.001 M	—	—
Polypeptides (MW <3000)	—	✓	—	—
Potassium chloride	2 M	1 M	—	—
Potassium phosphate	0.5 M	—	—	—
Pyrophosphate	—	0.2 M	—	—
rRNA	—	0.25 mg/ml	—	—
SB 3–10	0.1%	—	—	—
SDS	0.03%	0.1%	10%	—
Sodium acetate, pH 4.8	0.2 M	—	—	—
Sodium azide	0.5%	—	0.05%	—
Sodium bicarbonate	0.2 M	—	—	—
Sodium carbonate	0.1 M	—	—	—
Sodium chloride	2.5 M	5 M	—	—
Sodium citrate, pH 4.8 or 6.4	0.2 M	—	—	—
Sodium hydroxide	0.1 M	—	0.5 M	2.5 M
Sodium phosphate	0.5 M	✓	—	—
Streptomycin sulfate	—	20%	—	—
Sucrose	10%	—	—	—
TBP	0.005 M	—	—	0.002 M
TBS	0.5x	—	—	—
TCEP	0.02 M	—	—	—
Thesit	—	—	1%	—
Thiourea	1 M	—	—	—
Thymidine	—	0.001 M	—	—
Total RNA	—	0.30 mg/ml	—	—
Tricine	0.05 M	✓	—	—
Triethanolamine, pH 7.8	0.05 M	—	—	—
Tris, pH 8	1 M	2 M	0.1 M	0.5 M
Tris/glycine	1x	—	—	—
Tris/glycine/SDS	0.5x	—	—	—
Triton X-100	0.05%	0.1%	1%	2%
tRNA	—	0.4 mg/ml	—	—
Tween 20	0.01%	—	1%	2%
Tyrosine	—	0.001 M	—	—
Urea	4 M	6 M	4 M	—
Vitamins	—	✓	—	—

Reagents tested for compatibility. Concentrations represent maximum concentrations for standard assay. ✓ = compatible; — = not tested; X = not compatible.

*Adapted from the method of Bradford, Bradford MM (1976). A rapid and sensitive method for the quantitation of microgram quantities of protein utilizing the principle of protein-dye binding. Anal Biochem 72, 248-254. **Adapted from the method of Lowry, Lowry OH et al. (1951). Protein measurement with the Folin phenol reagent. J Biol Chem 193, 265-275. Coomassie is a trademark of BASF Aktiengesellschaft. Triton is a trademark of Union Carbide Corporation. Tween and Brij are trademarks of ICI Americas Inc.