

NanoOrange[®] Protein Quantitation Kit (N-6666)

Quick Facts

Storage upon receipt:

Component A:

- 4°C or –20°C
- Desiccate
- Protect from light

Component B:

- Room temperature

Component C:

- 4°C or –20°C

Abs/Em: 470/570 nm

Introduction

Molecular Probes' NanoOrange[®] Protein Quantitation Kit represents a technological breakthrough for the quantitation of proteins in solution, combining ease of use and superb sensitivity. The unique chemistry of the NanoOrange reagent allows for accurate detection of proteins in solution at concentrations between 10 ng/mL and 10 µg/mL (Figure 1). This level of sensitivity is considerably better than that achieved with the BCA




Figure 1. Quantitative analysis of bovine serum albumin (BSA) using the NanoOrange Protein Quantitation Kit. The inset corresponds to the shaded area in the lower left corner of the plot (0 to 500 ng protein per mL) and illustrates the detection limit of 10 ng/mL. Fluorescence measurements were carried out on a fluorometer using excitation/emission wavelengths of 485/590 nm.

method, the Bradford assay, the Lowry assay or absorption at 280 nm (Table 1).¹⁻⁴ The NanoOrange assay also shows less protein-to-protein variability than the Bradford assay.

To perform a protein assay, the protein sample is simply added to the diluted NanoOrange reagent, and this mixture is heated at 95°C for ten minutes. Fluorescence can be measured as soon as the mixture has cooled to room temperature. Alternatively, samples can be read up to six hours after preparation with no loss in sensitivity, as long as samples are protected from light.

The NanoOrange reagent is virtually nonfluorescent in aqueous solution, but upon interaction with proteins it undergoes a dramatic fluorescence enhancement. When bound to proteins in the diluent provided, the NanoOrange quantitation reagent has a broad excitation peak centered at about 470 nm and a broad emission peak centered at about 570 nm (Figure 2); it is suitable for use with a variety of instruments, including standard fluorometers, minifluorometers or fluorescence microplate readers. With a fluorometer or minifluorometer, the NanoOrange assay can detect from 10 ng protein per mL to 10 µg/mL using a single dye concentration; with fluorescence-based microplate readers, the assay is reliable from 100 ng/mL to 10 µg/mL.

Table 1. Comparison of methods for quantitating proteins in solution.

Assay Method	Useful Range	Comments
NanoOrange assay	10 ng/mL to 10 µg/mL	<ul style="list-style-type: none"> • Samples can be read up to six hours later without any loss in sensitivity • Low protein-to-protein signal variability • Detection not influenced by reducing agents or nucleic acids
BCA method ¹	0.5 µg/mL to 1.5 mg/mL	<ul style="list-style-type: none"> • Samples must be read within 10 minutes • Not compatible with reducing agents
Bradford assay ²	1 µg/mL to 1.5 mg/mL	<ul style="list-style-type: none"> • Proteins precipitate over time • High protein-to-protein signal variability • Not compatible with detergents
Lowry assay ³	1 µg/mL to 1.5 mg/mL	<ul style="list-style-type: none"> • Lengthy, multistep procedure • Not compatible with detergents, carbohydrates or reducing agents
Absorbance at 280 nm ⁴	50 µg/mL to 2 mg/mL	<ul style="list-style-type: none"> • High protein-to-protein signal variability • Detection influenced by nucleic acids and other contaminants




Figure 2. The normalized fluorescence excitation and emission spectra of the NanoOrange reagent in the presence of protein quantitation diluent and 150 $\mu\text{g}/\text{mL}$ BSA.

Materials

- **NanoOrange protein quantitation reagent** (Component A), 1.0 mL of a 500X concentrate in DMSO
- **NanoOrange protein quantitation diluent** (Component B), 50 mL of a 10X concentrate containing 2 mM sodium azide
- **Bovine serum albumin (BSA) standard** (Component C), 0.5 mL of a 2 mg/mL solution in water, containing 2 mM sodium azide

The kit supplies sufficient material for performing 200 assays using the protocol described. This assay is performed in 2.5 mL volumes, from which 2.0 mL is transferred to a standard cuvette for fluorescence determination in a fluorometer or minifluorometer. If the assay volume is reduced by 10-fold, there is sufficient material for 2000 assays using 200 μL volumes for detection in a fluorescence microplate reader.

Storage and Handling

Storage

For long-term storage, the NanoOrange protein quantitation reagent should be stored at either 4°C or -20°C, desiccated and protected from light. Before opening, the vial containing the NanoOrange reagent should be allowed to warm to room temperature and then briefly centrifuged in a microfuge to deposit the DMSO solution at the bottom of the vial. For convenient short-term storage (up to one week), the NanoOrange reagent may be kept at room temperature, protected from light.

The NanoOrange protein quantitation diluent should be stored at room temperature. At cooler temperatures, some of the components may precipitate; however, the precipitate will readily redissolve if the diluent is warmed to room temperature and gently mixed.

The BSA standard should be stored at 4°C; it may be frozen at -20°C for long-term storage.

When stored properly, all three stock solutions are stable for six months to one year.

Handling

We must caution that no data are available addressing the toxicity of the NanoOrange reagent. The DMSO stock solution

should be handled with particular caution as DMSO is known to facilitate the entry of organic molecules into tissues. We strongly recommend using double gloves when handling DMSO stock solutions.

Disposal

Solutions of the NanoOrange reagent should be poured through activated charcoal prior to disposal. The charcoal must then be incinerated to destroy the dye. The amount of NanoOrange reagent present in one vial can be fully adsorbed from normal staining solutions using 1 g of activated charcoal.

Experimental Protocol

This section describes the protocol for 2 mL assays. For 200 μL microplate assays, adjust all volumes accordingly.

Reagent Preparation

1.1 Prepare 1X protein quantitation diluent: Dilute the concentrated NanoOrange protein quantitation diluent (Component B) 10-fold in distilled water. For each assay, 2.5 mL of 1X protein quantitation diluent will be required.

1.2 Prepare 1X NanoOrange working solution: Dilute the NanoOrange protein quantitation reagent (Component A) 500-fold into the 1X protein quantitation diluent. For example, to prepare 50 mL of 1X NanoOrange working solution (enough for 20 assays), first prepare the 1X diluent by mixing 5 mL of the 10X diluent stock with 45 mL of distilled water; next add 100 μL of the NanoOrange reagent and mix thoroughly. Protect the 1X NanoOrange working solution from light to prevent photodegradation of the NanoOrange dye. For best results, the working solution should be used within a few hours of its preparation.

Notes About Heating Solutions in Microplates

This assay requires heating the sample to 90°C for 10 minutes to denature the protein. When samples are heated in a microplate, the plate must be covered with a plastic cover or heat-resistant film to prevent evaporation. During the heating and cooling steps, water may condense on the cover, unless a special heater cover is used. To achieve accurate results, these water droplets must be collected by spinning the microplate briefly in a centrifuge. Alternatively, samples may be processed in microcentrifuge tubes, spun in a microcentrifuge and then transferred to microplate wells for reading. If transferring samples, we recommend transferring only 200 μL of the 250 μL sample, as it is difficult to retrieve all 250 μL using a pipette.

Protein Standard Curve

The reference standard curve is used not only to convert fluorescence to $\mu\text{g}/\text{mL}$ protein, but also to control for variation between fluorometers and for day-to-day variation in the performance of a single fluorometer.

Ideally, the protein type used for the standard curve should be the same as the protein that is used in the experiment; however, as with other protein assays, bovine serum albumin (BSA) serves as a convenient reference standard. The NanoOrange Kit includes a 2 mg/mL sample of BSA (Component C) that can be used to prepare a standard curve.

The standard curve may be generated to cover the full assay range, 0–10 $\mu\text{g}/\text{mL}$, or to cover a selected range. This section

Table 2. Protocol for preparing a standard curve using BSA.

Volume (µL) of BSA Solution*	Volume of 1X NanoOrange Working Solution	Final BSA Concentration
0 mL	2.50 mL	0 µg/mL
2.50 mL of 10 µg/mL	0 mL	10 µg/mL
1.50 mL of 10 µg/mL	1.00 mL	6 µg/mL
0.75 mL of 10 µg/mL	1.75 mL	3 µg/mL
0.25 mL of 10 µg/mL	2.25 mL	1 µg/mL
1.50 mL of 1 µg/mL	1.00 mL	0.6 µg/mL
0.75 mL of 1 µg/mL	1.75 mL	0.3 µg/mL
0.25 mL of 1 µg/mL	2.25 mL	0.1 µg/mL
1.50 mL of 0.1 µg/mL	1.00 mL	0.06 µg/mL
0.75 mL of 0.1 µg/mL	1.75 mL	0.03 µg/mL
0.25 mL of 0.1 µg/mL	2.25 mL	0.01 µg/mL

* The BSA solutions must be made up in 1X NanoOrange working solution, as described in the text.

describes how to generate a simple standard curve with points corresponding to 0, 1, 3, 6 and 10 µg BSA per mL. If desired, serial dilutions can be made to create additional standards ranging from 0.01 to 0.6 µg/mL, to fill out the standard curve in the low range.

2.1 Prepare a 10 µg/mL solution of BSA: Dilute the 2 mg/mL BSA standard (Component C) 1:200 into the 1X NanoOrange working solution. For instance, dilute 30 µL of BSA standard into 5.97 mL of 1X NanoOrange working solution (prepared in step 1.2).

2.2 Dilute the 10 µg/mL BSA solution to make 0, 1, 3, 6 and 10 µg/mL standards, as described in Table 2.

2.3 If desired, prepare 0.1, 0.3 and 0.6 µg/mL standards, as described in Table 2, by diluting a 1 µg/mL BSA solution. Prepare the 1 µg/mL BSA solution by diluting 300 µL of 10 µg/mL BSA (made in step 2.2) into 2.70 mL of 1X NanoOrange working solution.

2.4 If desired, prepare 0.01, 0.03 and 0.06 µg/mL standards, as described in Table 2, by diluting a 0.1 µg/mL BSA solution. Prepare the 0.1 µg/mL BSA solution by diluting 300 µL of 1 µg/mL BSA (made in step 2.3) into 2.70 mL of 1X NanoOrange working solution.

2.5 Incubate samples at 90°C to 96°C for 10 minutes, protected from light.

2.6 Cool to room temperature for 20 minutes, protected from light.

2.7 Measure the fluorescence: Transfer 2.0 mL of the sample to a standard disposable fluorescence cuvette and measure the fluorescence using a fluorometer equipped with filters or settings capable of allowing excitation at about 485 nm and capturing the emission at about 590 nm. To ensure that the sample readings

remain in the detection range of the fluorometer, set the instrument's gain so that the 10 µg/mL sample yields a fluorescence intensity near the maximum. To minimize photobleaching effects, keep the time for fluorescence measurement as short as possible and constant for all samples.

2.8 Generate a standard curve: Subtract the fluorescence value of the reagent blank (0 µg/mL) from that of each sample. Use these corrected values to generate a standard curve of fluorescence versus protein concentration (for example, see Figure 1).

Sample Analysis

3.1 Dilute the experimental protein solution in 1X NanoOrange working solution (prepared in step 1.2) to achieve a final volume of 2.5 mL. It is best if the sample volume is no more than ~4% of the total volume. If higher volumes must be used, we recommend

Table 3. Tolerance levels for contaminants in the NanoOrange protein quantitation assay.

Contaminating Compound	Maximum Tolerable Concentration*
DTT	100 mM
β-mercaptoethanol	100 mM
urea	1 M
sodium chloride	20 mM
potassium chloride	20 mM
magnesium chloride	1 mM
calcium chloride	1 mM
zinc chloride	0.5 mM
sodium acetate	20 mM
sodium phosphate	20 mM
ammonium sulfate	10 mM
HEPES	10 mM
sodium azide	10 mM
EDTA	5 mM
sodium hydroxide	10 mM
hydrochloric acid	10 mM
ascorbic acid	10 mM
sucrose	10 mM
glycerol	10%
PEG	1%
DNA	100 ng/mL
amino acids	10 µg/mL
Tween® 20	below 0.001%
SDS	0.01%
Triton® X-100	0.001%

* Compounds present in the final assay solution at or below the indicated concentration do not appreciably interfere with the NanoOrange protein quantitation assay. Whenever feasible, the blank and protein standards should be prepared in a solution closely matching that of the experimental samples.

performing a standard curve using similar volumes. You may wish to use two or three different dilution factors for a given sample. Higher dilution factors may dilute contaminants to acceptable levels (see Table 3 for contaminant tolerance limits in the final assay solution); however, extremely small sample volumes should be avoided as they are difficult to pipet accurately.

3.2 Incubate samples at 90°C to 96°C for 10 minutes, protected from light.

3.3 Cool to room temperature for at least 20 minutes, protected from light.

3.4 Measure the fluorescence: Transfer 2.0 mL of the sample to a standard disposable fluorescence cuvette and measure the fluorescence using instrument parameters identical to those used in generating the standard curve (step 2.6). To minimize photo-bleaching effects, keep the time for fluorescence measurement as short as possible and the same as that used for the standard curve.

3.5 Determine the protein concentration: Subtract the fluorescence value of the reagent blank from that of the sample and use the standard curve generated in step 2.7 to determine the protein concentration of the sample.

References

1. Anal Biochem 150, 76 (1985); 2. Anal Biochem 72, 248 (1976); 3. J Biol Chem 193, 265 (1951); 4. Scopes, R.K., *Protein Purification, Principles and Practice, 2nd Edition*, Springer-Verlag (1987).

Product List

Current prices may be obtained from our Web site or from our Customer Service Department.

Cat #	Product Name	Unit Size
N-6666	NanoOrange® Protein Quantitation Kit *200-2000 assays*	1 kit

Contact Information

Further information on Molecular Probes' products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Leiden, the Netherlands. All others should contact our Technical Assistance Department in Eugene, Oregon.

Please visit our Web site — www.probes.com — for the most up-to-date information

Molecular Probes, Inc.

PO Box 22010, Eugene, OR 97402-0469
Phone: (541) 465-8300 • Fax: (541) 344-6504

Customer Service: 7:00 am to 5:00 pm (Pacific Time)
Phone: (541) 465-8338 • Fax: (541) 344-6504 • order@probes.com

Toll-Free Ordering for USA and Canada:
Order Phone: (800) 438-2209 • Order Fax: (800) 438-0228

Technical Assistance: 8:00 am to 4:00 pm (Pacific Time)
Phone: (541) 465-8353 • Fax: (541) 465-4593 • tech@probes.com

Molecular Probes Europe BV

PoortGebouw, Rijnsburgerweg 10
2333 AA Leiden, The Netherlands
Phone: +31-71-5233378 • Fax: +31-71-5233419

Customer Service: 9:00 to 16:30 (Central European Time)
Phone: +31-71-5236850 • Fax: +31-71-5233419
eurorder@probes.nl

Technical Assistance: 9:00 to 16:30 (Central European Time)
Phone: +31-71-5233431 • Fax: +31-71-5241883
eurotech@probes.nl

Molecular Probes' products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Several of Molecular Probes' products and product applications are covered by U.S. and foreign patents and patents pending. Our products are not available for resale or other commercial uses without a specific agreement from Molecular Probes, Inc. We welcome inquiries about licensing the use of our dyes, trademarks or technologies. Please submit inquiries by e-mail to busdev@probes.com. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2001, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.