For life science research only. Not for use in diagnostic procedures. FOR *IN VITRO* USE ONLY.

Rapid Translation System RTS Amino Acid Sampler

Amino acid solutions for use in RTS 500 cell-free protein expression

Cat. No. 3 262 154

Store the kit at -15 to $-25^{\circ}C$

For 5 RTS 500 reactions

Instruction Manual

Version 1, August 2001

The Rapid Translation System RTS website: www.proteinexpression.com



1. Preface

1.1 Table of contents

1.	Preface	2
1.1	Table of contents	
1.2	Kit Contents	3
2.	Introduction	5
2.1	Product overview	5
2.2	Background Information	6
3.	Procedures and required materials	7
3.1	Before you begin	7
3.2	Preparation of Amino Acid Mixture	8
3.3	Preparation of reactions components	
3.3.1	Reconstitution of reaction components from the RTS 500 E.coli HY kit	10
3.3.2	Reaction Components from the RTS Amino Acid Sampler	11
3.3.3	Preparation of working solutions	
3.3.4	Running an experiment	11
4.	Typical results	12
5.	Appendix	14
5.1	Trouble shooting	14
5.2	References	
5.3	How to contact Roche Molecular Biochemicals	
5.4	Related products	15

1. Preface, continued

1.2 Kit Contents

Vial	Label	Contents and function
1	Alanine Ala, A	
2	Arginine Arg, R	
3	Asparagine Asn, N	
4	Aspartic Acid Asp, D	
5	Cysteine Cys, C	
6	Glutamine Gln, Q	
7	Glutamic Acid Glu, E	
8	Glycine Gly, G	
9	Histidine His, H	
10	Isoleucine Ile, I	One vial each, containing 1.0 ml 42 mM buffered
11	Leucine Leu, L	L-amino acid solution.
12	Lysine Lys, K	
13	Methionine Met, M	
14	Phenylalanine Phe, F	
15	Proline Pro, P	
16	Serine Ser, S	
17	Threonine Thr, T	
18	Tryptophan Trp, W	
19	Tyrosine Tyr, Y	
20	Valine Val, V	
21	DTT 1,4-Dithiothreitol	Five vials, each containing 0.5 ml 40 mM DTT

1. Preface, continued

Safety Information None of the bottles contains hazardous substances in reportable quantities. The usual precautions taken when handling chemicals should be observed. Used reagent can be disposed off in waste water in accordance with local regulations. In case of eye contact flush eyes with water. In case of skin contact wash off with water. In case of indestion seek medical advice. Laboratory To avoid contamination with RNases we recommend using RNase-free materials and Requirements wearing gloves. The Rapid Translation System RTS Amino Acid Sampler is designed for use with the Additional ٠ reagents required Rapid Translation System RTS 500 E. coli HY Kit (Cat. No. 3 246 817; 2 reactions and 3246 949: 5 reactions). · Modified amino acids for combination with the amino acids provided in the sampler must be obtained separately.

2. Introduction

2.1 Product overview

System components and application	The Rapid Translation System RTS 500 consists of the RTS 500 Instrument and the RTS 500 <i>E.coli</i> HY Kit which contains reagents and reaction devices for 2 (Cat. No. 3 246 817) or 5 (Cat. No. 3 246 949) coupled transcription/translation reactions.		
	The amino acids in the RTS 500 <i>E.coli</i> HY kit are provided in a separate tube apart from other reaction components and can therefore be exchanged for customized amino acid mixtures, e.g. for production of labeled proteins for NMR-spectroscopy. The RTS Amino Acid Sampler provides the basis for preparing these customized mixes, in which single amino acids are exchanged. By providing each amino acid separately, any mixture can be prepared.		
Stability	This product is stable at -15° C to -25° C until the expiration date printed on the label.		
	All components have been tested for their stability after freezing and thawing. No decrease in expression level was observed after 5 freeze-thaw cycles.		
	DTT (Vial 21, five vials) should be used only once after opening a vial. A separate vial is used for each RTS 500 reaction.		

2.2 Background Information

Coupled <i>in vitro</i> transcription/ translation using the Rapid Transla- tion System	NMR spectroscopy plays an increasingly important role in structural biology since it not only allows to solve the 3D-structure of the protein of interest but also to investi- gate dynamics and intermolecular interactions (1, 2, 3). Cell-free protein synthesis has been shown to have exceptional advantages when proteins have to be labeled for NMR:
	• the protein of interest is expressed from an exogenous plasmid template and is the only protein produced and labeled,
	 the translation reaction is performed in a small volume, requiring relatively low amounts of labeled amino acids
	• the metabolism of amino acids in a cell free system is lower compared to cellular systems, leading to more amino acid-specific labeling
	 no expensive media are required since no cells need to be cultivated
	Using the RTS 500 <i>E.coli</i> HY Kit in combination with the RTS Amino Acid Sampler and selected labeled amino acid(s), it is possible to label any amino acid of choice in the protein of interest. The complexity of NMR-spectra can be greatly reduced by amino acid-specific labeling (2, 3).
Testing protein expression in the Rapid Translation System	The RTS 500 <i>E.coli</i> HY Kit is based on an <i>E.coli</i> lysate and has restrictions for the expression of certain classes of proteins. To test whether a protein is successfully expressed in RTS, we recommend to start with a small-scale test expression in the RTS 100 <i>E.coli</i> HY Kit (Cat. No. 3 186 148). This kit allows the fast and economic production of protein in a 50 μ l reaction.
	After positive reactions we recommend to proceed with the expression on a larger scale using the RTS 500 <i>E.coli</i> HY Kit.
	For a detailed description of the Rapid Translation System please refer to the pack insert of the RTS 500 <i>E.coli</i> HY Kit or visit www.proteinexpression.com

3. Procedures and required materials

3.1 Before you begin

Reagents	 First prepare the amino acid mixture (section 3.2), then reconstitute the required reaction components of the RTS 500 <i>E. coli</i> HY Kit (Section 3.3).
	 The contents of the RTS Amino Acid Sampler can be thawed at room temperature. Some amino acids show some precipitation after thawing, which disappears after warming to 37–42°C.
	 Reconstitute only the bottles from the RTS 500 <i>E.coli</i> HY Kit needed for the experiment, using only the reconstitution buffer provided in the RTS 500 <i>E.coli</i> HY Kit.
	Keep reconstituted reagents and working solutions on ice before use.
Equipment	 Use the Rapid Translation System RTS 500 Instrument from Roche Molecular Biochemicals for optimal reaction performance and convenience (Cat. No. 3 064 859)
	 Use only calibrated pipettes.
	 Use RNase-free plastic and glassware.

3.2 Preparation of Amino Acid Mixture

Preparation of
labeled amino acidA 42 mM stock solution of a labeled amino acid is recommended for labeling during
labeled amino acid RTS 500 *E.coli* HY protein expression. The amino acid solutions provided in the
sampler are all adjusted to 42 mM, leading to a final concentration of 2.1 mM for each
amino acid in the amino acid mix and 0.5 mM for each amino acid in the final RTS 500
reaction.

For adjustment of concentration the use of the Reconstitution Buffer (Bottle 6) from the RTS 500 E. coli HY Kit is recommended.

Some amino acids require additional treatment for dissolving and stability. When dissolving

• Tryptophan, adjust pH to 7.1 with HCl 0.1 M (approx. 50 μl required for 1 ml Trp),

The amino acids in the RTS 500 Amino Acid Sampler are arranged in a grid (see Fig.

- Tyrosine, adjust pH to 10.0 with KOH 2M (approx. 165 μl required for 1 ml Tyr),
- · Cysteine, add DTT to a final concentration of 8 mM,
- Methionine, add DTT to a final concentration of 4 mM.

Preparation of amino acid labeling mix



Fig. 1: Before preparing the mix all amino acids are in their basic position.

We recommend the removal of the amino acid(s) to be labeled from the grid and the replacement by the labeled amino acid stock solution(s) on the open position(s) (as indicated by a black capped tube for Gly in Fig. 1 and 2). The mix can now be prepared by pipetting 150 μ l per amino acid, using a fresh pipette tip for each amino acid. Use an RNase free tube (e.g. 15 ml Falcon) for preparation of the mix. The final volume of the mix is 3.0 ml.

By moving the vial to the next row after pipetting there is a visual check for which amino acids already have been combined. The risk of leaving amino acids out or double pipetting of amino acids is thus minimized.

3.2 Preparation of Amino Acid Mixture, continued



Fig. 2: After being pipetted from its tube, each amino acid is moved to the next row (result after 5 amino acids shown boxed). The progress of mixing can easily be monitored this way.

After pipetting all 20 amino acids the solution is mixed thoroughly by vortexing. The amino acid labeling mix is now ready for addition to reaction and feeding solution (See Section 3.3.3.).

3.3 Preparation of reactions components

3.3.1 Reconstitution of reaction components from the RTS 500 E.coli HY kit

Using the RTS 500 The Rapid Translation Kit RTS 500 *E.coli* HY Kit is available in 2 pack sizes: *E.coli* HY Kit Cat. No. 3 246 817 for two 1 ml reactions Cat. No. 3 246 949 for five 1 ml reactions

When combined with the RTS Amino Acid Sampler, the Amino Acid Mix without Methionine (Bottle 4, brown cap) and Methionine (Bottle 5, yellow cap) contained in the Rapid Translation Kit RTS 500 *E.coli* HY Kit are not required.

Component	Reconstitution/Preparation of working solution	For use in
<i>E. coli</i> Lysate Bottle 1, red cap	Reconstitute the lyophilizate with 0.525 ml of Reconstitution Buffer (bottle 6). Mix gently. DO NOT VORTEX!	Section 3.3.3 Reaction Solution
Reaction Mix Bottle 2, green cap	Reconstitute the lyophilizate with 0.25 ml of Reconstitution Buffer (bottle 6), mix by rolling or shaking.	Section 3.3.3Reaction Solution
Feeding Mix Bottle 3, blue cap	Reconstitute the lyophilizate with 8.1 ml of Reconstitution Buffer (bottle 6), mix by rolling or shaking.	Section 3.3.3Feeding Solution
Reconstitution Buffer Bottle 6, white cap	 Ready-to-use solution The solution is stable at 4°C but can also be stored at -20°C 	

3.3.2 Reaction Components from the RTS Amino Acid Sampler

Component	Reconstitution/Preparation of working solution	For use in
Amino acids labeling mix Vial 1-20	Prepare a 2.1 mM stock solution of all amino acids in 3 ml of Reconstitution Buffer as described under Section 3.2.	Section 3.3.3. Reaction Solution and Feeding Solution
DTT Vial 21	Ready-to-use solution Open each vial only once Repeated freezing and thawing is allowed	Section 3.3.3. Reaction Solution and Feeding Solu- tion

3.3.3 Preparation of working solutions

Content	Reconstitution/Preparation of working solution
Feeding Solution	To the content of bottle 3 (Feeding Mix) add 2.65 ml of the amino acids labeling mix (see section 3.2.) and 0.3 ml of DTT-solution (Vial 21). Mix by rolling or shaking.
Reaction Solution	To the content of bottle 1 (<i>E.coli</i> Lysate) add 0.225 ml of the reconstituted Reaction Mix (bottle 2), 0.27 ml of the amino acid labeling mix (see Section 3.2.) and 30 μ l of DTT solution (Vial 21). Add 10–15 μ g of DNA template in a maximum volume of 50 μ l. Mix carefully by rolling or gentle shaking. Do not vortex!

3.3.4 Running an experiment

Running the experiment After preparation of the working solutions proceed with the procedure for loading the RTS 500 device as described in the pack insert of the RTS 500 *E.coli* 500 HY Kit under "Running an experiment".

4. **Typical results**

Using a mixture of amino acids from the RTS Amino Acid Sampler, the protein synthesis rate was compared to the synthesis rate obtained with the Amino Acid Mix without Methionine (bottle 4) and Methionine (bottle 5) from the RTS 500 E.coli HY Kit. The amount of protein produced in the 1 ml standard reaction is similar (see Fig. 3)



0.1 µl sample

standard



Fig. 3: Comparision of RTS reactions performed with the ready-to-use amino acid mixture contained in the RTS 500 E.coli HY kit or with amino acids provided separately in the Amino Acid Sampler. After running the RTS reaction, various amounts of sample were loaded in duplicates on a Coomassie-stained SDS-Gel. A GFP standard was loaded in parallel for subsequent densitometric analysis.

continued on next page

To test the incorporation of a labeled amino acid, an amino acid mixture containing ¹⁵N-Glycine (Isotech, Miamisburg, OH) was prepared according to the procedures described in this pack insert. The final concentration of labeled amino acid in the RTS reaction device Solution was 0.5mM. After synthesis the produced protein was analyzed on SDS-PAGE (See Fig. 4), showing the production of protein of expected length.



Fig. 4: Expression of phosphoserine phosphatase (PSP) in RTS (A) Unlabeled protein from a total RTS 500 *E.coli* HY reaction mix (lane 1) and after purification (lane 2). M: 10 kD molecular weight marker. (B) ¹⁵N-Glycine labeled protein in a total reaction mix (lane 4) compared to a reaction without expression template (lane 3) Note: this experiment was performed with a slightly modified proteol and a with a prototype of the amino acid sampler .

5. Appendix

5.1 Trouble shooting

Problem	Possible Cause	Recommendation	
No protein produced	The protein is difficult to express in the Rapid Trans- lation System	Test expression first in the RTS 100 <i>E.coli</i> HY Kit. Follow the recommenda- tions for increasing success rate provided with the kit	
No protein produced although RTS 100 protein production was positive	The amino acid labeling mix is missing one or more amino acids	Repeat the experiment using freshly prepared mixtures. Compare expres- sion with the expression using the amino acids provided in the RTS 500 <i>E.coli</i> HY Kit.	
Some amino acids show precipitation after thawing		Warm the tubes to 37– 42°C. Precipitates must dissolve.	

5.2 References

- 1 Riek, R., Pervushin, K., Wuethrich, K. (2000). TROSY and CRINEPT: NMR with large molecular and supramolecular structures in solution *TIBS* **25** (2000), p. 462.
- 2 Gardner, K.H., Kay, L. E. (1998): The use of H, C, N multidimensional NMR to study the structure and dynamics of proteins. *Annu Rev Biophys Biomol Struct* 27 (1998), p. 357.
- 3 Kigawa T, Yabuki T, Yoshida Y, Tsutsui M, Ito Y, Shibata T, and Yokoyama S (1999): Cell-free production and stable isotope labeling of milligram quantities of proteins. *FEBS Lett* **442** (1999), pp.15-19.

5.3 How to contact Roche Molecular Biochemicals

Ways to contact usTo contact Roche Molecular Biochemicals for technical assistance, please choose one of the following:

IF you are using	THEN
Internet	access our web-site at: http://www.proteinexpression.com, or http://biochem.roche.com.
E-mail	Please refer to the address that corre- sponds to your particular location, on the last page of this instruction manual.
Telephone	Please refer to the address that corre- sponds to your particular location, on the last page of this instruction manual.

5.4 Related products

Product	Pack Size	Cat. No.
Rapid Translation System RTS 500 Instrument	1 instrument	3064 859
Rapid Translation System RTS 100 <i>E.coli</i> HY Kit	24 reactions (50 μl each)	3186 148
Rapid Translation System RTS 500 <i>E.coli</i> HY Kit	2 reactions 5 reactions	3246 817 3246 949
Rapid Translation System RTS pIVEX His-tag Vector Set	5 vectors (10 μg each)	3253 538
Rapid Translation System RTS GroE Supplement	5 × 125 μl	3263 690
Restriction Protease Factor Xa	3 × 30 μg 3 x 100 μg 3 x 250 μg	1179 888 1585 924 1179 896
Restriction Protease Factor Xa Cleavage and Removal Kit	Kit I : 3 × 30 μg Kit II: 3 × 100 μg	1644 777 1644 785
rGFP	50 µg	1814 524
Anti-His ₆	100 µg	1922 416
Anti-His ₆ -Peroxidase	50 U	1965 085
Complete Protease Inhibitor Cocktail Tablets, mini, EDTA-free	25 tablets (each tablet for 10 ml extract)	1836 170

	E-mail Address	Country	E-mail Address	Country
	argentina.biochem@roche.com	Argentina	Raitis@invitros.lv	Latvia
	biochem.au@roche.com	Australia	Sakkijha@rdleb.com	Lebanon
	Gerhard.Muehlbauer@roche.com	Austria	Gintaras@eksma.lt	Lithuania
	biochem.be@roche.com	Belgium	diagnostics@prophac.lu	Luxembourg
	Valent@mbox.cit.bg	Bulgaria	Vccl@vol.net.mt	Malta
	africhem@camnet.cm	Cameroon	Aiouche.echo@dounia.net.ma	Morocco
	biochem.ca@roche.com	Canada	biocheminfo.nl@roche.com	Netherlands
	biochem.cn@roche.com	China	biochem.nz@roche.com	New Zealand
	Info@medisell.com.cy	Cyprus	bofungwu@linkserve.com.ng	Nigeria
	Bm-comp@bm-comp.cz	Czech Republic	biochem.se@roche.com	Norway
<	dk.biochem@roche.com	Denmark	biochem.pt@roche.com	Portugal
E	ou.melestrum@neti.ee	Estonia	Topdiag@fx.ro	Romania
9	pharsc.et@telecom.net.et	Ethiopia	biochem.sg@roche.com	Singapore
99	helsinki.biochem_diagnostics@roche.com	Finland	roche.diagnostics@siol.net	Slovenia
2	biochem.fr@roche.com	France	south_africa.bioboffin@roche.com	South Africa
E .	mannheim.biocheminfo@roche.com	Germany	biochem.es@roche.com	Spain
ñ	Bm_roche@hotmail.com	India	biochem.se@roche.com	Sweden
12	h.hajian@tebtech.com	Iran	BiochemInfo.CH@roche.com	Switzerland
æ	tubanegin@istn.irost.com	Iran	Jean-Marie.kindbeiter@roche.com	Tunisia
5	Dyn@netvision.net.il	Israel	bmuae@emirates.net.ae	United Arab Emirates
8	it.biochem@roche.com	Italy	uk.biochem@roche.com	United Kingdom
	biochemicals@rdj.co.jp	Japan	biochemts.us@roche.com	USA
69	pharmakp@net2000ke.com	Kenya	Mvalentiner@telcel.net.ve	Venezuela
2	Bmskorea@chollian.net	Korea	dusica@eunet.yu	Yugoslavia
<u><u> </u></u>	react@ncc.moc.kw	Kuwait	biochemts.row@roche.com	All other countries
.02		_		
ö	http://hiochom	rocho con	n/nack_incort/276716/19	ndt

http://biochem.roche.com/pack-insert/3262154a.pdf

Argentina 541 954 5555; Australia (02) 9899 7999; Austria (01) 277 87; Belgium (02) 247 4930; Brazil +55 (11) 3666 3565; Bulgaria +35929625408; Cameroon 237-370269; Canada (450) 086 7050; (800) 361 2070; Chile 00 56 (2) 22 33 737 (central) 00 56 (2) 22 32 099 (Exec); China 86 21 6427 5586; Colombia 0057-1-3412797; Cyprus +357-2-311362; Czech Republic (0324) 45 54, 58 71-2; Denmark +45 363 999 58; Egypt 20-2-3619047; Estonia 372-7-447600; Ethiopia 251-1-552799; Finland +358 9 525 333 66; France 04 76 76 30 87; Germany (0621) 759 8568; Greece 3 (01) 67 40 238; Hong Kong (852) 2485 7596; India +91-22-8379906; Indonesia 62 (021) 252 3820 ext. 755; Iran +98-21-8072374 / +98-21-8797027; Israel 972-6-6380569; Italy 039 247 4109-4181; Japan 03 5443 5284; Kenya +254-2-750112; Korea 82-2-3471-6500; Kuwait +965-4837859; Latvia 371-787828309; Lebanon Fax: 00961-1-399667; Lithuania 370-2-729715; Luxembourg +352-496098; Malta Fax: +356-341087; Morocco Fax: +212-2-944040; Malaysia 60 (03) 755 5039; Mexico (5) 227 8967; Netherlands (036) 539 4911; New Zealand (09) 276 4157; Nigeria +234-1-521767; Norway (47) 23 373300; Philippines (632) 810 7246; Poland +48 (22) 22 66 84 305; Portugal (01) 4171717; Republic of Ireland 1 800 40 90 41; Romania +40-1-2123763; Russia (49) 621 759 8636 Fax: (49) 621 759 8611; Saudia Arabia +966-1-4010364; Singapore 0065 272 9200; Slovenia +386 61 1363528; South Africa (011) 886 2400; South Korea 02 569 6902; Spain (93) 201 4411; Sweden (08) 404 8800; Switzerland +41 (41) 799 6161; Taiwan (02) 736 7125; Thailand 66 (2) 274 07 08 (12 line); Turkey 0090 212 216 32 80: United Arab Emirates +971-4-694351: United Kingdom (0800) 521578: USA (800) 428 5433. Venezuela Fax: +0058-4810697; Yugoslavia +381 11 137163.



Roche Diagnostics GmbH Roche Molecular Biochemicals Sandhofer Strasse 116 D-68305 Mannheim Germany