

# Technical Review

## • Gel Stains



FEATURING ...

GelCode<sup>®</sup> Blue  
Stain Reagent

GelCode<sup>®</sup>  
SilverSNAP<sup>™</sup> Stain


GelCode<sup>®</sup> E-Zinc<sup>™</sup>  
Reversible Stain Kits

PIERCE

**AL**  
PERBIO

# Protein Staining

## FEATURE STORY



### Water Wash-Enhanced Protein Staining With GelCode® Blue Stain Reagent

#### Introduction

Coomassie® Brilliant Blue is the most common dye used to stain proteins on polyacrylamide gels. The dye was originally developed as an acid wool dye, but found favor among life scientists because of its ease of use and relatively good sensitivity. Although alternative methods, such as silver staining,<sup>1</sup> heavy metal staining<sup>2</sup> and fluorescent dye stains<sup>3</sup> have been developed, they are not as easy to use as Coomassie® stains. These methods have only been used for some specific applications.

Coomassie® staining requires an acidic environment to enhance the ionic interaction between the dyes and basic amino acids. Classically, the R-250 form of the dye is dissolved in a solution containing 20%-45% methanol and 5%-10% acetic acid at concentrations between 0.1%-0.25%.<sup>4</sup> After the gel is stained, protein bands are not visible before destaining with methanol/acetic acid to remove the characteristic background staining. The main drawback to traditional Coomassie® Brilliant Blue staining methods is the long time period required for staining and destaining. Overnight destaining periods are usually necessary to achieve a clear background. It is common that the staining intensity of protein bands is decreased by prolonged

destaining periods, resulting in reduced sensitivity. Using this method, a gel typically cannot be analyzed until the day after it is run.

The G-250 form of Coomassie® Brilliant Blue is typically used as a colloidal suspension. The colloidal particles do not penetrate the gel, but the dye molecules are extracted from the colloid by the proteins in the gels. The proteins become tinted by the dye molecules and, because there is little coloration in the bulk of the gel, destaining is not required. A method developed by Neuhoff, V., *et al.*<sup>5</sup> has improved the sensitivity of Coomassie® Brilliant Blue G-250 stain to the nanogram level. The main drawbacks of this method is that it requires long (12-hour) staining times to reach a maximum sensitivity, and the band signal-to-noise ratio is low.

Pierce recently developed a Coomassie® Brilliant Blue G-250-based protein stain, GelCode® Blue Stain Reagent, which uses the colloidal properties of Coomassie® G-250 for protein staining in polyacrylamide gels. This exclusive Pierce formulation yields staining properties unique to a Coomassie® Brilliant Blue G-based protein stain. For example, only the protein is stained, not the gel, allowing protein bands to be viewed directly within the staining reagent during the staining process. In addition, the staining sensitivity reaches its maximum within one hour. Longer (overnight) staining does not increase the background. In addition, a simple water wash after staining significantly increases the staining sensitivity, and the gel background is crystal-clear without having to destain with methanol/acetic acid.

#### Materials and Methods

##### 1) Protein Electrophoresis

Protein molecular weight standards consisting of myosin H-chain (200 kDa), phosphorylase b (97.4 kDa), bovine serum albumin (68 kDa), ovalbumin

*Only the protein is stained, not the gel, allowing protein bands to be viewed directly within the staining reagent during the staining process.*



(43 kDa), carbonic anhydrase (29 kDa),  $\beta$ -lactoglobulin (18.4 kDa) and lysozyme (14.3 kDa) were reconstituted with deionized water at an approximate concentration of 1 mg/ml of each protein. The reconstituted proteins were further diluted serially in 1X SDS-PAGE buffer (50 mM Tris•HCl, pH 6.8, 100 mM dithiothreitol, 2% SDS, 0.1% bromophenol blue, 10% glycerol) and boiled for 5 minutes to denature the protein samples. Samples were loaded identically onto two 4%-20% gradient, 8 cm x 7.5 cm mini gels, and run with constant current at 20 mA for each gel until the bromophenol blue dye reaches the bottom.

### II) Protein Staining

After electrophoresis, one of the gels was rinsed with 200 ml of deionized water, 3 times for 5 minutes and then stained with 20 ml of GelCode® Blue Stain Reagent for 1 hour. After staining, the staining reagent was replaced with 200 ml of deionized water. A gentle agitation was used during the entire process to ensure that the gel was treated uniformly. The other gel was stained with a commercial staining reagent containing 0.25% Coomassie® Blue G-250, 25% methanol and 5% acetic acid for one hour and destained with several changes of 25% methanol and 5% acetic acid.

### III) Staining of Proteins After Western Blot Transfer

Approximately 50  $\mu$ g of bacterial *E. coli* cell lysates expressing polyhistidine-tagged (PHT) mouse dihydrofolate reductase and baculovirus polyhedrin protein respectively, were lysed with SDS-PAGE loading buffer (1X) and boiled for 5 minutes before loading onto two 4%-20% gradient mini gels. After electrophoresis, one of the gels was stained with GelCode® Blue Stain Reagent as previously described, and the other gel was transferred with a semi-dry transfer unit for 1 hour at a current of 0.8 mA/cm<sup>2</sup>. After transfer, the membrane was

processed for detection of PHT proteins using INDIA™ HisProbe™-HRP\* (Product # 15165) following the manufacturer's instructions. The transferred gel was rinsed once with deionized water for 20 seconds and stained with GelCode® Blue Stain Reagent for one hour, followed by equilibration in water.

## Results

### Protein Bands Visible Without Destaining


Colloidal staining of proteins with GelCode® Blue Stain Reagent following PAGE allows highly sensitive and background-free staining. Standard proteins were serially diluted, run in 4%-20% SDS-PAGE gels and washed three times with deionized water before staining with GelCode® Blue Stain Reagent. The colloidal dye stains protein rapidly; protein bands at concentration above 1  $\mu$ g were visible within only a few minutes in the staining reagent. With one-hour staining, the staining intensity of protein bands reaches a maximum and standard proteins at concentrations above 250 ng were clearly visible (Figure 1A). Unlike other staining systems, the GelCode® Blue Stain Reagent only stains protein; the bulk of the gel remains virtually clear. The background staining did not increase even with an overnight staining (data not shown).

### Water Wash Enhancement™ Step

After a one-hour staining, the staining reagent was replaced with deionized water. The staining sensitivity was significantly increased after placing the gel in deionized water for a period of 1-2 hours. On the water-equilibrated gel, standard protein bands at a concentration of 31 ng were clearly visible. Some of these standard proteins were even visible at concentrations as low as 8 ng (Figure 1B). Variations were observed among different standard proteins. This result could be due to discrepancy in

*The staining sensitivity was significantly increased after placing the gel in deionized water for a period of 1-2 hours.*



the amount of each individual protein in the standard or because individual proteins interact with the dye differently. In contrast, no protein bands were visible on the gel stained with the commercial staining reagent for one hour (Figure 1C) before destaining, and the gel retained very high background even after destaining with 25% methanol/5% acetic acid for 1 hour (Figure 1D).



**Figure 1.** Comparison of GelCode® Blue Stain Reagent With a Classical Coomassie® Stain Formulation. Standard proteins were serially diluted and run on two 4%-20% gradient gels. One gel was stained with GelCode® Blue Stain Reagent for one hour (A), and then washed with deionized water for one hour (B). The other gel was stained with 0.25% Coomassie® Blue G-250 in 25% methanol and 5% acetic acid for one hour (C), and then destained with 25% methanol and 5% acetic acid for one hour (D). The amount of each standard protein, as described in the text, was approximately 3,000 ng, 2,000 ng, 1,000 ng, 500 ng, 250 ng, 125 ng, 62 ng, 31 ng, 16 ng and 8 ng, in each lane (left to right).

### Linear Range Staining

Scanning the visible spectrum of the stained gel using a Bio-Rad Molecular Imager® GS-700 indicated that most of the standard proteins, such as bovine serum albumin, ovalbumin, carbonic anhydrase, lactoglobulin and lysozyme, responded to the staining reagent at a linear range between 61 ng to 2,000 ng per band, while myosin H-chain and phosphorylase b exhibited a linear range between 500 ng to 3,000 ng (Figure 2).




**Figure 2.** The Linear Range of Standard Proteins Stained With GelCode® Blue Stain Reagent. The gels shown in Figure 1B were scanned within the visible spectrum using a Bio-Rad Molecular Imager® GS-700. The data was processed with Microsoft Excel and converted to graphs. Figure 2B is the inset of Figure 2A at concentrations between 8 ng to 250 ng.

### Convenient Protein Staining After Western Blot Transfer

It has been observed that high concentrations of SDS in the PAGE gel interfere with the GelCode® Blue Stain Reagent. A simple water wash is usually sufficient to remove SDS from the gel, achieving excellent results. For Western Blot-transferred gels, most of the SDS has been removed during transfer. Therefore, this washing step can be simplified to a

*Following a one-hour staining and a one-hour water equilibration, trace amounts of proteins remaining on the gel can be easily detected.*



**Figure 3.** Staining of SDS-PAGE Gel With GelCode® Blue Stain Reagent After Western Blot Transfer. *E. coli* lysates were run on two 4%-20% SDS-PAGE gels. After electrophoresis, one gel was washed three times with deionized water and stained with GelCode® Blue Stain Reagent as described in the text (A). The other gel was subjected to Western Blot transfer. The transferred membrane was probed with INDIA™ HisProbe™-HRP\*

(Product # 15165) to detect six histidine-tagged recombinant proteins (B) and the transferred gel was rinsed with deionized water and stained with GelCode® Blue Stain Reagent as described in the text (C). Lanes 1 and 4 were the uninduced *E. coli* lysates. Lanes 2, 3, 5 and 6 were *E. coli* lysates expressing PHT mouse dihydrofolate reductase and baculovirus polyhedrin protein, respectively. Lanes M were standard protein markers.

single rinse with deionized water. Protein bands on a transferred gel are developed more rapidly than protein bands on untransferred gel during the staining process. Following a one-hour staining and a one-hour water equilibration, trace amounts of proteins remaining on the gel can be easily detected (Figure 3). This feature is particularly useful for evaluating protein loading and for checking if proteins on SDS-PAGE gels are uniformly and completely transferred.

## Conclusion

As demonstrated by our results, the GelCode® Blue Stain Reagent represents a new generation of Coomassie® Blue Stain with several unique properties: (1) easily viewed protein bands during the staining process, (2) elimination of the methanol/acetic acid destaining step and (3) enhancement of staining sensitivity by a simple water wash step. Other benefits of GelCode® Blue Stain Reagent include a crystal-clear gel background, no gel shrinkage, significantly shortened timeline to a result, timesavings and reduced costs because there is no need to stock and formulate methanol and acetic acid for destaining purposes. Although high concentrations of SDS present in the gel interfered with the reagent, a simple water wash usually is sufficient to remove SDS from the gel. Even without prewashing, the background can

be easily removed by the final water wash. It was found that different prefixing and washing steps only affect the band visibility during the staining process. The sensitivity remains virtually the same after the gel has been equilibrated in water.

## REFERENCES

- 1.] Sammons, D.W., Adams, L.D. and Nishizawa, E.E. (1981). Ultrasensitive silver-based color staining of polypeptide in polyacrylamide gels. *Electrophoresis* **2**, 135-141.
- 2.] Fernandez-Patron, C., Castellanos-Serra, L. and Rodriguez, P. (1982). Reverse staining of sodium dodecyl sulfate polyacrylamide gels by imidazole-zinc salts: sensitive detection of unmodified proteins. *BioTechniques* **12**, 564-573.
- 3.] Javier Alba, F., Bermudez, A., Bartolome, S. and Daban J.-R. (1996). Detection of five nanograms of protein by two minute Nile red staining of unfixed SDS gels. *BioTechniques* **21**, 625-626.
- 4.] Bollag, D.M., Rozycki, M.D. and Edelstein S.J. (1996). *Protein Methods*, 126-134, Wiley-Liss, Inc., NY.
- 5.] Neuhoﬀ, V., Arold, N., Taube, D. and Ehrhardt, W. (1988). Improved staining of proteins in polyacrylamide gels including isoelectric focusing gels with clear background at nanogram sensitivity using Coomassie® Brilliant Blue G-250 and R-250. *Electrophoresis* **9**, 255-262.

## Ordering Information

Product #	Description	Pkg. Size
24590	<b>GelCode® Blue Stain Reagent</b> Sufficient for staining up to 25 SDS-PAGE 8 cm x 10 cm mini gels	500 ml
24592	<b>GelCode® Blue Stain Reagent</b> Sufficient for staining up to 175 SDS-PAGE 8 cm x 10 cm mini gels	3.5 liters
72300	<b>Pump</b> (fits 3.5 liter package only) <b>FREE, upon request, with purchase of Product # 24592GX</b>	1 pump

Coomassie® is a registered trademark of ICI Americas.

\*Patent Pending

*Benefits of  
GelCode® Blue Stain  
Reagent include a  
crystal-clear gel  
background,  
no gel shrinkage,  
significantly short-  
ened timeline to a  
result, time-savings  
and reduced costs.*



# GelCode® Stain Kits

## FEATURE STORY

### Rapid and Ultrasensitive PAGE Gel Stains:

#### Introducing GelCode® SilverSNAP™ Stain and GelCode® E-Zinc™ Reversible Stain Kits

Coomassie® stain,<sup>1</sup> silver stain and heavy metal (Zn) stain are often used for protein detection in polyacrylamide gels following electrophoresis. Although Coomassie® is the dye most frequently used for protein staining in polyacrylamide gels, silver staining has a tendency to be much more sensitive. The main disadvantage associated with common Coomassie® stains is that they require a time-consuming and foul-smelling methanol-acetic acid destaining procedure. GelCode® Blue Stain Reagent, a new Coomassie® stain developed by Pierce Chemical Company, does not require cumbersome destaining steps.<sup>1</sup> The reagent stains only protein,

not the gel, resulting in a crystal-clear background after a Water Wash Enhancement™ Step. Thus GelCode® Blue Stain Reagent is the ideal stain for many applications. In addition to GelCode® Blue Stain Reagent, Pierce recently developed two new related kits: GelCode® SilverSNAP™ and GelCode® E-Zinc™ Reversible Stains.

#### GelCode® SilverSNAP™ Stain Kit Provides Superior Sensitivity

The lower detection limit of GelCode® Blue Stain Reagent is 8-10 ng, which is sufficient for most applications. In some applications such as two-dimensional (2-D) PAGE gel staining, however, a higher sensitivity is required. The GelCode® SilverSNAP™ Stain Kit (Product # 24602) provides greater sensitivity than Coomassie® stain kits. The kit consists of three ready-to-use components and a simple procedure that can be completed within 90 minutes. As shown in Figure 1, as little as 0.25 ng of protein molecular weight standards separated by a 4%-20% SDS-PAGE gel are clearly visible. The yellowish background associated with common

*GelCode® Blue Stain Reagent does not require cumbersome destaining steps.*

### Three GelCode® Stain Products are Complementary

GelCode® Blue Stain Reagent, GelCode® E-Zinc™ Reversible Stain and GelCode® SilverSNAP™ Stain are designed for different applications. Table 1 shows


each stain's characteristics. No matter what your application, one of these stains will deliver satisfactory results.

Table 1. Comparison of Three GelCode® Stain Products

Stains	GelCode® E-Zinc™ Reversible Stain	GelCode® SilverSNAP™ Stain	GelCode® Blue Stain Reagent
Number of Components	3	3	1
Number of Steps	2	3	2
Staining Time (minutes)	15	35	60
Sensitivity (ng)	0.25	0.25	8
Applications	<ul style="list-style-type: none"><li>• Protein elution from PAGE gels</li><li>• Staining before Western blotting</li><li>• Diluted protein samples</li><li>• Purified protein samples</li></ul>	<ul style="list-style-type: none"><li>• Two-dimensional (2-D) gels</li><li>• Diluted protein samples</li><li>• Purified protein samples</li><li>• DNA and RNA</li></ul>	<ul style="list-style-type: none"><li>• SDS-PAGE gels</li><li>• Native PAGE gels</li><li>• PhastGel® and IEF gels</li><li>• PVDF membranes</li></ul>



silver stain is eliminated with the GelCode® SilverSNAP™ Stain. In comparative experiments, some proteins that were not detectable by the Silver Stain from Vendor B are clearly visible with GelCode® SilverSNAP™ Stain. Thus the GelCode® SilverSNAP™ Stain Kit is more reliable. This kit has also been tested for DNA and RNA staining on PAGE gels.




**Figure 1.** Comparison of GelCode® SilverSNAP™ Stain With a Silver Stain From Vendor B. Protein molecular weight standards, as indicated at the top of the gel, were separated on 15-well pre-cast 4%-20% SDS-PAGE gels. For GelCode® SilverSNAP™ Stain (top), following fixing and washing, the gel was stained with SilverSNAP™ Stain for 30 minutes, developed with SilverSNAP™ Developer for approximately 5 minutes, and stopped with 5% acetic acid. For Silver Stain (Vendor B) (bottom), the vendor's procedure was followed. Note that GelCode® SilverSNAP™ Stain not only provides higher sensitivity, but also stains all proteins. Because all GelCode® SilverSNAP™ Stain Kit components are ready-to-use, the background caused by mixing multiple components – as seen when the vendor B stain is used – is eliminated.

## GelCode® E-Zinc™ Reversible Stain Kit is Ideal for Recovering Proteins From PAGE Gels and Western Blots After Staining


The length of time required to use Coomassie® stains and silver stains is sometimes cited as a drawback of both techniques. In addition, once the protein is stained, it is difficult to remove the dye from the protein-dye complex. The dye-contaminated protein is not appropriate for downstream applications, such as antibody generation and/or immunological detection. As shown in Figure 2, GelCode® E-Zinc™ Reversible Stain Kit (Product # 24582) is ideal

*E-Zinc™ Reversible Stain Kit is ideal for rapid detection of proteins on a PAGE gel with a sensitivity equivalent to silver stain.*



**Figure 2.** Comparison of GelCode® E-Zinc™ Reversible Stain With Zn Stain (Vendor B), GelCode® Blue Stain Reagent and GelCode® SilverSNAP™ Stain. Protein molecular weight standards, the same as those used in Figure 1, were separated on 15-well pre-cast 4%-20% SDS-PAGE gels. For Zn stain, the gels were stained with the stain solution for 10 minutes and developed with the developer solution for approximately 2 minutes. GelCode® SilverSNAP™ Stain was used as indicated in Figure 1. GelCode® Blue Stain Reagent was used as described earlier.<sup>1</sup>

for rapid detection of proteins on a PAGE gel with a sensitivity equivalent to silver stain. The two-step procedure takes less than 15 minutes and the stain can be removed completely by the E-Zinc™ Eraser. For example, Figure 3 demonstrates the use of the GelCode® E-Zinc™ Reversible Stain Kit for the purification of a recombinant protein from PAGE gel.



**Figure 3.** Purification of Inclusion Body Protein From PAGE Gel After GelCode® E-Zinc™ Reversible Stain. The inclusion body of a GST fusion protein expressed in *E. coli* was extracted and purified by B-PER® Bacterial Protein Extraction Reagent (Product # 78248). As some of the bacterial proteins are associated with the GST-TSP recombinant protein within the inclusion body, the inclusion body preparation was further purified through an SDS-PAGE gel to obtain homogeneous recombinant protein. On a preparative SDS-PAGE gel, 4 mg of crude GST-TSP inclusion body was dissolved in 1X SDS-PAGE sample buffer, loaded and separated from contaminants. The preparative gel was stained with GelCode® E-Zinc™ Reversible Stain Kit and the fusion protein band was excised. After removing Zn by destaining with the EZ-Zinc™ Stain Eraser, the GST fusion protein was electro-eluted from the gel and the purity was analyzed by loading various amounts of the eluted sample on an SDS-PAGE mini gel stained with GelCode® Blue Stain Reagent. **Lane 1:** uninduced total bacterial protein; **Lane 2:** induced bacterial protein; **Lane 3:** the soluble protein extracted by B-PER® Reagent; **Lanes 4-5:** two washes with 1:10 diluted B-PER® Reagent; **Lane 6:** purified inclusion body protein; and **Lanes 7-9:** 1, 2 and 4 µg of GST fusion protein electro-eluted from PAGE gel after Zn stain. Note that the insoluble GST fusion protein seen in Lane 2 is not present in Lane 3.

By using the reversible stain, the final protein preparation is free of dye contamination and is ideal for antibody generation. In addition to reversibility, GelCode® E-Zinc™ Reversible Stain may also be used for rapid analysis of purified protein and the composition of protein components in a diluted sample on PAGE gels.

#### REFERENCE

1.J. Chu, R. and Vigna, R.A. (1997). Water wash-enhanced protein staining with GelCode® Coomassie® Blue Stain Reagent. *Previews* **1**(4), 18-21.

#### Ordering Information

Product #	Description	Pkg. Size
24590	<b>GelCode® Blue Stain Reagent</b> Sufficient for staining up to 25 SDS-PAGE 8 cm x 10 cm mini gels.	500 ml
24592	<b>GelCode® Blue Stain Reagent</b> Sufficient for staining up to 175 SDS-PAGE 8 cm x 10 cm mini gels.	3.5 liters
24602	<b>GelCode® SilverSNAP™ Stain Kit</b> Sufficient for staining up to 20 SDS-PAGE 8 cm x 10 cm mini gels. Includes: SilverSNAP™ Stain SilverSNAP™ Developer SilverSNAP™ Enhancer	Kit 500 ml 500 ml 25 ml
24582	<b>GelCode® E-Zinc™ Reversible Stain Kit</b> Sufficient for staining up to 80 SDS-PAGE 8 cm x 10 cm mini gels. Includes: E-Zinc™ Stain E-Zinc™ Developer E-Zinc™ Eraser	Kit 500 ml 500 ml 500 ml

Coomassie® is a registered trademark of ICI Americas.

PhastGel® is a registered trademark of Pharmacia LKB Biotechnology AB.

#### Recommended reading or information on the use of GelCode® Blue Stain Reagent for zymogram staining:

1.J. Meade-Tollin, L.C. (1998). Rapid staining of gelatin zymograms with GelCode® Blue Stain Reagent. *Previews* **2**(2), 2-4.


No matter what  
your application,  
one of these  
stains will deliver  
satisfactory results.



## Gel Forming Reagents, Markers and Staining Kits

### Ammonium Persulfate

*A useful catalyst for acrylamide gel polymerization.*



#### Ordering Information

Product #	Description	Pkg. Size
17874	Ammonium Persulfate	4 x 25 gm

### ImmunoPure® Lane Marker Sample Buffers

*Ready-to-use lane markers that can be used for Western blots!*

#### Features/Benefits:


- Bright pink hydrophobic tracking dye for SDS-PAGE
- Ability to see proper transfer of proteins from gel to nitrocellulose membrane
- Easy cutting of sample lanes when different antibody systems are used in the “shaker” method
- Dye front is easily shown on both the gel and the nitrocellulose membrane for determination of molecular weight (Rf values)
- Easy alignment of sample lanes containing transferred protein in the Pierce Fast Blot-Developer™ Systems
- DTT is used in place of β-Mercaptoethanol in the Reducing Sample Buffer, so strong odors and auto-oxidation are avoided

#### Ordering Information

Product #	Description	Pkg. Size
39000	<b>ImmunoPure® Lane Marker Reducing Sample Buffer</b> 0.3 M Tris•HCl, pH 6.8; 5% SDS; 50% Glycerol; 100 mM Dithiothreitol (DTT); Lane Marker Tracking Dye	5 ml
39001	<b>ImmunoPure® Lane Marker Non-Reducing Sample Buffer</b> 0.3 M Tris•HCl, pH 6.8; 5% SDS; 50% Glycerol; Lane Marker Tracking Dye	5 ml

### 2-Mercaptoethanol (2-ME)

*Reduces dithiol bonds in peptides and proteins.*




#### Ordering Information

Product #	Description	Pkg. Size
35600	2-Mercaptoethanol (2-ME)	500 gm

### SDS

*Useful for resolving viral proteins during gel electrophoresis and for when renaturation after SDS-PAGE is required.*



#### Typical Alkyl Chain Length Distribution for Pierce SDS Lauryl Grade

C <sub>10</sub>	C <sub>12</sub>	C <sub>14</sub>	C <sub>16</sub>
—	63.5%	29.5%	7.0%

#### Ordering Information

Product #	Description	Pkg. Size
28364	SDS (Sodium Dodecyl Sulfate, Lauryl)	100 gm
28365	SDS (Sodium Dodecyl Sulfate, Lauryl)	1 kg

## SDS C<sub>12</sub>, Sequanal Grade

Contains very low levels of C<sub>16</sub> alkyl sulfate.

### Specifications:

- White powder
- Greater than 99% alkyl sulfate
- Greater than 98% C<sub>12</sub> alkyl sulfate
- 10% solution is clear, colorless and free of particulate matter

### Typical Alkyl Chain Length Distribution for Pierce C<sub>12</sub> Sequanal Grade

C <sub>10</sub>	C <sub>12</sub>	C <sub>14</sub>	C <sub>16</sub>
0.4%	99.4%	0.2%	0.0%

### Ordering Information

Product #	Description	Pkg. Size
28312	SDS C <sub>12</sub> , Sequanal Grade	500 gm

See the Detergents section for more information on SDS.

## TEMED

Greater than 99% pure!



TEMED  
M.W. 116.21

### Specifications:


- Purity: > 99.9%
- Refractive Index: 1.417-1.419
- Boiling Range: 119-121°C

### Ordering Information

Product #	Description	Pkg. Size
17919	TEMED (N,N,N',N'-Tetramethylethylenediamine)	25 ml

## Urea (Sequanal Grade)

A low UV-absorbing protein denaturant.



Urea  
M.W. 60.06

### Specification:

- A<sub>280</sub> nm: < 0.100

### General Electrophoresis References

1. An der Lan, B.C. (1988). *Methods for Protein Analysis*. American Oil Chemists' Society, Champaign, IL, Chapter 5. Options for Gel Electrophoresis, pages 52-69.
2. Hames, B.D. and Rickwood, D., eds. (1984). *Gel Electrophoresis of Proteins: A Practical Approach*. Washington D.C.: IRL Press.
3. Lacks, S.A., et al. (1979). Effect of the composition of sodium dodecyl sulfate preparations on the renaturation of enzymes after polyacrylamide gel electrophoresis. *Anal. Biochem.* **100**, 357-363.

### Ordering Information

Product #	Description	Pkg. Size
29700	Urea (Sequanal Grade)	1 kg

Tel: 800-874-3723 or 815-968-0747 • Fax: 815-968-7316 • E-mail: [TA@piercenet.com](mailto:TA@piercenet.com) • Internet: [www.piercenet.com](http://www.piercenet.com)