

# Pierce<sup>®</sup> Protein L Agarose

20510 20512

0778.5

Number	Description
20510	<b>Pierce Protein L Agarose</b> , 2 ml settled resin
20512	<b>Pierce Protein L Agarose</b> , 10 ml settled resin Support: Crosslinked 6% beaded agarose supplied as 50% slurry (e.g., 2 ml of settled resin is equivalent to 4 ml of 50% slurry) containing 0.02% sodium azide Binding Capacity: ~5-10 mg human IgG/ml of resin

**Storage:** Upon receipt store product at 4-8°C. Product is shipped at ambient temperature.

## Table of Contents

Introduction .....	1
Important Product Information .....	1
Gravity-flow Column Procedure for Antibody Purification .....	2
A. Additional Materials Required .....	2
B. Immunoglobulin Purification Procedure .....	2
Troubleshooting .....	3
Additional Information .....	3
Related Thermo Scientific Products .....	3
General References .....	4
Product References .....	4

## Introduction

Protein L is an immunoglobulin-binding protein that was isolated from the bacteria *Peptostreptococcus magnus* and is now produced recombinantly. Protein L binds to immunoglobulin kappa light chains without interfering with the antigen-binding site and binds a wider range of Ig classes and subclasses than other antibody-binding proteins such as Protein A or Protein G. Protein L binds to all classes of Ig (i.e., IgG, IgM, IgA, IgE and IgD). Protein L also binds single chain variable fragments (Scfv) and Fab fragments. The Thermo Scientific Pierce Protein L Agarose is prepared using a leak-resistant coupling method that ensures excellent resin stability and binding characteristics.

## Important Product Information

- Protein L **only** binds to immunoglobulins containing light chains of type kappa I, III, and IV in human and kappa I in mouse. Protein L also may be specific for certain kappa subgroups in other species. Protein L binds scfv without interfering with antigen binding.
- Protein L binds weakly to rabbit immunoglobulins and does not bind immunoglobulins from bovine, goat or sheep; nor does it bind to lambda light chains.
- The total IgG content of serum is approximately 10-15 mg/ml and, therefore, 2 ml of settled resin has only enough binding capacity to purify antibody from ~1 ml of serum. The concentration of antibody in tissue culture supernatant varies considerably among hybridoma clones.
- The described purification procedure is for human IgG. For optimal recovery, use a sample size such that the expected Ig load on the column is less than 80% of the maximum binding capacity.
- Serum samples, ascites fluid, plasma or tissue culture supernatant may be used with this product.

---

## Gravity-flow Column Procedure for Antibody Purification

**Note:** The following protocol is for using a column packed with 2 ml of Pierce Protein L Agarose. When using columns containing other resin volumes, reagent amounts must be adjusted accordingly.

### A. Additional Materials Required

- Disposable column capable of containing at least 2 ml resin bed volume such as the Disposable Polypropylene Columns (Product No. 29922) or the Column Trial Pack (Product No. 29925), which contains two each of three column sizes.
- Binding Buffer: 0.1 M phosphate, 0.15 M sodium chloride; pH 7.2 (Product No. 28372)
- Elution Buffer: 0.1 M glycine, pH 2-3 or Pierce IgG Elution Buffer (Product No. 21004 and 21009)
- Neutralization Buffer: 1 ml of high-ionic strength alkaline buffer such as 1 M phosphate or 1 M Tris, (pH 7.5-9)
- Thermo Scientific Slide-A-Lyzer Dialysis Cassette or Zeba Desalt Spin Column (for example, Product No. 89894) for buffer exchange

### B. Immunoglobulin Purification Procedure

1. Equilibrate the Protein L agarose and reagents to room temperature before use.
2. Carefully pack the Protein L resin slurry (4 ml) into the column (see Additional Information section for the Tech Tip Protocol). Equilibrate the Protein L column by adding 10 ml of the Binding Buffer and allowing the solution to drain through the column.
3. Dilute sample at least 1:1 with Binding Buffer before applying onto the Protein L column to maintain optimal ionic strength and pH for binding.

**Note:** If plasma is used, the sample may appear hazy after adding the Binding Buffer caused by lipoprotein precipitation. For optimal Ig recoveries, centrifuge the diluted sample at  $10,000 \times g$  for 20 minutes and apply the supernatant to the equilibrated Protein L column.

4. Apply up to 4 ml of the diluted sample to the column and allow it to flow completely into the resin. The column will stop flowing automatically when the liquid level reaches the top disc. Larger volumes may be applied provided the total amount of Ig is less than 80% of column capacity.

**Note:** If the sample contains more Ig than can bind to the Protein L column (or is an antibody type that does not bind to Protein L), the flow-through will contain excess antibody. By saving the flow-through, non-bound antibody can be recovered and examined by antibody-specific assays.

5. Wash the Protein L column with 10-15 ml of Binding Buffer.

**Note:** If desired, verify that all non-bound proteins are removed from the column by collecting separate 2 ml fractions as the solution drains and measuring their absorbance at 280 nm. The last fractions should have absorbances similar to Binding Buffer alone.

6. Elute antibodies with 6-10 ml of Elution Buffer and collect 0.5-1 ml fractions. Immediately adjust eluted fractions to physiologic pH by adding 100  $\mu$ l of the Neutralization Buffer to 1 ml of eluate. Monitor the elution by measuring the absorbance at 280 nm or by protein assay such as BCA™ Protein Assay Kit (Product No. 23225).
7. Pool the eluted Ig fractions that contain the highest absorbance. The purified antibodies may be used directly for SDS-PAGE, or the buffer may be exchanged to a system compatible with the specific downstream application (see optional procedure that follows).
8. Regenerate column by washing with 12 ml of Elution Buffer.
9. For storage, wash column with 5 ml of water containing 0.02% sodium azide. When approximately 3 ml of solution remains, replace the bottom cap followed by the top cap on the column. Columns may be regenerated a minimum of 10 times without significant loss of binding capacity.

**Note:** If required, perform a buffer exchange using a Slide-A-Lyzer Dialysis Cassette or a Zeba Desalt Spin Column.

## Troubleshooting

Problem	Possible Cause	Solution
Flow of the column is exceedingly slow (i.e., < 0.5 ml/minute)	Outgassing of buffers or sample on the column, which results in blockage of resin pores with microscopic air bubbles	Degas buffers and remove air bubbles from column (see Additional Information section for suggested Tech Tip protocol)
No protein detected in any elution fractions	Sample devoid of antibody species or isotype that binds to Protein L	Use Protein A, Protein G or Protein A/G resin
Considerable antibody purified, but no specific antibody of interest detected	Antibody of interest is at low concentration	Use serum-free medium for cell supernatant samples
		Affinity purify the antibody using the specific antigen coupled to an affinity support (Product No. 44894)
Antibody of interest purified, but it is degraded (as determined by lack of function in downstream assay)	Antibody is sensitive to low-pH Elution Buffer	Try Gentle Ag/Ab Elution Buffer (see Related Products)
	Downstream application is sensitive to neutralized Elution Buffer	Desalt or dialyze eluted sample into suitable buffer

## Additional Information

Visit our website for additional information including the following items:

- Tech Tip Protocol: Packing Gel (resin) into Polypropylene Columns
- Tech Tip Protocol: Remove Air Bubbles from Columns
- Tech Tip Protocol: Degas Solutions for use in Affinity Columns
- Tech Tip: Binding Characteristics for Immunoglobulin Proteins and Proteins L, A, G, and A/G
- Tech Tip: Protein Stability and Storage

## Related Thermo Scientific Products

<b>66382</b>	<b>Slide-A-Lyzer Dialysis Cassette Kit</b> , 10 dialysis cassettes, each appropriate for 0.5-3.0 ml samples
<b>66528</b>	<b>Slide-A-Lyzer Concentrating Solution</b> , 200 ml
<b>21027</b>	<b>Gentle Ag/Ab Elution Buffer</b> , 500 ml
<b>44894</b>	<b>AminoLink® Plus Immobilization Kit</b>
<b>37501</b>	<b>Monoclonal Antibody Isotyping Kit I (HRP/ABTS)</b>
<b>20421</b>	<b>Pierce Protein A/G Agarose</b> , 3 ml settled resin with column sample kit
<b>20422</b>	<b>Pierce Protein A/G Agarose</b> , 15 ml settled resin with column sample kit
<b>20398</b>	<b>Pierce Protein G Agarose</b> , 2 ml settled resin
<b>20399</b>	<b>Pierce Protein G Agarose</b> , 10 ml settled resin
<b>20333</b>	<b>Pierce Protein A Agarose</b> , 5 ml settled resin
<b>20334</b>	<b>Pierce Protein A Agarose</b> , 25 ml settled resin

---

## General References

- Åkerström, B. and Björck, L. (1989). Protein L - an immunoglobulin light chain-binding bacterial protein. Characterization of binding and physicochemical properties. *J Biol Chem* **264**:19740-6.
- Beckingham, J.A., *et al.* (1999). Interactions between a single immunoglobulin-binding domain of protein L from *Peptostreptococcus magnus* and a human  $\kappa$  light chain. *Biochem J* **340**:193-9.
- Björck, L. (1988). Protein L - a novel bacterial cell wall protein with affinity for immunoglobulin light chains. *J Immunol* **140**:1194-7.
- Enokizono, J., *et al.* (1997). NMR Analysis of the interaction between protein L and Ig light chains. *J Mol Biol* **270**:8-13.
- Kastern, W., *et al.* (1992). Structure of *peptostreptococcal* protein L and identification of a repeated immunoglobulin light chain-binding domain. *J Biol Chem* **267**:12820-5.
- Nilson, B.H.K., *et al.* (1993). Purification of antibodies using protein L-binding framework structures in the light chain variable domain. *J Immunol Meth* **164**:33-40.

## Product References

- Kinoshita, Y., *et al.* (2001). Pescadillo, a novel cell cycle regulatory protein abnormally expressed in malignant cells. *J. Biol. Chem.* **276**(9):6656-65.
- Nozawa, K., *et al.* (2001). Preferential blockade of CD8+ T cell responses by administration of Anti-CD137 ligand monoclonal antibody results in differential effect on development of murine acute and chronic graft-versus-host diseases. *J. Immunol.* **167**:4981-6.

Slide-A-Lyzer<sup>®</sup> Dialysis Cassette Technology is protected by U.S. Patent # 5,503,741 and 7,056,440.

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

**No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).**

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current versions of product instructions are available at [www.thermo.com/pierce](http://www.thermo.com/pierce). For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2008 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.