

Immobilized Jacalin

3747 N. Meridian Road
P.O. Box 117
Rockford, IL 61105

20395

0735

Product Description

Number

20395

Description

Immobilized Jacalin, 5 ml gel

See product label for lot-specific binding capacity.

Loading: Approximately 4.5 mg Jacalin per ml of gel.

Store at 4°C. Do not freeze.

Introduction

Jacalin is an α -D-galactose binding lectin extracted from jack-fruit seeds (*Artocarpus integrifolia*).¹ This lectin is a glycoprotein of approximately 40,000 MW composed of four identical subunits. Jacalin immobilized on supports such as agarose has been useful for the purification of human serum or secretory IgA.²⁻¹² IgA can be separated from human IgG and IgM in human serum or colostrum.² IgD is reported to bind to jacalin.⁴ This support is also useful for removing contaminating IgA from IgG samples. Immobilized jacalin is preferable to protein A because protein A will bind some IgA along with IgG and is not effective for IgA removal.³

Materials

- A. Immobilized Jacalin, 5 ml
- B. Human Serum
- C. Binding Buffer: Phosphate buffered saline (PBS), 0.1 M sodium phosphate, 0.15 M NaCl, pH 7.4
- D. Elution Buffer: 0.1 M melibiose or 0.1 M α -D-galactose in PBS
- E. Column

Packing a Column

1. Allow the gel slurry and deionized water to come to room temperature.
2. Place the bottom cap on the end of a column.
3. Place the column in a test tube and add 2.75 ml of degassed water to a 2.0 ml column, 6.0 ml to the 5.0 ml column or 12.0 ml to the 10.0 ml column. Tap the end of the column sharply to dislodge any air bubbles. A 0.035% solution of detergent such as Surfact-Amps™ X-100 (Product No. 28314) can be substituted for degassed water to make it easier to remove entrapped air.
4. Float one of the porous Polyethylene discs on top of the liquid within the column.
5. Use the reverse end of a serum separator or pipette tip to push the disc to the bottom of the column. As the disc passes to the bottom, the air will be purged out through the disc.
6. Empty the column of water and add gel slurry to the column.
7. Allow 30 minutes for the gel to settle in the column.

Telephone: 800-8-PIERCE or 815-968-0747 • Fax: 815-968-7316 or 800-842-5007

Internet: <http://www.piercenet.com>

8. Place one of the discs on top of the liquid within the column and depress it to just above the settled gel level. Leave approximately 1-2 mm of space between the top of the gel bed and the bottom of the top disc.
9. Wash the inside top part of the column to clean out any gel that may have remained along the sides during packaging. The column is now ready for use. It may be stored in 0.02% sodium azide at 4°C.

Method

1. Equilibrate the column with at least 5 column volumes of binding buffer.
2. Mix human serum 1:1 with binding buffer, centrifuge if necessary and add to column.
3. Wash the column with 5 column volumes of binding buffer or until baseline is reached using absorbance at 280 nm to measure effluent.
4. Elute the bound IgA1 using 2 ml aliquots of elution buffer until baseline absorbance is reached. Collect the eluate from the aliquots in separate test tubes and monitor absorbance at 280 nm to determine baseline.
5. Buffer exchange the sample into binding buffer using a gel filtration/desalting column to remove the melibiose or galactose.
6. Regenerate the column by washing with 20 column volumes of PBS.

References

1. Kumar, G.S., Appukuttan, P.S. and Basu, D. (1982). α -D-Galactose-specific lectin from jack fruit (*Artocarpus integra*) seed. *J. Biosci.* **4**(3), 257-261.
2. Roque-Barreira, M.C. and Campos-Neto, A. (1985). Jacalin, an IgA-binding lectin. *J. Immunol.* **134**(30), 1740-1743.
3. Van Kamp, G.J. (1979). IgA contamination of IgG prepared on a protein A column. *J. Immunol. Meth.* **27**, 301-305.
4. Aucouturier, P., Mihaesco, E., Mihaesco, C. and Preud'Homme, J.-L. (1987). Characterization of jacalin, the human IgA and IgD binding lectin from jackfruit. *Mol. Immunol.* **24**(5), 503-511.
5. Biewenga, J., Steneker, I. and Hameleers, D.M.H. (1988). Effect of serum albumin on the recovery of human IgA1 from immobilized jacalin. *J. Immunol. Meth.* **115**, 199-207.
6. Bunn-Moreno, M.M. and Campos-Neto, A. (1981). Lectins(s) extracted from seeds of artocarpus integrifolia (Jackfruit): potent and selective stimulators(s) of distinct human T and B cell functions. *J. Immunol.* **127**(2), 427-429.
7. Gregory, R.L., Rundegren, J. and Arnold, R.R. (1987). Separation of human IgA1 and IgA2 using jacalin-agarose chromatography. *J. Immunol. Meth.* **99**, 101-106.
8. Jarvis, G.A. and Griffiss, J.M. (1989). Human IgA1 initiates complement-mediated killing of *Neisseria Meningitidis*. *J. Immunol.* **143**(5), 1703-1709.
9. MacDermott, R.P., Nash, G.S., Bertovich, M.J., Mohrman, R.F., Kodner, I.J., Delacroix, D.L. and Vaerman, J.-P. (1986). Altered patterns of secretion of monomeric IgA and IgA subclass 1 by intestinal mononuclear cells in inflammatory bowel disease. *Gastroent.* **91**(2), 379-385.
10. Saxon, A., Tsui, F. and Martinez-Maza, O. (1987). Jacalin, an IgA-binding lectin, inhibits differentiation of human B cells by both a direct effect and by activating T-suppressor cells. *Cell. Immunol.* **104**, 134-141.
11. Kondoh, H., Kobayashi, K., Hagiwara, K. and Kajii, T. (1986). Jacalin, a jackfruit lectin, precipitates IgA1 but not IgA2 subclass on gel diffusion reaction. *J. Immunol. Meth.* **88**, 171-173.
12. Wood, G.M., Trejdosiewicz, L.K. and Losowsky, M.S. (1987). ELISA for measurement of secretory IgA distinct from monomeric IgA. *J. Immunol. Meth.* **97**, 269-274.
13. Mestecky, J., Kulhavy R., Kraus, F.W. (1971). Method of Serum IgA Isolation. *J. Immunol.* **107** (2), 605-607.
14. Coppo, R., Arico, S., Piccoli, G., Basolo, S., Roccatello, D., Amore, A., Tabone, M., De LA Pierre, M., Sessa, A., Delacroix, D.L., Vaerman, J.P. (1985). Presence and Origin of IgA1 and IgA2-Containing Circulating Immune Complexes in chronic Alcoholic Liver diseases with and without glomerulonephritis. *Clin. Immunol. Immunopath.*, **35**, 1-8.
15. Hammarstrom, L., Holm, G., Palmblad, J., Persson, M.A.A., Smith, C.I.E. (1984). Lack of IgG in a Healthy Adult: A Rare Case of Dysgammaglobulinemia with Undetectable Serum IgG, IgA2, and IgE. *Clin. Immunol. Immunopath.*, **30**, 1-10.
16. Despont, J.P., Abel, C.A. (1974). Glycopeptides of heavy chains from Human IgA myeloma Proteins. *J. Immunol.*, **112** (5), 1623-1627.

Triton is a registered trademark of Rohm & Haas.

© Pierce Chemical Company, 11/1997. Printed in U.S.A.