

# Basics about CIM® technology and key applications

Aleš Štrancar

March, 2011



# Leaders in Monolith Chromatography

- BIA Separations was founded in September 1998 as a spin-off from BIA d.o.o. founded in 1989. Headquarters in Austria, R&D and Production in Slovenia.
- BIA Separations USA established in September 2007 - sales and tech support office.
- BIA Separations China established in January 2011 - sales and tech support office.
- Main focus: To develop and sell methacrylate monolithic columns & develop methods and processes for large biomolecules separation and purification.
- Pioneers and leaders in proprietary monolithic technology (CIM<sup>®</sup>). 4 USA patents granted including their foreign equivalents, more pending.



# Important Milestones

- 2002: **First Drug Master File (DMF) for CIM<sup>®</sup> DEAE supports.**
- 2002: Pass first FDA audit for one of the projects.
- 2004: First monolith used [for the industrial cGMP purification for plasmid DNA](#) at Boehringer Ingelheim provide [15-fold increase in productivity](#)
- 2006: **Drug Master File (DMF) for CIM<sup>®</sup> QA supports.**
- 2006: First cGMP production of a vaccine (influenza) using CIM<sup>®</sup>.
- 2008: Partnership with Agilent Technologies – develop and produce analytical monolithic columns
- 2009: Pass second FDA audit for one of the projects.
- 2010: **Drug Master File (DMF) for CIM<sup>®</sup> SO3 supports.**
- 2001 - 2010: Pass many audits by Novartis, Boehringer Ingelheim, Octapharma, ...



# Moving to the new facility in summer



# BIA Separations CIM® Monolithic Columns are Becoming Industry Standard for Production of Complex Biomolecules

- Drug Master Files (DMF) for CIM® DEAE, QA and SO3 columns in place, HIC in preparation.
- First drug purified using CIM monoliths passed CPIII trial (pDNA for gene therapy).
- More than 15 projects in CPI – CPIII trials (various Influenza, various Adenovirus, bacteriophages, various IgMs, Inter-alpha-inhibitors).
- More than 200 projects in pre-clinical trials (Influenza A and B virus (eggs, Vero and MDCK cells), Rabies virus, Rotavirus, AAV, various Adenovirus subtypes, Hepatitis A, Vaccinia, Mulv, MVM, Feline calicivirus, Japanese encephalitis, Crimean-Congo hemorrhagic fever, Hantaan virus, VLP (Hepatitis B, HPV, Influenza, Adenovirus), bacteriophages (Lambda, T4, VDX10, Pseudomonas phage), Tomato and Pepino Mosaic virus, pDNA, IgM, various proteins).



# Short Monolithic Columns Technology






# Many novel drug targets are large and complex, and in some cases used live

- These include different viral particles, pDNA, protein complex, IgM.
- "Whilst highly effective for the purification of proteins and smaller molecules, chromatographic techniques are not necessarily well suited to purification of these newer, larger targets." (N. Willoughby, J Chem Tech & Biotech, 84, 2008, 145).
- Why?



# Conventional Liquid Chromatography Media

**Packed bed of porous particles - two types of void volume exist!**



**Interparticle void volume**  
(preferential flow path)


**Intraparticle void volume**  
(contains majority of  
binding sites: > 90 %)





# Diffusion Limitations

Binding capacity at high flow rate:  $f_3 > f_2 > f_1 > f_0$




**Larger the molecule faster the capacity drop**



# Diffusion Limitations


Resolution @ linear gradient elution at high flow rate:  $f_2 > f_1$



**Larger the molecule wider the peak – lower the resolution**



# Another Challenge – the Size of the Molecule of Interest



**Pores too small for large solutes!**

(Binding mostly on outer surface)

results in


**Very low binding capacity for large solutes**

**(behave like nonporous particles)**



# Working with big molecules

- Molecule size: surface accessibility




Molecule	nm
Proteins	1-3
IgM	25
Plasmids	150-250
Rotavirus	130
Poxvirus	200 x 500
T4	220 x 85

Courtesy P. Gagnon [www.validated.com](http://www.validated.com)



# Yet Another – Shear Forces



Gray areas indicate particles.

The white area indicates the void space between particles.

Black arrowheads indicate primary flow.

Red arrowheads indicate countercurrent flow.

*The frictional differential between particle surfaces and the deep void space creates eddies — areas of persistent countercurrent flow.*



**Eddies create shear forces that damage labile biomolecules (similar to the effect when using ultracentrifugation).**

**Eddy-generated shear is proportional to flow rate.**


Courtesy of Pete Gagnon, Validated Biosystems, USA – details at [www.validated.com](http://www.validated.com)



# What are the Alternatives?



## Membranes – Stack of very thin Monoliths




Stacks of thin polymeric layers – supplied in single piece but in fact they are discontinuous unit.

### Problems with:

- resolution due to void volumes
- share forces due to eddies.

## Monoliths





Single piece continuous units with a homogeneous open pore structure in all 3 directions (flow through channels).





# Convective Transport: Consequences

- Flow independent properties











*Podgornik et al., Anal. Chem. 72 (2000) 5693*




# Monoliths

Monoliths are chromatography media that are cast as a single block and inserted into a chromatography housing. They are characterized by a highly inter-connected network of channels, sometimes compared to a sponge.




# Sum up: CIM<sup>®</sup> Monolithic Columns are purpose designed for the chromatography of big biomolecules



## Traditional approach - Porous particle:

1. Diffusive mass transport – slow process or **lower resolution**
2. Pores too small – very **low capacity**
3. Countercurrent flow - shear forces – **lower yields**



## Novel approach – CIM monoliths:



1. Convective mass transport – **flow independent resolution and capacity**, very fast process
2. Big channels – **high capacity**
3. Laminar flow - No shear forces – **better yields**



# CIM® Monoliths - the Only Material Engineered to Address the Needs of Large Molecule Separation


CIM® monolithic supports are highly porous rigid polymers with:

- High porosity (over 60 %)
- Flow-through channels (“pores”) having large diameter (1.5  $\mu\text{m}$ ), for Vaccinia special monolith (3-4  $\mu\text{m}$ )
- Biocompatible with uniform channel connectivity in 3D (homogeneous structure)
- Ligands (active groups) for **AEX, CEX, HIC, RPC, Affinity, Activated, Bioreactor.**




# *Comparison of support performance*

## ❖ For proteins & peptides




# ***Comparison of support performance***

❖ **For large proteins, DNA & viruses**






# Effect of the Molecule Size on Dynamic Binding Capacity



Yamamoto S. and Kita A., *Trans IChemE, Part C, Food and Bioproducts Processing*, 84 (2006) 72-77.



# Plasmid DNA Binding Capacity



**CIM DEAE binding capacity  
= ~10 mg/ml**

**Currently used for CP III trials**

## 15-fold increase in productivity

- High binding capacity at relevant flow rates
- High elution concentration - pDNA eluted in lower volume (important for SEC!)
- Fast process (no product loss due to oxidative degradation or enzymatic attack)

*Urthaler et al., J.Chrom. A, 1065 (2005), 93-106*



# Surface accessibility for CIM<sup>®</sup> Monoliths

High capacity for IgM, viruses and DNA

Molecule	Column	Capacity
IgM	CIM QA, SO3	25-50 mg/ml
Plasmid DNA	CIM DEAE	8 mg/ml
Genomic DNA	CIM DEAE	15 mg/ml
Endotoxins	CIM QA	>115 mg/ml
ToMV	CIM QA	2.0E+14 vp/ml
Influenza virus	CIM QA	2.0E+10 vp/ml
Adenovirus	CIM QA	3.0E+12 vp/ml
Ad3 VLPs	CIM QA	7.3E+16 VLP/mL





# CIM<sup>®</sup> Chemical Structure


Made of  
**highly cross-linked**  
**porous rigid monolithic**  
*poly(glycidyl methacrylate-  
co-ethyleneglycol dimethacrylate)*  
or  
*poly(styrene-divinylbenzene)*  
polymers

## Well proven and biocompatible:

- Toyopearl<sup>®</sup> from Tosoh
- Fractogel<sup>®</sup> from Merck / EMD



# Available Chemistries




\* on request



# CIM<sup>®</sup> Columns Design to Allow High Volumetric Flow Rates - High Productivity


Designed as Short Chromatographic Layers

Disk dimensions




$$V = 0.34 \text{ ml}$$

Disks used with housing




Typical flow rates





# What About the Separation on Short Chromatographic Layers?



Due to an almost **rectangular adsorption isotherm**, macromolecules remain adsorbed on the column almost **irreversibly**.


They are eluted by changing the mobile phase composition commonly applying linear or step **gradients**.

## “Theory of short chromatographic layers”:

Protein remains adsorbed at the top of the column until the eluting power of the mobile phase reaches the point at which a small change in the composition of the mobile phase causes the movement of the protein without any retention (Yamamoto, 1988). As a result, **even very short columns can provide very good separations**.




# Very Short CIM<sup>®</sup> Monolithic Columns Offer Outstanding Resolution



Anion Exchange Semi-Preparative Purification of a 16-mer Oligodeoxynucleotide on a 0.34 ml CIM<sup>®</sup> DEAE **Disk** Monolithic Column (**3 mm long** x 12 mm ID column)



# CIM<sup>®</sup> CLC - Multidimensional chromatography



Separation of IgG (8 mg), Transferrin (within range of 1mg/mL) and Albumin (within range of 1mg/mL) on monolithic column consisting of **two Protein G** and **one CIM<sup>®</sup> QA** Disks.

**Buffer A:** 20 mM Tris-HCl, pH 7.4

**Buffer B:** 20 mM Tris-HCl, 1 M NaCl, pH 7.4

**Buffer C:** 0.1 M Gly-HCl, pH 2.6

**Flow:** 4 ml/min

**Injection volume:** 250  $\mu$ l

Order of elution: 1- Transferrin  
2- Albumin  
3- IgG



# Fast Method Development

	Short Layer Monolith (3x12 mm ID)	Porous Particles
Column volume	0.34 ml	1 ml
Flow rate applied	4 ml/min	1 ml/min
Flow rate applied	12 CV/min	1 CV/min
Time – loading (5 CV)	0.4 min	5 min
Time – elution (10 CV)	0.9 min	10 min
Time – equilibration (5 CV)	0.4 min	5 min
Time – total per run	1.7 min	20 min
Time for 20 runs	0.6 h	6.7 h
Time for 100 runs	2.8 h	33.3 h



# Lab Scale Columns – CIM® Disks

- Smaller units in disk format - intended for media screening, method development, laboratory purification

## CIM Disk Monoliths with housing



## CIM Disk housing parts



**Color of ring denotes the disk chemistry**

**Blue** Strong AEx (QA)  
**Green** Weak AEx (DEAE)


**Red** Strong CEx (SO3)  
**Black** Weak CEx (CM)



# Industrial Scale Units – CIM<sup>®</sup> Tubes



80, 800, and 8.000 ml CIM Monoliths





# CIM<sup>®</sup> Tube Column Structure



1. Piston – Collector  
with flow-out

2. Seal

3. Frit

4. Monolith

5. Housing - Distributor

6. Upper plate  
with flow-in



# Currently Available Sizes

**Performance and process time maintained!**




**0.34 ml disk**

**3-8 ml/min**



**8 ml column**

**10-40 ml/min**



**80 ml column**

**40-250 ml/min**






**800 ml column**

**400-2000 ml/min**



**8000 ml column**

**2000-10000 ml/min**



# Applications



# Clotting Factor VIII/von Willebrand Factor Complex




# Separation of Clotting Factor VIII on a CIM<sup>®</sup> QA Disk Monolithic Column

Separation of clotting factor VIII/von Willebrand factor complex from contaminating (model) proteins

Peaks 1 & 2 - Transferrin and IgG

Peak 3 – Human Serum Albumin

Peak 4 – FVIII/vWF




(Štrancar A. et al. J. Chromatogr. A 1997; 760: 117-123)



# Semi-Preparative Isolation of Clotting Factor VIII/von Willebrand Factor Complex

Separation of sample No. 3 from FVIII production resolved cryoprecipitate after  $\text{Al}(\text{OH})_3$  precipitation and S/D virus inactivation, obtained on a QA anion-exchange compact porous tube (53 mm long, 23 mm diameter and with a 1-mm inner hole). Conditions: Buffer A, 10 mM sodium citrate, 120 mM glycine, 1 mM  $\text{CaCl}_2$ , pH 7.0; buffer B, buffer A containing 1 M NaCl; injection volume, 20 ml.



(Štrancar A. et al. J. Chromatogr. A 1997; 760: 117-123)






# Clotting Factor IX



# FIX – Production Process




# Anion exchange chromatography: Purification of FIX out of the Sephadex eluate

Column	Buffer	Step gradient	Dynamic capacity	Recovery (%)	Specific activity
<b>CIM® DEAE tube</b>	<b>Citrate</b>	<b>10-50-100</b>	<b>109.4</b>	<b>70.1</b>	<b>12.5</b>
	<b>Phosphate</b>	<b>10-50-100</b>	<b>187.5</b>	<b>76.2</b>	<b>10.6</b>
<b>DEAE conventional column</b>	<b>Citrate</b>	<b>10-50-100</b>	<b>93.4</b>	<b>100</b>	<b>3.4</b>
<b>CIM® QA tube</b>	<b>Citrate</b>	<b>5-30-100</b>	<b>156.1</b>	<b>86.3</b>	<b>14.5</b>
	<b>Phosphate</b>	<b>5-30-100</b>	<b>177.1</b>	<b>97.5</b>	<b>12.8</b>

*Courtesy of Prof. Dj. Josic, Octapharma, Vienna, Austria*



# Process Design: Use of Step Gradient Disk-Shaped Monolithic Column (340 $\mu\text{l}$ )



# Process Design

## Disk-Shaped Monolithic Column (340 $\mu$ l)


<b>FIX-yield</b>	<b>–</b>	<b>77.1 %</b>
<b>FIX-specific activity Protein</b>	<b>–</b>	<b>44.4 IU/mg</b>

*Courtesy of Prof. Dj. Josic, Octapharma, Vienna, Austria*



# Scaling-up

## 8 ml Tube Monolithic Column






# Optimization

## Scaling-up, 8 ml tube

	FVII (IU)	FIX (IU)	FX (IU)	PA x10 <sup>3</sup> (U)	VN (mg)	Protein (g)	Vol. (ml)	Spec. act. (IU FIX/mg)
Load	240	500	485	19.85	12.50	0.423	10	1.02
F/T	4.32	0.24	12.7	0.47	0.48	0.157	12	-
W	27.50	0.9	25.2	0.12	2.30	0.167	45	-
E <sub>1</sub>	292	80	372	0.08	1.98	0.124	40	-
E <sub>2</sub>	1.5	487.5	15	0.13	1.28	0.018	25	27.05
E <sub>3</sub>	0.25	15	1	25.2	9.0	0.165	25	-
Σ	560.75	682.4	893	25.25	26.06	0.428	140	-
%	1.25	1.18	27.6	1.21	1.26	1.28	-	-



# Optimization Separation with 800 ml CIM<sup>®</sup> DEAE Tube Monolithic Column



# Optimization Separation with 800 ml CIM<sup>®</sup> DEAE Tube Monolithic Column

## Scaling-up, 800 ml tube sephadex eluate

	FVII (IU)	FIX (IU)	FX (IU)	PA x10 <sup>3</sup> (U)	VN (mg)	Protein (g)	Vol. (ml)	Spec. act. (IU FIX/mg)
Load	28800	40680	58200	2382	1500	39	1200	1.02
F/T	853	1114	1651	132	152.7	8.05	1376	-
W	6173	-	4079	21.2	747.8	22.5	4249	-
E <sub>1</sub>	23233	3319	33927	17.5	545.8	8.85	1843.9	-
E <sub>2</sub>	40	40172	331	9.33	86.4	1.2	1004.3	33.3
E <sub>3</sub>	30	1719	69	1806	288	2.0	955.9	-
Σ	30329	46324	40057	1986	1819	42.6	-	-
%	105	113	68.8	83.3	121	109	-	-



# Major advantages using CIM<sup>®</sup> columns for FLX purification

- Very fast method development.
- Very fast scale-up to industrial scale.
- Highly reproducible results regarding speed, specific activity and yield from the disk to industrial columns.
- About 20 times faster purification than using column packed with bulk support.
- Specific activity in the eluate 5 times higher than in corresponding purification step using column packed with bulk support.


*Josic, PPB 03, Curacao*



# Pegilated proteins



# Faster Separation of PEGylated Proteins



Q-Sepharose HiTrap  
25mm x 7mm (1ml)  
Amersham/GE




QA-CIM®  
3 disks (1ml)  
BIA Separations

**- 5-fold faster run!**  
**- Sharper peaks!**





# Separation of PEGylated Proteins on cation exchange CIMac™ columns



CIMac™ SO3 and CM columns, 15 mm length; Gradient: 0- 0,5 min 100 % Buffer A, than 0-100 % Buffer B in 52 Column Volumes; Buffer A: 20 mM phosphate, pH 6,5, Buffer B: Buffer A + 0,3 M NaCl




# Platform IgM purification processes




# Immunoglobulin - IgM

Pentameric, 0.96 Md



Hexameric, 1.15 Md



# Are IgMs really difficult to purify?

IgMs have some characteristics that can limit the application of standard purification tools:

- They tend to be less soluble than IgGs and more susceptible to denaturation at extremes of pH. This can limit application of affinity chromatography.
- Low solubility is compounded by low conductivity. This can limit ion exchange chromatography.
- They are generally tolerant of high salt concentrations, but susceptible to denaturation from exposure to strongly hydrophobic surfaces. This can limit hydrophobic interaction chromatography.
- Large size corresponds with slow diffusion constants.
- Porous particle based chromatography media depend on diffusion for mass transport.
- Slow diffusion constants translate into lower capacity and lower resolution, and/or lower flow rates.
- This is a particular limitation for size exclusion chromatography because it already suffers from low capacity and low flow rate.



# Are IgMs really difficult to purify?

On the other hand:

- IgMs are typically more charged than IgGs. They bind more strongly than IgG to anion exchangers or cation exchangers.
- They also bind more strongly than IgG to hydroxyapatite, and much more strongly than most contaminants.
- HIC on moderately hydrophobic supports usually elutes IgM in a well defined peak at reasonably low salt concentration.
- A new generation of industrial ion exchangers is available that does not rely on diffusion.
- Convection is independent of size and flow rate, so capacity and resolution are not affected by the large size of IgM, nor does flow rate need to be reduced.

For details visit [validated.com](http://validated.com)



# Purification of Clinical Grade Human IgM from Cell Culture Supernatant

STEP 1

Cation exchange chromatography  
on a **CIM® SO3 Tube Monolithic Column**

STEP 2

Anion Exchange Chromatography  
on a **CIM® QA Tube Monolithic Column**

STEP 3


Polishing on a hydroxyapatite column

For details visit [validated.com](http://validated.com)



# STEP 1:

## Cation Exchange Chromatography



**Column:** CIM<sup>®</sup> SO<sub>3</sub> Tube  
Monolithic Column (V=8 mL)  
**Flow rate:** 20 mL/min


Equilibrate: 50 mM MES, pH 6.0  
Titrate dilute load: 1500 mL  
Wash: 50 mM MES, pH 6.0  
Wash: 25 mM NaPO<sub>4</sub> pH 7.0  
Elute: LG to 225 mM NaPO<sub>4</sub>  
Clean: 500 mM NaPO<sub>4</sub> pH 7.0  
Sanitize/store: 1.0/0.01 M NaOH





# STEP 2:

## Anion Exchange Chromatography



Column: CIM® QA Tube Monolithic  
Column (V=8 mL)  
Flow rate: 20 mL/min


Equilibrate: 10 mM NaPO<sub>4</sub> pH 7.0  
Load 1: pH 7.0 **Load 2: pH 8.0**  
Wash: 75 mM NaPO<sub>4</sub> pH 7.0  
Elute: LG to 225 mM NaPO<sub>4</sub>  
Clean: 500 mM NaPO<sub>4</sub> pH 7.0  
Sanitize/store: 1.0/0.01 M NaOH

For details visit [validated.com](http://validated.com)



# STEP 3:

## Hydrophobic Interaction using Hydroxyapatite



Media: 10 mL CHT type II 40 $\mu$ m  
Column: Hydroxyapatite column  
Flow rate: 3.4 mL/min

Equilibrate: 10 mM NaPO<sub>4</sub> pH 7  
Load: Adjusted eluate  
from STEP 2

Wash: 10 mM NaPO<sub>4</sub> pH 7.0  
Wash: 75 mM NaPO<sub>4</sub> pH 7.0  
Elute: LG to 225 mM NaPO<sub>4</sub> pH 7  
Clean: 500 mM NaPO<sub>4</sub> pH 7.0  
Sanitize/store: 1.0/0.1 M NaOH

For details visit [validated.com](http://validated.com)



# IgM Purification Process Summary

## Process summary

	<i>Cation exchange 8 mL monolith</i>	<i>Anion exchange 8 mL monolith</i>	<i>Hydroxyapatite 10 mL column</i>
<i>Sample volume, mL</i>	250	25	25
<i>Diluted Sample, mL</i>	1250	250	250
<i>Diluted Sample, CV</i>	156	31	25
<i>Flow rate, mL/min</i>	20	20	3.34
<i>Flow rate, CV/min</i>	2.5	2.5	0.67
<i>Application time</i>	62.5	12.5	75.0
<i>Total volume<sup>1</sup></i>	1950	950	650
<i>Total time, min</i>	98	48	195
<i>Recovery %</i>	78(86) <sup>2</sup>	84	88
<i>Purity %</i>	~90	~95	~99


<sup>1</sup>Includes equilibration, sample application, wash, elution, cleaning.

<sup>2</sup>Includes the IgM that eluted prematurely in the wash.

Courtesy of Pete Gagnon, Validated Biosystems, USA – details at [www.validated.com](http://www.validated.com)



# Separation of IgG, IgA and IgM using CIM® QA



Column: **CIM® QA (strong anion exchanger)** column (5,2 mm I.D. x 4,95 mm L; V = 100 µl)

Sample: A mixture of human IgG (Octapharma, Md = 150 kDa), IgA (Sigma, xxxx, Md = 160 kDa) and IgM (Sigma, yyy, Md = 950 kDa) dissolved in 20 mM Tris-HCl buffer, pH 7,4

Mobile phase A: 20 mM Tris-HCl buffer, pH 7,4      Mobile phase B: 20 mM Tris-HCl buffer + 1.0 M NaCl, pH 7,4

Flow rate: 1,0 ml/min      Gradient: A linear gradient from 0 % buffer B to 35 % buffer B in 4 min (40 column volumes).

Detection: UV at 280 nm      Column pressure: 15 bar (1,5 MPa)



# IgG impurity removal processes




# IgG Impurity Removal

- Sample:
  - Chimera
    - 1 mg/mL purified monoclonal IgG + 350 ml of 0.1 mg/mL DNA
- Comparison of polishing applications of 3 different resins:
  - Microparticulate
    - Q Sepharose™ Fast Flow (1mL HiTrap™ column)
  - Membranes
    - Sartobind™ Q nano (1mL)
  - Monoliths
    - CIM® QA Disk Monolithic Columns (1 mL - 3 disks)



# Dynamic Binding Capacity for DNA




Note 50 times higher dynamic binding capacity than particle based resin while operating at 4-fold higher flow rate!

Courtesy of Pete Gagnon, for details visit [validated.com](http://validated.com)



# DNA Binding Efficiency



**Membrane:** earlier breakthrough, shallower slope

**CIM<sup>®</sup> QA:** later breakthrough, steeper slope



Lower Binding Capacity

4.8 mg/mL no-bt capacity = Membrane

Higher Binding Capacity

14.3 mg/mL no-bt capacity = CIM<sup>®</sup> QA

**Important implications for manufacturing of therapeutic antibodies**

Courtesy of Pete Gagnon, for details visit [validated.com](http://validated.com)






# Platform plasmid DNA purification proces



# Plasmid DNA Process

Speed + Capacity = Productivity = < Costs



**CIM DEAE binding capacity  
~10 mg pDNA/ml**

**Passed CP III trials**


15-fold increase in productivity reported by Boehringer Ingelheim. Used for gene therapy, DNA vaccines.

- High binding capacity at relevant flow rates
- High elution concentration - pDNA eluted in lower volume (important for SEC!)
- Fast process (no product loss due to oxidative degradation or enzymatic attack)

*Urthaler et al., Boehringer Ingelheim, J.Chrom. A, 1065 (2005), 93-106*



# RNA-se Free pDNA Purification Process



*Urthaler et al., Chem.Eng.Technol., 28 (2005), 1408-1420*

HIC is a bottleneck due to low capacity and slow process – larger columns are needed, more buffers are consumed.



# Novel plasmid DNA purification process

The goal:


to introduce a resin with higher capacity  
and better flow properties as a capturing

step –



smaller columns and less buffer  
consumption



# Novel pDNA purification process design



# RNA Removal - Precipitation with $\text{CaCl}_2$




# Selective Precipitation with $\text{CaCl}_2$

- Alkaline lysis
  - 50 mM Tris pH 8.0/ 10 mM EDTA
  - 0.2 M NaOH/ 1% SDS
  - 3 M potassium acetate, pH 5.0
- Adjustment to proper  $\text{CaCl}_2$  concentration (0,3-0,8 M)
- Incubation for 15 minutes at 4°C
- Centrifugation/Filtration



# Anion Exchange Step



Equilibration buffer: 50 mM Tris, 10 mM EDTA, pH 7.2

Washing buffer: 50 mM Tris, 10 mM EDTA, 0.6 M NaCl, pH 7.2

Elution buffer: 50 mM Tris, 10 mM EDTA, 1 M NaCl, pH 7.2


Regeneration buffer: 50 mM Tris, 10 mM EDTA, 2 M NaCl, pH 7.2

**Working capacity: 6 mg/ml DEAE monolith**





# Hydrophobic Interaction Step



Equilibration and washing buffer: 50 mM Tris, 10 mM EDTA, 1.7 M  $(\text{NH}_4)_2\text{SO}_4$ , pH 7.2

Elution buffer: 50 mM Tris, 10 mM EDTA, 0.4 M  $(\text{NH}_4)_2\text{SO}_4$ , pH 7.2

Regeneration buffer: 50 mM Tris, 10 mM EDTA, pH 7.2

**Working capacity: 3 mg/ml C4 HLD monolith**




# Yield and Purity

	Alkaline lysate	CIM® DEAE-8	CIM® C4-8
pDNA (µg/ml)	28	630	300
pDNA (mg)	40	38	34
Homogeneity (%SC)	94	95	98
Endotoxins (EU/mg pDNA)	80000	12,4	1,1
Host cell proteins (µg/ml)	5000	20	1,1
gDNA (µg/mg pDNA)	30	74,3	3,4
RNA (µg/ml)	N.D.	0	0
Step yield (%)	100%	95%	90%



# Novel Plasmid DNA Purification Process



# Consistent Scale-up

	Monolith Bed Volume (AEX and HIC)	pDNA produced per batch
<b>CIM-1</b>	1 ml	6 mg
<b>CIM-8</b>	8 ml	48 mg
<b>CIM-80</b>	80 ml	480 mg
<b>CIM-800</b>	800 ml	4,8 g
<b>CIM-8000</b>	8.000 ml	48 g



# Costs Comparison of Monolith and Particle Based Plasmid DNA Purification Processes

## 1 ml CIM® monolith

Calculations	
Buffer	76,3 ml buffer/mg pDNA
Time	23,6 min/mg pDNA
Recovery	85%
Purity	cGMP grade

### Costs using columns for 1 Run

Quantity of purified pDNA	5,10 mg pDNA
€ (Column costs)	114 €/mg pDNA
€ (column+buffer)	114 €/mg pDNA
€(column+buffer+work)	123 €/mg pDNA

### Costs using columns for 10 Runs

Quantity of purified pDNA	51 mg pDNA
€ (Column costs)	11,4 €/mg pDNA
€ (column+buffer)	11,8 €/mg pDNA
€(column+buffer+work)	21,1 €/mg pDNA

### Costs using columns for 20 Runs\*

Quantity of purified pDNA	102 mg pDNA
€ (Column costs)	5,7 €/mg pDNA
€ (column+buffer)	6,1 €/mg pDNA
€(column+buffer+work)	15,4 €/mg pDNA

## Particle based

Calculations	
Buffer	108,0 ml buffer/mg pDNA
Time	70,0 min/mg pDNA
Recovery	79%
Purity	cGMP grade

### Costs using columns for 1 Run

Quantity of purified pDNA	4 mg pDNA
€ (Column costs)	227 €/mg pDNA
€ (column+buffer)	228 €/mg pDNA
€(column+buffer+work)	257 €/mg pDNA

### Costs using columns for 10 Runs

Quantity of purified pDNA	40 mg pDNA
€ (Column costs)	23 €/mg pDNA
€ (column+buffer)	24 €/mg pDNA
€(column+buffer+work)	53 €/mg pDNA

### Costs using columns for 20 Runs

Quantity of purified pDNA	79 mg pDNA
€ (Column costs)	11 €/mg pDNA
€ (column+buffer)	12 €/mg pDNA
€(column+buffer+work)	42 €/mg pDNA



# Platform Flu Vaccine Purification Proces




# Uhhhhh, how do we Scale-up our vaccine production process?



# Influenza Viruses

- *Orthomyxoviridae*: enveloped, negative strand RNA virus
- Segmented genome: eight segments encode 11 viral proteins
- Host range
  - Influenza A: humans, swine, horses, domestic and wild avian species
  - Influenza B: humans
- Size: 80 – 120 nm





# Conventional purification based on Continuous Flow Ultracentrifugation

- 40 years of reliable operation carrying out runs every day for much of the year
- The units of choice for 85% of current Influenza vaccine market
- Already used with cell culture-based virus production
- Combines Concentration and Purification replacing multiple steps by alternative technologies
- Nevertheless:
  - Long processing time (10-18hours per run)
  - Expensive equipment / potential safety hazards
  - Multiple runs may be needed for impurity removal

*Dr. Charles Lutsch, Process Development, Sanofi Pasteur Inc., USA, MSS2008, Portorož*



# Conventional Vaccine Purification is Based on Continuous Flow Ultracentrifugation

- In addition:
  - Expensive service contracts and maintenance of centrifugation equipment.
  - Footprint of the equipment (size of the facility) might be an issue, as well the completely closed processing loop.
  - How to address in-process control (PAT)?
  - Can centrifugation address requests for “Disposable factory”?
  - What about the yield when life virus particle is in question (Conventional methods of virus purification using ultracentrifugation frequently result in distorted particles with low levels of biological activity)?
  - Can centrifugation address purification of non-viral particle candidates (proteins, pDNA)?



# Evaluation of Different Supports for Purification of Influenza A

Average values	CIM QA	Mustang <sup>®</sup> Coin Q	Q Sepharose <sup>™</sup> XL	Celufine Sulfate
Virus Recovery	54%	35%	35%	27%
DNA Depletion	96%	95%	95%	91%
Protein Depletion	95%	94%	98%	99%
Dynamic Binding Capacity	<b>10.3</b> log <sub>10</sub> TCID50/mL Support	<b>10.3</b> log <sub>10</sub> TCID50/mL Support	<b>9.0</b> log <sub>10</sub> TCID50/mL Support	<b>8.4</b> log <sub>10</sub> TCID50/mL Support



# To Sum Up: Novel Chromatography Supports Shall be Used for Vaccine Purification Processes

- Monolithic and membrane adsorbers offer an order or two higher capacity to particle based supports (1 L can replace 10 to 100 L column).
- Much smaller footprint of the facility.
- Much lower buffer consumption when using monoliths or membranes.
- Monolithic and membrane adsorbers offer much shorter process time.
- As a result; lower production costs, lower COGs of the vaccine.
- Ease of use (easier to handle smaller units).
- No column packing needed.
- Purity of the product is the same or better.
- Membranes (Stack of very thin Monoliths) may offer lower yield and/or lower purity than single piece Monoliths.




# IEX Monoliths Screening

Virus	IEX	HA Yield (%)	
		Flowthrough	Elution
<b>H1N1</b>	CIM QA	0	43.7
	CIM DEAE	0	47.3
	CIM SO3	0	49.6
<b>H3N2</b>	CIM QA	0	50.0
	CIM DEAE	0	37.1
	CIM SO3	0	62.1
<b>H5N1</b>	CIM QA	0	84.0
	CIM DEAE	-	-
	CIM SO3	<b>3.6</b>	114.6
<b>FLUB</b>	CIM QA	0	51.2
	CIM DEAE	0	35.4
	CIM SO3	<b>37.1</b>	30.2



# Mobile Phase Screening using CIM QA Columns




Zwitterionic buffers stand out for two important reasons:

- first that they have no inherent conductivity,
- are immune from binding to charged groups which results in more robust process.



# Purification of Clinical Grade Flu Vaccine Using CIM QA Columns




CIM QA (strong anion exchange column), 8 mL tube

Flow rate: 45 ml/min (150 cm/h)

Load: 140 ml of H1N1 TUF concentrate



# Final Polishing and Buffer Exchange using SEC





# Mobile Phase Optimisation for SEC Step

	SPG	SPG <sub>modified</sub>
Virus recovery*	59%	100%
DNA depletion	54%	48%
Protein depletion	79%	81%

\*TCID50

SPG


Sucrose-phosphate-glutamate buffer

SPG<sub>modified</sub>

Sucrose-phosphate-glutamate buffer modified



# Chart of the Clinical Grade Flu Vaccine Manufacturing Process



**Currently used for CP I and CP II trials**



# Process Yields and Purity Obtained

Step	Aim	Virus yield*
Clarification	DNA depletion, cell debris clearance	50-100%
TUF	Protein and DNA depletion, virus concentration	80-100%
CIM QA	Protein and DNA depletion, virus concentration	50-100%
SEC	Protein and DNA depletion, buffer exchange	100%
Overall	From harvest to purified vaccine bulk	≥ 25%


\*TCID50

	H1N1	H5N1
DNA depletion	≥ 99.9%	≥ 99.9%
Protein depletion	≥ 99.0%	≥ 99.5%




# Centrifugation versus Chromatography Based Flu Purification Process

## CENTRIFUGATION BASED PURIFICATION PLATFORM



Infectious virus yield	11.4 %
DNA removal	99.50 %
Protein removal	97.4 %

## MONOLITH BASED PURIFICATION PLATFORM



Infectious virus yield	47.3 %
DNA removal	99.96 %
Protein removal	97.8 %



# Centrifugation versus Chromatography Based Phage Purification Process

## CsCl-ULTRACENTRIFUGATION

- None
- 5 days
- costs & time
- $10^{10}$ - $10^{11}$  phages
- not for CsCl-sensitive phages
- 0,01 – 0,1 %



OPTIMIZATION

TIME

UP-SCALING

CAPACITY

APPLICABILITY

RECOVERY

(P. Putida phage  $\phi$ 15)

## CIM MONOLITHIC COLUMNS

- 10 days (pH-stability & cond.)
- 3 days
- Multiple runs on the same column
- $10^{11}$ - $10^{14}$  phages
- all phages
- about 71%




# To Sum Up: Chromatography Using Novel Supports can Better Address the Needs of Novel Vaccine Purification Platforms than Ultracentrifugation

- Yield of live virus particle is (much) higher (lower product degradation).
- Much smaller footprint of the facility and more flexible equipment.
- Lower production and equipment maintenance costs.
- Purity of the product is the same or better.
- Faster and easier scale-up, scale-down.



# Multiuse Disposable Units - “Plug and Play”



CIMmultus™ from BIA Separations (1 mL – 8 L)

Carbon fibre reinforcement embedded into epoxy thermoset resin (carbon fibre composite); tough, light material; 5-times lower density than stainless-steel; operate at 20 bar (291 psi).

**Replace the SS housings.**

**More inert than polypropylene, less leachables, less unspecific binding.**





# Introduction of Continuous LC (BioSMB™)



**Disposable monolithic or membrane columns**





# Disposable and Continuous Liquid Chromatography Systems Fit to “Single use” Vaccine Production Facility




*Courtesy of Xcellerex*



**Do we know enough about  
our complex molecule  
samples to develop safe  
product and/or run the  
production process?!**





# In-process control (PAT) using CIMac™ monolithic HPLC columns



# Monolithic Analytical Columns for In-process Control (PAT)



## CIMac™ HPLC Columns



10 ml/min = 4500 cm/h = 360 CV/min (res. time: 0,1 s) = faster than biosensor



# CIM<sup>®</sup> ImmunoDisk allows rapid quantification of biomarkers



**MAb 69.31 ImmunoDisk standard curve**  
**Fresh Frozen Human Plasma - 250 mg/L of INTER-ALPHA INHIBITOR**


*Courtesy of Prof. Yow-Pin Lim, ProThera Biologics, Providence, USA*



# Results on progres of Sepsis within few minutes CIM<sup>®</sup> ImmunoDisk


## Inter-alpha Inhibitors

Normal range 600-800 mg/mL



## Granzyme K

Normal range 50-100 aU/mL



Septic Patient (#242)

- @ Admission, 0 hr
- 12 hr later
- 30 hr later @ ICU

## SERIAL PLASMA STUDY OF SEVERE SEPTIC PATIENT


Courtesy of Prof. Yow-Pin Lim, ProThera Biologics, Providence, USA





# Rapid IgM Content Analysis

*Rapid analysis, IgM content of cell culture supernatant*





Reproduced from A high speed monolithic assay for IgM quantitation in cell culture production and purification process monitoring: Pete Gagnon, Richard Richieri, Simin Zaidi, Roy Sevilla, Alexander Brinkman, 3<sup>rd</sup> Wilbio Conference on Purification of Biological Products, September 24-26, 2007, Waltham, MA USA





# Separation of Plasmid DNA Isoforms Using Monolith Columns – PAT of pDNA Production



**AGE of different isoforms**

- 1 Initial pDNA sample
- 2 Linear form
- 3 OC form


Column: CIMas DEAE, 5.2 mm x 4.95 mm, max. pressure: 150 bars

Buffer A: 20 mM Tris-HCl; pH 8.5, Buffer B: Buffer A + 1 M NaCl, Flow rate: 1 ml/min, Gradient: 60 to 75 % buffer B within 100 CV, Injection volume: 2  $\mu$ l, Detection: UV at 260 nm.



# PAT - Phage and Impurity Growth Monitoring in Pilot Scale Bioreactor

CIMac QA disk  
Take 50  $\mu$ l from bioreactor  
Filter using 1  $\mu$ m filter  
Inject in the HPLC




Growth monitoring

Monitoring impurity profile



# PAT of Ad5 Production Using Monolithic HPLC




## Eden Biodesign Platform Ad5 Process



C. Sims et al., Eden Biodesign, MSS 2010, Portoroz, Slovenia, 2010



# Monolithic HPLC used for Ad5 Production Process Development – Basis for the PAT




*P. Ball et al., Eden Biodesign, MSS 2008, Portoroz, Slovenia, 2008*




# Monolithic HPLC used for Ad5 Production Process Development – Basis for the PAT


DNase Pre-treatment



DNase Post-treatment







- Species eluting at 2 minutes are proteinaceous (280 nm > 260 nm). Likely to be HCP or Ad5 proteins.
- Species eluting at 4 minutes a mixture of Ad5 particles, possibly free DNA.



*P. Ball et al., Eden Biodesign, MSS 2008, Portoroz, Slovenia, 2008*




# Monolithic HPLC used for the Ad5 Production Process PAT



*C. Sims et al., Eden Biodesign, MSS 2010, Portoroz, Slovenia, 2010*




# Monolithic HPLC Used for the Ad5 QC/QA of the Final Product



*P. Ball et al., Eden Biodesign, MSS 2008, Portoroz, Slovenia, 2008*



# Monolithic RP HPLC used for the Recognition of Different Flu Strains



**Virus disruption with detergent (Zwittergent) followed by Trypsin treatment and Injection to CIMac SDVB.** Conditions: mobile phase A = 5% AcCN, 0.1 % TFA, mobile phase B = 90% AcCN, 0.1% TFA, Injection volume = 30  $\mu$ L,  $\lambda$  = 215 nm; flow rate = 0.8 mL/min.

*L. Urbas et al., MSS 2010, Portoroz, Slovenia, 2010*





# Take Home Message

Vaccine development and manufacturing can be supported by many new innovative technologies and materials but chemistry and biochemistry are staying the same.

Without understanding the basic of the technology and investment in R&D one cannot expect successful adoption. Miracles rarely happen  
- Good analytical methods are the key.

Unspecific binding, column clogging, product degradation and/or inactivation are NOT caused by innovation or new technology but materials used.

Even perfect DSP cannot solve problems of badly designed, not reproducible USP. PAT is a key to understand and design robust USP.



# BIA Separations CIM® Monolithic Columns are Becoming Industry Standard for Production of Complex Biomolecules

- Drug Master Files (DMF) for CIM® DEAE, QA and SO3 columns in place, HIC in preparation.
- First drug purified using CIM monoliths passed CPIII trial (pDNA for gene therapy).
- More than 15 projects in CPI – CPIII trials (various Influenza, various Adenovirus, bacteriophages, various IgMs, Inter-alpha-inhibitors).
- More than 200 projects in pre-clinical trials (Influenza A and B virus (eggs, Vero and MDCK cells), Rabies virus, Rotavirus, AAV, various Adenovirus subtypes, Hepatitis A, Vaccinia, Mulv, MVM, Feline calicivirus, Japanese encephalitis, Crimean-Congo hemorrhagic fever, Hantaan virus, VLP (Hepatitis B, HPV, Influenza, Adenovirus), bacteriophages (Lambda, T4, VDX10, Pseudomonas phage), Tomato and Pepino Mosaic virus, pDNA, IgM, various proteins).



