

# Affinity Chromatography Media

## Matrex<sup>®</sup> Cellufine<sup>™</sup> Sulfate

### ■ For Concentration, Purification and Depyrogenation of Virus, Viral/Microbial Antigens, Heparin Binding Proteins

Advances in vaccines and clinical diagnostics have created an increasing demand for large volumes of highly purified and concentrated virus and viral or microbial antigens. Matrex Cellufine Sulfate affinity media is a simple, rapid and effective means for concentration, purification and depyrogenation of these important products.

Matrex Cellufine Sulfate eliminates cumbersome, time-consuming and potentially unsafe classical ultracentrifugation and density gradient methods. It can also provide a significant improvement in concentration and purity. Matrex Cellufine Sulfate offers flow properties superior to hydroxyapatite and can reduce or eliminate the expense, ligand leakage and reproducibility problems associated with immobilized dextran sulfate, chondroitin sulfate or heparin.

Matrex Cellufine Sulfate consists of a rigid spherical cellulose matrix of 3,000 Dalton exclusion limit, with a low concentration of sulfate ester functionality on the 6-position of cellobiose (see Figure 1). This low density of cation exchange groups provides affinity binding properties for a wide range of viruses and antigens. In normal physiological buffers, virus particles or viral/microbial antigens bind to the surface of the beads due to an affinity interaction. Pyrogens and most contaminating proteins, as well as nucleic acids and endotoxins, pass through the column unbound.


Elution of the bound product is effected through simple stepwise or gradient increases in ionic strength. Additionally, Matrex Cellufine Sulfate is similar enough to heparin and dextran sulfate to simulate these materials in many protein affinity purification applications.

A wide variety of lipid-enveloped viruses have been found to bind to the gel, including Influenza, Hepatitis B, Encephalitis (Japanese), Herpes, Rabies.

### FEATURES


- Affinity for a wide range of live, killed or disrupted viruses, viral or microbial antigens and heparin-binding proteins.
- More effective than ultracentrifugation at removing contaminants from culture media and host cells.
- Closed column operation assures safety and product sterility.
- Avoids excessive product handling and safety concerns, particularly with viral preparations.
- Endotoxins do not bind, allowing a rapid and contaminant free depyrogenation.
- Simultaneous concentration and purification improve yield, reduce processing steps, time and costs.
- Gentle binding and elution conditions provide high capacity and product yield.
- Rigid, high-strength beads resist compression, providing rapid flow for high-speed processing, even in large columns, making it easily scaleable.
- Autoclavable, resistant to chemical depyrogenation with base and chemically sterilizable with formalin.

MILLIPORE



**Figure 1**  
**Partial Structure of Cellufine Sulfate**

## FLOW PROPERTIES



**Figure 2**  
**Pressure/Flow Curves**

The nearly rigid properties of the spherical cellulose support matrix allow outstanding flow properties, particularly in large production columns.

Column A: 90 x 200mm  
Column B: 350 x 200mm

## CHARACTERISTICS

Support Matrix:	Cellulose
Particle Size:	44 - 105µm
Particle Shape:	Spherical
Gel Exclusion Limit:	3kD
Activated Group:	Sulfate Ester
Total Sulfur:	>700µg/g dry
Protein Binding Capacity	
Lysozyme:	>3mg/ml
Hepatitis B Surface Antigen:	7mg/ml
Environmental Resistance:	Resistant to 0.1M NaOH, 0.1% of 37% Formalin
Operating Pressure:	<2 bar (30 psi)
Autoclavable:	In suspension at neutral pH: 30 min at 121°C
Supplied:	Suspension in 20% Ethanol

## APPLICATIONS

### Virus, Viral/Microbial Antigens


There are many applications of Matrex Cellufine Sulfate in the concentration or purification of viral and microbial antigens, proteins and viruses - see Table 1.

<b>Viruses</b>
Rabies*
Influenza*
Japanese Encephalitis*
Feline Leukemia
Feline Herpes
Feline Calicivirus
Respiratory Syncytial Virus
Human Herpes Simplex
Human Measles
Human Parainfluenza
<b>Viral/Microbial Antigens</b>
Herpes Simplex gA and gB Glycoprotein Subunits*
Hepatitis B Surface Antigen
Filamentous Hemagglutinin from <i>B. pertussis</i> *
Leucocytosis Promoting Factor Hemagglutinin*
* These applications are covered by US and foreign process patents. Please inquire regarding details and licensing arrangements.

**Table 1**

## Purification of Rabies Virus

The following example illustrates the high degree of concentration, purification and yields obtained with Matrex Cellufine Sulfate on typical viral preparations.



**Figure 3**  
**Purification of Rabies virus from chick embryo tissue culture fluid**

Column: 50 x 70mm (140ml)  
Buffer: 0.01M Phosphate (pH 7.2)  
Eluant: 1M NaCl/0.01M Phosphate (pH 7.2)

	Load	Eluate
Volume (ml)	4,200	50
Virus Titer	32	4,096
Protein (ug/ml)	8.5	14
Yield (%)	100	152
Purification Factor	-	79x
Concentration Factor	-	126x

Concentration and Purification of Virus with Cellufine Sulfate

**Table 2**

## Purification of Influenza Virus

Hen's egg allantoic fluid was loaded directly onto a 333ml gel bed and 94.5% virus was recovered in the eluate fraction.


	Volume ml	Virus titer	TCA-N ug/ml	Recovery %	Fold Purification
Allantoic Fluid	4200	77	337.1	100	1
Wash	6700	1	209.2	2.1	-
Eluate	170	1797	448.0	94.5	20.1

**Table 3**

Column: 50 x 170mm  
Buffer: 0.01M Phosphate pH 7.4  
Wash: 0.01M Phosphate pH 7.2 + 0.2M NaCl  
Elution: 0.01M Phosphate pH 7.0 + 1.5M NaCl

## Antigenic Protein Purification and Depyrogenation

Matrex Cellufine Sulfate is ideal for depyrogenating virus and other microbial extracts because it does not bind endotoxins. Figure 4 shows the purification of filamentous hemagglutinin (FHA) from the whooping cough bacterium *Bordetella pertussis*.



**Figure 4**

### Purification of filamentous hemagglutinin from *B. pertussis*

Column: 16 x 70mm (20ml)  
Sample: 800ml *B. pertussis* culture fluid (endotoxin liter > 1015 by *Limulus* lysate test)  
Buffer: 0.01M Phosphate (pH 7.6)  
Eluant: 1M NaCl/0.01M Phosphate (pH 7.6)  
FHA Yield: 94%  
Purification Factor: 20 x  
Concentration Factor: 28x (30ml product)  
Endotoxin: Below standard level by *Limulus* lysate, rabbit pyrogen and mouse toxicity tests

## Protein Purification

Matrex Cellufine Sulfate mimics the affinity of heparin or dextran sulfate for many proteins. It can function as an affinity support for selected plasma proteins, cellular growth factors and lipases. Its capacity is comparable to conventional heparin gels.


Binding Proteins	Non-Binding Proteins
Antithrombin III β-Lipoprotein Complement C5, C6, C8 Complement C3 activator Trypsin Trypsin Inhibitor Chymotrypsinogen lysozyme Urease Catalase Factor IX	Albumin α-Lipoprotein Complement C3, C9 Complement C1, C3b inactivators IgG Ceruloplasmin α2-Macroglobulin RNase Bacitracin Glucose Oxidase
Binding and elution are extremely rapid and very fine separations can be generated in gradient mode	

**Table 4**

## ORDERING INFORMATION

Type	Quantity	Catalogue No.
Matrex Cellufine Sulfate	50ml	19845
Matrex Cellufine Sulfate	500ml	19846
Matrex Cellufine Sulfate	5 liters	19847
Matrex Cellufine Sulfate	10 liters	19849
Matrex Cellufine Sulfate	20 liters	19848

## Purification of Partially Purified Casein Kinase II from Calf Thymus



**Figure 5**

Column: 10 x 20mm  
 Sample: 7ml  
 Buffer: 50mM Tris-HCl (pH 7.9)  
       + 50mM MgCl<sub>2</sub> + 0.1mM EDTA  
       + 0.1mM PMS + 0.5mM DTT  
       + 25% glycerol  
 Eluant: 0.05 - 1.0M NaCl in buffer

## To Place an Order or Receive Technical Assistance

For additional information call the nearest number listed below:

In the U.S. and Canada, call toll-free 1-800-MILLIPORE (1-800-645-5476)

In the U.S., Canada and Puerto Rico, fax orders toll-free 1-800-MILLIFX (1-800-645-5439)

On the Internet <http://www.millipore.com>

E-mail: [tech\\_service@millipore.com](mailto:tech_service@millipore.com)

## Millipore Worldwide

### AUSTRALIA

Tel. 1 800 222 111  
or (02) 9888 8999  
Fax (02) 9878 0788

### AUSTRIA, CENTRAL EUROPE, C.I.S., AFRICA, MIDDLE EAST AND GULF

Tel. (43) 1 877-8926  
Fax (43) 1 877-1654

### BALTIC COUNTRIES

Tel. 358-9-804 5110  
Fax 358-9-256 5660

### BELGIUM AND LUXEMBOURG

Tel. (02) 726 88 40  
Fax (02) 726 98 84

### BRAZIL

Tel. (011) 548-7011  
Fax (011) 548-7923

### CANADA

Tel. 1-800-645-5476  
Fax 1-800-645-5439

### CHINA, PEOPLE S REPUBLIC OF

Beijing:  
Tel. (8610) 6500-8063  
Fax (8610) 6500-7372  
Guangzhou:  
Tel. (8620) 8755-4021  
Fax (8620) 8755-4350  
Hong Kong:  
Tel. (852) 2803-9111  
Fax (852) 2513-0313  
Shanghai:  
Tel. (8621) 5306-9100  
Fax (8621) 5306-0838

### CZECH REPUBLIC

Tel. 02-205 138 41  
02-205 138 42  
Fax 02-205 14294

### DENMARK

Tel. 70 10 00 23  
Fax 70 10 13 14

### FINLAND

Tel. (09) 804 5110  
Fax (09) 256 5660

### FRANCE

Tel. (01) 30.12.70.00  
Fax (01) 30.12.71.80

### GERMANY

Tel. (06196) 494-0  
Fax (06196) 43901

### HUNGARY

Tel. 1-205 9784  
Fax 1-205 9792

### INDIA

Tel. (91) 80-839 46 57  
Fax (91) 80-839 63 45

### ITALY

Vimodrone (Milano):  
Tel. (02) 25.07.81  
Fax (02) 26.50.324  
Roma:  
Tel. (06) 52.03.600  
Fax (06) 52.95.735  
JAPAN  
Tel. (03) 5442-9711  
Fax (03) 5442-  
9736 Analytical  
9737 BioProcess  
9734 Lab Water

### KOREA

Tel. (822) 551-0990  
Fax (822) 551-0228

### MALAYSIA

Tel. 603-757 1322  
Fax 603-757 1711

### MEXICO

Tel. (525) 576-9688  
Fax (525) 576-8706

### THE NETHERLANDS

Tel. 076-5022000  
Fax 076-5022436

### NORWAY

Tel. 22 67 82 53  
Fax 22 66 04 60

### POLAND

Tel. 22-669 12 25  
22-663 70 31  
Fax 22-663 70 33

### PUERTO RICO

Tel. (787) 273-8495  
Fax (787) 747-6553

### SINGAPORE

Tel. (65) 842 1822  
Fax (65) 842 4988

### SPAIN AND PORTUGAL

Madrid:  
Tel. 917 283 960  
Fax 917 292 909  
Barcelona:  
Tel. 934 525 530  
Fax 934 516 048

### SWEDEN

Tel. 08-628 69 60  
Fax 08-628 64 57

### SWITZERLAND

Volketswil:  
Tel. (01) 908-30-60  
Fax (01) 908-30-80  
Lausanne:  
Tel. (021) 648 51 18  
Fax (021) 648 51 37

### TAIWAN

Taipei:  
Tel. (886-2) 700-1742  
Fax (886-2) 755-3267  
Hsin Chu City:  
Tel. (886-3) 571-0178  
Fax (886-3) 572-9520

### U.K. AND IRELAND

Tel. (01923) 816375  
Fax (01923) 818297

### U.S.A.

Tel. (781) 533-6000  
Fax (781) 533-3110

### OTHER COUNTRIES

Millipore Intertech (U.S.A.)  
Tel. +1 (781) 533-8622  
Fax +1 (781) 533-8630

# MILLIPORE

Millipore and Matrex are registered trademarks of Millipore Corporation or an affiliated company. Cellufine is a trademark of Chisso, Corp. Lit. No. PF1007EN00 Rev A 7/99. © 1999 Millipore Corporation or an affiliated company. All rights reserved.